The Lighted Ear Curette™ with Magnification

INSTRUCTIONS FOR USE

MD	CE
Ⓢ	EC REP
R_x ONLY	Obelis s.a Boulevard Général Wahnis 33 1030 Brussels, BELGIUM Tel: +(32) 2.732.59.54 Fax: +(32) 2.732.60.03 E-Mail: mail@obelis.net
MADE IN THE U.S.A.	Bionix LLC 1670 Indian Wood Circle Maumee, Ohio 43537, United States Phone: +1 419.727.8421 Fax: +1 419.727.8426 www.Bionix.com

Importer, Importeur, Importör, Importador, Importateur, Importeur, Importer, Dovezce, Maahantuaja, Importatore, Εισαγωγέας, Importuotojas:		MedEnvoy Prinses Margrietplantsoen 33 - Suite 123 2595 AM The Hague The Netherlands
---	---	---

	©2023 Bionix LLC 1670 Indian Wood Circle Maumee, Ohio 43537 United States Phone: 800.551.7096 Fax: 800.455.5678 Email: bionix@bionix.com *Patents Pending www.bionix.com RM95-0006 2023/03/23 Rev. 04
--	---

ENGLISH

Intended use: The intended use of the Lighted Ear Curette is to remove cerumen and/or foreign bodies from and illuminate a patient's ear canal without causing injury.

Indications for use: The Lighted Ear Curette is indicated for use in cases of impacted cerumen and/or foreign bodies occluding the ear canal.

Assembly/Disassembly:

- To assemble the lighted ear curette, insert the end of the curette into the light source by aligning the pegs on the curette end with the grooves inside the light source. Push down and twist **FIRMLY** clockwise until rotation stops. This will lock the curette into place and activate the light.
- To attach the magnification lens, snap the lens onto the neck of the light source so that the top of the lens angles towards the ear curette tip.
- To disassemble, push down on the curette and rotate counter clockwise to remove the curette from the light source. The light will go off when the curette is removed.

Instructions For Cerumen Removal

- Use an otoscope to examine the ear canal and determine cerumen location.
- Select the appropriate style lighted ear curette for the type of patient and cerumen occlusion. Assemble according to instructions.
- Use lighted ear curette to remove cerumen from the ear canal.
- Repeat otoscope exam. Determine location of any remaining cerumen.
- Repeat curettage procedure until sufficient cerumen is removed to allow adequate visualization of the tympanic membrane.
- Disassemble the lighted ear curette, and discard the used curette in the usual fashion. Save the light source for the next curettage procedure. LEC Tips are intended for single patient use only.

CAUTION: Discontinue curettage immediately if bleeding, irritation, or other trauma to the ear canal or tympanic membrane occurs.

CAUTION: This product is a spring-loaded device; improper assembly may cause ejection of the ear curette.

CAUTION: DO NOT hold device by light source as accidental ejection of curette may occur if not assembled properly.

WARNING: Cross-contamination risk. Do not reuse disposable curette tips as this may spread contamination from one patient to another patient.

Medical Device Reporting: Notice to Users and/or Patients in EU: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Clinical Benefits: Indicated in Intended Use.

Contraindications: There are no absolute contraindications for cerumen removal.

Residual Risk: Risk associated with the use of this product has been reduced as far as possible, but the product cannot completely eliminate potential patient or user harm arising from the following:

- Harm from mechanical hazards
- Harm from misuse, or use error
- Harm from unanticipated origins

MD	Medical Device		Do Not Use if Package is Damaged
	Manufacturer		Do Not Reuse Single-Use Device
	Date of manufacture		Single patient-multiple use
	Use-By Date		Consult Instructions For Use (IFU)
LOT	Lot Code		Caution
REF	Reorder Number		Warning
SN	Serial number	EC REP	Authorized representative in the European Community
STERILE R	Sterilized using irradiation	CE	European Conformity
	Do Not Resterilize	R_x ONLY	Prescription Only or "For Use by or on the order of a licensed medical professional"
	Non-sterile		Temperature limitation

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