

QUALITY PRODUCTS | VALUE PRICED | TIMELY DELIVERY

# 2025/2026 PRODUCTS CATALOGUE



[www.canadamedical.net](http://www.canadamedical.net)

# About Canada Medical Ltd.

Founded over 35 years ago Canada Medical Ltd. is proudly Canadian-owned and operated, with roots firmly planted in the Canadian healthcare community since 1986. For nearly four decades, we have been serving clinics and healthcare professionals across the country. Being Canadian-owned means we understand the unique needs of our local medical practices and are deeply committed to supporting them. Our long history and homegrown values translate into a personal, reliable approach that clinics nationwide have come to trust.

**Our mission:** 100% dedicated to Canada's Medical supply needs.

**We are proud to specialize in:**

- Competitive pricing with exceptional service
- Specializing in clinic supplies and equipment
- Canada-wide inventory and fast delivery
- Sterilization consultations and high standards

Whatever your needs, our experienced team is always ready to share insights and help healthcare providers find the right solutions.

Connect with us:



**Canada Medical Ltd.**

Head Office:

16503 116 Ave NW  
Edmonton, AB T5M 3V1

Local: 780.465.2020

Local Fax: 780.468.0991

Toll-Free: 1.800.232.7286

Fax Toll-Free: 1.800.278.8503

[orders@canadamedical.net](mailto:orders@canadamedical.net)

[www.canadamedical.net](http://www.canadamedical.net)



\*Visit [canadamedical.net](http://canadamedical.net) for more details

# PRODUCT CATEGORIES

**EXAM ROOM EQUIPMENT & FURNITURE**

**PAGE 4 - 18**

**DIAGNOSTICS**

**PAGE 19 - 38**

**INFECTION PREVENTION**

**PAGE 39 - 43**

**HOUSEKEEPING / JANITORIAL**

**PAGE 44 - 46**

**STERILIZATION**

**PAGE 47 - 51**

**INSTRUMENTS REUSABLE**

**PAGE 52 - 57**

**INSTRUMENTS DISPOSABLE**

**PAGE 58 - 62**

**SURGERY AND MINOR PROCEDURES**

**PAGE 63 - 77**

**INFUSION & BLOOD COLLECTION**

**PAGE 78 - 84**

**LABORATORY**

**PAGE 85 - 90**

**NEEDLES & SYRINGES**

**PAGE 91 - 97**

**PATIENT CARE & EXAM ROOM SUPPLIES**

**PAGE 98 - 107**

**PROTECTIVE APPAREL**

**PAGE 108 - 111**

**RESPIRATORY & OXYGEN**

**PAGE 112 - 115**

**SKIN & WOUND CARE**

**PAGE 116 - 125**

**DENTAL**

**PAGE 126 - 129**

**TATTOO**

**PAGE 130**

**SPA AND BEAUTY**

**PAGE 131 - 133**



## RITTER 224/225 BARRIER-FREE® EXAMINATION CHAIR

Ritter 224 & 225 Barrier-Free chairs are designed to offer low and high positioning heights so you can provide the best care to all patients, whether treating the elderly, expectant mothers or patients with disabilities. Ritter Barrier-Free Examination Chairs by Midmark adjust in height from 18-37", improving accessibility to the exam chair and reducing the likelihood of staff injury.



- Power Base
- 224 Model Air Spring Assist Back
- 225 Model Power Back
- 500lb Patient Weight Capacity
- Foot Control

RIT224-001 .....	Base with Manual Back
RIT224-002 .....	Base with Manual Back, Pelvic Tilt and Drawer Heater
RIT225-002 .....	Base with Power Back
RIT225-003 .....	Base with Power Back, Clean Assist Roller System
RIT002-2009-(Colour Code) .....	Seamless Top

## RITTER® STANDARD COLOURS

Mist-813	Stone-814	Restful Path-858	Iron Ore-851	Lunar Gray-845	Obsidian-857	Oasis-843	Soothing Blue-855	Cranberry-859	Dark Linen-856
Latte-860	Dream-847	Mineral-844	Healing Waters-848	Shaded Garden-853	Robust Brown-852	Deep Earth-854	Citrus-849	Sandy Retreat-850	Curative Copper-846



## RITTER® 204 MANUAL EXAMINATION TABLE

What makes Ritter the most purchased brand of exam tables in the industry? It's our willingness to research, study and listen so we can meet your needs—and the needs of your patients. If you need an exam table with easy access for patients to make them feel secure, three side drawers for optimal supply management, modern styling to coordinate with your office and an additional weight capacity up to 500lb, the Ritter 204 Manual Exam Table is the examination table for you.



- One piece seamless drawers
- Stylish and easy to clean
- Removable top system
- Enlarged footstep
- Storage access and volume
- Seamless upholstered top
- Manual back with air spring assist



RIT204-011 .....	Base with Manual Back
RIT002-10138-(Colour Code) .....	Seamless Top



## 626 BARRIER-FREE® EXAMINATION CHAIR

The Midmark 626 is the first chair of its kind designed to promote proper patient positioning for a more accurate blood pressure measurement. Connect the IQvitals® device to the Midmark 626 Barrier-Free Exam Chair with the integrated IQscale® to weigh patients quickly and discreetly from the comfort of the exam chair. Then, send data to the EMR. The Midmark IQconnect® framework provides a seamless connection from your devices to your EMR system. In other words, it allows you to operate all of your devices and transfer data to the EMR from one, single pane of glass—your computer screen.

### Standard Features:

- Clean Assist™ Roller System
- The retractable roller base allows you to easily and safely move the chair for cleaning
- Hand and Foot Control
- Stirrups and Pelvic Tilt
- Drawer Heater
- Active Sensing Technology pauses table movement when disrupting impact is detected
- Quick Exam® and Return to Home Quick Exam elevates the chair to the exam position with one touch of the control, while Return to Home returns the chair to its home position
- Patient Weight Capacity: 650lb
- 15.5" L - 37" H

RIT626-001..... Base Only  
RIT002-2000-(Colour Code ..... Seamless Top



## RITTER 244 BARIATRIC BARRIER-FREE® EXAM CHAIR

Barrier-Free® Accessibility - Safe and Efficient Care for your Patient.

- 850lb Weight Capacity
- 18" Low Height. By lowering the seat to 18" from the floor, the Ritter 244 performs this sometimes difficult task. It becomes a very simple routine, requiring minimal assistance from you and your staff. The patient's dignity is maintained and the likelihood of injury is greatly reduced, providing efficient, effective patient care
- 34" High Height - The Ritter 244 raises the seat to a height of 34", designed to bring the patient to you while helping to minimize the stress and strain on your back. By having the patient in the right position, the exam can be completed properly and efficiently

RIT244-001 ..... Base Only  
RIT002-0861-(Colour Code ..... Seamless Top



## MIDMARK® 631 PROCEDURE CHAIR

The right procedure room equipment can make a crucial difference in both patient outcomes and caregiver safety. Particularly for patients with mobility limitations, accessible in-office procedure environments are essential. Unfortunately, many procedure chairs on the market today can create barriers to patient access. That's why we designed the Midmark 631 Procedure Chair.

### Accessibility

- The Midmark 631 has a US Access Board compliant 17-inch low seat height (while maintaining a 39-inch high seat height.)
- Transfer patients onto the 28 or 32-inch wide transfer surface with assistance from compliant Patient Support Rails.
- The chair is compatible with patient lifts and helps ensure safe patient transfer with no interference.
- Utilize US Access Board compliant Knee Crutches to support the patient's thigh, knee, and calf during lower body exams and procedures.

RIT631-001..... Base Only  
RIT002-10609-(Colour Code)..... Seamless Top





## HI-LOW POWER TWO SECTION EXAM TABLE

### Two Section Table

32" .....	Maximum Height
17.5" .....	Minimum Height
27" .....	Width
78" .....	Length
39.5" .....	Head Section
24.5" .....	Seat Section



SOLAYR1008

High-Low Exam Table design provides an 17.5" profile in low position to allow level transfer from standard wheelchair heights, complete with upright seat back.

- Adjustable pneumatic backrest from either side, raising up to 80° from flat
- Chair-to-Table design facilitates seamless transition from consultation to examination level
- Hands-free operation with a low voltage foot control
- Sturdy, all welded powder coated steel frame for heavy-duty use
- Welded steel support frames in the top sections
- 2" firm density foam tops for optimal comfort
- Accommodates paper rolls of 18" to 21" with a convenient paper roll holder
- Impressive weight capacity of 500lbs

## HI-LOW POWER THREE SECTION EXAM TABLE

### Three Section Table

32" .....	Maximum Height
17.5" .....	Minimum Height
27" .....	Width
78" .....	Length
39.5" .....	Head Section
24.5" .....	Seat Section
13" .....	Foot Section



SOLAYR1009

SOLAYR1009E.....with stirrups

## EXAM ROOM BASIC TREATMENT TABLE

- Table Dimensions (24.5" W x 75" L x 28.5" H)
- 350lb weight capacity
- All steel construction with strong 4-leg base
- Nylon floor levelling glide helps prevent scratches on floor

SOL2001



## EXAM ROOM BASIC TREATMENT TABLE WITH DRAWERS

- Table Dimensions 24.5" W x 75" L x 30" H
- 350lb weight capacity
- Mechanical 3-position headrest with adjustable degrees of 25°, 40° & 60°
- 2" thick soft bedded black vinyl top
- Paper roll holder
- Bottom storage shelf
- Two front drawers 20" W x 18.5" D x 3.5" H

SOL22401



## RITTER SIDE CHAIRS

The Midmark 280 is soft, comfortable seating for your patients in the exam room, waiting room and front office. The wall-saver kickback leg design allows the chair to be close to the wall without damaging it.

### COMFORT

Contoured, upholstered seat and back section

### STABILITY

Sturdy, powder-coated black steel tube frame with leveling feet

### OPTIONS

Available with or without high-impact plastic chair arms

RIT280-001-(Colour Code)\* .....Basic chair without arms

RIT280-002-(Colour Code)\* .....Basic chair with arms

\*See Page 4 for colour codes.



## RITTER 282 BARIATRIC SIDE CHAIR

The Midmark 282 Side Chair can enhance your clinical workflow while complementing the style of your space with a range of upholstery color options.

### COMFORT

Contoured, upholstered seat and back section

### STABILITY

Sturdy, powder-coated black steel tube frame with levelling feet, a 30" wide seat and a 650lb weight capacity

RIT282-001-(Colour Code)\* .....Basic chair with arms

\*See Page 4 for colour codes.





## RITTER® AIR LIFT HEIGHT STOOL, 5 CASTER

These stools offer a black, plastic composite five leg caster base and hood and a single lever height control. The seat height adjusts pneumatically. Black vinyl upholstery is stocked for quick delivery. Other colours available with lead time.

RIT273-001-(Colour Code) .....Stool with back

RIT272-001-(Colour Code) .....Stool without back

\*See Page 4 for colour codes.



## STEP ON STOOLS

AMG116-460 .....Foot Stool

DUK4349-2 .....Foot Stool with Handrail

AMG116-465 .....Replacement rubber tips

DUK4350 .....Bariatric Foot Stool

DUK4348 .....Bariatric Foot Stool with Handrail

**Health o meter®**  
Professional



## DIGITAL PEDIATRIC TRAY SCALE

Scale function that allows user to select and use only pounds or kilograms, reducing the risk of incorrect data entry.

- Capacity: 44lb / 20kg
- Resolution: 0.5oz / 10g
- Display: 0.6" LCD
- Power Source: 9V battery included
- Warranty: 2 Year limited warranty
- Carrying case: optional carrying case (order # 64769)
- Tray Size: 19 1/4" W x 13" D x 2 5/8" H

CON549KL .....Pediatric scale

**Health o meter®**  
Professional



## DIGITAL PEDIATRIC SEAT SCALE

Plastic, easy-to-clean seat limits infant movement while accurately calculating weight measurement.

- Capacity: 50lb / 23kg
- Resolution: 0-20lb / 0.2oz; 20-50lb / 0.5oz; 0-9kg / 5g; 9-23kg / 10g
- Power Source: 6 AA batteries included, optional 120V adapter (order # ADPT31)
- Connectivity: USB
- Seat Size: 24 1/4" W x 14 1/4" D x 18 1/2" H
- Warranty: 2 Year limited warranty

CON524KL .....Digital scale with seat



## REMOTE DISPLAY DIGITAL FLOOR SCALE

Interfaces with Welch Allyn® Spot Vital Signs® LXi Monitor. The remote display provides versatility with the ability to be mounted on a wall or placed on a table.

- Capacity: 400lb / 180kg
- Resolution: 0.2lb / 0.1kg
- Platform Size: 12 1/2" W x 12" D x 1 7/8" H
- Display: 1/2" LCD
- Power Source: 9V battery included, optional 120V power adapter (Order#: ADPT40)
- Connectivity: serial port
- Warranty: 2 Year limited warranty

CON349KLX .....Digital floor scale



## LARGE PLATFORM FLOOR SCALE

- Capacity: 400lb / 180kg
- Platform Size: 11" W x 12 1/2" D x 3" H
- Graduation: 1lb / 1kg
- Warranty: 1 Year limited warranty

CON160KL .....Platform floor scale



## MECHANICAL BEAM SCALE

The easy-to-use height rod provides height measurements from 23" to 84".

- Capacity: 500lb / 200kg - with included counterweights
- Platform Size: 10 1/2" W x 14" D x 3 1/4" H / 267mm W x 356mm D x 83mm H
- Height Rod: 23 5/8" - 84" / 60cm - 213cm
- Graduation: 1/4lb / 100g
- Warranty: 10 Year limited warranty

CON450KL .....Beam scale



## EYE LEVEL DIGITAL SCALE

Scale has the ability to interface with a computer, EMR software or other electronic device via various Health-o-meter® Professional Connectivity Solutions. The easy-to-use height rod provides height measurements from 30" to 84"/76cm - 214cm.

- Capacity: 500lb / 220kg
- Resolution: 0.2lb / 0.1kg
- Platform Size: 13 3/4" W x 16 1/2" D x 2 3/8" H
- Wheels: 2 Wheels for ease of movement
- Power Source: 6 AA batteries included, optional 120V adapter (order # ADPT31)
- Warranty: 2 Year limited warranty - eligible for ScaleSurance extended warranty

CON500KL .....Digital scale

CON500HB.....Handlebars for 500 series scale

CON500KLHB .....Digital Scale with Handlebars

**Health o meter®**  
Professional





## DIGITAL PHYSICIAN SCALE WITH INTEGRATED DIGITAL HEIGHT ROD

The newest member of the number-one-selling family of physician scales, the all-digital Health o meter® 502KL features a fully integrated digital height rod and EMR capabilities, providing the optimum in quality, reliability, and ease of use. The 502KL, created to cater to various healthcare settings like doctor's offices, hospitals, clinics, and acute and long-term care facilities, showcases a market-leading built-in digital height rod, a spacious non-skid platform, and a reinforced base that ensures enhanced stability. The scale offers providers and patients an eye-level weight indicator, accommodating patients weighing up to 660 pounds.

- Capacity: 660 lb / 300 kg
- Resolution: 0.2 lb / 0.1 kg
- Display: 1" / 25 mm, LCD display
- EMR Connectivity: Via optional Pelstar® wireless technology
- Power Source: 6 AA batteries included, optional adapter (order# ADPT31)
- Warranty: 2-year limited manufacturer warranty

CON502KL .....Scale with height rod

**Health o meter®**  
Professional



## DIGITAL PLATFORM SCALE WITH HEIGHT ROD

EMR Connectivity Scale has the ability to interface with a computer, EMR software or other electronic device via various Health-o-meter® Professional Connectivity Solutions. Calculates Body Mass Index (BMI).

- Capacity: 1000lb / 454kg
- Resolution: 0.2lb / 0.1kg
- Platform Size: 22" W x 15 ¾" D x 2 1/2" H
- Height Rod: Included
- Display: 3/4" LCD, 180° Swivel
- Wheels: 2 Wheels for ease of movement
- Power Source: 110-240V adapter included, 6 AA batteries (not included)
- Connectivity: USB
- Warranty: 2 Year limited warranty

CON1100KLHR .....Scale with height rod

**Health o meter®**  
Professional



## HEIGHT MEASURING TAPE, WALL MOUNT

- The space-saving solution with roll-up mechanism
- Easy fixing on wall with just one screw
- Display window in the head piece, it's easy to read the height
- The durable tape is made of metal

SOLHX400C .....Height Measuring Tape



## RETRACTABLE TAPE MEASURE

- The space-saving solution with roll-up mechanism
- Easy fixing on wall with just one screw
- Display window in the head piece so it's easy to read the height
- The durable tape is made of metal

AMG116-860 .....Measuring Tape

## WALL MOUNTED MECHANICAL MEASURING ROD

The Seca 216 has a variable scale with a range of 1.4" to 90" that is inserted into a 55.3" frame.

- 2 different measurement ranges
- Easy read out
- Adjustable head positioner
- Choice of two units of measurement
- Simple assembly

SEC2161814009 .....Measuring rod

**seca**  
Precision for health

## TELESCOPIC HEIGHT MEASURING ROD

Mechanical telescopic measuring rod with large measuring range 60-200cm.

- Designed specifically for wall mounting
- The telescopic measuring rod ensures that the result can be easily read at eye level
- When not in use, the measuring slice can be folded down

SOLHX200A.....Measuring Rod

## INFANT HEIGHT MEASURING DEVICE

The Seca 210 is a lightweight and space-saving solution for precise length measurement of babies and toddlers while lying down. The Seca 210 is washable and foldable. The fixed head piece and the sliding foot positioner make it simple to use.

SEC2101821004 .....Measuring device



## Ritter® 253 LED EXAM LIGHT

Designed to be easy to use, easy to clean, durable and aesthetically pleasing, the Ritter 253 exam light is the simple lighting solution. Because the controls are located on the light head within easy reach, you'll never have to reach around the patient to adjust the lighting, making this the perfect light for lengthy procedures. With table mount and mobile options available, the 253 is a flexible, spacesaving addition perfect for any exam room.

- Life: greater than 50,000 hours
- Flexible arm
- Power cord: 8ft (2.5m) hospital grade plug

RIT253-017 .....Light with wall mount

RIT253-012 .....Light with table mount

RIT253-016 .....Light with mobile stand

**midmark**



## Ritter® 255 LED PROCEDURE LIGHT

Whether you practice in an office environment or specialized procedure area, your work requires superior, white, shadowfree illumination that is safe and cool to the touch. The Ritter 255 LED Procedure Light by Midmark is designed to deliver both quality and affordability with features that can help you provide the best care to your patients.

- Illumination at 39" (1m)
- LED Light Module: 16 LEDs arranged in a 4x4 array, 32W rating
- Fully autoclavable handle
- Available in single or dual ceiling mount versions

RIT255-001 .....Single head (8' ceiling mount)

RIT255-004 .....Dual head (9' ceiling mount)

**midmark**





**WelchAllyn®**



## THE GREEN SERIES EXAM LIGHT IV

Perfect for the OB/GYN office or for specialties that require very focused light in specific areas.

- One LED with 3X the typical light output than halogen
- Highly flexible light pipe that is twice as long as conventional pipes with a compact, unobtrusive head that ensures pinpoint accuracy
- Touchless On/Off
- No bulb replacement

WEL48810 ..... GS Exam Light IV with wall mount

WEL48950 ..... Mobile stand for Green Series Light



**WelchAllyn®**



## GREENSERIES 300 GENERAL EXAM LIGHT

Perfect for the OB/GYN office or for specialties that require very focused light in specific areas.

- LED with 3X the typical light output than halogen
- Highly flexible light pipe that is twice as long as conventional pipes with a compact, unobtrusive head that ensures pinpoint accuracy
- Touchless on/off
- No bulb replacement
- 5 Year limited warranty

WEL44410 ..... GS Exam Light 300 with wall mount

WEL44400 ..... GS Exam Light 300 with mobile stand



**WelchAllyn®**



## 600 MINOR PROCEDURE LIGHT

Ideal for minor procedures and exams in any healthcare environment.

- Touchless intensity controls provide the output you need, when you need it
- Features 3 LEDs and a wider head specifically for demanding viewing conditions that require good light
- 5 Year limited warranty

WEL44610 ..... GS600 Light with table/wall mount

WEL44600 ..... GS600 Light with mobile stand



**WelchAllyn®**



## GS900 PROCEDURE LIGHT

Perfect for minor procedures/procedures in the hospital, ambulatory care center, outpatient facility or physician's office. Features 6 LEDs with intensity controls, outstanding durability, extraordinary maneuverability and unparalleled spot quality. Ceiling mount options gives ability to maximize floor space and direct light wherever it is needed during a procedure. 5 Year limited warranty.

WEL44900 ..... with mobile stand

WEL44900-C ..... with ceiling mount

WEL44900-W ..... with wall mount



## AARON® A942 HYFRECATOR HIGH FREQUENCY DESICCATOR

- 40W of power with higher resolution in Bi-Polar mode
- Increase and decrease in 1/10th watt increments from 0 to 10W
- Universal Power with 100-240 VAC
- More User Friendly with features that can save you time and money
- Less to do with fewer steps to accomplish the same procedures
- Less to buy as there are fewer accessories needed for additional functionality
- 3-Button handpiece

AARA942 ..... Aaron® 942 Desiccator



## BANTAM|PRO A952 HF DESICCATOR WITH CUT

The Bovie Bantam | PRO, A952 is designed and engineered to be the most reliable and durable generator available today.

- Offers you up to 50W in Cut mode, Blended Cut mode, and Pinpoint Coagulation mode. 40W – Fulguration, Bipolar and Micropolar
- The unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat
- Unit cleanup is simplified by its non-absorbing, non-staining, high-impact ABS enclosure
- Standard package includes the unit, power cord, reusable handpiece, non-sterile drapes, 2 sterile drapes, 5 sterile disposable electrodes, 6 disposable sterile dermal tips, 12 disposable non-sterile dermal tips, 2 reusable grounding cord, 5 disposable split grounding pad, wall mount kit and users guide

AARA952 ..... Aaron® 952 Desiccator



## DISPOSABLE DERMAL TIPS

### Non-Sterile

AARA804-IMC ..... Sharp Dermal Tip..... Box/100

AARA806-IMC ..... Blunt Dermal Tip ..... Box/100

### Sterile (packaged individually)

AARA805 ..... Sharp Dermal Tip..... Box/50

AARA807 ..... Blunt Dermal Tip ..... Box/50



## ACCESSORIES

AARA812 ..... Mobile stand for A941, A942 & A952

AARA808T ..... Mayo tray for stand

AARA803 ..... Footswitch for A941, A942 & A952

AARA827V ..... Bipolar forcep cord for A942 & A952

AARA825 ..... Bipolar forcep reusable 4.75" straight

AARESRE ..... Disposable split grounding pads Box/50

AARA902 ..... Reusable 3-Button handpiece for A941, A942 & A952

AARA901 ..... Reusable 3-Button handpiece for A940 & A950







Bovie

## STERILE HIGH TEMPERATURE CAUTERIES

High temperature cauteries are popular for many types of procedures including pinpoint hemostasis, evacuation of a subungual hematoma, controlling diffuse bleeders, vasectomies and sculpting woven grafts. Individually packaged sterile.

AARAA01 .....	Fine tip 2200° F (1204° C)
AARAA03 .....	Loop tip 2200° F (1204° C)
AARAA21 .....	Elongated tip 1800° F (982° C) vasectomy



Bovie

## DISPOSABLE ELECTRODES

AARES01 .....	Standard blade .....	Box/50
AARES37 .....	Modified blade .....	Box/25
AARES20 .....	Standard 3/16" ball .....	Box/25
AARES02 .....	Standard needle .....	Box/25
AARES18 .....	Angled blade .....	Box/5
AARES55 .....	Extended blade 4" .....	Box/25
AARES03 .....	Extended needle 6" .....	Box/25
AARES07 .....	Extended 5mm ball .....	Box/5
AARES60 .....	SuperCut fine needle 3cm .....	Box/5
AARES62 .....	SuperCut fine needle 3cm .....	Box/5
AARES24 .....	Short shaft loop 3/8" .....	Box/5



Bovie

## SMOKE SHARK™ II EVACUATOR

A modern, light-weight, compact design allowing for simplified use, storage and stackability. Never again order the Smoke Shark™ and still need to order a filter upon receipt, the SE02 is packaged with both the SF35 Extended-Life Filter offers up to 35 hours of performance reducing procedure cost by 38% (+/-) from the previous model, when used at the lowest setting. Covered 3-port filter design accommodates a variety of tubing sizes [1/4" (6.4mm), 3/8" (9.5mm), and 7/8" (22mm)] providing versatility in procedure and accessory setup. The ultra-quiet design offers the same excellent performance with less noise. The SE02 Smoke Shark Smoke Evacuator is -3 and RoHS compliant. Package includes: Unit, Power Cord, User's Guide, Pneumatic Footswitch (FSSEP) and Long Life 35-hour Filter (SF35).

AARSE02 .....	Smoke Shark™
---------------	--------------



Bovie

## ACCESSORIES

AAR786T	.....7/8" x 6' tubing, non-sterile	.....Box/24
AAR786TS	.....7/8" x 6' tubing, sterile	.....Box/24
AARSEPAT	.....Pencil adapter 6' tubing, sterile	.....Case/10
AARSF35	.....Replacement 35hr filter for SE02	
AARSF18	.....Replacement 18hr filter for SE01	
AARSEAS	.....Arm stand to hold tubing	

## PREMIER PHARMACY/VACCINE UNDERCOUNTER REFRIGERATOR BUILT-IN

- 4.6 Cu. Ft. Capacity
- One solid swing door, right hinged
- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- 2 probes; 1 in air and 1 in glycol bottle
- Adjustable operating temperature range: 36°F - 46°F (2° - 8° C)
- Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable.
- Three shelves (two adjustable/one fixed) with guard rail on back to prevent products from being pushed against the back wall
- 23 <sup>3</sup>/<sub>4</sub>" W x 24 <sup>1</sup>/<sub>2</sub>" D x 33 <sup>3</sup>/<sub>8</sub>" H

ABSPH-ABT-HC-UCBI-0404



## PREMIER PHARMACY/VACCINE UNDERCOUNTER REFRIGERATOR FREESTANDING

- 2.5 Cu. Ft. Capacity
- One solid swing door, right hinged, not reversible
- Microprocessor temperature controller
- Digital temperature display
- Adjustable operating temperature range: 36°F - 46°F (2°C - 8°C)
- Refrigerator/Freezer freestanding thermometer with 3 Year certificate of calibration, with min/max memory, °F/°C switchable
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- 2 probes; 1 in air and 1 in sample bottle
- Four adjustable shelves (two large, two small) with guard rail on back to prevent product from being pushed against the back wall
- 17 <sup>3</sup>/<sub>4</sub>" W x 19 <sup>1</sup>/<sub>2</sub>" D x 28 <sup>7</sup>/<sub>8</sub>" H

ABSPH-ABT-HC-UCFS-0204



## SOLID DOOR PHARMACY PREMIER LABORATORY REFRIGERATOR FREESTANDING

- 10.5 Cu. Ft. Capacity
- Its ample room and six adjustable shelves provide ample space for a variety of contents, while a microprocessor temperature controller makes sure that it's always at a correct, safe temperature conveniently shown on a digital display
- Its average temperature cycle is 41°F, while its adjustable operating range is between 36 and 46°F
- High-density, urethane foam insulates the low-temp refrigerator and works in tandem with HFC-free refrigerant (R600a)
- Its door features a magnetic gasket for a positive seal, as well as a keyed door lock to keep contents secure
- Audible and visual alarms are always on hand to alert nearby users of high or low temperature fluctuations
- 59" H x 23 <sup>3</sup>/<sub>4</sub>" W x 26 <sup>1</sup>/<sub>2</sub>" D

ABSPH-ABT-HC-10PS





## SINGLE PROBE TEMPERATURE MONITORING DEVICE

- Accuracy of  $\pm 0.5^{\circ}\text{C}$  ( $\pm 1^{\circ}\text{F}$ )
- High and Low Temp Visual/Audio Alarm
- 3 years certification of calibration
- "buffered" probe in the product simulated solution
- Min/Max historical temp point
- $^{\circ}\text{F}/^{\circ}\text{C}$  Switchable Readout

ABSTMD-18



## Plastic Utility Cart

- Plastic composition
- Multi-swivel locking castors
- Easy to assemble
- With 3 Shelves
- 21.5"L x 14"W x 37.6"H

SOL101B



## STAINLESS STEEL UTILITY CART, 3 SHELVES

- Stainless steel cart
- Strong and elegant
- 4 noiseless smooth casters
- 3 shelves
- 26" W x 34.75" H x 17.3" D

SOLF20



## STAINLESS STEEL CART WITH 1 DRAWER

- 25.5"L x 25.5"W x 39.2"H
- 4" noiseless, smooth casters with breaks
- Includes 2 shelves and 1 drawer

SOLF22

## CRASH CART

- High-impact ABS material frame with aluminum alloy plastic plate, cold rolled steel frame, durable outer double-layer frame
- Handles on both sides for easy maneuvering
- Two small drawers, two medium drawers, and one large drawer, all equipped with removable dividers for easy cleaning
- Drawer with transparent medicine card
- Defibrillator shelf (1 pc)
- Dust bin (2 pcs)
- ABS Side basket
- Sharp container
- Adjustable I.V. pole, 4 hooks on top
- CPR Board at back
- Oxygen tank holder on side
- Sliding side shelf
- Two keys for centralized locking
- 5" Casters, two in front with brakes
- 26.3" × 20.8" × 40"



SOLNJ801

## CRASH CART

- High-impact ABS material frame with four aluminum alloy columns
- Stainless steel guard rails
- Handles on both sides for convenient handling
- Two small drawers, two medium drawers, and one large drawer, all equipped with removable dividers for easy cleaning
- Dust bin (1 pc)
- ABS Side basket (1 pc)
- Sharp container (1 pc)
- Adjustable I.V. pole, 4 hooks (1 pc)
- Oxygen tank holder (1 pc)
- Sliding side shelf (1 pc)
- Transparent file box (1 pc)
- Two keys for centralized locking (1 set)
- Casters (5 inches), two with brakes
- 33.4" × 20.4" × 38"



SOLJ101

## 5-LEG MOBILE INSTRUMENT STAND

- Adjustable height (28 - 48") for comfortable use while sitting or standing
- 1" (2.54cm) chrome plated steel tubing
- Stable five-leg 20" base
- Removable stainless steel tray: 12.5" x 19"
- 2" (5.08cm) dual wheel nylon



DUK4365

## MAYO STAND WITH TRAY

Mayo Stands with adjustable height, and stainless steel tray.

DUK4368 .....Double Post

DUK4352 .....Single Post





## DYNARIDE SERIES 3 LITE WHEELCHAIR

- Powder-coated steel frame with silver vein finish
- Black, padded, nylon upholstery
- Detachable, flip-back, padded armrests (desk length)
- Includes built-in chart pocket
- Tool-free, adjustable, elevating leg rest
- 24" solid rubber rear wheels feature hand rims and push-lock brakes
- Dual-axle provides adjustment from standard height to hemi level
- Weight capacity 350 lb
- Limited Lifetime Warranty

DYN10231



## DYNARIDE TRANSPORT WHEELCHAIR

The DynaRide Transport Wheelchair is a compact and dependable chair that will last through wear and tear. The carbon steel construction is inclusive of a foldable back with release hinge allowing it to fold down for easy storage and transport.

- 19" x 16" with Fixed Full Arm
- Powder-coated steel frame with silver vein finish
- Equipped with tool-free, adjustable, swing-away footrests with loops
- Includes built-in chart pocket
- Weight capacity 250 lb.
- Limited Lifetime Warranty

DYN10243



## BARI+MAX BARIATRIC WHEELCHAIR

The Bari+Max Wheelchair with Elevating Leg Rest provides reliable support and improved mobility for larger users

- 24" x 18"
- Weight capacity 500 lb
- Padded arm rests provide improved comfort
- Black, padded, double-embossed vinyl upholstery
- Dual-axle provides adjustment from standard height to hemi level
- Tool-free, adjustable, elevating leg rests with padded calf pads and extra-wide foot plates
- 24" composite real wheels with solid rubber tires
- Push-to-lock brakes
- Includes built-in chart pocket for added convenience
- Padded arm rests provide improved comfort
- Black, padded, double-embossed vinyl upholstery
- Dual-axle provides adjustment from standard height to hemi level
- Tool-free, adjustable, elevating leg rests with padded calf pads and extra-wide foot plates
- 24" composite real wheels with solid rubber tires
- Push-to-lock brakes
- Includes built-in chart pocket for added convenience

DYN10234





## STRETCHERS

- Adjustable H 76"x22"x16"-28.8"
- Made of high strength aluminum alloy thick wall pipes
- Adjustable back
- Foldaway legs are controlled by the handles on both ends of the stretcher
- The legs can be folded under the stretcher

SOLAG4H



## RECOVERY STRETCHER

- Hydraulic stretcher for ICU and Emergency
- Steel frame with Aluminum alloy retracting side rails, I/V pole included
- 2 gas cylinders for easy back lift control
- Trendelenburg: 18° & Reverse Trendelenburg: 18°
- Back rest max angle: 0°-90° & Footrest max angle: 0°-40°
- Weight capacity: 220Kg/485 Lbs
- 7.8" central lock silent casters
- Fifth wheel in middle for easy maneuvering
- Storage rack at the bottom with 22lbs weight capacity
- Washable blue zippered mattress
- Height Adjustment: 25" x 38"
- Length: 83.5"
- Width: 30.5"

SOLFF05



# Working through a list?

Let us help you check it off.

# Welch Allyn Connex® Vital Signs Monitor 6000 Series

Experience accurate, automatic spot-check vitals and continuous monitoring sent to your EMR from the bedside to help support improved patient safety and clinical decisions.

The Connex Vital Signs Monitor represents the next generation in vital signs monitoring with an easy-to-use color touchscreen and advanced vital signs technology to monitor vitals, respiration and motion. Plus, EMR integration gives clinicians immediate access to the information they need—where and when they need it.



## Connect to the future of vital signs today

### Designed specifically for busy floors like yours to increase efficiency

- > Captures both spot-check vitals and continuous monitoring
- > EMR integration for accurate electronic vitals documentation
- > Vivid color touchscreen simplifies training and operation
- > Comprehensive point-of-care documentation including customizable data fields
- > Advanced Patient ID including two-way name confirmation
- > Patient data protection for HIPAA compliance
  - Protect patient identifiable information with customizable options like abbreviations
  - Allow automated deletion of patient data after sending to EMR
  - Apply rules to check record completeness before saving
- > Easy-to-understand error messages for simple troubleshooting

### Industry-leading parameters for accurate patient information

- > SureBP® blood pressure technology
  - Adult, pediatric and neonatal blood pressure modes
  - Customizable blood pressure monitoring programs help nurses follow hospital protocol
  - 15-second readings for greater patient comfort
- > SureTemp® Plus and/or Braun ThermoScan® PRO 6000 ear thermometry (optional—choose one or both)
- > Masimo® or Nellcor® SpO<sub>2</sub> algorithms (optional)—both sensors and signal processing
- > Masimo noninvasive continuous monitoring of total hemoglobin (optional)
- > Manual measurements entry (weight, height, respiration rate, pain) and documentation of modifiers like body position and O<sub>2</sub> therapy
- > Customizable documentation fields to fully document at the point of care

## Configurations Available

CONNEX VSM	6700	6800
<b>Parameters</b>		
SureBP NIBP	●	●
SureTemp Plus Thermometry	●	●
Braun PRO 6000 Ear Thermometry	●	●
Masimo or Nellcor SpO <sub>2</sub>	●	●
EarlySense Contact-Free Motion	⬆	⬆
Covidien Microstream EtCO <sub>2</sub>	⬆	⬆
Masimo RRA	⬆	⬆
Masimo SpHb	⬆	⬆
<b>Communications</b>		
4 USB Ports	●	●
Ethernet	●	●
802.11 a/b/g Wireless Radio	⬆	●
<b>Accessories</b>		
Barcode Scanner	●	●
Weight Scale	●	●
EMR Connectivity	●	●
Printer	●	●

● Included    ● Optional    ⬆ Upgradeable

#### Features:

- SureBP technology can provide BP in only 15 seconds BP averaging to help ensure accurate hypertension detection and diagnosis
- SureTemp Plus Oral / Axillary Thermometry
- Option for Nonin® SpO2 technology
- 7" color touchscreen for a simplified workflow
- Options for documenting patient observations (i.e. Respiration Rate, BMI)
- EMR connectivity helps minimize documentation errors
- 8 hour Lithium ion battery to maximize time before recharging
- Each device ships with one adult blood pressure cuff



WEL44XT .....	Spot 4400 with NIBP and SureTemp® Thermometer
WEL44WT .....	Spot 4400 with NIBP, SureTemp Thermometer and Nonin® SpO2
WEL4400-DST .....	Spot 4400 Desk Stand
WEL4400-MBS .....	Spot 4400 Mobile Stand
WEL77794-M4400 .....	Green Series™ 777 Integrated Wall System for Spot 4400
WEL77794-2M4400 .....	Green Series 777 Integrated Wall System with PanOptic™ Ophthalmoscope for Spot 4400
WELBATT22 .....	Connex® Spot, 2-Cell LI-Ion Battery
WEL106275 .....	Micro USB Cable
WELREUSE-11 .....	Cuff Reusable Adult
WELREUSE-12 .....	Cuff Reusable Large Adult
WEL4500-35 .....	Fast BP Hose with FlexiPort®, 10 ft

The Welch Allyn ProBP 2400 is not a typical electronic blood pressure device. It is an affordable, professional device that you can trust to accurately measure your patients' blood pressures despite many common underlying circumstances. The ProBP 2400 is specifically designed for practices like yours. It gives you the quality, durability and features you don't get from home blood pressure devices that are mainly designed for general household use, which some practices resort to because of cost. It provides you consistent technique that's not assured with manual blood pressure gauges. And it does so at a budget-friendly price that other professional blood pressure monitors with the same features don't offer.

#### Features

- More affordable than similarly featured professional blood pressure monitors
- 3-Reading Average Mode can give a broader view of patient's true blood pressure
- Irregular heartbeat detection provides data to help you make testing decisions
- Manual Mode allows it to serve as a digital gauge when auscultation is indicated
- Microprocessor-controlled deflation rate provides consistent technique
- Built-in folding handle allows you to carry it safely and comfortably
- Optional mobile stand makes transporting easy for maximum location flexibility
- Optional wall mount bracket ensures device is exactly where and when you need it



WEL2400 .....	ProBP 2400 BP Unit with Adult & Large Adult Cuffs
WEL2400-MS .....	Mobile Stand for ProBP 2400
WEL2400-WMB .....	Wall Bracket for ProBP 2400
WELREUSE-09-2400 ....	Blood Pressure Cuff, Child (14 - 22cm)
WELREUSE-11-2400 ....	Blood Pressure Cuff, Adult (22 - 32cm)
WELREUSE-12-2400 ....	Blood Pressure Cuff, Large Adult (32 - 42cm)



The Connex® Spot Monitor features an easy-to-use, vivid touchscreen display and provides accurate vital signs measurement including blood pressure averaging, spot checking, interval monitoring and custom scoring across patient populations. The device connects wirelessly to your EMR and sends vitals to the patient chart from the point of care to help improve efficiency. Its upgradeable design and custom configurations include optional wall unit integration, industry-leading thermometry and pulse oximetry choices, connectivity and other mounting/mobility options. In addition, the Accessory Power Management Stand with WhisperDrive™ technology features additional power backup, providing up to 17 hours of on-time, alleviating the need to charge the monitor as frequently. Ships standard with one Adult and one Large Adult cuff.

### Features

- SureBP® technology can provide BP in only 15 seconds
- Connex Scoring App calculator for custom protocols (e.g., NEWS, MEWS)
- BP averaging to help ensure accurate hypertension diagnosis
- Supports resident list workflows for long-term care settings
- SureTemp® Plus or Braun ThermoScan® PRO 6000
- Thermometry Partnerships with more than 90 leading EMRs including Epic® and Cerner®
- Choice of leading SpO2 technologies: Masimo®, Nellcor® or Nonin®
- Add to integrated wall system for complete patient exam workflow
- Scale connectivity to Detecto, Health-o-meter® and SECA scales
- Mobile stand option includes work surface with integrated battery for extended usage

WEL71XX-B .....	Digital BP
WEL71WX-B .....	Digital BP, Oximeter
WEL71WT-B .....	Digital BP, Oximeter & SureTemp Thermometer
WEL7000-MS3 .....	Mobile Stand
WEL7000-GCX .....	Wall Mount
WEL7000-DST .....	Desk Stand



Affordable and easy to use, Welch Allyn FlexiPort Reusable Blood Pressure Cuffs offer all the advantages of FlexiPort technology—including simplified connectivity—in a dependable, durable design.

## Features

- FlexiPort technology enables every cuff to be used as a one- or two-tube cuff
- Tested to withstand demanding clinical environments
- Latex-free material minimizes risk of allergic reactions
- Rotatable port, improves patient comfort and cuff durability
- Folded edge reduces risk of cuts and scrapes, maximizing patient comfort
- Meets all latest clinical guidelines for proper fit from AAMI and AHA
- Primary sizes are color-coded to make it easy to find the proper cuff size
- 3 Year warranty

## FlexiPort™ Cuffs (No Tubing)

WELREUSE-06 ..... Small Infant  
 WELREUSE-07 ..... Infant  
 WELREUSE-08 ..... Small Child  
 WELREUSE-09 ..... Child  
 WELREUSE-10 ..... Small Adult  
 WELREUSE-11 ..... Adult  
 WELREUSE-12 ..... Large Adult  
 WELREUSE-13 ..... Thigh  
 WELREUSE-MLT ..... 8-Cuff Kit

## Complete Inflation System (includes Cuff, Tubing, Bulb and Valve)

WELREUSE-09-2BV ..... Child  
 WELREUSE-11-2BV ..... Adult  
 WELREUSE-12-2BV ..... Large Adult  
 WELREUSE-13-2BV ..... Thigh



## WelchAllyn® INFLATION BULB & VALVE COMBINATION

WEL5088-01 ..... Bulb & Valve Combination  
 WEL5087-01 ..... Air Valve Only  
 WEL5086-01 ..... Bulb Only



## WelchAllyn® COIL TUBING with TWIST FIT CONNECTOR

WEL5089-13 ..... 8"







**WelchAllyn®**

## 767 SERIES ANEROIDS

The versatile 767 Aneroid is the smart choice for accurate and reliable readings throughout your facility. Improved ergonomics and a new modern design make for an easy-to-use blood pressure aneroid, developed with staff and facility workflow efficiency in mind. Recessed, laser-engraved dial face with high-contrast pointer provides certified accuracy to reduce errors and improve patient care.

WEL7670-01 ..... Aneroid with adult cuff, wall mount

WEL7670-03 ..... Aneroid with adult cuff, mobile stand



**WelchAllyn®**

## DS58 HAND ANEROID

With gear-free DuraShock™ technology, you can count on the DS58 aneroid to remain in calibration longer than traditional aneroids, providing more accurate readings over time and a lower lifetime cost, helping to improve practice efficiencies and your bottom line.

- Gear-free DuraShock™ technology
- Certified accuracy to  $\pm 3$  mmHg
- Laser-engraved dial face for accuracy
- Jeweled movement for long life
- Lifetime calibration warranty
- Latex-free for safety
- Thumbscrew bleed valve and ergonomic inflation bulb
- Unsurpassed reliability in a traditional, handheld design
- Designed to enhance comfort
- Withstands a 60" drop, exceeds AAMI's shock resistant requirement

WELDS58-11 ..... With adult cuff and zippered case

WELDS58-MC ..... With 4 sized cuffs and zippered case



**WelchAllyn®**

## DIAGNOSTIC SET, COAXIAL

Welch Allyn 3.5V Diagnostic Set with PanOptic Basic LED Ophthalmoscope, MacroView Basic LED Otoscope, one Lithium Ion Plus Rechargeable Power Handle, and Hard Case.

WEL71-PM2LXU .....Diagnostic set

## WelchAllyn® 3.5V NICKEL-CADMIUM RECHARGEABLE HANDLE

Includes #72300 battery. A choice of convenient rechargeable power sources that fit all the WelchAllyn 3.5V instrument heads. Smooth, non abrasive knurled finish ensures a no-slip grip. Chrome-plated brass construction.

WEL71000-A



## WelchAllyn® 3.5V LITHIUM ION RECHARGEABLE HANDLE

Twice the on-time, half the weight and size. Fits all Welch Allyn 3.5V Instrument Heads. 120 minutes of on-time (compared to 60 minutes for standard handles). Low battery indicator, no "memory" effect. Includes AC charging module.

WEL71900



## WelchAllyn® COAXIAL OPHTHALMOSCOPE HEAD

- Halogen HPX lamp illumination for true tissue color
- 18 unique aperture/filter combinations for greater versatility
- Polarizing filter virtually eliminates corneal blue filter
- Compatible with all existing Welch Allyn 3.5V power sources
- Look for the traditional "flat-top" design you've come to trust

WEL11720-L ..... Ophthalmoscope head only



## WelchAllyn® PANOPTIC™ BASIC LED OPHTHALMOSCOPE

- Welch Allyn PanOptic Basic LED Ophthalmoscope; power handle not included
- 5X larger view of the fundus vs. standard ophthalmoscopes in an undilated eye
- Digitally capture, store and share fundus images by combining with iExaminer
- Greater working distance improves comfort for both practitioner and patient
- Dynamic focusing wheel allows continuous, smooth action and more precise control
- Hallogen HPX® lamp provides bright, white light
- Compatible with all existing Welch Allyn 3.5V power sources

WEL118-2 ..... Head only





## WelchAllyn® MACROVIEW™ OTOSCOPE HEAD

- A nearly complete view of the tympanic membrane
- Better clarity and definition of landmarks
- The ability to adjust focus for variable ear canal length or farsighted eyes
- Lip grip for secure fastening and ejection of ear specula
- Fiber optics produce cool light with no reflections, no obstruction
- Twice the field of view and 30% greater magnification than traditional otoscopes

WEL238-2 ..... MacroView Basic LED Otoscope

WEL238-3 ..... MacroView Plus LED Otoscope



## WelchAllyn® DISPOSABLE OTOSCOPE SPECULUM

WEL52432-U ..... Pediatric (2.75mm) Bag / 850

WEL52434-U ..... Adult (4.25mm) Bag / 850




## WelchAllyn® SPOT™ VISION SCREENER

The Spot Vision Screener is a handheld, portable device designed to help users quickly and easily detect vision issues on patients from 6 months of age through adult. Spot screens both eyes at once from a nonthreatening 3' distance. The touch-screen display allows for one-touch activation, simple management of patient data entry, and easy configuration for both vision screening and autorefraction applications.

### Features:

- Supports AAP screening guidelines for early detection of amblyopic risk factors
- Automated screening provides thorough and easy-to-understand results
- Easy to use, easy to implement with minimal user training required
- Adequately powered for mass screening events
- Captures readings 97% of the time helping to screen otherwise difficult patients
- Wireless printing and WiFi-enabled for easy export of data
- Manual exports available via USB port
- Uses lights and sounds to help engage children

WELVS100S-B



**Our shelves are deeper  
than these pages.**

Ask us.

## WelchAllyn® REPLACEMENT BULBS

WEL03000-U6 .....For Ophthalmoscope (11710)  
 WEL03100-U6 .....For Otoscope  
 WEL04200-U6 .....For Exam Light (48600 & 48700)  
 WEL04900-U6 .....For Ophthalmoscope (11720)  
 WEL06300-U6 .....For Exam Light (44100 & 44110)  
 WEL06500-U6 .....For MacroView™ (23810 & 23820)  
 WEL07800-U6 .....For Vaginal Illuminator (78000)  
 WEL08800-U6 .....For Vaginal Illuminator (78810)



## WelchAllyn® RECHARGEABLE BATTERIES

WEL72200 .....Fits 71000 & 71670 Handles (Black)  
 WEL72300 .....Fits 71000-A Handle (Orange)



## 3M LITTMANN™ CARDIOLOGY IV STETHOSCOPE

- Outstanding acoustics with better audibility of high-frequency sounds (adult side, firm pressure)
- Tunable diaphragm; hear high or low frequency sounds by slightly adjusting pressure on the chestpiece
- Dual-lumen tubing: two sound baths in one tube eliminates the rubbing noise of traditional double tubes

3M 6152.....Black	3M 6153 ..... Burgundy
3M 6154.....Navy Blue	3M 6155..... Hunter Green
3M 6156.....Plum	3M 6157..... Carribean Blue
3M 6158.....Raspberry	3M 6161 ..... Rose Pink



## DOUBLE HEAD STETHOSCOPE WITH ALUMINUM CHESTPIECE

ALM10-1510..... Black  
 ALM10-1530 .... Dark Blue

## DOUBLE HEAD STETHOSCOPE ECONOMY MODEL

ALM10-1700 ...Black

## LITTMANN™ CLASSIC III STETHOSCOPE

The 3M™ Littmann® Classic III™ Stethoscope offers high acoustic sensitivity for exceptional performance when doing general physical assessments. It features dual tunable diaphragms and an updated design that is easier to clean and maintain.

The 3M™ Littmann® Classic III™ Stethoscope is the latest version of the stethoscope that helps millions of medical professionals achieve their best. The Classic III stethoscope offers high acoustic sensitivity for exceptional performance, plus a versatile two-sided chestpiece with tunable diaphragms. The Classic III stethoscope is used by students and medical professionals alike to identify, listen to, and study heart, lung, and other body sounds for physical assessment, patient monitoring, and diagnosis.

Tunable diaphragm technology—a 3M invention—lets clinicians hear different frequency sounds by simply adjusting the pressure on the chestpiece. Hold the chestpiece with light pressure to hear low frequency sounds; press a little more to hear higher frequency sounds. This time-saving feature allows the healthcare provider to focus on the patient, while the non-chill diaphragm provides patient comfort.

Headset tension is easily adjusted for head size and comfort by squeezing together or pulling apart the eartubes. Snap-tight, soft-sealing eartips conform to individual ears for an excellent acoustic seal and comfortable fit. Eartips snap firmly onto the ends of the eartubes and, for safety, require a significant effort to remove.

Resilient next-generation tubing retains its shape and flexibility even after folding tightly into a pocket. It provides longer stethoscope life due to improved resistance to skin oils and alcohol, and is less likely to pick up stains. Next-generation tubing is not made with natural rubber latex, a benefit to sensitive medical professionals and their patients. And, to play its part in the protection of human health and the environment, the next-generation tubing does not contain phthalate plasticizers.

The Classic III stethoscope is made in the USA for consistent quality. It comes complete with large and small soft-sealing eartips, non-chill bell sleeve, and instructions.

### Features:

- Two-sided chestpiece with tunable diaphragms on both the adult and pediatric sides
- Single-piece tunable diaphragm is easy to attach, and easier to clean because its surface is smooth without crevices
- Pediatric side converts to a traditional open bell by replacing the single-piece diaphragm with a non-chill rim
- Next-generation tubing provides longer life due to improved resistance to skin oils and alcohol; less likely to pick up stains
- Small tunable diaphragm is useful for pediatric, small, or thin patients; around bandages; and for carotid assessment
- Stainless steel chestpiece is precision-machined into an aesthetically pleasing, less angular shape. The stem features open side indicator
- Snap-tight eartips have a soft, smooth surface providing a comfortable acoustic seal and comfortable fit

3M 5620..... Black  
 3M 5626..... Raspberry  
 3M 5832..... Lavender  
 3M 5839..... Lemon Lime  
 3M 5639..... Rose Pink  
 3M 5621..... Grey  
 3M 5629..... Orange  
 3M 5829..... Lime Green

3M 5622..... Navy Blue  
 3M 5627..... Burgundy  
 3M 5835..... Turquoise  
 3M 5840..... Emerald  
 3M 5624..... Hunter Green  
 3M 5623..... Caribbean Blue  
 3M 5630..... Ceil Blue  
 3M 5831..... Plum





## DIGITAL FINGER TIP PULSE OXIMETER

Moisture proof, one-button operation, requires two AAA batteries, 1 Year limited warranty.

- Checks heart rate & respiratory function
- Suitable for finger thickness: .37" - .86" (9.5mm - 21.9mm)
- Not made with natural rubber

DYN7090 .....Adult

DYN7090 .....Pediatric



## PHYSIO LOGIC DIGIPULSE OXIMETER

- One-Touch Operation
- Accurate Readings (Oxygen Level, Pulse Rate)
- OLED Multi-Directional Display
- Compact and Portable
- Low Power Consumption
- Auto Shut-Off
- Adjustable Alarm, Sound, and Brightness

AMG107-100



## PULSE OXIMETER, HAND HELD

The Edan H100B pulse oximeter with powerful memory (300 hours) and configurable alarms provides an easy to use tool to measure oxygen levels (SpO2) and heart rate information. The Edan H100B has a wide variety of sensors options (neonate, pediatric, adult) to fit any patient sizes. Memory from the Edan H100B can be downloaded and printed using the optional software and USB cable.

### Edan H100B Hand Held Pulse Oximeter Includes:

- Operator's Manual & Quick Reference Card
- Adult Finger Clip Sensor
- 4 AA alkaline batteries
- Operator's Manual
- Deluxe Carrying Case

SOLH100B .....Pulse Oximeter

SOL02.01.210121.....Pediatric Sensor

SOL01.57.040199.....Infant Disposable Sensor

SOL03.25.110260.....Data Management Kit (Data from the H100B can be transferred to PC for storage, review and printing)



## EDAN F6 FETAL MONITOR

- 10.1" foldable, color TFT-LCD screen
- Display waveforms and digitals simultaneously
- Signals Overlap Verification to separate twins FHR
- Professional CTG analysis system
- Compatible with 150/152 mm universal recording paper
- 60 hours waveform storage and playback
- Wired & wireless central monitoring solution
- Telemetry Transducer System compatibility

SOLF6..... Fetal Monitor

SOL01.57.471869..... Recording Paper





**Welch Allyn**

**Hillrom**

## THERMOSCAN® PRO 6000 TYMPANIC EAR THERMOMETER

- Innovative PerfectTemp™ technology adjusts for variability in probe placement
- ExacTemp™ technology detects stability of the probe during measurement
- Designed to take quick and accurate readings, to improve patient satisfaction
- 60 second pulse timer assists with manual measurement of pulse rate, respiration
- Electronic security features help prevent theft, loss
- Small or large cradles for probe cover storage that pop-up for easy attachment
- Optional charging station stores 200 probe covers, with rechargeable battery
- 6'-9' security tether options keep device attached to cradle

WEL06000-200 ..... Braun PRO 6000 with Small Cradle (holds 1-pc box)

WEL06000-005 ..... Braun PRO 6000 Probe Covers Box/200



**Welch Allyn**

**Hillrom**

## SURE TEMP® PLUS 690 ORAL THERMOMETER

- Fast: quickly captures accurate readings on all patient types
- Large LCD displays temperatures in Fahrenheit or Celsius
- Convenient storage for 25 probe covers
- Last temperature recall button
- Optional security wall mount
- Takes approximately 6,000 readings on 3 AA batteries
- Interchangeable, removable wells to help reduce the risk of cross-contamination

WEL01690-200 ..... Thermometer

WEL05031-105 ..... Probe Covers ..... Box/250



## PROSCAN NON-CONTACT INFRARED THERMOMETER

- 3-in-1 functionality for body, room, or surface/liquid temperature checks
- Non-contact readings at 3–5 cm (1–2 inches)
- Easy-to-read, color-coded fever alert system
- Switchable temperature display in Celsius (°C) or Fahrenheit (°F)
- Comes with 2 AA batteries for immediate use
- 2 year limited warranty

AMG016-660

## EDAN® SONOTRAX DOPPLER

- A high performance Fetal Doppler with heart rate
- Compact design
- Built-in audio recorder
- Large LCD display & backlit
- Accurate FHR detection with clear sound
- Automatic / Average / Manual fetal heart rate counting
- Highly sensitive interchangeable probes
- 9V alkaline battery



SOLMS9104661 .....3mHz OB Probe with LCD Display

## SONOTRAX VASCULAR DOPPLER

- SonoTrax Vascular has been widely adopted as standard equipment for the assessment of peripheral vascular disease.
- The versatile, compact and portable Sonotrax Vascular allows fast and simple but effective assessment of both arterial and venous disease
- 4 or 8 MHz waterproof probe
- Built-in audio recorder
- No digital display
- Accurate detection with clear sound
- Highly sensitive interchangeable probes



SOLMS314321..... SonoTrax Vascular Doppler with 8 MHz probe

SOLMS314322..... SonoTrax Vascular Doppler with 4 MHz probe

## EDAN U50 PORTABLE ULTRASOUND SYSTEM

- 12.1" high resolution LCD display with intuitive control panel
- Two transducer connectors
- Lithium battery for extensive remote use
- Multiple peripheral ports including VGA, composite video and USB
- Large capacity integrated hard disk for extended data storage
- Dicom supported
- eSRI: Speckle Reduction Imaging, Spatial Compound Imaging, Quick Imaging tools for B Mode and Color Mode, Multi-Beam Imaging
- Needle guide, PW, CW, PDI/DPDI, HPRF
- For multiple applications, including gastroenterology, obstetrics, pediatrics, gynaecology, cardiology, urology, MSK
- Probe Options available



SOLU50



## AQUA-SONIC® 100 ULTRASOUND GEL

### Blue Gel

PAR01-01 ..... 20ml Foil Package - Sterile ..... Box/48  
 PAR01-08 ..... 250ml  
 PAR01-50 ..... 5L

### Clear Gel

PAR03-08 ..... 250ml  
 PAR03-50 ..... 5L



## AQUAFLEX® GEL PAD

PAR04-02 ..... 2cm x 9cm ..... Box/6



## ECLIPSE® ULTRASOUND PROBE COVERS

PAR38-01 ..... Eclipse Pre-lubed Probe Cover Box/100  
 PAR38-03 ..... Eclipse 3D Pre-lubed Probe Cover Box/100



## THERMASONIC ULTRASOUND GEL WARMER

PAR82-01 ..... Single Bottle Warmer  
 PAR82-03 ..... Three Bottle Warmer



## TRANSEPTIC® CLEANSING SOLUTION

Multi-purpose spray for professional use in cleansing ultrasound transducer/probe surfaces between procedures.

PAR09-25 ..... 250ml Spray Bottle



## POLYSONIC® ULTRASOUND LOTION

- Multi-purpose ultrasound lotion with superior coupling efficiency in a pleasing cosmetic quality base. Recommended for all ultrasound procedures when a lotion is preferred
- Bacteriostatic
- No formaldehyde

PAR21-08 ..... 250 ml Bottle



## DERMLITE®

The versatile Dermlite has a multitude of applications. Viewing pigmented lesions, varicose veins, hair transplants, skin irritations, or skin lesions is just the beginning. Together with the Dermlite's compact size and lightweight design (only 3.5oz/99g), the possibilities are virtually infinite. To make sure you will not run out of power, the Dermlite's endurance measures over 1200 lesion examinations per battery, due to its diode lighting technology and long-life lithium power.

DERDL100 .....DL-100 Unit

## NEUROPEN® SENSATION SCREENING DEVICE

**Includes:** Neuropen, 1 x 10g Monofilament, 1 x Neurotip

The Neuropen® is a reusable device that features a calibrated spring mechanism to deliver a 40g sharpness test using Neurotips® single use neurological examination pins, to assess reduced sensation to sharpness/ pain in small nerve fibers and a 10g pressure test using Neuropen Monofilaments® to assess protective touch/pressure sensation in large nerve fibres.

MDPNT0100 .....Neuropen Device

MDPNT0104 .....10g Monofilaments Box/5

MDPNT5405 .....40g Neurotips Box/100



SE-1200 Express



## SE-1200 Express 12-channel ECG with Interpretation

- Large high resolution color touch screen
- Alphanumeric keyboard
- One-touch operation
- Enhanced weak signal detection
- Comprehensive filters & anti-noise technology
- Real-time ECG waveforms freezing and review
- Advanced Glasgow algorithm for automatic measurements and interpretation
- Built-in high resolution thermal printer with external USB printer (optional) support
- Extended printing
- Large internal memory which can be enlarged by USB flash disk
- Built-in data management with reports review
- LAN/Wi-Fi (optional) connections
- Bi-directional communication with SE-1515 Data Management System
- Multi-format reports export (optional)
- Barcode scanner support (optional)

SOLSE1200-EXB .....SE-1200 ECG Basic

SOLSE1200-EXPRESS .....SE-1200 ECG with touch screen



## SPIROMETER DISPOSABLES

WEL703419 .....Disposable Flow Transducer Pkg/100

WEL100680 .....Spirometry Noseclips Pkg/25







WelchAllyn®

Hillrom.

## CP 150™ ELECTROCARDIOGRAPH

The CP 150 combines the best traditional box ECG features (e.g. full-size printouts and clinical decision support with MEANS interpretation) with connectivity enhancements designed especially for today's primary care practices.

- ECG print previews can help save time and reduce paper costs
- 7" color, touchscreen display for easy entry of patient data
- Optional wireless connectivity helps increase device mobility and utilization
- Optional spirometry lets you perform two tests with one device
- Optional DICOM compatibility enhances interoperability
- Optional adult and pediatric MEANS ECG interpretation software
- 3, 6 or full 12-lead ECG preview options
- Battery operation lets you easily transport your device to the patient
- Simplified EMR connectivity without a complicated interface
- Instant-on powers the device up quickly to take an ECG at a moment's notice

WELCP150-1ENB ..... ECG Non-interpretive

WELCP150A-1ENB..... ECG Interpretive

WEL105353 ..... CP100/150/200 Recording paper.....Box/200

WEL105341 ..... CP150 - Office Cart

## WELCH ALLYN® DIAGNOSTIC CARDIOLOGY SUITE

WelchAllyn®

Hillrom.

The Welch Allyn Diagnostic Cardiology Suite now includes ECG and SPIROMETRY capabilities, helping you shorten the path to diagnosis for improved patient outcomes.

### Features:

- Captures best 10 seconds of ECG data from memory, which helps reduce need for repeat ECGs
- ECG complies with ACC, AHA and HRS recommendations for adult/pediatric patients
- Easily access prior spirometry exams to monitor treatment plans over time
- Utilize one device for spirometry and ECG exams for a common workflow and setup
- ECG and spirometry applications launch directly from the EMR
- Instantly verify spirometry exam quality and test performance



**WELCC-RXX-WAXX**..... Connex Cardio PC based 12-Lead Multi-Channel Resting ECG - software, Connex Cardio WAM wireless acquisition module, AHA 10 wire, 12-lead lead set with banana plugs, CardioSens Ultra II resting tab electrodes (pouch of 100), Clear Choice adapter clips (set of 10), Connex Cardio UTK

**WELCC-RXX-AAXX**.....Connex Cardio PC based 12-Lead Multi-Channel Resting ECG - software, Connex Cardio AM12 wired acquisition module, AHA 10 wire, 12-lead lead set with banana plugs, CardioSens Ultra II resting tab electrodes (pouch of 100), Clear Choice adapter clips (set of 10)

**WELCC-RXX-AADX** .....Connex Cardio PC based 12-Lead Multi-Channel Resting ECG - software, Connex Cardio AM12 wired acquisition module, AHA 10 wire, 12-lead lead set with banana plugs, CardioSens Ultra II resting tab electrodes (pouch of 100), Clear Choice adapter clips (set of 10), DICOM

**WELDCSS-CXX**.....Diagnostic Cardiology Suite Spirometry with Calibration Syringe

**WELDCSS-NXX**.....Diagnostic Cardiology Suite Spirometry without Calibration syringe

**WEL105341** .....Office Cart

## DISPOSABLE HOLTER POUCHES

Designed to secure the Digital Holter device to the patient's waist or around the neck during testing

- 4 x 6"
- 25 per Bag

NIOP100



## EDAN 3-CHANNEL HOLTER ECG WITH 5 LEADS

- Lightweight Edan 3-channel Holter ECG with 5 lead option
- Up to 7 days of recording with one AAA battery
- IP27 waterproof and dustproof
- Real-time waveform display
- Industrial-grade Micro SD card
- Pacemaker detection

SOLSE2003-7D      Edan 3-channel Holter ECG  
with 5 leads, color LCD screen, 7 Days recording  
SOLE2012SW..... Software for the Holter SE-2003-7D



## Need something specific?

A quick call can save you time.



## 3M RED DOT™ RESTING ECG TAB ELECTRODE

3M2330E ..... Pkg/100



## SKINTACT® MONITORING SOFT CLOTH ELECTRODES

- More comfort for patients due to breathable
- backing material, ideal for long-term monitoring
- Solid Gel

CMLW-601 ..... 30 Electrodes per Pack, 50 Packs per Box



## 3M RED DOT™ PEDIATRIC MONITORING ELECTRODE WITH 3M™ MICROPORE™ TAPE BACKING

- Smaller size for increased conformability and comfort, suitable for pediatric patients
- Lift tab for easy removal
- Adhesive is gentle on skin for minimal irritation

3M 2248-50 ..... 50 Electrodes per Bag



## 3M RED DOT™ MONITORING ELECTRODE

3M 2670-5 ..... Flexible Non-Woven and Sticky Gel..... Case/1000



## 3M RED DOT™ RESTING EKG DIAGNOSTIC ELECTRODE

Easy to apply and remove, 3M™ Red Dot™ Resting EKG Electrodes are ideal for the unique needs of resting EKG monitoring. The conductive adhesive provides intimate skin contact to ensure accurate EKG traces

- Little to no residue or adhesive transfer
- Can be repositioned
- Aggressive adhesive

3M 2360 ..... 100 per Bag



## ME MEDICAL SURGICAL CLIPPER

The ME Surgical Clipper is the ultimate in hair removal. Skin preparation has never been faster, easier, safer or more cost effective.

- Revolutionary head twists 45 degrees to contour to even the most difficult-to-reach areas.
- Easy, no touch insertion and ejection of disposable clipper blades.
- Easy-to-hold non-slip grip handle
- Rechargeable lithium battery fully charges in just one hour
- Two types of disposable clipper blades: the Universal designed for most hair and clipping needs, and the Neuro for coarse or dense hair

NIOSC001 ..... Rechargeable Surgical Clipper (with charging station and 2 blades)

NIOSC002 ..... Universal disposable blades 50 per Box

NioSC003 ..... Neuro disposable blades 50/ per Box

## ELECTRODE SKIN PREP PAD

Formulated to remove excess skin oil enhancing ECG electrode to skin contact.

DYN1508 .....Skin Prep Pad.....Box/100



## 3M RED DOT™ ECG TRACE PREP TAPE

Designed to remove the non-conductive skin layer.

3M 2236.....18mm x 5m.....Roll



## GALLANT DISPOSABLE PREP RAZORS

The Dynarex Gallant Disposable Prep Razors feature a stainless steel blade and open design allowing an optimal view of the shaving area. Gallant prep razors are easy to use, requiring only light gentle pulling instead of pressure.

- Stainless Steel Blade
- Disposable
- Ideal for prepping prior to procedures such as ECGs

DYN4251 ..... 50 per Box



## PERSONNA PREP RAZOR (PREMIUM)

- Specially designed comb (rake)
- Wrapped blade reduces cuts and nicks
- Polymer Coated Stainless Steel Blade
- Heavy weight handle design
- Fixed head
- Blades precision ground for comfort

PER75-4004 .....Single Edge.....Box/100



## TWIN BLADE RAZORS

- Long-lasting, superior quality blade
- Clear plastic guard
- Disposable
- Unisex

DYN4250 .....Twin Blade.....Box/50





**ZOLL**  
Advancing Resuscitation. Today.™

**AED PLUS**

## SEMI-AUTOMATIC DEFIBRILLATOR



The AED Pro is designed for both basic and advanced life support. While it offers guidance to the BLS rescuer, it also provides more advanced capabilities to the professional rescuer. The AED Pro not only relays advanced CPR feedback, it's also a vital signs patient monitor that can be configured to the way you work.

- A three-lead cable allows basic monitoring using ECG electrodes
- Flexible operating modes let the professional rescuer apply manual override to access more advanced capabilities, including heart rate displays, and to decide if and when to deliver treatment
- A high-resolution auto-gaining LCD display keeps the ECG trace within a specific area, always displaying the gain, so the rescuer can easily see the patient's rhythm

ZOL22600010102011060



## AED PLUS METAL WALL CABINET WITH ALARM

ZOL8000-0855



## AED PLUS ELECTRODES CPR-D PADZ®

ZOL8900-0810-01..... Pediatric

ZOL8900-0800-01..... Adult



## AED PLUS LITHIUM BATTERIES

ZOL8000-0807..... Pkg/10

**ZOLL**  
Advancing Resuscitation. Today.™



## ZOLL AED 3 DEFIBRILLATOR

Enhanced Real CPR Help, featuring a full-color display with vivid rescue images, CPR cycle timer, and a large color bar gauge that shows CPR compression depth. RapidShock Analysis™ enables the industry's shortest rhythm analysis for more continuous life-saving CPR. CPR Uni-padz™ electrodes can be used on both adult and pediatric patients, so there is no need for two separate sets of pads. Cloud connectivity that can automatically report AED readiness. 5-year shelf life for pads and batteries increases readiness and decreases frequency of maintenance. Available in semi-automatic and fully automatic models.

ZOL8501-001102-06



## CAVICIDE1® SURFACE DISINFECTANT

CaviCide1, a multi-purpose disinfectant/decontaminant cleaner, can be used on hard, non-porous surfaces. CaviCide1 is effective in killing bacteria and fungi in 1 minute, including TB, Norovirus, HIV-1, HBV, HCV and many more. When used as directed, fragrance-free CaviCide1 will effectively clean and disinfect surfaces and can help reduce the risk of cross-contamination.

MET11-5024 ..... CaviCide1, 24oz spray  
MET11-5001 ..... CaviCide1, 1 Gallon  
MET11-5100 ..... CaviWipes1, Small, 6" x 7"  
MET11-5150 ..... CaviWipes1 X-Large, 9" x 12"

**Metrex™**  
INFECTION PREVENTION



## CAVIWIPES™ BLEACH

Disinfecting towelettes that demonstrates effectiveness against over 40 relevant microorganisms. 1:10 Bleach dilution. No PPE required\*, leave the gloves behind and reach for CaviWipes Bleach

MET11-1200 ..... CaviWipes Bleach (6"x10.5") 90 Wipes per Canister

**Metrex™**  
INFECTION PREVENTION



## CAVIWIPES™ 2.0

2-minute universal contact time for bacteria, viruses, and pathogenic fungi (Candida auris)

MET11-1410 ..... Caviwipe 2.0 (6"x6.75") 160 Wipes per Canister  
MET11-1450 ..... Caviwipe 2.0 X-Large (9"x12") 65 Wipes per Canister

**Metrex™**  
INFECTION PREVENTION



## CAVIWIPES™ HP

CaviWipes HP is the newest member of the trusted CaviWipes family providing an effective 1-minute universal contact time from an alcohol-free hydrogen peroxide formulation. The High Performance disinfection protocol offered by CaviWipes HP protects both patients and staff.

MET11-1610 ..... CaviWipes HP (6"x6.75") 160 wipes per Canister  
MET11-1650 ..... Caviwipes HP X-Large (9"x12") 65 Wipes per Canister

**Metrex™**  
INFECTION PREVENTION



## SANI-CLOTH® PLUS GERMICIDAL DISPOSABLE CLOTH

- Effective against 16 microorganisms in 3 minutes
- Quaternary ammonium

PDIQ90172 ..... Sani-Cloth® Plus 160 Wipes per Canister





## GLUTE-OUT NEUTRALIZER FOR GLUTARALDEHYDE & OPA

Disposing of used disinfectants can be extremely dangerous as glutaraldehyde and OPA products remain highly toxic and biocidal even after 14 days of re-use. Improper disposal of disinfectants by pouring down a drain has proven to cause high levels of exposure for the clinician and is considered illegal in many regions. Glute-Out Neutralizer offers a solution to dangerous disinfectant disposal problems. Glute-Out deactivates glutaraldehydes and OPA by using a glycine powder base that neutralizes solutions in just five minutes.

CV-610-1045 ..... 2oz bottle (Neutralizes 1 Quart)

CV-610-1099 ..... 5oz bottle (Neutralizes 1 Gallon)



## CIDEX OPA SOLUTION, DEVICE DISINFECTANT SOLUTION

- 14-Day High Level Medical Device Disinfectant Solution
- Sterilization in 10 hours at 25°C
- Activation required

JOH20394 ..... 3.8L

JOH20393 ..... Test Strips..... Bottle/15 Box/2



## METRICIDE 28® DISINFECTANT SOLUTION

MetriCide 28 must maintain the minimum effective concentration (MEC) level of 1.8% glutaraldehyde or greater during the 28-day use life to effectively sterilize and high-level disinfect. It is recommended to test the glutaraldehyde solution prior to each use to help prevent healthcare-associated infections (HAIs).

- Recommended for immersible and heat-sensitive equipment and instruments
- Contains surfactants

MET11-2801 ..... 3.8L

MET10-304..... Test Strips (1.8%)..... Bottle/60



## EASYTASK CLEANING SYSTEM 200WIPES/BOX

10"x12" wipe, add any cleaning or disinfecting solutions to make a durable wipe. Resealable bag keeps wipes moist. Holds 1/2 gallon of liquid, low lint.

CFSN-ETA110CZGW.....Dry Wipe Box/200 wipes

## EMPOWER® DUAL ENZYMATIC DETERGENT

EmPower® is used to pre-clean critical or semi-critical medical devices prior to terminal sterilization or high level disinfection. It contains two protease enzymes that offer the broadest cleaning action on a variety of protein soils. Inclusion of surfactants provides additional cleaning on carbohydrates, lipid and protein soils. Works fast to remove blood, tissue, mucous and other protein-rich bodily fluids from instruments.

- Dual-enzymatic detergent cleaner
- Gentle on instruments
- Less corrosive than distilled water
- Low-foaming formula is ideal for use in automated machines
- Fragrance free option for respiratory-sensitive users

MET10-4100..... 3.8L



## EMPOWER® FOAM SPRAY

EmPower® Foam is used to pre-clean critical or semi-critical medical devices prior to terminal sterilization or high level disinfection. It is a ready-to-use, foaming dual-enzymatic spray that starts the pre-cleaning process by breaking down blood, tissue, mucus and other protein-rich bodily fluids from instruments. The foaming spray delivery system is ideal for use in the operating room, sterile processing, and other departments where instrument cleaning is not immediately available.

- Ideal for instruments with the presence of moderate to heavy organic debris
- Covers instruments, reducing splashes and spills
- Ready-to-use in a convenient sprayer
- Covers instruments to rapidly begin the cleaning process

MET10-4224..... 24oz



## ENZYCLEAN PROTEASE INSTRUMENT DETERGENT

Enzyclean Protease Enzyme Detergent is specially formulated for use as a presoak or manual detergent.

The protease and amylase enzymes break down organic and inorganic soils, residues, films and deposits and ensures complete cleaning.

WEIEBL1HC.....3.78L



## CUTS-IT GEL INSTRUMENT PRE-SOAK

Clings to and dissolves blood, fat, tissue and other body fluids and keeps bioburden from drying on instrumentation.

Neutral pH, Non-corrosive, Phosphate-free, Biodegradable. Perfect for a transport gel.

WEIM60005HC.....710ml





## ES4 MANUAL HAND SANITIZER DISPENSER

### ES4 Hand Sanitizer Dispensers

GOJ5020-01..... Wall Dispenser, White

GOJ5024-01..... Wall Dispenser, Black/Chrome

### Purell Advanced Hand Sanitizer

GOJ5060-02.....Gel for ES4 Manual Dispenser, refill



## ES6 TOUCH FREE HAND SANITIZER SYSTEM

### ES6 Touch Free Hand Sanitizer Dispensers

GOJ6420-01..... Wall Dispenser, White

GOJ6424-01..... Wall Dispenser, Black

### Purell Advanced Hand Sanitizer

GOJ6460-02..... Gel for ES6 Touch Free Dispenser, 1200ml Refill



## GEL HAND SANITIZER

Kills 99.9% of most common germs that may cause illness. Specially formulated with moisturizers to leave hands feeling soft and refreshed.

### Gel

GOJ3770-12..... 354ml Pump Bottle

GOJ9674-12..... 8oz Bottle with Aloe

GOJ5015-04..... 1.5 Litre Pump Bottle



## HAND SANITIZING GEL 540ML



Ecolab's waterless hand sanitizer formulated to promote skin health and reduce the spread of infections. Contains moisturizers to condition and protect the skin and contains 70% ethyl alcohol.

ECO6000325 ..... Quikcare Gel Hand Sanitizer 540ml

## Purell ES4 MANUAL SOAP DISPENSER

### ES4 Hand Soap Dispensers

GOJ5030-01.....Wall Dispenser, White/White  
GOJ5034-01.....Wall Dispenser, Black

### Healthcare Foam Soap

GOJ5072-02.....1200ml Refill



## Purell ES6 TOUCH FREE SOAP DISPENSER

### ES6 Touch Free Soap Dispensers

GOJ6430-01.....Wall Dispenser, White  
GOJ6434-01.....Wall Dispenser, Black

### Health Care Foam Soap

GOJ6472-02.....1200ml Refill



Refill Sold Separately



Refill Sold Separately

## ECOLAB® BACTI-STAT® LIQUID ANTIMICROBIAL HANDWASH

Ecolab's Bacti-Stat is formulated to provide outstanding efficiency and exceptional mildness. Ecolab's hand soaps are designed specifically for the handwashing needs of healthcare workers and help reduce the potential for transmission of healthcare acquired infections, as well as promote skin health.

ECO6017285..... 4L



## ECOLAB® BACTISTAT HAND SOAP 540ML

- 0.3% Triclosan formula meets FDA requirements for Healthcare Personnel Handwash
- Lathers quickly, rinses freely which minimizes residual lather than can cause irritation
- Contains Aloe Vera which helps protect against irritation and dry, cracked skin

ECO6048722 .....Bactistat Hand Soap 540ml



## GREEN FOAM HAND SOAP 222ML

A mild, biodegradable foam soap formula certified by Green Seal™ and EcoLogo® to meet their joint standard for hand cleaners and hand soaps. The certifications mean that the product has been thoroughly evaluated by these unbiased, third-parties. Unscented.

GOJ5715-06.....Green Foam Hand Soap 222ml







## SCOTT® ESSENTIAL FOLDED PAPER TOWELS

UNI01700 ..... 250 Sheets per Pack,  
16 Packs per Case (4,000 Sheets/Case)



## TORK® PREMIUM HAND TOWELS

CFS101298 ..... Pack/100 Case/21 Packs  
CFS553120 ..... Dispenser



## TORK PREMIUM HAND TOWEL ROLLS

CFS774820 ..... Dispenser  
CFS8634550 ..... Roll Towel Case/12 Rolls



## WHITE SWAN® 2 PLY HOUSEHOLD PAPER TOWELS

CFS01890 ..... 90 Sheets/Roll Case/21 Rolls



## WHITE SWAN® 2 PLY QUICK WIPES

CFS08500 ..... 80 Wipes/Box Case/135 Boxes



## KLEENEX® PROFESSIONAL FACIAL TISSUE

UNI21606 ..... 125 Tissues per Box, 48 Boxes per Case



## TORK ADVANCED BATH TISSUE ROLL

- 2 Ply
- Individually wrapped rolls for hygiene and cleanliness

CFSTM6130S ..... 48 Rolls per Case



## TORK® ADVANCED MINI-JUMBO BATH TISSUE

CFS12024402.....2 Ply, 751' Roll Case/12  
CFS555020A.....Dispenser



## GARBAGE CAN, WHITE

- Steel square outer can and lid with removable plastic inner liner
- Enamel Finish Size
- 32 quart (8 gallon) capacity
- 11.5" D x 12.25" W x 18.25" H (29.21cm x 31.12cm x 46.36cm)

BRE35266



## GARBAGE/RECYCLING BAGS

CFSTT-20W..... Garbage Bag, Regular Strength, White (20" x 22") ..... Cs/500  
CFSTT-24..... Garbage Bag, Regular Strength, Black (22" x 24") ..... Cs/500  
CFSTT-24W..... Garbage Bag, Regular Strength, White (22" x 24") ..... Cs/500  
CFSTT-32..... Garbage Bag, Regular Strength, Black (26" x 36") ..... Cs/250  
CFSTT-52..... Garbage Bag, Regular Strength, Black (30" x 38") ..... Cs/250  
CFSTT-32TB..... Recycling Bag, Regular Strength, Blue (26" x 36") ..... Cs/250



## METRIMIST® ODOUR ELIMINATOR

MET10-1158.....240ml, 8oz



## LYSOL® DISINFECTANT SPRAY

Lysol® Disinfectant Spray kills 99.9% of viruses and bacteria on hard surfaces and eliminates odours by killing the bacteria that cause them. Also effective on mold and mildew.

UNI34052 ..... Crisp Linen Scent





## CLOROX® TOTAL 360® DISINFECTANT CLEANER

Kill COVID-19 (SARS-CoV-2, on hard nonporous surfaces), cold and flu viruses, norovirus, MRSA and more than 40 additional pathogens in two minutes

- Can be used with a Sprayer Device

CLX31650 ..... 3.78L bottle



## CLOROX® GERMICIDAL BLEACH 3.5L 3/CS

- Concentrated 7.4% bleach
- Health Canada registered to kill a total of 35 microorganisms, including C. difficile spores in 5 minutes
- DIN# 02408961 Concentrated 7.4% bleach

CLX01383..... Clorox® Germicidal Bleach 3.5L 3/CS



## DAWN ULTRA DISH SOAP 443ML

UNI439001..... Dawn Ultra Dish soap 443mL



## TIDE LIQUID HE LAUNDRY DETERGENT 1.36L

UNI41825 .....32 Loads per Jug

## midmark™ M11 STEAM STERILIZER

The 11" x 18" chamber makes this one of the largest of any standard countertop sterilizers on the market.

RITM11-050 ..... Midmark® M11  
- Automatic large tabletop with automatic door



## midmark™ M9 STEAM STERILIZER

Pack all the reliable sterilizing power you need into a compact unit perfect for areas with limited space

RITM9-050 ..... Midmark® M9  
- Fully automatic sterilizer with automatic door



## Tuttnauer EZ11PLUS STERILIZER

Features:

- 11" diameter X 19.8" deep chamber and flexible loading options
- Printer Included
- USB Port enables users to easily transfer cycle data via a flash drive to a PC
- EZGlide door features a sleek modern look and effortless opening and closing
- Closed door drying, meaning the door remains closed through the dry cycle to maintain sterility
- Comprehensive 2 Year parts & labor warranty and 10 Year chamber warranty

TutEZ11Plus-P



## MIDMARK™ QUICKCLEAN™ ULTRASONIC CLEANERS

Table top unit, basket is included. QuickClean™ washes instruments for sterilization, allowing caregivers and staff to spend their time with patients, rather than washing instruments.

**POWERFUL CLEANING. CONSISTENT RESULTS. SAFER FOR EVERYONE.**

RITQC1-01 ..... 1.2 Gallon  
RITQC3-01 ..... 3.3 Gallon  
RITQC6-01 ..... 6.6 Gallon



## MINIBIO AUTO READER

Tuttnauer MiniBio will allow 20 Minute & 1 hour process times simultaneously, holds 3 indicators with different incubation times, at the same temperature, USB connection with free software download, thermal printer.

TUTWTL198-0057



## ULTRA-RAPID BIOLOGICAL INDICATORS 50/BOX

Tuttnauer Ultra-Rapid BI's are to be used on the Bionova or MiniBio units. 20 minute read out, Can be used in gravity or pre/post vacuum autoclaves.

TUTWTL198-0072



## EZTEST® STERILIZATION MONITORING SYSTEM

Complete in office sterilization monitoring system - incubator, EZTest (25), English log book. Everything you need for effective monitoring of steriliation cycles in your office. Visual results in 24 hours.

SSP4145-SK/25B



## DELTA STEAM BIOLOGICAL INDICATOR

- EZTest is a self-contained biological indicator (BI) for monitoring sterilization
- Steam, ethylene oxide, propylene oxide, hydrogen peroxide
- Bacterial spores on a carrier material packaged within a small, thermoplastic culture tube
- Inside is a sealed-glass ampoule of specially formulated soybean casein digest culture medium containing a color indicator which turns a dramatic yellow when growth occurs
- Easy-to use in-office; no laboratory transfer needed
- Dramatic visual results within 24 hours (48 hours for ethylene oxide)

SSP4145-BV5-100 Box/100





## 3M COMPLY™ STERIGAGE™ CHEMICAL INTEGRATOR

Chemical integrators are used for pack control which validates that the sterilant has penetrated to the point of placement in the pack and confirms that sufficient exposure conditions have been met. 3M™ Comply™ Steam Chemical Integrators consist of a paper wick and steam and temperature sensitive chemical pellet contained in a paper/film and foil laminate envelope of laminate materials.

- Moving front ink technology provides easy to read Accept or Reject results
- Can be used in all 250-275°F / 121-135°C steam sterilization cycles
- Response correlates with a Biological Indicator at 3 time/temperature relationships under ideal steam sterilization conditions
- Class 5



3M 1243B ..... 100/Bag

## 3M COMPLY™ STEAM CHEMICAL INDICATOR STRIPS

The 3M™ Comply™ Indicator Strip enables pack control by validating that the sterilant has penetrated to the point of placement in the pack and confirming that sufficient exposure conditions have been met. Use in steam sterilization cycles.

3M 1250 ..... Class 4 Box/480



## STERILIZATION WRAP SHEETS (CSR)

MAXGEM2112.....	12"x12"	Cs/1000
MAXGEM2115.....	15"x15"	Cs/1000
MAXGEM2118.....	18"x18"	Cs/1000
MAXGEM2120.....	20"x20"	Cs/500
MAXGEM2124.....	24"x24"	Cs/500
MAXGEM2130.....	30"x30"	Cs/500



## STERILIZATION POUCHES, SELF SEALING

SPMSP22540 .....	2.25" x 4"	Box/200
GSM442100 .....	3.5" x 6.25"	Box/200
SPMSP35900.....	3.5" x 9"	Box/200
SPMSP52510.....	5.25" x 10"	Box/200
SPMSP71300.....	7.5" x 13"	Box/200
GSM222555 .....	12" x 15"	Box/200
SPMSP27590.....	2.75"x9"	Box/200



## INSTRUMENT AUTOCLAVABLE BRUSHES

SRG2872550-50 ..... Nylon  
SRG2972550-50 ..... Stainless Steel



## 3M COMPLY™ AUTOCLAVE TAPE

3M 1322-18MM..... 3/4" x 60 Yards  
3M 1322-24MM..... 1" x 60 Yards



## DISTILLED WATER

RW 050405 ..... Per Case of 4



## SPEED CLEAN® AUTOCLAVE CLEANER

Speed Clean is an autoclave cleaner solution that can be used in all brands of autoclaves, to clean and remove chlorides.

RIT002-0396-05..... 16oz



## CHAMBER BRITE® AUTOCLAVE CLEANER

Chamber Brite powder is simply the easiest way to keep autoclaves clean and free of water deposit build-up. Designed specifically for steam sterilizers, Chamber Brite effectively cleans and descales the reservoir, chamber and tubing in minutes.

TUTCB0010..... Box/10



## INSTRUMENT TRAYS, STAINLESS STEEL

AMG020-142 .....	10.5" W x 6.5" H x 2.5" D/L
AMG020-140 .....	Cover For 020-142
AMG020-162 .....	12.75" W x 7" H x 2.5" D/L
AMG020-160 .....	Cover For 020-162




## FLAT TRAY, STAINLESS STEEL

AMG020-527 .....	10" W x 6.5" H x 3/4" D/L
AMG020-525 .....	15" W x 10.5" H x 3/4" D/L



## STERILIZER POUCH RACK

RIT002-2108-00..... Rack




## INSTRUMENT PROTECTORS

Instrument Protectors will protect instrument tips and edges from damage during handling and sterilization.

### Variations Include:

- perforated protectors (holes) which reduce vacuum seal & moisture build-up
- steotome protectors
- double protectors
- retractor protectors
- validated for pre-vacuum sterilization
- available in a wide array of sizes for almost every application
- steam and gas permeable
- non sterile
- latex free

**Assorted Sizes and Colours Available. Call for details.**



## STERILIZATION MARKERS

- Approved to be used to write on self seal pouches and wraps

CML 13601 ..... Non-Sterile, withstands temperatures up to 500°F, Black



## THUMB DRESSING FORCEP

ALMM06-0120..... 5"



## TISSUE FORCEP

ALMM08-0120 ..... 5", 1 x 2 teeth  
 ALMM08-0140 ..... 5.5", 1 x 2 teeth  
 ALMM08-0160..... 6", 1 x 2 teeth




## ADSON TISSUE FORCEP

ALMM08-0420 ..... 4.75", 1 x 2 teeth



## ADSON DRESSING FORCEP

ALMM06-0520..... 4.75", Serrated



## SPLINTER FORCEP

ALMM10-0120 ..... 4"



## WILDE NASAL & EAR FORCEP

ALMM06-0710..... 6"  
 ALMM06-0720..... 7"



## HARTMANN ALLIGATOR FORCEP

ALMM06-0800..... 3.5"  
 ALMM06-0830..... 6.5"



## KELLY HEMOSTAT FORCEP

ALMM18-0320 ..... 5.5" Straight  
ALMM18-0340 ..... 5.5" Curved



## HALSTED MOSQUITO FORCEP

ALMM18-0200 ..... 5" Straight  
ALMM18-0220 ..... 5" Curved




## BOZEMAN UTERINE DRESSING FORCEPS

ALMM14-0300 ..... 10" Straight  
ALMM14-0320 ..... 10" Curved



## FOERSTER SPONGE FORCEP

ALMM20-0180 ..... 7" Straight Serrated  
ALMM20-0200 ..... 7" Curved Serrated  
ALMM20-0220 ..... 7" Straight Smooth  
ALMM20-0240 ..... 7" Curved Smooth  
  
ALMM20-0260 ..... 9.5" Straight Serrated  
ALMM20-0280 ..... 9.5" Curved Serrated  
ALMM20-0300 ..... 9.5" Straight Smooth  
ALMM20-0320 ..... 9.5" Curved Smooth



## SIMS UTERINE SOUND

ALMM14-0500 ..... 12.5"



## SCHROEDER TENACULUM FORCEP

ALMM14-0340 ..... 9.5" Straight  
ALMM14-0360 ..... 9.5" Curved






## OPERATING SCISSORS

ALMM04-0220 ..... 5.5" Straight, Blunt/Blunt  
 ALMM04-0240 ..... 5.5" Straight, Sharp/Blunt  
 ALMM04-0260 ..... 5.5" Straight, Sharp/Sharp



ALMM04-0500 ..... 5.5" Curved, Blunt/Blunt  
 ALMM04-0520 ..... 5.5" Curved, Sharp/Blunt  
 ALMM04-0540 ..... 5.5" Curved, Sharp/Sharp



## IRIS SCISSOR

ALMM04-1200 ..... 3.5" Straight  
 ALMM04-1220 ..... 3.5" Curved  
 ALMM04-1240 ..... 4.5" Straight  
 ALMM04-1260 ..... 4.5" Curved



## MAYO SCISSORS

ALMM04-0800 ..... 5.5" Straight  
 ALMM04-0840 ..... 6.75" Straight  
 ALMM04-0870 ..... 9.0" Straight

ALMM04-0820 ..... 5.5" Curved  
 ALMM04-0830 ..... 6.0" Curved  
 ALMM04-0872 ..... 9.0" Curved



## DOUBLE ENDED PROBE

ALMM34-0130 ..... 6.0"



## GROOVED DIRECTOR

ALMM34-0240 ..... 5.5"



## METZENBAUM SCISSORS

ALMM04-0910 ..... 5.5" Straight  
 ALMM04-0920 ..... 7.0" Straight  
 ALMM04-0960 ..... 9.0" Straight

ALMM04-0912 ..... 5.5" Curved  
 ALMM04-0940 ..... 7.0" Curved  
 ALMM04-0980 ..... 9.0" Curved



## SPENCER STITCH SCISSORS

ALMM04-1100 ..... 3.5"  
 ALMM04-1120 ..... 4.5"



## UNIVERSAL BANDAGE SCISSORS

Available in Black, Blue, Orange or Red

ALMM02-0500 ..... 6.0"



## LISTER BANDAGE SCISSORS

ALMM02-0140 ..... 5.5"  
 ALMM02-0160 ..... 7.0"



**CRILE WOOD NEEDLE HOLDER**

ALMM12-0420 ..... 6"

**DERF NEEDLE HOLDER**

ALMM12-0120 ..... 4.75"

**HALSEY NEEDLE HOLDER, SERRATED**

ALMM12-0160 ..... 5"

**MAYO-HEGAR NEEDLE HOLDER**

ALMM12-0320 ..... 5"

ALMM12-0340 ..... 5.5"

**OLSEN-HEGAR NEEDLE HOLDER**

ALMM12-0280 ..... 5.5"

ALMM12-0300 ..... 6.5"

ALMM12-0310 ..... 7.0"

**BAUMGARTNER NEEDLE HOLDER**

ALMM12-0260 ..... 5.5"

**NAIL NIPPER (WITH BARREL SPRING)**

ALMM29-0720 ..... 5.25"


**NAIL NIPPER (WITH DOUBLE SPRING)**

ALMM29-0744 ..... 5.5"



## TAYLOR PERCUSSION HAMMER

ALMM24-0140 ..... 8"



## QUEEN SQUARE HAMMER

ALMM24-0170 ..... 9"

ALMM24-0180 ..... 15.5"



## FINGER RING SAW

ALMM30-0920 ..... with Blade Protector

ALMM30-0940 ..... Blade Only



## CURETTES - STAINLESS STEEL - AUTOCLAVABLE

### Straight Buck Currettes

ALMM22-0300 ..... Size 00

ALMM22-0320 ..... Size 0

ALMM22-0340 ..... Size 1

ALMM22-0360 ..... Size 2

ALMM22-0380 ..... Size 3

ALMM22-0400 ..... Size 4

ALMM22-0420 ..... Size 5

### Angled Buck Currettes

ALMM22-0500 ..... Size 00

ALMM22-0510 ..... Size 0



ALMM22-0520 ..... Size 1

ALMM22-0530 ..... Size 2

ALMM22-0540 ..... Size 3

ALMM22-0550 ..... Size 4

ALMM22-0560 ..... Size 5



### Billeau Ear Currettes

ALMM22-0600 ..... Size 1

ALMM22-0620 ..... Size 1

ALMM22-0640 ..... Size 1 Each

## FOX DERMAL CURETTE - AUTOCLAVABLE

ALMM30-0300 ..... 3mm

ALMM30-0320 ..... 4mm

ALMM30-0340 ..... 5mm

ALMM30-0360 ..... 6mm



## OTOCLEAR® EAR WASH SYSTEM

Otoclear spray wash kit, manual, intermittent lavage

- Delivers greater volume compared to syringes or other manual devices
- Efficient design allows for two-handed procedure
- Built-in temperature gauge on bottle for water temperature

BNX7290..... Ear Wash System w/ ear basin

BNX7200 ..... Otoclear Irrigation Tips Box/40 - Disposable



## EAR SYRINGE, METAL

MILV919-385 ..... 4oz



## EAR BASIN, PLASTIC

GFI1280 ..... Autoclavable



## TUNING FORKS

MIL19-102..... C128

MIL19-104..... C256

MIL19-106..... C512



## ELEPHANT EAR WASH SYSTEM

Designed by a physician, the spraybottle provides enough water pressure to quickly and thoroughly clean ears, and holds enough solution to clean both ears without refilling. The unique hose design allows for easy positioning of the solution stream.

STA555-EW ..... Ear Washer with 3 Tips

STA555-TW..... Tips Only Pkg/20 - Disposable

STA555-SHIELD ..... Shield Only



# Have a question?

Call us — no robots, just people.



## STERILE SINGLE-USE INSTRUMENTS


Our single-use sterile instrument trays are designed to give you a good floor quality stainless steel instrument in a peel top plastic package that is sterile and ready for use.



### IRIS SCISSOR, STERILE

Stainless Steel instrument in sterile peel top plastic package


SUPSS.390.11..... 4.5" Straight Scissor  
SUPSS.393.11..... 4.5" Curved Scissor



### LITTAUER STITCH SCISSOR, STERILE

Stainless Steel instrument in sterile peel top plastic package


SUPSS.210.14 ..... 5.5" Straight Scissor



### O.R. SCISSORS, STERILE

Stainless Steel instrument in sterile peel top plastic package


SUPSS.030.14..... 5.5" Sharp/Sharp Scissor  
SUPSS.020.14..... 5.5" Sharp/Blunt Scissor



### METZENBAUM SCISSORS, STERILE

Stainless Steel instrument in sterile peel top plastic package


SUPSS.120.18..... 7" Straight Scissor  
SUPSS.121.18..... 7" Curved Scissor



### LISTER BANDAGE SCISSORS, STERILE

Stainless Steel instrument in sterile peel top plastic package


SUPPL.014.18 ..... 5.5" Scissor



## ADSON FORCEP, STERILE

Stainless Steel instrument in sterile peel top plastic package

SUPDT.318.12 ..... 4.75" Serrated  
 SUPDT.320.12 ..... 4.75" 1x2 Teeth



## HALSTEAD MOSQUITO FORCEP, STERILE

Stainless Steel instrument in sterile peel top plastic package

SUPAF.060.13 ..... 5" Straight  
 SUPAF.061.13 ..... 5" Curved



## KELLY FORCEP, STERILE

Stainless Steel instrument in sterile peel top plastic package


SUPAF.090.14 ..... 5.5" Straight  
 SUPDT.320.12 ..... 5.5" Curved



## BAUMGARTNER NEEDLE HOLDER, STERILE

Stainless Steel instrument in sterile peel top plastic package


SUPND.120.14 ..... 5.25"



## TENACULUM FORCEP, STERILE

Stainless Steel instrument in sterile peel top plastic package

GOJ5015-04 ..... 1.5 Liter Pump Bottle



## PENLIGHT

AAR66RN ..... Disposable Pkg/6



## LIGHTED EAR CURETTES WITH MAGNIFICATION

### Standard Box (50 per box & 1 Light Source & 1 Magnifying Lens)

BNX2210 .....Lighted AngleLoop® Box  
 BNX2220 .....Lighted MicroLoop® Box  
 BNX2230 .....Lighted WaveCurette™ Box  
 BNX2240 .....Lighted FlexLoop® Box  
 BNX2250 .....Lighted CeraSpoon® Box  
 BNX2260 .....Lighted InfantScoop® Box  
 BNX2270.....Lighted VersaLoop® Box  
 BNX2245 .....Variety Pack (Includes 7 tip styles) Box

### Clinic Packs (200 per box & 4 Light Sources & 4 Magnifying Lens)

BNX2209 ..... Variety Pack (Includes 7 tip styles) Box  
 BNX2200 ..... Light Source Only Box



## LIGHTED PLACEMENT TOOL

Gives Audiologists the ability to safely and accurately place foam blocks in the ear.

BNX2850 .....Lighted Placement Tool Box/15



## LIGHTED SUCTION FOR CERUMEN REMOVAL

This combination provides the ability to quickly remove cerumen or foreign bodies with the utmost ease and safety. Suction unit sold separately or works with any existing suction unit.

BNX2625 .....Lighted Suction Box/25



## LIGHTED FORCEPS FOR FOREIGN BODY REMOVAL

BNX2750 ..... Forceps with Light Source & Magnifying Lens Box/10  
 BNX2775 ..... Forceps Only Box/10



## LIGHTED ARTICULATING EAR CURETTE

The easy to maneuver tip articulates up to 90 degrees to assist in the removal of cerumen from the ear canal and foreign bodies from the ear and nose. The Lighted Articulating Ear Curette is the only lighted curette that articulates to assist in the removal of cerumen from the ear canal and works because of the unique hinged design.

BNX2511..... Lighted Articulating Curette Box/50



## MILTEX DISPOSABLE DERMAL CURETTES

MIL33-52 .....	2mm.....	Box/50
MIL33-53 .....	3mm.....	Box/50
MIL33-54 .....	4mm.....	Box/50
MIL33-55 .....	5mm.....	Box/50
MIL33-57 .....	7mm.....	Box/50



BNX9555.....	Safe Ear Currettes™ - White FlexLoop® 4mm .....	Box/50
BNX4888.....	Safe Ear Currettes™ - Blue InfantScoop® 2mm.....	Box/50
BNX1222.....	Safe Ear Currettes™ - Green MicroLoop® 3mm .....	Box/50
BNX6333.....	Safe Ear Currettes™ - Yellow CeraSpoon® 4mm .....	Box/50
BNX3444.....	Safe Ear Currettes™ - Red AngleLoop® 4mm.....	Box/50
BNX4111.....	Safe Ear Currettes™ - Purple VersaLoop® 3mm .....	Box/50
BNX2999.....	Safe Ear Currettes™ - Orange ControLoop® 4mm.....	Box/50
BNX5777.....	Safe Ear Currettes™ - Variety Pack (10/each style & 15/white) .....	Box/75



## INSTRUMENT CLEANING SUPPLIES

MIL3-740 .....	Surgical Instrument Stain Remover, 3oz
MIL3-720 .....	Surgical Instrument Cleaner, 8oz
SRG2872550-50 .....	Surgical Instrument Brush, Nylon
SRG2972550-50 .....	Surgical Instrument Brush, Stainless Steel



## METRILUBE® INSTRUMENT LUBRICANT

MetriLube is an oil-based instrument lubricant that coats metal instruments to prevent spotting, staining, and rusting of surgical instruments and helps keep hinges working freely.

MET-10-3425 .....	24oz
-------------------	------



**Metrex**  
INFECTION PREVENTION

## INSTRUMENT LUBRICANT MILK

Instrument Lubricants inhibit rust and prolong the life of surgical instruments. Weiman® lubricants contain dual-purpose corrosion inhibitors that protect all types of metal instruments, including carbon steel. Do not interfere with steam sterilizing.

WEIT2 .....	1 Gallon
-------------	----------





## PRI-MED SOLACE NITRILE EXAM GLOVES

- Chemo Drug Tested and Approved
- Micro Textured Fingertips
- Powder Free
- Blue
- 250 Per Box

PRIPM6-2905 ..... X-Small  
 PRIPM6-2906 ..... Small  
 PRIPM6-2907 ..... Medium  
 PRIPM6-2908 ..... Large  
 PRIPM6-2909 ..... X-Large (230 per Box)



## PRI-MED BLUEBERRY NITRILE EXAM GLOVES (ACCELERATOR FREE)

Specially formulated without the use of any chemical accelerators or sulfur, the PrimaTouch blueberry glove protects both the wearer and patient from Type I (natural rubber latex protein) and Type IV hypersensitivities (Chemical Accelerators: Thiurams, Thiazoles & Carbamates).

- Chemotherapy Drug Tested and Approved
- Textured Fingertips
- Powder-Free
- Berry Blue
- 250 per Box

PRIPM6-2805 ..... X-Small  
 PRIPM6-2806 ..... Small  
 PRIPM6-2807 ..... Medium  
 PRIPM6-2808 ..... Large  
 PRIPM6-2809 ..... X-Large



## PRI-MED SUSTAIN BIODEGRADABLE NITRILE EXAM GLOVES

PRIMED Sustain™ nitrile gloves are specially formulated to attract microbes found in landfills to enhance the enzymatic biodegradation of the material. Sustain™ achieves 90% biodegradability in 490 days versus 0.8% in 626 days for standard nitrile formulations.

- Chemo Drug Tested and Approved
- Textured Fingertips
- Purple
- 200 Per Box

PRIPM6-2025 ..... X-Small  
 PRIPM6-2026 ..... Small  
 PRIPM6-2027 ..... Medium  
 PRIPM6-2028 ..... Large  
 PRIPM6-202 ..... X-Large





## COPOLYMER GLOVE

AMG014-200 .....One Size Fits All Box/100



## GLOVE BOX HOLDER WIRE

BOWGL011-0613 .....Single Glove Wire  
 BOWGL022-0613 .....Double Glove  
 BOWGL033-0613 .....Triple Glove Wire  
 BoWGL044-0613.....Quad Glove Wire



## GLOVE BOX HOLDER PLASTIC

BOWGP-013 .....Single Glove Plastic  
 BOWGP-014 .....Double Glove Plastic  
 BOWGP-015 .....Triple Glove Plastic  
 BOWGP-016 .....Quad Glove Plastic



## TRIPLE GLOVE BOX HOLDER - WIRE

- Holds three boxes of standard examination gloves.
- Dimensions: 17.1 x 4.25 x 7.5 inch (43.7 x 10.8 x 19 cm) White Wire
- No hardware with product

AMG011-910 ..... Triple Glove Box Holder - Wire



# Working through a list?

Let us help you check it off.

## VIBRANT™ POWDER-FREE LATEX EXAM GLOVES

- Powder free
- Beaded cuff eliminates roll-up
- Micro Textured
- Non-Sterile
- 4.5ml thickness

SPM98225 .....X-Small  
 SPM98226 .....Small  
 SPM98227 .....Medium  
 SPM98228 .....Large  
 SPM98229 .....X-Large

**AURELIA**  
 Sign of Comfort™



## DELIGHT PF VINYL EXAM GLOVES

- Powder free
- Smooth finish
- Latex-free
- Hypoallergenic
- Protein-free
- Non-sterile
- Clear white

SPM38225 .....X-Small  
 SPM38226 .....Small  
 SPM38227 .....Medium  
 SPM38228 .....Large  
 SPM38229 .....X-Large



## BLUSH® NITRILE PINK EXAMINATION GLOVES - BOX OF 200

- Nitrile (100% Latex-Free)
- Non Sterile
- Glove Design & Feature: Ambidextrous and Finger Textured
- Cuff: Beaded
- Color: Blush Pink
- 2.5ml thickness
- Packaging: 200 pieces per box / 10 boxes per case

SPM78885 ..... X-Small  
 SPM78886 ..... Small  
 SPM78887 ..... Medium  
 SPM78888 ..... Large



## ABSOLUTE 100 BLACK NITRILE EXAM GLOVES

SPM9899A5 .....X-Small  
 SPM9899A6-S .....Small  
 SPM9899A7-M .....Medium  
 SPM9899A8-L .....Large  
 SPM9899A9-XL .....X-Large



## NITRIDERM® NITRILE STERILE EXAM GLOVES SERIES 108

Nitrile synthetic contains no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort. Pair packaged, inner wrapper in ozone resistant multi-vac packaging.

IHC108100.....	Small
IHC108200.....	Medium
IHC108300.....	Large
IHC108350.....	X-Large
50 Pairs/Box	



## PRESTIGE® DHD LATEX SURGICAL GLOVES

Damp don gel interior provides effortless donning. Curved finger design increases comfort. Powder-Free to reduce risk of powder induced allergies, adhesion, and granulomas.

IHC139550.....	Size 5.5
IHC139600.....	Size 6.0
IHC139650.....	Size 6.5
IHC139700.....	Size 7.0
IHC139750.....	Size 7.5
IHC139800.....	Size 8.0
IHC139850.....	Size 8.5
IHC139900.....	Size 9.0
50 Pairs/Box	



## NITRIDERM® NITRILE STERILE EXAM GLOVES SINGLE

Nitrile synthetic contains no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort. Pair packaged, inner wrapper in ozone resistant multi-vac packaging.

IHC109100.....	Small
IHC109200.....	Medium
IHC109300.....	Large
IHC109350.....	X-Large
100/Box	



## NITRIDERM® NITRILE STERILE EXAM GLOVE

NitriDerm Nitrile Surgical Gloves, These nitrile synthetic surgical gloves contain no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort, fit, and protection, while offering exceptional tactile sensitivity.

IHC135260.....	Size 6.0
IHC135265.....	Size 6.5
IHC135270.....	Size 7.0
IHC135275.....	Size 7.5
IHC135280.....	Size 8.0
IHC135285.....	Size 8.5
IHC135290.....	Size 9.0
50 Pairs/Box	



## PRESTIGE® CHLOROPRENE SURGICAL GLOVES

Made from non-latex, polychloroprene synthetic polymer, these new surgical gloves offer all of the benefits of natural rubber latex, with none of the allergen problems associated with latex. The unique formulation allows them to be worn for extended periods without hand fatigue. Slightly roughened finish provides the ideal combination of gripping ability and tactility.

IHC134550.....	Size 5.5
IHC134600.....	Size 6.0
IHC134650.....	Size 6.5
IHC134700.....	Size 7.0
IHC134750.....	Size 7.5
IHC134800.....	Size 8.0
IHC134850.....	Size 8.5
IHC134900.....	Size 9.0
25 Pairs/Box	



## DERMASSIST® LATEX SURGICAL GLOVES

IHC133550.....	Size 5.5
IHC133600.....	Size 6.0
IHC133650.....	Size 6.5
IHC133700.....	Size 7.0
IHC133750.....	Size 7.5
IHC133800.....	Size 8.0
IHC133850.....	Size 8.5
IHC133900.....	Size 9.0
50 Pairs/Box	



## ANSELL ENCORE® LATEX TEXTURED SURGICAL GLOVES

Latex, powder-free surgical gloves ensuring secure grip for confident instrument handling

- Chemo-tested
- Sterile
- 50 pairs per box

Size 5.5.....	ANS5785000
Size 6.0.....	ANS5785001
Size 6.5.....	ANS5785002
Size 7.0.....	ANS5785003
Size 7.5.....	ANS5785004
Size 8.0.....	ANS5785005
Size 8.5.....	ANS5785006
Size 9.0.....	ANS5785007



## MAGNA DISPOSABLE SCALPELS

These scalpels are individually packaged and are attached to a plastic handle with a plastic blade cover for protection.

ALMM92-10 .....	#10
ALMM92-11 .....	#11
ALMM92-12 .....	#12
ALMM92-15 .....	#15
ALMM92-20 .....	#20
ALMM92-21 .....	#21
ALMM92-22 .....	#22
ALMM92-23 .....	#23

Box/10



## MEDI-CUT® DISPOSABLE SAFETY SCALPELS

Blade retracts safely with one simple motion, and can be permanently locked and covered for safe disposal.

DYN4160.....	#10
DYN4161.....	#11
DYN4165.....	#15

Box/10



## MAGNA STAINLESS STEEL SURGICAL BLADES


ALMM90-10 .....	#10
ALMM90-11 .....	#11
ALMM90-12 .....	#12
ALMM90-15 .....	#15
ALMM90-20 .....	#20
ALMM90-21 .....	#21
ALMM90-22 .....	#22
ALMM90-23 .....	#23
ALMM90-24 .....	#24
ALMM90-25 .....	#25

Box/100



## SCALPEL HANDLE STAINLESS STEEL


ALMM36-0100.....	#3 (Fits 10, 11, 12, 15 blades)
ALMM36-0120.....	#4 (Fits 20 - 25 blades)





MILTEX® DISPOSABLE BIOPSY PUNCHES

MIL33-31AA	1.0mm
MIL33-31A	1.5mm
MIL33-31	2.0mm
MIL33-31B	2.5mm
MIL33-32	3.0mm
MIL33-33	3.5mm
MIL33-34	4.0mm
MIL33-35	5.0mm
MIL33-36	6.0mm
MIL33-37	7.0mm



DYNAREX DISPOSABLE BIOPSY PUNCHES

DYN4086	1.5mm
DYN4087	2.0mm
DYN4088	2.5mm
DYN4089	3.0mm
DYN4090	3.5mm
DYN4091	4.0mm
DYN4093	5.0mm
DYN4094	6.0mm
DYN4096	8.0mm
DYN4097	Assorted Includes: Sizes 2mm, 3mm, 4mm, 5mm, 6 mm (5 of each size)



BIOPBLADE

- Used for skin/tissue saucerization and shave biopsy, sterile, disposable
- Flexible design allows the blade to be positioned at various angles

MIL33-100	50 per box
-----------	------------



SUTUREOUT KIT

- SutureOut is a sterile, disposable product that is intended for removing multiple stitches from a single patient
- Kit includes SutureOut tool, plastic forcep and 1ea 4x4 gauze

CHS85-4337



## BRYMILL LIQUID NITROGEN GUN

Units come with four spray tips with apertures of various sizes (.04 - .016") plus a 20G straight spray and 20G bent spray for more accurate, controlled treatments.

CYGCERY-AC.....	500ml
CYGCERY-AC-3.....	350ml
CYGCERY-BABY .....	150ml "BABY"



## LIQUID NITROGEN TANKS AND ACCESSORIES

CYGTANK XC20 .....	20L Dewar Tank
CYGD12 .....	Dipper Cup for Liquid Nitrogen
CYGCERY GLOVES .....	Liquid Nitrogen Safety Gloves

## HISTOFREEZER® PORTABLE CRYOSURGICAL SYSTEM

- On the spot treatment for common warts, skin tags, benign lesions and plantar warts
- 2 Year shelf life
- Freezing Temp: -55°C
- Each kit includes 2 80ml cannister. 24ea 2mm foam tips for small lesions, 36ea 5mm foam tips for large lesions
- 40 treatments per kit
- Instruction manual



PAL2879100580

## HISTOFREEZER® STAT PORTABLE CRYOSURGICAL SYSTEM

Histofreezer STAT is a revolutionary advance in portable cryosurgery that delivers a colder

- 50 Treatments, includes: 1 Canister, 30 buds, 40 Cones (3,5,7,9 and 12mm), 4 cone adapters
- NEW handle design allows for single hand delivery of gas via a cone
- NEW cone design allows for deeper freezing with less gas used

CCL1500-1008



## CANTHARONE® WART TREATMENT SOLUTIONS

- For topical use of removal of warts

DOR9001-975M.....	Cantharone Wart Treatment 7.5ml
DOR9002-975M.....	Catharone Plus Wart Treatment 7.5ml





## LOOK™ SUTURES

### Nylon

LOO917B ..... 5-0 Nylon with C17 Needle ..... Box/12  
LOO921B ..... 4-0 Nylon with C17 Needle ..... Box/12  
LOO916B ..... 6-0 Nylon with C17 Needle ..... Box/12



### Reverse Cutting Needles

	Look	Ethicon
11mm	3/8 Circle C1	
12mm	C17	C-2
14mm	•	•
16mm	C22	FS-3
18mm	C6	FS-2
24mm	C7	FS-1
26mm	C8	FS
1/2 Circle	C31	X1

### Precision Reverse Cutting Needles

	Look	Ethicon
3mm	3/8 Circle •	G-6
6mm	•	P-6
11mm	C1	P-1
13mm	C3	P-3
16mm	PC34	PS-3
18mm	PC31	PS-2
24mm	PC33	PS-1

Canada Medical carries the complete line of Look® and Ethicon Sutures.  
Contact us for item codes and pricing.

# ETHICON

a Johnson & Johnson company

## ETHICON™ SUTURES

### Chromic Gut

ETH634G.....	5-0 Chromic Gut with FS-2 Needle.....	Box/12
ETH635H.....	4-0 Chromic Gut with FS-2 Needle.....	Box/36
ETH636H.....	3-0 Chromic Gut with FS-2 Needle.....	Box/36

### Plain Gut

ETHH820G.....	5-0 Plain Gut with FS-2 Needle .....	Box/12
ETHH821H.....	4-0 Plain Gut with FS-2 Needle .....	Box/36
ETHH822H.....	3-0 Plain Gut with FS-2 Needle .....	Box/36

### Ethilon

ETH668G.....	5-0 Nylon with C-2 Needle.....	Box/12
ETH661G.....	5-0 Nylon with FS-2 Needle.....	Box/12
ETH662G.....	4-0 Nylon with FS-2 Needle.....	Box/12

### Silk

ETH682G.....	5-0 Silk with FS-2 Needle.....	Box/12
ETH683G.....	4-0 Silk with FS-2 Needle.....	Box/12
ETH684G.....	3-0 Silk with FS-1 Needle.....	Box/12

### Prolene

ETH8661G.....	5-0 Prolene with FS-2 Needle.....	Box/12
ETH8683G.....	4-0 Prolene with FS-2 Needle.....	Box/12
ETH8665G.....	3-0 Prolene with FS-2 Needle.....	Box/12



### Reverse Cutting Needles

	Look	Ethicon	Syneture
11mm	3/8 Circle	C1	C-1
12mm		C17	
14mm		•	•
16mm		C22	FS-3
18mm		C6	FS-2
24mm		C7	FS-1
26mm		C8	FS
1/2 Circle	C31	X1	C23

### Precision Reverse Cutting Needles

	Look	Ethicon	Syneture
7mm	3/8 Circle	•	G-6
8mm		•	P-6
11mm		C1	P-1
13mm		C3	P-3
16mm		PC34	PS-3
18mm		PC31	PS-2
23mm		PC33	PS-1

Canada Medical carries the complete line of Ethicon™ Sutures.

Contact us for item codes and pricing.

## 3M STERI-STRIP™ SKIN CLOSURES

3M R1540 .....	1/8" x 3"	.....	5/Env Box/50 Env
3M R1541 .....	1/4" x 3"	.....	3/Env Box/50 Env
3M R1542 .....	1/4" x 1-1/2"	.....	6/Env Box/50 Env
3M R1547 .....	1/2" x 4"	.....	6/Env Box/50 Env



## PRI-MED PRIMAFLEX SKIN CLOSURES

PRIPM3-71831.....	1/8" x 3"
PRIPM3-71431 .....	1/4" x 3"
PRIPM3-71432 .....	1/4"x1.5"
PRIPM3-71241 .....	1/2"x4"
PRIPM3-71221 .....	1/2"x2"
PRIPM3-71151 .....	1"x5"



## GLUSTITCH TISSUE ADHESIVE

BLK200.....	Violet	.....	0.2ml	Box/12
BLK201 .....	Clear	.....	0.2ml	Box/12



## TWIN BLADE RAZOR

- Long-lasting, superior quality blade
- Clear plastic guard
- Disposable
- Unisex

DYN4250 Box/50



## PREP RAZOR

- Double Edge Premium Razor

PER75-4004 ..... Personna Prep Razor Box/100



## SAFETY PREP RAZOR

DYN4251 ..... Gallant Safety Prep Razor Box/50



## SURGICAL SKIN MARKERS WITH RULER, STERILE

SRGS-3000-R0 ..... Multi Tip Box/50





## PROXIMATE™ PLUS MD SKIN STAPLER

The PROXIMATE Plus MD skin stapler is designed to deliver rectangular, stainless steel staples for routine wound closure. It permits multi-directional release with an ergonomic design.

- Sterile, single-patient use
- 0.58 diameter


ETHPMW35 ..... 6 Per Box



## SUTURE REMOVAL TRAY, STERILE

Suture Removal Trays help you run your procedure more efficiently and help reduce the risk of cross-contamination.

DYN4521-IMC.....Suture Removal Tray



DYN4521-IMC

## SKIN STAPLE REMOVER, DISPOSABLE, STERILE

BUS716-IMC.....Skin Staple Remover



## SKIN STAPLE REMOVER WITH PLASTIC HANDLES STERILE, DISPOSABLE

BUS772



## BASIC DRESSING TRAY, STERILE

**This tray contains:**

- |                              |                            |                    |
|------------------------------|----------------------------|--------------------|
| • 1 Three Compartment Tray   | • 4 10x10 cm Gauze Sponges | • 1 Safety Pin     |
| • 1 Plastic Transfer Forceps | • 4 5x5 cm Gauze Sponges   | • 2 Towel/Underpad |
| • 2 Metal Forceps            | • 1 Waste Bag with Tie     |                    |

TRA85-2000TY



## IUD INSERTION KIT, STERILE

**Kit Includes:**

- |                                    |                       |
|------------------------------------|-----------------------|
| • 1 Scissors Metal, 7" Curved      | • 10 4x4" Gauze       |
| • 1 Tenaculum Metal, 9.5" Straight | • 30x30" Sterile Wrap |
| • 1 Uterine Sound Metal            |                       |

TRASMA30IUOPB



## IUD INSERTION KIT, STERILE, DISPOSABLE

**Kit Includes:**

- |                          |                         |                      |
|--------------------------|-------------------------|----------------------|
| • Tenaculum Forceps, 10" | • Littauer Scissors, 8" | • Uterine Sound, 13" |
|--------------------------|-------------------------|----------------------|

SUPK-161312



## LACERATION TRAYS, DISPOSABLE, STERILE

- |  |                                 |                                   |
|--|---------------------------------|-----------------------------------|
| • 1 Tray   | • 1 4-1/2" sharp/sharp scissors | • 5 4x4" gauze sponges            |
| • 1 5" floor-grade, Baumgartner (serrated jaw) needle holder | • 1 5" metal tissue forceps     | • One polylined sterile field     |
|  | • 2 absorbent towels            | • One polylined fenestrated drape |

BUS751




## XYLOCAINE®

APNVR89896CA1 .....Xylocaine 1% Plain 10ml 50/BX  
APNVR44571CA1 .....Xylocaine 1% Plain 20ml  
APNVR89892CA1 .....Xylocaine 1% Plain 2ml Poly Ampoule 50/BX  
APNVR44569CA1 .....Xylocaine 1% Plain 50ml  
APNVR89893CA1 .....Xylocaine 1% Plain 5ml Poly Ampoule 50/BX Box/50 vials  
APNVR44553CA1 .....Xylocaine 1% with Epinephrine 20ml  
APNVR89897CA1 .....Xylocaine 2% Plain 10ml 50/BX  
APNVR44567CA1 .....Xylocaine 2% Plain 20ml  
APNVR89894CA1 .....Xylocaine 2% Plain 2ml Poly Ampoule 50/BX  
APNVR44565CA1 .....Xylocaine 2% Plain 50ml  
APNVR89895CA1 .....Xylocaine 2% Plain 5ml 50/BX  
APNVR44551CA1 .....Xylocaine 2% with Epinephrine 20ml  
APNVR89868CA1 .....Xylocaine Endotracheal Spray 10mg/50ml  
APNVR44900CA1 .....Xylocaine Jelly 2% 30ml Tube  
APNVR44927CA1 .....Xylocaine Ointment 5% 35gm Tube  
APNVR44677CA1 .....Xylocard 100mg 5ml 10amp/box



## MARCAINE®

WIN02073020.....Marcaine 20ml, 0.25% Plain Single Use  
WIN02083020.....Marcaine 20ml, 0.5% with Epinephrine Single Use Box/10  
WIN02080020.....Marcaine 20ml, 0.25% with Epinephrine Multi-Dose Box/10  
WIN1587050 .....Marcaine 50ml, 0.25% Plain Multi-Dose  
WIN1610050 .....Marcaine 50ml, 0.5% Plain Multi-Dose



## LIDOCAINE

TEL0122AG01.....Lidocaine 1% Plain 20ml  
TEL0122AJ01.....Lidocaine 1% Plain 50ml  
TEL0126AF01 .....Lidocaine 2% Plain 10ml Tray/20  
TEL0127AG01.....Lidocaine 2% Plain 20ml  
TEL0127AJ01.....Lidocaine 2% Plain 50ml  
TEL0126AD01.....Lidocaine 2% Plain 5ml 20/bx Box/20  
TEL0128AG01.....Lidocaine 2% with Epinephrine 20ml



**Canada Medical also distributes  
pharmaceuticals including  
prescription drugs.  
Don't see what you need? Call us.**

## PAIN EASE® INSTANT TOPICAL ANESTHETIC

Gebauer's Pain Ease® instant topical anesthetic skin refrigerant (vapocoolant) is FDA cleared to temporarily control the pain associated with needle procedures and minor surgical procedures such as immunizations, allergy injections, foreign body removal, suturing and suture removal, IV starts and blood draws. Pain Ease can be used prior to, or in some cases in place of lidocaine injections.

SEK0008-02 .....3.5oz Mist Spray  
SEK0008-03 .....3.5oz Medium Spray



## EMLA® CREAM

- Topical Anesthetic

APNVR44812CA1 .....30g tube



## MINIMS® TETRACAINE HYDROCHLORIDE DROPS

BL 10666CA .....0.5% Pkg/20  
BL 10667CA .....1.0% Pkg/20



## MINIMS® FLOURESCIN SODIUM DROPS

BL 10692CA .....0.5ml Drop, 2% Pkg/20



## DISPOSABLE TOWEL DRAPES, STERILE

- For use in minor surgical or therapeutic procedures where a sterile field or drape is required
- Poly-interlined, 3 ply construction - tissue/poly/tissue - provides absorbency while preventing possible moisture or bacteria migration
- Eliminates all labor costs associated with cloth towels or drapes
- Available in both fenestrated and plain

DYN4409-IMC .....Fenestrated, 18"x26" Box/50  
DYN4410-IMC .....Plain, 18"x26" Box/50



## PROVIDINE IODINE SOLUTION 10%

Providine Iodine solution is an antiseptic cleanser for wound cleansing. Effective in destroying certain bacteria and removing impurities to provide antiseptic cleansing. Air dries quickly requiring no wiping. For external use only.

LER109-08 .....500ml  
LER109-09 .....115ml



## BAXEDIN® PRE-OP ANTISEPTIC SOLUTION

- 0.5% Chlorhexidine Gluconate (CHG) with 70% Isopropyl Alcohol
- For skin asepsis prior to surgery
- Disinfectant

OMGL0000003.....500ml Bottle, Tinted



## STANHEXIDINE® LOW SUDSING AQUEOUS SOLUTION

- 2% Chlorhexidine Gluconate (CHG) with 4% Isopropyl Alcohol
- Anti-bacterial for the skin
- Treatment for superficial skin infections

OMGL0000009 .....450ml Bottle



## BENZALKONIUM CHLORIDE WIPES

For use wherever antiseptic cleansing/refreshing is needed for sting free sensitive skin areas. Saturated with benzalkonium chloride 0.13%. Skin antiseptic.

LER126-01 .....Wipes..... Box/100



## PVP PREPS & SWABS

PDIB51200 .....Povidone Iodine 10% Prep Pad ..... Box/100  
PDIS42050.....Povidone Iodine 10% Swabstix.....Box/50



## NYLON SCRUB BRUSH

AMG118-140



## E-Z SCRUB SPONGE

Proprietary polyethylene soft bristles enhance mechanical action while reducing skin irritation. Soft sponge provides gentle cleaning on irritated areas or when firm mechanical action is not required.

BD 377479 .....4% Chlorhexidine.....Box/30



## **Baxter** NORMAL SALINE (0.9% SODIUM CHLORIDE) I.V. SOLUTIONS

BAXJB1322P .....	250ml
BAXJB1323 .....	500ml
BAXJB1324 .....	1000ml



## **Baxter** LACTATED RINGERS, I.V. SOLUTIONS

BAXJB2322P .....	250ml
BAXJB2323 .....	500ml
BAXJB2324 .....	1000ml



## **Baxter** DEXTROSE 5%, I.V. SOLUTIONS

BAXJB0062P .....	250ml
BAXJB0063 .....	500ml
BAXJB0064 .....	1000ml



## **Baxter** DEXTROSE 3.3% & 0.3% SALINE, I.V. SOLUTIONS

BAXJB1033 .....	500ml
BAXJB1034 .....	1000ml



## **B. BRAUN** SHARING EXPERTISE STERILE WATER, I.V. SOLUTIONS

TRABBL8502 .....	250ml
TRABBL8501-01 .....	500ml



## **Baxter** MINI BAG® I.V. SOLUTIONS

### Dextrose 5% Water Solution

BAXJB0081P .....	50ml per Bag, 4 bags/pack
BAXJB0082P .....	100ml Cs/96

### Sodium Chloride 0.9% Solution

BAXJB1301PK .....	50ml per Bag, 4 bags per pack
BAXJB1302PK .....	100ml per Bag, 4 bags per pack





## Baxter 4 WAY STOPCOCKS

Stopcock 4-way large bore, with rotating male luer lock adapter.

BAX2C6204 ..... Cs/50



## Baxter CLEARLINK® NEEDLE-FREE I.V. SETS

IV Needleless administration set with Clearlink® luer activated valve and male luer lock adapter. 76" Length, 10 Drops/ml

BAX2C8401 ..... Cs/48



IV Needleless Continu-Flo solution set with 2 Clearlink luer activated valve and male luer lock adapter. 104" Length, 10 Drops/ml

BAXJC8519 ..... Cs/48



## Baxter SOLUTION SET (10 DROPS / ML)


BAXJC8109 ..... with Luer Lock Adapter (101")



## Baxter I.V. EXTENSION SETS

IV Microbore catheter, 8" extension set with Interlink® injection site and luer lock adapter.

BAX2N3374 ..... Cs/200



## STERILE WATER FOR IRRIGATION

BAXJF7623P ..... 500ml Pour Bottle  
BAXJF7624 ..... 1000ml Pour Bottle



## STERILE WATER FOR INJECTION

**Bacteriostatic**  
WIN03977030 ..... 30ml vial

**No Preservative**  
OMGL0010005 ..... 10ml vial  
WIN04887020 ..... 20ml vial Cs/25



## NORMAL SALINE 0.9% (NaCl) FOR IRRIGATION

BAXJF7633P ..... 500ml Pour Bottle  
BAXJF7634 ..... 1000ml Pour Bottle



## NORMAL SALINE 0.9% (SODIUM CHLORIDE) FOR INJECTION

**No Preservative**  
WIN04888010 ..... 10ml vial  
WIN04888020 ..... 20ml vial Cs/25

**Bacteriostatic**  
WIN01966010 ..... 10ml vial Cs/25  
WIN01966020 ..... 20ml vial  
WIN01966030 ..... 30ml vial



## **Baxter** MINI BAG® SOLUTION SETS

### Dextrose 5% Water Solution

BAXJB0081P ..... 50ml per Bag, 4 bags/pack  
BAXJB0082P ..... 100ml Cs/96

### Sodium Chloride 0.9% Solution

BAXJB1301PK ..... 50ml per Bag, 4 bags per pack  
BAXJB1302PK ..... 100ml per Bag, 4 bags per pack



## SALJET® SINGLE DOSE STERILE SALINE 30ML

WLB60188151000000 ..... 20 Amps per Box



BD Insyte™ Autoguard™ BC Pro Shielded IV Catheter with Blood Control

Technology is designed to provide improved control with a tall push tab and three ridges for successful placement and advancement of the catheter the first time. Instaflash™ Needle Technology (included on all 18-24-gauge catheters) is designed to confirm immediate vessel entry at the point of insertion and improve first-stick success, reducing painful hit-and-miss insertions.

BD Insyte™ Autoguard™ BC Pro Shielded IV Catheter with Blood Control Technology

- 50 per Box

381012 .....	24 G x 0.75"
381023 .....	22 G x 1.00"
381024 .....	22 G x 1.25"
381033 .....	20 G x 1.00"
381034 .....	20 G x 1.16"
381037 .....	20 G x 1.88"
381044 .....	18 G x 1.16"
381047 .....	18 G x 1.88"
381054 .....	16 G x 1.16"
381057 .....	16 G x 1.77"

**With Wings:**

382912 .....	24 G x 0.75"
382923 .....	22 G x 1.00"
382924 .....	22 G x 1.25"
382933 .....	20 G x 1.00"
382934 .....	20 G x 1.16"
382937 .....	20 G x 1.88"
382944 .....	18 G x 1.16"
382947 .....	18 G x 1.88"
382954 .....	16 G x 1.16"
382957 .....	16 G x 1.77"



**Don't see what you're looking for?**

We'll help you track it down. Let's talk!

## BD INSYTE™ I.V. CATHETER

New Standard in IV therapy. When faced with the everyday challenges of infusion therapy, Insyte™ I.V. catheters deliver.

- Increases the ease of penetration
- Reduces the risks of infiltration
- Exhibits greater kink resistance
- Reduces the risks of phlebitis up to 50%
- Reduces IV restarts due to reduced complications
- Reduces the total IV therapy cost to use due to longer indwelling times

BD 381212 .....	24G x 0.75" .....	Yellow
BD 381223 .....	22G x 1.00" .....	Blue
BD 381233 .....	20G x 1.00" .....	Pink
BD 381234 .....	20G x 1.16" .....	Pink
BD 381237 .....	20G x 1.88" .....	Pink
BD 381244 .....	18G x 1.16" .....	Green
BD 381247 .....	18G x 1.88" .....	Green
BD 381254 .....	16G x 1.16" .....	Gray
BD 381257 .....	16G x 1.77" .....	Gray
BD 381267 .....	14G x 1.75" .....	Orange

Box/50



## BD NEXIVA™ DIFFUSICS™ CLOSED IV CANNULA SYSTEM WITH BD VIALON™ BIOMATERIAL FOR RADIOGRAPHIC POWER INJECTION

- 20 Per Box

383691.....	24 G X 0.75"
383692.....	22 G X 1.00"
383693.....	20 G X 1.00"
383694.....	20 G X 1.25"
383695.....	18 G X 1.25"

## BD VACUTAINER® SAFETY-LOK™ BLOOD COLLECTION SETS

The BD Vacutainer® Safety-Lok™ Blood Collection Set is simple, easy to use, and is safety-engineered. Single-handed safety activation with irreversible locking mechanism, confirmed with an audible click.

BD 367296.....	21G x 0.75", 12" tubing
BD 367297 .....	23G x 0.75", 12" tubing
BD 367298 .....	25G x 0.75", 12" tubing

Box/50



## BD NEXIVA™ CLOSED IV CANNULA SYSTEM WITH BD VIALON™ BIOMATERIAL, SINGLE PORT

383511.....	24 G x 0.75"
383510.....	24 G x 0.56"
383512.....	22 G x 1.00"
383518.....	20 G x 1.75"
383517.....	20 G x 1.25"
383516.....	20 G x 1.00"
383520.....	18 G x 1.75"
383519.....	18 G x 1.25"

Box/20




## TERUMO SURSHIELD® SAFETY WINGED INFUSION SET

Terumo's Surshield Safety Winged Needle Sets with One-Handed Activation provides a solution with a cost-effective, easy to use design. Press shield over the needle until it locks into place.

- Unique needle locking mechanism helps reduce the risk of needlestick injuries
- One handed activation for safety and ease of use
- No change in needle insertion technique required

TERSV.S19BL .....	19G x 3/4", 12" tubing
TERSV.S19BLS .....	19G x 3/4", 3.5" tubing
TERSV.S21BL .....	21G x 3/4", 12" tubing
TERSV.S21BLS .....	21G x 3/4", 3.5" tubing
TERSV.S23BL .....	23G x 3/4", 12" tubing
TERSV.S23BLS .....	23G x 3/4", 3.5" tubing
TERSV.S25BL .....	25G x 3/4", 12" tubing
TERSV.S25BLS .....	25G x 3/4", 3.5" tubing

Box/50



## SOL-CARE™ SAFETY BLOOD COLLECTION NEEDLE

The SOL-CARE™ Safety Blood Collection Needle is a safety-engineered device with a sliding safety shield that clicks when activated, telling you the needle is safely locked

With Luer Adapter, 50 Per Box:

SOM110201030030.....	21G x 3/4" 7" Tubing
SOM110201030031.....	21G x 3/4" 12" Tubing
SOM110201030032.....	23G x 3/4" 7" Tubing
SOM110201030033.....	23G x 3/4" 12" Tubing
SOM110201030034.....	25G x 3/4" 7" Tubing
SOM110201030035.....	25G x 3/4" 12" Tubing

with Pre-attached Holder, 25 Per Box:

SOM110201030004.....	21G x 3/4" 7" Tubing
SOM110201030002.....	21G x 3/4" 12" Tubing
SOM110201030008.....	23G x 3/4" 7" Tubing
SOM110201030006.....	23G x 3/4" 12" Tubing
SOM110201030012.....	25G x 3/4" 7" Tubing
SOM110201030010.....	25G x 3/4" 12" Tubing



## TOURNIQUETS, DISPOSABLE

Latex Free, 1x18"

PRIPM5-1018 Box/25



## 3M TEGADERM® TRANSPARENT FILM DRESSING

3M 1622W.....4.4cm x 4.4cm Box/100

3M 1624W.....6cm x 7cm Box/100



## 3M SoluPrep™ ANTISEPTIC SWABS

Chlorhexidine gluconate and 70% v/v alcohol

3M 10106 .....2% CHG & Alcohol Pad Box/200

3M 10202 .....0.5% CHG & Alcohol Swabstick Box/50

3M 10203 .....2% CHG & Alcohol Swabstick Box/50



## BD I.V. ACCESS CANNULA

The BD Plastic Cannula is a truly universal cannula and can be used with all types of split-septum IV injection sites and vials designed for needleless access. It provides simplified access to the Interlink injection site with smoother insertion, reduced penetration force and greater fluid dispersion.

BD 303345 .....Blunt Plastic Cannula Box/100



## BD PRN LUER LOCK IV ADAPTERS

Does not contain natural rubber latex. 0.1ml/min. Sterile.

BD 385111 .....Short 0.75" Injection site, sterile, latex free Box/50

BD 385112 .....Mid-Size 0.0625" Injection site, sterile, latex free Box/50



## BD Q-SYTE™ IV NEEDLELESS CONNECTOR, LUER ACCESS

The BD Q-Syte™ Luer Access Split Septum combines the convenience of a Luer access with the benefits of a split septum design to provide closed access for IV therapy. The benefits of this unique split septum design include reliable performance, high flow rates, and a visible fluid path.

BD 385100 .....Split Septum Box/50



## DELUXE I.V. POLE, 4 LEG

This IV has a four point base, twin wheeled swivel casters, and a circular "Ram's Horn" design. As well, the Economy 2 Hook IV Stand features a positive friction grip, chrome plated steel pole, stable base design and friction grip lock adjustment.

AMG775-751 .....IV Pole



## CLARITY AUTOMATIC URINE ANALYZER

- Large LCD Display
- Auto calibration and self-check
- Flagging of abnormal patient results
- Internal Memory of last 2000 Tests
- RS 232 port or USB for connectivity



CLACD-MCPROMO5 ..... Clarity Urocheck Urine Analyzerch

Includes: Analyzer, 4 Bottles of 10SG, and 1 Bottle of Microalbumin

CLACD-UCTL30 ..... Clarity Urocheck Urine Control 10ml 3/BXch

For use with Clarity Urine Strips and Analyzers only, REFRIGERATE UPON ARRIVAL.

Control solution must be used with each new box of strips

CLADTG-10SG..... Clarity Urocheck urine strips 10SG 100/Bt Bottle

CLADTG-UAPPR..... Clarity Analyzer paper rolls 4/PK Pack

## BIOSTRIP® URINE TEST STRIPS

URINE TEST STRIP	BLOOD	LEUKOCYTES	UROBILINOGEN	BILIRUBIN	NITRATE	SPECIFIC GRAVITY	GLUCOSE	PROTEIN	KETONES	pH
4192-2P										
4192-U7										
4192-U10										

SSP4192-2P ..... Urinalysis for protein and glucose Bottle/100

SSP4192-U7 ..... Urinalysis reagent strip, 7 parameters Bottle/100

SSP4192-U10 ..... Urinalysis reagent strip, 10 parameters Bottle/100



## SIEMENS® URINE TEST STRIPS

URINE TEST STRIP	BLOOD	LEUKOCYTES	UROBILINOGEN	BILIRUBIN	NITRATE	SPECIFIC GRAVITY	GLUCOSE	PROTEIN	KETONES	pH
Labstix										
Multistix 5										
Multistix 8 SG										
Multistix 10 SG										
Uristix										

SEM2300A ..... Multistix 10 SG Bottle/100

SEM2304A ..... Multistix 8 SG Bottle/100

SEM2309 ..... Multistix 5 Bottle/100

SEM2855U ..... Uristix Bottle/100



## CHEMSTRIP™ URINE TEST STRIPS

URINE TEST STRIP	SG	pH	Leukocytes	Nitrite	Protein	Glucose	Ketone	Urobilinogen	Bilirubin	Blood
Chemstrip 10										
Chemstrip 9										
Chemstrip 7										
Chemstrip 5										
Chemstrip GP										

BMC1120347911 ..... Chemstrip 10 Bottle/100

BMC11008552119 ..... Chemstrip 7 Bottle/100

BMC11893467119 ..... Chemstrip 5 Bottle/100

BMC1020074311 ..... Chemstrip GP Bottle/100



## BD VERITOR® STREP & FLU TEST KIT

### For Rapid Detection of Flu A+B & Strep

- Technology-verified Performance-Incorporates state of the art technologies that deliver results you can trust
- Multiple Sample Types-Validated to be used with nasopharyngeal washes/aspirates, and nasopharyngeal swab in transport
- Convenient Kit Configuration-Full customized kit to test and delivery results in minutes
- Objective Results-Easy to read digital display eliminates subjective interpretation



BD 256066.....	BD Veritor Plus System Analyzer
BD 256040 .....	Veritor Strep Test Kits Box/30
BD 256045 .....	Veritor Flu A+B Test Kits Box/30
BD 256088 .....	Veritor SARS, Cov-2 and Flu A+B Test Kits

## CLEARVIEW® INFECTIOUS MONONUCLEOSIS (MONO) RAPID TESTS

Detects Mononucleosis antibodies to aid in diagnosis of infectious mononucleosis. Whole venous blood, fingerprick blood, serum, or plasma.

- Cassette test
- Results in 5 minutes

SSP4010-CV506815N..... 20 Tests Per Box



## INNOVACON STREP A RAPID TEST

- Rapid detection of Group A Streptococcus bacterial antigen
- High sensitivity at low colony count
- Accurate results directly from throat swab
- Results in 5 minutes
- Easy to read control and test lines
- All materials included in the kit
- Available in strips and cassettes format

SSP4010-IST-501 ..... Strips 25 Per Box

SSP4010-IST-502 ..... Cassettes 20 Per Box



## BIOSTRIP® HCG URINE PREGNANCY TESTS

SSP4220-5 ..... 25 Strips Per Box



## QUICKSTICK® HCG URINE PREGNANCY TESTS

- Strips are individually Wrapped
- SSP4220-4K ..... 50 Strips Per Box



## QUICKSTEP® PLUS HCG URINE PREGNANCY TESTS

SSP4220-33 ..... 20 Cassettes Per Box



## PANBIO COVID-19 ANTIGEN SELF TEST

SSP4010-41FK701 ..... 10 Tests Per Box



## BINAXNOW™ INFLUENZA A&B CARD

- To aid in the rapid differential diagnosis of influenza A and B viral infections.
- Rapid detection of influenza A and B nucleoprotein antigens.
- Used for nasopharyngeal (NP) swab, nasal swab samples and nasal wash/aspirate specimens.
- Results in 15 minutes.

SSP4010-WL416-000 ..... 22 tests Per Box



## BINAXNOW™ RSV CARD

- To aid in the diagnosis of RSV infections in patients under the age of 5.
- Rapid detection of respiratory syncytial virus (RSV) fusion protein antigen.
- Used for nasopharyngeal (NP) swab and nasal wash specimens.
- Results in 15 minutes.

SSP4010-WL430-022 ..... 22 Tests Per Box



## BINAXNOW™ MALARIA

- To aid in the rapid diagnosis of human malaria infections.
- Rapid detection of Plasmodium antigens circulating in human venous and capillary EDTA whole blood.
- Targets the histidine-rich protein II (HRP II) antigen specific to Plasmodium falciparum (P.f.) and a pan-malarial antigen common to all four malaria species capable of infecting humans - P. falciparum, P. vivax (P.v.), P. ovale (P.o.), and P. malariae (P.m.).
- Results in 15 minutes.

SSP4010-WL660-000 ..... 22 Tests Per Box





## BIOSTRIP® LH OVULATION PREDICTION TEST

SSP4220-13.....Box/6 Tests



## DRUGS OF ABUSE TEST KIT

### Cups

SSP4215-DUA-1147A-125.....QuickCheck 14 Cup Urine Drug Screen Test 25 Test Cups per Box

### Panel

SSP4215-DUA-1144-125 .....QuickCheck 14 Panel Urine Drug Screen Tests 25 Test Panels Per Box

**Tests for:** Amphetamine 500 ng/mL, Barbiturates 300 ng/mL, Benzodiazepine 300 ng/mL, Buprenorphine 10 ng/mL, Cocaine 150 ng/mL, Methadone metabolite 100 ng/mL, Ethyl-  $\beta$ -D-Glucuronide 500 ng/mL, Fentanyl 20 ng/mL, Hydromorphone 300 ng/mL, Methamphetamine 500 ng/mL, Methadone 300 ng/mL, Opiate 2000 ng/mL, Oxycodone 100 ng/mL, Marijuana 50 ng/mL, Creatinine test



## AMNICATOR® AMNIOTIC FLUID INDICATOR SWABS

- A rapid and simple method of detecting ruptured amniotic membranes during pregnancy
- Highly sensitive indicator
- Non threatening to patients
- Each to read-color change visible in seconds

SSP4080-MWAM1-10.....Box/10



## ROVERS® CERVEX-BRUSH® COMBI WITH REMOVER TUBE

SSP4008-38010100RT .....Bag/25 Brushes



## CYTOLOGY BRUSH

SSP4201-CB-8B .....8" Non-Sterile..... Bag/100




## ENDOMETRIAL BIOPSY SAMPLER

Bioteque's 3.0mm Preferred Curette is used to obtain samples of mucosal lining for histological study in cancer screening, endometrial dating, or bacterial culturing.

MDG022721 .....MEDGYN 3mm Curette Sampler

BTQM0014 .....BIOTEQUE 3mm Curette Sampler



## MEDGYN ENDO SAMPLER STERILE KIT

MDG022720 .....3mm with 10cc Syringe





## URINE SPECIMEN CONTAINERS, STERILE

- Individually Wrapped

DYN4253-IMC..... 4oz Bag/100



## URINE COLLECTOR

- Fits securely on standard toilets

DYN4232



## CONDOMS

ANB5711 .....Non-Lubricated Box/144

ANB30070907058009-007 .....Lubricated Box/144



## BIOHAZARD SPECIMEN BAGS

- Outer Document Pouch
- WHMIS Standards Biohazard Symbol
- English and French Instructions Printed on Bag

SSP4180-SB69-5R 6" x 9" .....1000 Bags Per Box

SP4180-SB912-5R 9" x 12" .....500 Bags Per Box



## CYTOSPRAY FIXATIVE SPRAY

Ensures cells adhere to the cytology slide during staining and that cell morphology is maintained.

SSP4240-CYT-017..... 100ml



## DRIMOP® LIQUID ABSORBENT PACKETS

- Efficient containment of spills
- Use in specimen shipping containers and specimen bags
- Swiftly absorbs and solidifies up to 300 times its own weight of liquid
- 0.5g packet

SSP4090-02-01083CG17 ..... 1500 Per Box



## ZIP-IQ TT TEST TUBE CENTRIFUGE

The LW Scientific Zip-IQ TT centrifuge is an impressive 6-place test tube centrifuge for spinning blood, urine, and other fluids in 1.0ml to 15ml test tubes and microtubes, motor spins at 5,000rpm (2,500g) producing serum in as little as 5 minute.

LWS06AD-15T3



## BLOOD DRAWING CHAIR (PHLEBOTOMY CHAIR)

- Contoured padded seat and padded, arched backrest
- Padded flip arm and armrests
- Height adjustable armrest
- Comfortable, 2" padding throughout (5cm)
- Easy-clean, durable vinyl
- 1-1/4" square, heavy duty, all-welded, tubular steel frame
- Adjustable height and depth flip arm
- Flip arm positions in place when in use, then flips out of the way for easy exit
- Neutral Gray, powder-coated frame
- Adjustable foot levelers
- 400lb load capacity under normal use

SOLS573



## LORIS SAFETY LANCETS, STERILE

- Tri-bevel tip for virtually painless sampling
- Needle automatically retracts
- For capillary blood samples
- No lancing device required
- Pressure activated lancet ready to use
  - no loading needed
- Latex free
- 1.4mm puncture depth



LER201-01 .....Cs/100

## MEDIPOINT BLOOD LANCETS

- 36 Gauge
- Stainless Steel
- Single Use
- Tri Belvel Tip
- 200 Per Box

CML19570.....200 Per box




### Needle Only

BD 305757	30G x 1/2"
BD 305758	27G x 1/2"
BD 305767	25G x 1-1/2"
BD 305761	25G x 1"
BD 305759	25G x 5/8"
BD 305762	23G x 1"
BD 305763	22G x 1-1/2"
BD 305768	22G x 1"
BD 305765	21G x 1-1/2"
BD 305764	21G x 1"
BD 305766	18G x 1-1/2"

### Syringe & Needle Combo

BD 305778	1ml 30G x 1/2"	Box/50
BD 305789	1ml 27G x 1/2"	Box/50
BD 305780	1ml 25G x 5/8"	Box/50
BD 305787	3ml 25G x 1"	Box/50
BD 305781	3ml 25G x 5/8"	Box/50
BD 305782	3ml 23G x 1"	Box/50
BD 305783	3ml 22G x 1-1/2"	Box/50
BD 305779	3ml 21G x 1"	Box/50
BD 305785	5ml 22G x 1-1/2"	Box/50
BD 305786	10ml 22G x 1-1/2"	Box/50

### Safety Engineered Needles



### BD SYRINGES ONLY

BD 309659	1ml Slip tip	Box/200
BD 309628	1ml Luer Lock	Box/100
BD 309657	3ml Luer Lock	Box/200
BD 309646	5ml Luer Lock	Box/125
BD 302995	10ml Luer Lock	Box/200
BD 302830	20ml Luer Lock	Box/48
BD 302832	30ml Luer Lock	Box/58
BD 309653	60ml Luer Lock	Box/40
BD 305819	Sterile Caps	Box/200



## **PRECISION GLIDE™ NEEDLES ONLY**

BD 305128 .....	30G x 1" .....	Box/100
BD 305106 .....	30G x 1/2" .....	Box/100
BD 305136 .....	27G x 1-1/4" .....	Box/100
BD 305109 .....	27G x 1/2" .....	Box/100
BD 305111 .....	26G x 1/2" .....	Box/100
BD 305122 .....	25G x 5/8" .....	Box/100
BD 305125 .....	25G x 1" .....	Box/100
BD 305127 .....	25G x 1-1/2" .....	Box/100
BD 305145 .....	23G x 1" .....	Box/100
BD 305155 .....	22G x 1" .....	Box/100
BD 305156 .....	22G x 1-1/2" .....	Box/100
BD 305167 .....	21G x 1-1/2" .....	Box/100
BD 305195 .....	18G x 1" .....	Box/100
BD 305196 .....	18G x 1-1/2" .....	Box/100
BD 305198 .....	16G x 1-1/2" .....	Box/100



## **ALLERGY TRAYS**

BD 305535 .....	1/2cc 27G x 1/2"
BD 305540 .....	1cc 27G x 1/2"
BD 305541 .....	1cc 27G x 3/8"



## **SYRINGE & NEEDLE COMBINATIONS**

BD 309623 .....	1ml 27G x 1/2" .....	Box/100
BD 309625 .....	1ml 26G x 3/8" .....	Box/100
BD 309626 .....	1ml 26G x 3/8" .....	Box/100
BD 309570 .....	3ml 25G x 5/8" .....	Box/100
BD 309571 .....	3ml 23G x 1" .....	Box/100
BD 309572 .....	3ml 22G x 1" .....	Box/100
BD 309574 .....	3ml 22G x 1-1/2" ....	Box/100
BD 309575 .....	3ml 21G x 1" .....	Box/100
BD 309577 .....	3ml 21G x 1-1/2" ....	Box/100
BD 309578 .....	3ml 20G x 1" .....	Box/100
BD 309579 .....	3ml 20G x 1-1/2" ....	Box/100
BD 309581 .....	3ml 25G x 1" .....	Box/100
BD 309582 .....	3ml 25G x 1-1/2" ....	Box/100



## SOL-CARE™ SAFETY NEEDLE

- Easy, one-handed technique of activating the protective arm
- Audible click indicates safety mechanism activation
- Patented locking mechanism ensures needle is secure and completely locked in protective arm
- Two options for activation of the protective arm
- Compatible with all standard luer lock and slip tip syringes
- 100/Box



SN3012.....	SOL-CARE™ Safety Needle 30G x .5"	SN2215.....	SOL-CARE™ Safety Needle 22G x 1.5"
SN2712.....	SOL-CARE™ Safety Needle 27G x .5"	SN2210.....	SOL-CARE™ Safety Needle 22G x 1"
SN2515.....	SOL-CARE™ Safety Needle 25G x 1.5"	SN2115.....	SOL-CARE™ Safety Needle 21G x 1.5"
SN2510.....	SOL-CARE™ Safety Needle 25G x 1"	SN2110.....	SOL-CARE™ Safety Needle 21G x 1"
SN2558.....	SOL-CARE™ Safety Needle 25G x 5/8"	SN2015.....	SOL-CARE™ Safety Needle 20G x 1.5"
SN2315.....	SOL-CARE™ Safety Needle 23G x 1.5"	SN2010.....	SOL-CARE™ Safety Needle 20G x 1"
SN2310.....	SOL-CARE™ Safety Needle 23G x 1"	SN1815.....	SOL-CARE™ Safety Needle 18G x 1.5"

## SOL-CARE™ SAFETY NEEDLE WITH SYRINGE

12558SN .....	SOL-CARE 1ml Luer Lock Syringe w/Safety Needle 25G x 5/8"
12705SN .....	SOL-CARE 1ml Luer Lock Syringe w/Safety Needle 27G x .5"
13005SN .....	SOL-CARE 1ml Luer Lock Syringe w/Safety Needle 30G x .5"
32110SN .....	SOL-CARE 3ml Luer Lock Syringe w/Safety Needle 21G x 1"
32115SN .....	SOL-CARE 3ml Luer Lock Syringe w/Safety Needle 21G x 1.5"
32210SN .....	SOL-CARE 3ml Luer Lock Syringe w/Safety Needle 22G x 1"
32215SN .....	SOL-CARE 3ml Luer Lock Syringe w/Safety Needle 22G x 1.5"
32310SN .....	SOL-CARE 3ml Luer Lock Syringe w/Safety Needle 23G x 1"
32510SN .....	SOL-CARE 3ml Luer Lock Syringe w/Safety Needle 25G x 1"
32558SN .....	SOL-CARE 3ml Luer Lock Syringe w/Safety Needle 25G x 5/8"

## SOL-M LUER LOCK SYRINGE WITHOUT NEEDLE

SOM180001PP .....	SOL-CARE 1cc Luer Lock Syringe .....	Box/100
SOM180003 .....	SOL-CARE 3cc Luer Lock Syringe .....	Box/100
SOM180005 .....	SOL-CARE 5cc Luer Lock Syringe .....	Box/100
SOM180010 .....	SOL-CARE 10cc Luer Lock Syringe.....	Box/100
SOM180020 .....	SOL-CARE 20cc Luer Lock Syringe.....	Box/100
SOM180030 .....	SOL-CARE 30cc Luer Lock Syringe.....	Box/50
SOM180060 .....	SOL-CARE 60cc Luer Lock Syringe.....	Box/30

## SOL-M LUER SLIP SYRINGE WITHOUT NEEDLE

SOM180011.....	SOL-CARE 1cc slip tip syringe 100/bx
----------------	--------------------------------------

## SOL-M NEEDLE ONLY, 100/BOX

SOM11810 .....	18G x 1"	SOM12310 .....	23G x 1"
SOM11815 .....	18G x 1-1/2"	SOM12315 .....	23G x 1-1/2"
SOM11910 .....	19G x 1"	SOM12506 .....	25G x 5/8"
SOM11915 .....	19G x 1-1/2"	SOM12510 .....	25G x 1"
SOM12010 .....	20G x 1"	SOM12515 .....	25G x 1-1/2"
SOM12015 .....	20G x 1-1/2"	SOM12612 .....	26G x 1/2"
SOM12110 .....	21G x 1"	SOM12638 .....	26G x 3/8"
SOM12115 .....	21G x 1-1/2"	SOM12658 .....	26G x 5/8"
SOM12210 .....	22G x 1"	SOM12705 .....	27G x 1/2"
SOM12215 .....	22G x 1-1/2"	SOM13005 .....	30G x 1/2"

## SOL-M BLUNT FILL NEEDLE

SOM110022.....	18g x 1-1/2
----------------	-------------



## **TERUMO** SYRINGE ONLY

TERSS-01T .....	1cc TB Syringe only, Slip tip	Box/100
TERSS-03L .....	3cc TB Syringe only, Luer Lock tip	Box/100
TERSS-05L .....	5cc TB Syringe only, Luer Lock tip	Box/100
TERSS-10L .....	10cc TB Syringe only, Luer Lock tip	Box/100
TERSS-20L2 .....	20cc TB Syringe only, Luer Lock tip	Box/50
TERSS-30L .....	30cc TB Syringe only, Luer Lock tip	Box/25
TERSS-60L .....	60cc TB Syringe only, Luer Lock tip	Box/25



## **TERUMO** NEEDLE ONLY

TERNN-1825R .....	18G x 1"
TERNN-1838R .....	18G x 1-1/2"
TERNN-2025R .....	20G x 1"
TERNN-2038R .....	20G x 1-1/2"
TERNN-2125R .....	21G x 1"
TERNN-2138R .....	21G x 1-1/2"
TERNN-2225R .....	22G x 1"
TERNN-2238R .....	22G x 1-1/2"
TERNN-2325R .....	23G x 1"
TERNN-2516R .....	25G x 5/8"
TERNN-2525R .....	25G x 1"
TERNN-2713R .....	27G x 1/2"

100 per box



## **TERUMO** SYRINGE WITH NEEDLE

TERSS-01T2516 .....	1cc TB Syringe w/25G x 5/8" needle
TERSS-01T2713 .....	1cc TB Syringe w/27G x 1/2" needle
TERSS-03L2025 .....	3cc Syringe w/20G x 1" needle
TERSS-03L2038 .....	3cc Syringe w/20G x 1-1/2" needle
TERSS-03L2125 .....	3cc Syringe w/21G x 1" needle
TERSS-03L2138 .....	3cc Syringe w/21G x 1-1/2" needle
TERSS-03L2238 .....	3cc Syringe w/22G x 1-1/2" needle
TERSS-03L2325 .....	3cc Syringe w/23G x 1" needle
TERSS-03L2525 .....	3cc Syringe w/25G x 1" needle

100 per box



## BD BLUNT FILL NEEDLE

- BD 305180 ..... Blunt fill needle only, 18G x 1-1/2"  
 BD 305211 ..... Blunt fill filter needle only, 18G x 1-1/2"

100 per box



## BD VACUTAINER® PRODUCTS

BD provides a full offering of BD Vacutainer® Blood Collection Tubes in both glass and plastic. BD Vacutainer® Plus Plastic Tubes offer a safe method of blood collection. They also reduce the potential for tube breakage and specimen spoilage, thereby reducing the potential for exposure to bloodborne pathogens.



100 per box



- BD 367290 ..... Vacutainer Sterile multi sample luer adapter, latex free  
 BD 367899 ..... Vacutainer Plus collection plastic tube, 13mm x 100mm, 6ml, Pink closure cap, sterile  
 BD 367925 ..... Vacutainer Plus collection plastic tube, 13mm x 100mm, 6ml, Grey closure cap, sterile  
 BD 367815 ..... Vacutainer Plus collection plastic tube, 13mm x 100mm, 6ml, Red closure cap, sterile  
 BD 367861 ..... Vacutainer Plus collection plastic tube, 13mm x 75mm, 4ml, Lavender closure cap, sterile  
 BD 367986 ..... Vacutainer Plus collection plastic tube, 13mm x 100mm, 5ml, Gold closure cap, sterile  
 BD 363083 ..... Vacutainer Plus Citrate plastic tube, 13mm x 7mm, 2.7ml, Light Blue closure cap, sterile  
 BD 363080 ..... Vacutainer Plus Citrate plastic tube, 13mm x 7mm, 1.8ml, Light Blue closure cap, sterile

## BD POSIFLUSH™ EXTERNALLY STERILE (XS) SALINE FLUSH SYRINGE

PosiFlush™ XS Syringes, filled with Saline, are designed for flushing vascular access devices within a sterile field, or where patients are particularly susceptible to infection.

306571 ..... 5mL 30 Per Box  
306572 ..... 10mL 30 Per Box



## BD SMARTSITE™ 20 MM NON-VENTED VIAL ACCESS DEVICE

Standalone, single-vial use, disposable devices which permit access to a medication vial without the use of a needle during medication preparation.

BD 2205E ..... Box/25



## BD INSULIN SYRINGES

- Ease of penetrating the skin
- Ease of dispensing the total dosage
- Ease of measuring to the correct dose line
- Ease of reading the scale on the syringe
- Consistent quality from syringe to syringe
- Comfort of injection into the skin
- Smooth action of the plunger when drawing up insulin



BD 329420 ..... Micro-fine IV 1ml syringe with 28G x 1/2" regular bevel sterile needle Box/100  
BD 329424 ..... Micro-fine 1ml syringe with 28G x 1/2" regular bevel sterile needle Box/100  
BD 329461 ..... Micro-fine 0.5ml Lo-Dose syringe with 28G x 1/2" sterile needle Box/100  
BD 320440 ..... Micro-fine 0.3ml syringe with 31G x 8mm thin wall sterile needle Box/100

## SURE COMFORT INSULIN SYRINGE

BB 22-6503 ..... 3/10CC Syringe with 31G x 5/16 Box/100  
BB 22-6005 ..... 1/2cc syringe with 30G x 5/16 Box/100



## BD SHARPS CONTAINERS

BD 300460 ..... 1.4L Sharps Container  
 BD 305970 ..... Non Locking Tray Style Bracket  
 BD 300450 ..... 3L Sharps Container  
 BD 300974 ..... 5.1L Sharps Container  
 BD 305447 ..... Locking Wall Bracket for 5.1L  
 BD 300440 ..... 7.6L Sharps Container  
 BD 305485 ..... Locking Wire Bracket  
 BD 300452 ..... 10.3L Sharps Container



## COVIDIEN PHLEBOTOMY SHARPS CONTAINER

Portable phlebotomy containers that fit into present blood drawing trays for the efficient disposal of blood, hypodermic and IV needles. These containers facilitate temporary and final closures in a variety of sizes and colors.

COV8906..... 1 Quart, Yellow



## COVIDIEN MULTI-PURPOSE SHARPS CONTAINER

Covidien offers a complete line of brackets and holders to ensure your safety and compliance with all OSHA regulations.

COV8970 ..... 7.6L with clear rotor opening, red  
 COV8975 ..... Wall bracket for 7.6L container  
 COV8963..... Locking wall bracket for 7.6L container



## COVIDIEN SHARPS CONTAINER, HORIZONTAL DOOR

Covidien offers a complete line of brackets and holders to ensure your safety and compliance with all OSHA regulations.

COV8967..... 7.6L Container w/clear horizontal drop lid, Yellow



## COVIDIEN SHARPS CONTAINER, SHARPSTAR™ LID

Covidien offers a complete line of brackets and holders to ensure your safety and compliance with all OSHA regulations.

COV8506Y ..... 5L, Yellow



## TABLE PAPER

BAN170912IMC ..... Smooth, 18" x 225' Cs/12  
 BAN170914IMC ..... Smooth, 21" x 225' Cs/12  
 BAN171002IMC ..... Crepe, 18" x 125' Cs/12  
 BAN171004IMC ..... Crepe, 21" x 125' Cs/12



## EXAM DRAPES, DISPOSABLE

BAN1710824IMC..... 40" x 48" Cs/100



## EXAM GOWNS, DISPOSABLE

BAN1710847IMC..... 30" x 42" Cs/50  
 BAN910521..... 32" x 42" (Extra Large, Blue) Cs/50



## GOWN (HALF WITH TIES), DISPOSABLE

BAN1710857IMC..... 30" x 21" Cs/100



## MEDICAL TOWELS, 3-PLY TISSUE, 13" x 18"

BAN918101..... White Cs/500





## PILLOWCASE

BAN179365IMC...21" x 30", Disposable, White Cs/100  
GFI3857.....21" x 27", Reusable, White, Vinyl



## PILLOWS (FLUID PROOF)

These pillows feature a soft, vinyl ticking that is inherently flame-retardant, antistatic and fluid proof. Product has antimicrobial properties built in to protect the product. The product does not protect users or others against bacteria, viruses or germs.

MAXMDT219684 ..... 18" x 24"



## EXAM SHORTS, DISPOSABLE

Disposable exam shorts have multi-layer design for comfort and fluid resistancy.

GRM10000 ..... Small/Medium, Blue Cs/50  
GRM10001 ..... Large/X-Large, Blue Cs/50



## APRONS, DISPOSABLE

AMG018-400 ..... Plastic, White Box/100



## UNDERPAD, DISPOSABLE, FLUFF FILLED

DYN1340 ..... 17" x 24" Bag/100



## SANITARY PADS

CFSMT-4..... Individually Wrapped Box/250



## WelchAllyn® KLEENSPEC® CORDLESS VAGINAL ILLUMINATION SYSTEM COMPLETE WITH CHARGER

- Superior Visualization - advanced LED light source provides better visualization of the exam area by supplying brighter, whiter light, right where it is needed
- Cordless and Rechargeable - no cord to get in the way or break during procedures, continuous on time of 80 minutes, recharge time of six hours, lithium ion battery technology
- Automatic on/off - no more hassle with turning the unit on and off

WEL80000.....Illuminator Only  
WEL80010.....Complete Illumination System



## WelchAllyn® KLEENSPEC® CORDED LIGHTING SYSTEMS

- Unique integrated light source transmits superior illumination for a more accurate diagnosis
- Eliminates positioning of an external light source
- Integrated on/off switch close to the task area for operational convenience
- New ergonomic design allows for easy cleaning
- New lamp design with improved lamp life

WEL78900.....Speculum Lighting System



## WelchAllyn® DISPOSABLE VAGINAL SPECULUM

Designed to eliminate the time and expense of reprocessing metal specula and to reduce the risk of cross-contamination. Single-use speculum is manufactured to deliver durability and reliability. Use with the KleanSpec Illumination System or alone with an external light source.

WEL-590XS.....KleanSpec 590 Series, X-Small..... Box/24  
WEL-59000.....KleanSpec 590 Series, Small..... Box/24  
WEL-59001.....KleanSpec 590 Series, Medium ..... Box/24  
WEL-59004.....KleanSpec 590 Series, Large ..... Box/18



## DISPOSABLE VAGINAL SPECULUM

DYN4911 .....Small ..... Bag/10  
DYN4912 .....Medium ..... Bag/10  
DYN4913 .....Large ..... Bag/10



## CLEARSPEC SINGLE-USE VAGINAL SPECULUM AND LED ILLUMINATOR

ClearSpec single use speculums save time, reduce the risk of cross contamination and eliminate the expense of reprocessing metal speculums. Clearspec is recommended for all routine gynecological exams and procedures. Additionally, every box of ClearSpec comes with a light source in the box which is good for up to 100 exams.

- 1 re-usable LED illuminator Per Box

CSC102NS..... Small..... 24 Per Box  
CSC103NS..... Medium ..... 24 Per Box  
CSC104NS..... Large..... 18 Per Box



## PEDERSON STAINLESS STEEL VAGINAL SPECULUM

ALMM14-0180 .....Small  
ALMM14-0200 .....Medium  
ALMM14-0220 .....Large



## GRAVES STAINLESS STEEL VAGINAL SPECULUM

ALMM14-0100 .....Virgin  
ALMM14-0120 .....Small  
ALMM14-0140 .....Medium  
ALMM14-0160.....Large



## ANOSCOPES, DISPOSABLE

ORI42301-1C..... SoloSpec Small .....25 Per Box  
ORI42301-2C..... SoloSpec Medium .....25 Per Box  
ORI42301-3C..... SoloSpec Large .....25 Per Box



## SIGMOIDOSCOPES, DISPOSABLE

WEL53130..... KleenSpec Disposable Sigmoidoscope .....Box/25  
(individually wrapped)



## UTERINE SOUND, DISPOSABLE

MIL30-6000..... Uterine Sound Sterile 9 3/4" ..... 25 Per Box  
MDG020250 ..... Uterine Sound Sterile 10" ..... 25 Per Box



## Miltex UTERINE SOUND, STAINLESS STEEL

MIL30-650SS.....13" Stainless Steel



## IUD EXTRACTOR/HOOK

MIL30-940.....10.25"



## IUD THREAD RETRIEVER DISPOSABLE

SMBRETRIEVER.....Sterile, Single Use Box/5



## CERVICAL SCRAPER DUE END

SSP51-200.....7" Ayres Cervical Scraper, Non-Sterile Box/200



## HEGAR UTERINE DILATOR

ALMM14-0610.....1 - 2mm, Double Ended

ALMM14-0620.....3 - 4mm, Double Ended

ALMM14-0630.....5 - 6mm, Double Ended



## IUD INSERTION KIT, STERILE

### Kit Includes:

- 1 Scissors Metal, 7" Curved
- 1 Tenaculum Metal, 9.5" Straight
- 1 Uterine Sound Metal
- 10 4x4" Gauze
- 30x30" Sterile Wrap

TRASMA30IUOPB




## CERVICAL DILATOR, SINGLE USE

Cervical Dilators are especially useful in locating and dilating the external OS, the cervical canal and the internal OS.

- Length: 215mm
- Taper: 23mm (1.75mm Tip to 3.8mm OD)
- Disposable
- Used on post-menopausal patients with cervical stenosis

BTQOS265BA-1



## BZK TOWELETTES

0.13% BZK Wipe, is saturated with 0.13% Benzalkonium Chloride and is used as an antiseptic topical skin cleanser.

- Will not sting
- Ideal for perineal and maternity care

CCPBZK100 .....Individually Packaged Box/100



## SILVER NITRATE APPLICATORS

AMG118-395 .....6", 75% Silver and 25% Nitrate Pack/100



## RAYON TOP PROCTOLOGY SWABS

DYN4321 .....8" Proctology Swab, Rayon Tipped Box/100

AMD57500 .....16" Proctology Swab, Rayon Tipped Box/50



## COTTON TIPPED APPLICATORS

SSP4202-40-100 .....3" Non-Sterile Pkg/100

SSP4202-40-200 .....6" Non-Sterile Pkg/100

DYN4305 .....6" Sterile (2/Pkg) Box/100 Pkg



## TONGUE DEPRESSORS

AMD50000 .....Junior Non-Sterile Box/500

SSP4202-20-500 .....Senior Non-Sterile Box/500

AMG118-421 .....Senior Sterile Box/100





## AQUAGEL® LUBRICATING JELLY

Convenient flip cap lid. Non-runny, viscous gel. Bacteriostatic. Clear, water-soluble, non-irritating lubricant. Compatible with natural and synthetic substances.

PAR57-05 .....5oz



## DYNALUBE LUBRICATING JELLY - STERILE

Water soluble and latex free lube is intended for medical purposes to lubricate body orifices and facilitate entry of diagnostic or therapeutic devices. Comes in 5g individual packets.

DYN1251 .....5g, Foil Packs Box/72



## SURGILUBE® STERILE SURGICAL LUBRICANT

- Carbomer-free formulation reduces interference with pap test results
- Health Canada recommended for use with Hologic ThinPrep Pap Test

HR 281020543 .....Foil Pack 144 Per Box

HR 281020537 .....Flip-Top Tubes



## SUNDRY JARS

DUK4010 .....Glass Sundry Jars, Labelled, Set of 5 Set

DUK4011 .....Plastic Sundry Jars, Labelled, Set of 5 Set

DUK4020 .....Wall Mount 5 Jar Rack

DUK4023 .....Wall Mount 3 Jar Rack



## RECTANGULAR WASH BASIN, 7L, NOT AUTOCLAVABLE

DYN4227 .....Plastic, Mauve



## RECTANGULAR WASH BASIN, 7L, AUTOCLAVABLE

AMG025-184 .....Plastic, Blue



## SPONGE BASIN, STAINLESS STEEL

AMG020-565 .....7.5" x 3"



## MEDICINE CUP, GRADUATED, STAINLESS STEEL

AMG90-6200 .....2oz



## IODINE CUP, STAINLESS STEEL

ALM90-6000 .....2oz (60cc)

ALM90-6005 .....3oz (90cc)

ALM90-6010 .....4oz (120cc)

ALM90-6020 .....6oz



## EMESIS BASIN, SINGLE USE

Made in plastic, this emesis basin has a 16oz (0.47L) capacity. Non-autoclavable.

AMG760-680 .....Plastic, Beige



## EMESIS BASIN, STAINLESS STEEL

AMG020-506 .....8", 13oz

AMG020-505 .....10", 26oz



## EYE WASH STATION AND SOLUTION

BIO MED WASH® is an all-natural, sterile first aid wash for cuts and lacerations, is the first product that replaces saline solution for washing abrasions and lacerations (saline is still preferred for deep intramuscular cuts). The product is invaluable in the initial washing of skin that has come in contact with irritating chemicals such as gasoline, pepper spray, tear gas or plant toxins; and the washing of eyes irritated by foreign particles, heavy smoke, dust or perspiration. BIO MED WASH contains 100% Bio-Logic Aqua® tissue-culture grade water, proven to have the ideal pH balance, hypotonicity and mild mineral nutrient content to maximize moisture absorption and promote healing.

RBW47085 .....Eyewash Station, 2 x 7oz Cannisters, Mirror

RBW40007 .....7oz Cannisters



## EYE WASH SOLUTION

Effective against redness, smarting and mild congestion of the eyes and eyelids due to light wind, dust, smoke or other pollutants.

SFA04077 .....Isotonic Solution Sterile, 32oz / 1L



## HYDROGEN PEROXIDE

RW 40153 .....500ml, 3%

RW 6040140162 .....250ml, 6%



## ISOPROPYL ALCOHOL

RW 040193 .....500ml, 70%

RW 040181 .....500ml, 99%



## ALCOHOL DISPENSER

DUK4024 .....9oz Dispenser



## COTTON BALLS

DYN3170 .....Medium Bag/2000



## PLASTIC MEDICINE CUPS, GRADUATED

DYN4258 .....1oz Pkg/100



## PAPER MEDICINE / PORTION CUPS

DYN4244 .....0.5oz Pkg/250



## PAPER DRINKING CUPS

SPMC00004 .....4oz Poly Coated Sleeve/100

SPMC00005 .....5oz Poly Coated Sleeve/100



## DRINKING CUPS, PLASTIC

DYN4236.....5oz, White, Plastic, Disposable Pkg/50

DYN4237 .....5oz, Blue, Plastic, Disposable Pkg/50

DYN4238 .....5oz, Mint Green, Plastic, Disposable Pkg/50

DYN4239 .....5oz, Mauve, Plastic, Disposable Pkg/50

DYN4240 .....5oz, Lavender, Plastic, Disposable Pkg/50



## FLEXIBLE DRINKING STRAW

DYN4260.....Individually Wrapped Box/400



## PROTECTIVE EYE GLASSES

DYN2198



## **PRI-MED** FACE SHIELD, FULL, DISPOSABLE

Pre-assembled, one-piece, face shield that includes an optical-grade lens, foam headband, and adjustable elastic strap to ensure a comfortable fit. Latex-free.

PRIPG4-9700



## EARLOOP FACE MASK

Earloop masks in Level 1, 2, and 3

SPM2110 ..... 3ply ASTM Level 1 Blue, 50 Per Box  
 SPM2120..... 3ply ASTM Level 2 Blue, 50 Per Box  
 SPM2122 ..... 3ply ASTM Level 2 Pink, 50 Per Box  
 SPM2130 ..... 3ply ASTM Level 3 Blue, 50 Per Box



## **PRI-MED** EARLOOP FACE MASK, WITH FACE SHIELD

Available with optically clear plastic wraparound face shield which provides protection for the face and eyes against pathogens contained in body fluids. Mask has fluid resistant shield. Latex-free.

PRIPG4-1573..... Primed ASTM Level 3, 25 Per Box



## **PRI-MED** CHILD SIZE FACE MASK

Child Size Face Mask with earloops is soft, non-irritating and easy to breathe through. The dye free mask is not made with natural rubber latex and is comfortable to wear for pediatric size patients.

PRIPG4-1680 ..... Child Size ASTM Level 1 Dinosaur Print 50 Per Box  
 PRIPG4-1693 ..... Small Adult ASTM Level 3 50 Per Box





## 3M QUALITATIVE FIT TEST KIT

3M Respirator Fit Test Kits can be used with any particulate respirator to test face mask to face seal fit. The 3M Fit Test Kit can help satisfy regulatory requirements for face mask fit testing. Available in either a bitter or sweet test solution kit. Simply spray the solution into the hood environment. If the wearer does not detect the test agent's taste, an acceptable respirator fit has been demonstrated. The 3M Respirator Fit Test Kit can fit test approximately 150 employees with the included supplies.

### Kit Includes:

- 55ml Sensitivity Solution
- 55ml Fit Test Solution
- Two nebulizers (FT-13)
- Laminated instruction book (English/Spanish/French/Portuguese)
- One collar assembly (FT-15)
- Two sets replacement nebulizer inserts
- One hood (FT-14)

3M FT-10.....Test Kit with Sweet Sensitivity Solution

3M FT-30.....Test Kit with Bitter Sensitivity Solution



### SENSITIVITY SOLUTION, SWEET, 55ML



3M FT-11.....Box/6

### TEST SOLUTION, SWEET, 55ML



3M FT-12.....Box/6

### SENSITIVITY SOLUTION, BITTER, 55ML



3M FT-31.....Box/6

### TEST SOLUTION, BITTER, 55ML



3M FT-32.....Box/6

## 3M N95 FACE MASK, FLAT FOLD



3M 9210PLUS ..... Box/20

## 3M N95 FACE MASK, MOLDED



Regular Size

3M 8210..... Box/20

Small Size

3M 8110S ..... Box/20

## 3M N95 FACE MASK, FLAT FOLD, FLUID RESISTANT



3M 1870PLUS ..... Box/20

## 3M N95 FACE MASK, MOLDED FLUID RESISTANT



Regular Size

3M 1860..... Box/20

Small Size

3M 1860S..... Box/20

## ISOLATION GOWN 50/CS

- Universal
- Offers economical, comfortable and reliable protection
- Full length for maximum coverage
- Elastic cuffs
- Sewn seams provide superior strength
- Latex free
- Extra length on waist ties allows gown to be easily secured



AMD8022 .....White, Impervious  
DYN2141 .....Yellow, Fluid Resistant

## **PRI-MED OVERHEAD PROTECTIVE FILM GOWN**

- Lightweight layer of protective polyethylene film with an open back design.
- AAMI Level 3

PRIPM4-88335 .....20 Per Box



## **PRI-MED CHEMOTHERAPY DRUG TESTED ISOLATION GOWN**

- ASTM F739-12 (using USP 800 Guidelines)
- Full Back Coverage
- Long Waist Ties
- Extra Soft Knit Cuff
- AAMI Level 2

PRIPM4-99450 .....chemotherapy approved ..... 10 per bag



## **PRI-MED AAMI 2 ISOLATION GOWN X-LARGE**

- Full Back
- Yellow
- Thumbloop with Elastic Cuff
- Ties

PRIPM4-57162XL ..... 100 Per Case



## SHOE COVERS

- Non-Skid, Non-Conductive
- Will cover shoe size from S to L shoes (SPM6007) and from M to XL (SPM6006) (gender non-specific)
- Manufactured of spunbond polypropylene
- Sewn seams, not glued - offers greater strength
- Blue in color
- Elastic opening provides a snug comfortable fit

SPM6007 ..... Universal ..... Bag/100  
SPM6006 ..... X-Large ..... Bag/100



## SHOE COVERS PLASTIC/WATERPROOF

Resistant to water, snow and slush, AMG plastic shoe covers are convenient, cost effective and safe. Made from embossed plastic to prevent slipping, these shoe covers are waterproof and feature an elastic band opening for universal fitting. Blue.

AMG018-285 ..... Universal Size ..... Cs/250 Pairs



## BOUFFANT CAP

- Lightweight and breathable spunbond

AMD8004.....24" Blue..... 100 per bag



- Adjustable with Ties

PRIPM4-801 ..... Surgeon's Cap, SMS, Blue, Universal, Fluid Resistant 100 Per Box

PRIPM4-804 ..... Surgeon's Cap, Kaycel, Blue, Large, Ultra Sweat Absorbent 100 Per Box



## COVER YOUR COUGH COMPLIANCE KIT

Cover Your Cough (CYC) Compliance Kits hold an assortment of personal protection equipment in a convenient and attractive dispenser. Kit 101 holds two boxes of earloop face masks, one - two boxes of facial tissue, and one hand sanitizer bottle or hand sanitizing wipe canister (consumables not included). The sign holder is vertical and a variety of signs can be downloaded and printed from [www.bowmandispensers.com](http://www.bowmandispensers.com). Proper and consistent use of the products in this dispenser reduces the chance of spreading infectious bacteria and viruses, such as those found in influenza and the common cold. CYC Compliance Kits are ideal for use in public waiting areas, offices, and lobbies. 18" W x 59" H

### Kit 101 includes the following products:

- BOWMAN® Respiratory Hygiene Station (RS001 0512)
- BOWMAN® Floor Stand (KS010-0412)
- BOWMAN® 8.5 x 11" Sign Holder - Vertical (MP-070)

BOWRS001-0212..... Hygiene Station Only

BOWKS010-0412 ..... Floor Stand Only, Steel

BOWMP-070..... 8.5" x 11" Sign Holder, Vertical, Clear Plastic



## INFECTION PROTECTION SYSTEM

Holds flat pack or launderable gowns, three boxes of gloves, one box of earloop face masks and one box of face shields. Keyholes for wall mounting or may be used with a wall or door hanger. Quartz ABS plastic.

BOWLD-070



## PPE CART II PROTECTION SYSTEM MOBILE

Easy roll 3" casters with toe-locking brakes, steel doors with magnetic closures, safety yellow powder-coated aluminum and steel.

BOWCT030-0000 .....Mobile PPE Cart



## PEAK FLOW METER

International, independent testing has proven that both the Mini-Wright and the Mini-Wright ATS are the most accurate peak flow meters available. All meters are individually calibrated by hand to ensure consistent accuracy and reproducibility.

MCA3103085 ..... Adult  
MCA3104050 ..... Pediatric



## DISPOSABLE MOUTHPIECES

Our disposable one-way-valved and single patient use mouthpieces are designed to provide protection against possible cross infection in common respiratory procedures.

MCA3122200 ..... Adult.....Bag/100  
MCA3123200 ..... Pediatric.....Bag/100



## OXYGEN MASKS & NASAL CANNULA

### Cannula

ZEN1012-SC..... Adult Cannula with Tubing  
ZEN1101 ..... Child Cannula with Tubing

### Masks-Medium Concentration

ZEN1041 ..... Adult Mask with Tubing  
ZEN1043 ..... Child Mask with Tubing

### Oxygen Tubing

ZEN1115-07..... 7'  
HUD1118..... 14'  
ZEN1115-25..... 25'



## NEBULIZER MASKS

ZEN1707 ..... Child Mask & Cup  
ZEN1805 ..... Adult Mask & Cup  
ZEN1809 ..... Mouthpiece & Cup



## CONNECTING TUBE, STERILE

CHS70-6072S..... 3/16" x 72"



## DAR® BARRIERBAC "S" VIRAL BREATHING FILTER

COV350/5879 ..... Box/25



## MANUAL RESUSCITATOR

Controlled ventilation. Disposable.

ZEN4050 ..... Adult, 66lb+  
ZEN4010 ..... Pediatric, 22-66lb  
ZEN4025 ..... Infant, 0-22lb



## GUEDEL AIRWAY KIT

Single use, colour coded airways. Eight sizes (40mm, 50mm, 60mm, 70mm, 80mm, 90mm, 100mm, 110mm). Comes in clear poly case with airway tray. Soft, disposable.

ZEN2200-1 ..... Guedel Airway Kit .....Set/8



## GUEDEL AIRWAY NON STERILE

Individual easy open packaging.

ZEN2104 ..... 40mm - Pink  
ZEN2105 ..... 50mm - Blue  
ZEN2106 ..... 60mm - Black  
ZEN2107 ..... 70mm - White  
ZEN2108 ..... 80mm - Green  
ZEN2109 ..... 90mm - Yellow  
ZEN2110 ..... 100mm - Red  
ZEN2111 ..... 110mm - Orange





## CPR FACESHIELD WITH KEYCHAIN CASE

CPR Protective Barrier with Belt/Keyring 11 x 11cm Pouch (Red), One Pair of Latex Free Gloves.

ZEN3052R



## CPR POCKET RESUSCITATOR

CPR Pocket Ventilator, Head Strap, Oxygen Inlet, Filter, Gloves and Complete with hard plastic clamshell carry case.

ZEN3100-10



## EMERGENCY RESPONDER KIT

- Includes Oxygen
- 1 Guedal Airway Kit
- 1 Universal Face Mask, BVM Reservoir, Oxygen Tubing
- 1 Cylinder "D" with 6LPM Regulator
- 1 Adult Mask 3 Valves with Tubing Non-Breathing
- 1 Pediatric Mask 3 Valves with Tubing Non-Breathing

ZEN911..... Complete Kit



## PORTEX® SOFT SEAL LMA MASK, DISPOSABLE

BCI100/222/200..... Size 2.0 Infant 10kg-20kg  
 BCI100/222/250..... Size 2.5 Infant 20kg-30kg  
 BCI100/222/300..... Size 3.0 Infant 30kg-60kg  
 BCI100/222/400..... Size 4.0 Infant 50kg-70kg  
 BCI100/222/500..... Size 5.0 Infant 70kg-100kg



## ASPIRATOR (MANUAL VACUUM/SUCTION PUMP)

State of the art, high vacuum, light, portable emergency suction device. Built in 3M HEPA Filter for bacterial and viral protection. Adult and Child settings, Built-in overflow protection seal, 100° swivel. Disposable collection jar and catheter.

ZEN7000 ..... Suction Machine  
ZEN7000-AD ..... Disposable Jar and Catheters



## OPTIVAC® AC/DC PORTABLE SUCTION UNIT

Long life rechargeable battery. Compact design. Heavy Duty. 3 Year warranty, 1 Year on battery.

AHPG180 ..... Suction Machine



## EPINEPHRINE INJECTION, SINGLE DOSE AMPOULE

WIN07241L0 ..... 1ml Ampule

\*DIN00721891



## EPIPEN® EPINEPHRINE AUTO INJECTOR

WIN00121 ..... EpiPen®  
WIN00111 ..... EpiPen® Jr.



## NORMAL SALINE, 0.9% SODIUM CHLORIDE

For respiratory therapy.

CHS116 ..... 15ml ..... Box/24



## MEDPRO PLASTIC ADHESIVE BANDAGES

AMG018-600 ..... 3/4" x 3" Box/100



## FABRIC STRIP ADHESIVE BANDAGES

AMDAF0334..... 3/4" x 3" Box/100



## SHEER SPOT ADHESIVE BANDAGES 7/8"

- Transparent on the skin

DUK1307033..... 7/8" Box/100



## FLEXIBLE FABRIC KNUCKLE BANDAGES

DUK99949..... 1.5" x 3" Box/100



## SHEER ADHESIVE BANDAGES 2" X 4" X-LARGE

DUK 11800..... 100 per Box



## KIDS BANDAGES

DIS912816..... Spiderman box/20

DIS912824..... Frozen box/20

DUK15606..... Emoji box/100



## Leukoplast® STRONG HEAVYWEIGHT FABRIC BANDAGES

BSN7646007 ..... 2.2 cm x 7.2 cm 100 Per Box  
 BSN7646005 ..... 2.2 cm x 3.8 cm 100 Per Box  
 BSN7646009 ..... 3.8 cm x 3.8 cm 100 Per Box  
 BSN7646013 ..... 5 cm x 7.2 cm 100 Per Box  
 BSN7646010 ..... 3.8 cm x 7.5 cm Knuckle 100 Per Box  
 BSN7646012 ..... 4.5 cm x 5 cm Fingertip 50 Per Box  
 BSN764601 ..... 4.5 cm x 7.5 cm Fingertip 100 Per Box  
 BSN7645908 ..... Doctor Bandages Set 101 Assorted per Box



**Leukoplast®****STRONG HEAVYWEIGHT FABRIC DRESSING**

BSN7645900..... 4cm x 5m 1 Roll Per Box  
 BSN7645900..... 6cm x 5m 1 Roll Per Box  
 BSN7645902..... 8cm x 5m 1 Roll Per Box

**DYNAREX PAPER SURGICAL TAPE**

DYN3552 ..... 1" x 10 Yards Box/12

**DYNAREX CLEAR SURGICAL TRANSPARENT TAPE**

DYN3572 ..... 1" x 10 Yards Box/12

**DYNAREX CLOTH SURGICAL TAPE**

DYN3562 ..... 1"  
 DYN3561 ..... 1/2"

**LEUKOPLAST HYPAFIX NON-WOVEN ADHESIVE**

A cross-elastic non-woven tape that provides light compression. Excellent for areas requiring flexibility and freedom from constraint while providing superior holding power. Air and exudate permeable. Hypoallergenic.

BSN7144300..... 2.5cm x 10m 1 roll per box  
 BSN7144301..... 5.0cm x 10m 1 roll per box  
 BSN7144302..... 10 cm x 10m 1 roll per box  
 BSN7144303..... 15 cm x 10m 1 roll per box  
 BSN7144304..... 20 cm x 10m 1 roll per box



## 3M MICROPORE™ PAPER TAPE

A latex free, hypoallergenic paper tape that is gentle to the skin yet adheres well and leaves minimal adhesive residue upon removal. An economical, general purpose, breathable surgical tape.

3M 1530-0.....	1/2" x 10 Yards.....	Box/24
3M 1530-1.....	1" x 10 Yards .....	Box/12
3M 1530-2.....	2" x 10 Yards .....	Box/6



## 3M TRANSPORE™ TRANSPARENT TAPE

A latex free, hypoallergenic, transparent and perforated plastic tape that offers strong adhesion. Transpore tape is easy to tear into very thin strips, and is easy to use even with gloves.

3M 1527-0.....	1/2" x 10 Yards.....	Box/24
3M 1527-1.....	1" x 10 Yards .....	Box/12
3M 1527-2.....	2" x 10 Yards .....	Box/6



## 3M BLENDERM™ TRANSPARENT TAPE

A latex free, hypoallergenic, occlusive and transparent plastic tape. This flexible medical tape is designed to provide a barrier from external fluids and contamination.

3M 1525-1.....	1" x 5 Yards.....	Box/12
3M 1525-2.....	2" x 5 Yards.....	Box/6



## TENOPLAST® ELASTIC ADHESIVE TAPE

BSN7206902.....	2.5cm x 4.5m	
BSN7205218.....	5.0cm x 4.5m .....	Box/6 RL
BSN7205219.....	7.5cm x 4.5m .....	Box/4 RL



## ADHESIVE REMOVER

- Effectively removes tape and residue from skin
- Reduces patient discomfort
- Contains C11-12 Isoparaffin, Isopropyl Alcohol
- Petroleum Naphtha, Isopropyl Palmitate, Isopropyl
- Myristate, Isopropyl Stearate

PDIB16400.....	Remover Pad .....	Box/100
----------------	-------------------	---------





## WOVEN GAUZE SPONGES, NON-STERILE

### 8-Ply Non-Sterile

AMDD3001 .....	2" x 2" .....	Pkg/200
DYN3232 .....	3" x 3" .....	Pkg/200
AMDD3005 .....	4" x 4" .....	Pkg/200

### 12-Ply Non-Sterile

DYN3223 .....	2" x 2" .....	Pkg/200
DYN3233 .....	3" x 3" .....	Pkg/200
DYN3243 .....	4" x 4" .....	Pkg/200



## WOVEN GAUZE SPONGES, STERILE 2'S

### 8-Ply Sterile 2's

DYN3322 .....	2" x 2" .....	Box/50 Pkg of 2
DYN3342 .....	4" x 4" .....	Box/25 Pkg of 2



## NON-WOVEN GAUZE SPONGES 4-PLY, NON-STERILE

AMDA2101-CH .....	2" x 2" .....	Pkg/200
AMD2100-IMC .....	4" x 4" .....	Pkg/200



## NON-WOVEN GAUZE SPONGES 4-PLY, STERILE 2'S

AMDA2240 .....	2" x 2" .....	Box/50 Pkg of 2
DYN3363 .....	3" x 3" .....	Box/50 Pkg of 2
AMD4440-IMC .....	4" x 4" .....	Box/50 Pkg of 2



## NON-ADHERENT PADS

AMDA1923 .....	2" x 3" .....	Box/100
COVZ6017 .....	2" x 3" w/Adhesive .....	Box/100
AMDA1934 .....	3" x 4" .....	Box/100
DYN3438 .....	3" x 8" .....	Box/50



## ABDOMINAL PADS

DYN3501 ..... Sterile, 5" x 9" .....Tray/20

DYN3503 ..... Sterile, 8" x 10" .....Tray/24



## EYE PADS, STERILE

DYN3360..... 1-5/8" x 2-5/8" ..... Box/50



## PACKING STRIPS, STERILE, PLAIN & WITH IODOFORM

### PLAIN, 1 per Bottle

DYN3401-IMC..... 1/4" x 5 Yards

DYN3402-IMC..... 1/2" x 5 Yards

DYN3403-IMC..... 1" x 5 Yards

### WITH IODOFORM, 1 per Bottle

DYN3411-IMC..... 1/4" x 5 Yards

DYN3412-IMC..... 1/2" x 5 Yards

DYN3413-IMC..... 1" x 5 Yards



## DRESSING TRAY, STERILE

- 1 Three Compartment Tray
- 1 Plastic Transfer Forceps
- 2 Metal Forceps
- 4 - 10 x 10cm Gauze Sponges
- 4 - 5 x 5cm Gauze Sponges
- 1 Waste Bag with Tie
- 1 Safety Pin
- 2 Towel/Underpad

TRA85-2000TY.....Each

TRA85-2000..... Cs/30



## FIRST AID KITS

Workplace, Type 2, Basic, Small, (2-25 Workers)

SFA51611..... Type 2 (AB, MB, NB, NL, NS, NT, ON, PEI, QC, SK & YT)



## CONFORM STRETCH GAUZE, NON-STERILE

DYN3102 .....	2" x 4.1 Yards .....	Box/12
DYN3103 .....	3" x 4.1 Yards .....	Box/12
DYN3104 .....	4" x 4.1 Yards .....	Box/12
DYN3106.....	6" x 4.1 Yards .....	Box/6



## CONFORM STRETCH GAUZE, STERILE

DYN3112 .....	2" x 4.1 Yards .....	Box/12
DYN3113 .....	3" x 4.1 Yards .....	Box/12
DYN3114 .....	4" x 4.1 Yards .....	Box/12
DYN3116.....	6" x 4.1 Yards .....	Box/6



## ELASTIC TENSOR BANDAGES

DYN3662.....	2" x 5 Yards
DYN3663.....	3" x 5 Yards
DYN3664.....	4" x 5 Yards
DYN3666.....	6" x 5 Yards



## STOCKINETTE TUBULAR 100% COTTON BANDAGE

DYN3692.....	2" x 25 Yards
DYN3693.....	3" x 25 Yards
DYN3694.....	4" x 25 Yards
DYN3696.....	6" x 25 Yards



## TRIANGULAR BANDAGE

DYN3681



## SURGILAST® ELASTIC NET BANDAGE, LATEX FREE, 25 Yards

ZENGL-701 .....	Size 1 Roll
ZENGL-702 .....	Size 2 Roll
ZENGL-703 .....	Size 3 Roll
ZENGL-704 .....	Size 4 Roll
ZENGL-705 .....	Size 5 Roll
ZENGL-706 .....	Size 5.5 Roll
ZENGL-707 .....	Size 6 Roll
ZENGL-708 .....	Size 7 Roll
ZENGL-709 .....	Size 8 Roll



## SURGITUBE® TUBULAR GAUZE BANDAGE FOR APPLICATOR, 50 Yards

ZENGL-219 .....	Size 1, Small Fingers & Toes Roll
ZENGL-220 .....	Size 2, Large Fingers & Toes Roll
ZENGL-221 .....	Size 3, Hands, Wrist & Feet Roll
ZENGL-222 .....	Size 4, Arms & Lower Legs Roll
ZENGL-223 .....	Size 5, Legs, Thighs & Head Roll



## SURGITUBE® APPLICATORS

ZENGL-230 .....	Size 0
ZENGL-231 .....	Size 1
ZENGL-232 .....	Size 2
ZENGL-234 .....	Size 3
ZENGL-236 .....	Size 4



## DYNAREX REUSABLE HOT AND COLD GEL PACKS

Dynarex Reusable Hot and Cold Gel Packs are convenient compresses designed to provide hot and cold therapy for temporary relief of muscle and joint soreness. These multi-purpose Packs can be easily frozen for conditions requiring cold therapy or microwaved for heat applications.

DYN4504 .....	4" x 6"
DYN4505 .....	5" x 11"
DYN4506 .....	6" x 9"



## HOT & COLD PACKS, INSTANT, DISPOSABLE

DYN4512 .....	Cold 5"x9 Cs/24
DYN4516 .....	Hot 5"x9



## CIRCULAR COOL GEL PACKS

Reusable.

RAP11331-48 .....	3" Cs/48
RAP11351-24 .....	5" Cs/24



## POLYSPORIN® OINTMENT

IDC75620..... 15gm Tube

IDC75647..... 30gm Tube



## BACITRACIN ZINC, FIRST AID ANTIBIOTIC OINTMENT

Helps prevent infection in minor cuts, scrapes, and burns.

SFA05796 ..... 0.9 g Foil Pack.....25/Box



## CUTICELL® DRESSING

Cuticell® is soothing and low-adherent minimizing pain at dressing change. The primary wound contact layer is made of smooth acetate fabric impregnated with non-medicated ointment on Eucerit® basis. The result: an ointment dressing that ensures that the wound is not disturbed and helps maintain a moist environment that supports wound healing.

BSN7253901..... 3" x 3" .....Box/10



## JELONET® DRESSING

SMN7403..... 5cm x 5cm.....Box/50



## TEGADERM™ TRANSPARENT FILM DRESSING

3M 1624W..... 6cm x 7cm..... Box/100

3M 1626W..... 10cm x 12cm.....Box/50



## OPSITE® TRANSPARENT FILM DRESSING

### Flexifix Film

SMN66000040 ..... 5cm x 10m

SMN66000041 ..... 10cm x 10m

### Flexigrid Dressing

SMN4628..... 6cm x 7cm..... Box/100

SMN4630..... 10cm x 12cm.....Box/50



## OPSITE® SPRAY

SMN66004978 ..... 100ml Spray Bottle





## SALJET® SINGLE DOSE NORMAL SALINE FOR IRRIGATION

Single use 30ml containers. Saves time, wastage and assures sterility.

WLB10-18815100009-08 ..... 30ml.....Box/40



## IRRIGATION SOLUTIONS

### Normal Saline

BAXJF7633P ..... 500ml

BAXJF7634 ..... 1000ml

### Sterile Water

BAXJF7623P ..... 500ml

BAXJF7624 ..... 1000ml



## STERILE WATER FOR IRRIGATION

MAXDYND40570 ..... 110ml cup.....Cs/48



## NORMAL SALINE 0.9% (SODIUM CHLORIDE) FOR IRRIGATION

MAXDYND40540 ..... 100ml



## DYNACAST® PRELUDE SPLINT SYSTEM (ROLL FORMAT)

Fiberglass.

BSN7141800.....2.5cm x 4.6m Roll

BSN7142000.....5.0cm x 4.6m Roll

BSN7142100.....7.5cm x 4.6m Roll

BSN7141500.....10cm x 4.6m Roll

BSN7141600.....12.5cm x 4.6m Roll

BSN7141700.....15cm x 4.6m Roll



## DELTA-LITE® PLUS CAST TAPE

Fiberglass.

BSN7345800.....	2.5cm x 1.8m, White .....	Box/10
BSN7345801.....	5.0cm x 1.8m, White .....	Box/10
BSN7345802.....	7.5cm x 3.6m, White .....	Box/10
BSN7345803.....	10cm x 3.6m, White .....	Box/10
BSN7345804.....	12.5cm x 3.6m, White.....	Box/10

\*Other colours available\*



## SOFFBAN CAST PADDING - SYNTHETIC

BSN7146401.....	5cm x 2.7m .....	Bag/12
BSN7146701.....	7.5cm x 2.7m .....	Bag/12
BSN7148605.....	10cm x 2.7m .....	Bag/12
BSN7146501.....	15cm x 2.7m .....	Bag/12



## GYPSONA® PLASTER OF PARIS CAST BANDAGE

BSNLPL23 .....	5cm x 2.7m .....	Box/24
BSNLPL33 .....	7.5cm x 2.7m.....	Box/12
BSNLPL43 .....	10cm x 2.7m.....	Box/12
BSNLPL63 .....	15cm x 2.7m.....	Box/12
BSNLPL65 .....	15cm x 4.5m.....	Box/12
BSNLPL83 .....	20cm x 2.7m.....	Box/12
BSNLPL85 .....	20cm x 4.5m.....	Box/12



## GYPSONA® PLASTER OF PARIS CAST SLABS

BSNL5330 .....	7.5cm x 75cm	Box/50
BSNL5430 .....	10cm x 10cm	Box/50
BSNL5530 .....	12.5cm x 75cm	Box/50



## ZIP STICK CAST REMOVAL AID

BSN7204629.....	2.5cm x 48cm
-----------------	--------------



## ALUMAFOAM FINGER SPLINTS

CCO6212 .....	1/2" x 18" .....	Box/12
CCO6234 .....	3/4" x 18" .....	Box/12
CCO6210.....	1" x 18" .....	Box/12



## DENTAL BIBS

DYN4491 .....	17-3/4" x 12-7/8", White
DYN4492 .....	17-3/4" x 12-7/8", Blue
DYN4493 .....	17-3/4" x 12-7/8", Mint Green
DYN4494 .....	17-3/4" x 12-7/8", Mauve
DYN4495 .....	17-3/4" x 12-7/8", Lavender



## DENTAL BARRIER FILM

Eliminates cross contamination while working with handheld equipment.  
Covers most surface areas.

DYN4488 .....	4" x 6", Clear .....	Box/1200
DYN4489 .....	4" x 6", Blue.....	Box/1200



## IMPRESSION TRAYS

- Disposable perforated plastic trays
- Consistent impressions
- Small, medium and large upper and lower trays, quadrant and anterior
- 12 pieces per bag, 10 bags per case
- Not made with natural rubber latex

DYN2220 .....	No. 1, Large Upper.....	Bag/12
DYN2221 .....	No. 2, Large Lower.....	Bag/12
DYN2222 .....	No. 3, Medium Upper.....	Bag/12
DYN2223 .....	No. 4, Medium Lower.....	Bag/12
DYN2224 .....	No. 5, Small Upper .....	Bag/12
DYN2225 .....	No. 6, Small Lower .....	Bag/12
DYN2234 .....	No. 7, Quadrant Tray, UL/LR .....	Bag/12
DYN2235 .....	No. 8, Quadrant Tray, UR/LL .....	Bag/12
DYN2238 .....	No. 9, Anterior Tray, Upper .....	Bag/12 (Bag Size 25x17)
DYN2239 .....	No. 10, Anterior Tray, Lower .....	Bag/12 (Bag Size 25x17)



## EARLOOP FACE MASK

Procedure Masks offer a flexible nose piece that helps form a closure for better protection against exposure. Pleated earloop procedure Masks with glass free filter. Not made with natural rubber latex.

Earloop Masks in level 1, 2, and 3 .....	Box/50
--	--------



## BOUFFANT CAP

- Lightweight and breathable spunbond

PRIPM4-6644.....Blue Small.....100 per bag



## SHOE COVERS

- Non-Skid, Non-Conductive
- Will cover shoe size from S to L shoes (SPM6007) and from M to XL (SPM6006) (gender non-specific)
- Manufactured of spunbond polypropylene
- Sewn seams, not glued - offers greater strength. Blue in color
- Elastic opening provides a snug comfortable fit

SPM6007..... Universal .....Bag/100

SPM6006..... X-Large .....Bag/100



## SHOE COVERS PLASTIC/WATERPROOF

Resistant to water, snow and slush, AMG plastic shoe covers are convenient, cost effective and safe. Made from embossed plastic to prevent slipping, these shoe covers are waterproof and feature an elastic band opening for universal fitting. Blue.

AMG018-285 ..... Universal Size .....Cs/250 Pairs



## SURGEONS CAPS, DISPOSABLE

PRIPM4-800 ..... Premium SMS..... Box/100

DYN2121 ..... Box/100



## ISOLATION GOWN TEAL UNIVERSAL 5/30bag/cs

- PriMed Overhead Fullback
- Thumbloop with Elastic Cuff
- 34gsm Teal SMS
- AAMI Level 3

PRIPM4-22430.....Bag/5



## NITRIDERM® NITRILE STERILE EXAM GLOVES

Nitrile synthetic contains no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort. Pair packaged, inner wrapper in ozone resistant multi-vac packaging.

IHC106100 .....Small  
IHC106200 .....Medium  
IHC106300 .....Large  
IHC106350 .....X-Large

50 Pairs/Box



## NOAH STRETCH NITRILE - BOX OF 100

- Finger Textured
- Nitrile
- Powder free
- 3.0ml thickness

SPM88896 .....Small  
SPM88897 .....Medium  
SPM88898 .....Large



## BLUSH® NITRILE PINK EXAMINATION GLOVES - BOX OF 200

- Nitrile (100% Latex-Free)
- Non Sterile
- Glove Design & Feature: Ambidextrous and Finger Textured
- Cuff: Beaded
- Color: Blush Pink
- 2.5ml thickness
- Packaging: 200 pieces per box / 10 boxes per case

SPM78885 .....X-Small  
SPM78886 .....Small  
SPM78887 .....Medium  
SPM78888 .....Large



## PRIMED PRIMATOUGH® FIT NITRILE GLOVES (NON-POWDERED) - BOX OF 200

PRIPM6-2005 .....X-Small  
PRIPM6-2006 .....Small  
PRIPM6-2007 .....Medium  
PRIPM6-2008 .....Large  
PRIPM6-2009 .....X-Large



See more glove options on pages 63-67,  
or check out the website for all available gloves.



## RUBBER BANDS

- Made of high quality rubber
- Fits tattoo machines that use size #12 & #13
- Helps keep your tattoo needle in place while tattooing
- A must-have for all standard tattoo machines
- One of the most commonly replaced parts on a tattoo machine



DYN2166..... 1L .....Bag

## GROMMET NIPPLES

- Gives the needle a snug fit and keeps it in place
- Easily slip your needle over the top of the nipple, secures needle tightly on the A-bar
- Rubber nipples designed specifically for tattoo machine use
- Packaged in bags of 100



DYN2179 ..... Black .....Bag/100

## INK CUPS

- Used to hold pigment or ink
- Enables you to place your ink cup directly on the work surface without the need for an ink cup stand
- Transparent to allow a better view of the ink color
- Disposable
- Non-Sterile
- Available in Small, Medium, Large & X-Large



DYN2175 ..... Medium, 14mm..... Bag/1000

## SOAKER PADS

- Absorbent center is wrapped and sealed on all four sides by a polyethylene film with perforations on one side
- Pads allow tattoo to breathe and not dry out
- Non-sterile
- Non-stick for easy removal
- Absorbs excess ink after tattoo application



DYN2168..... 4" x 7", Black..... Box/500

DYN2169..... 6.25" x 7", Black..... Box/500

## STERILIZATION POUCHES, SELF SEALING

- Self sealing, no tape or heat sealing needed
- Blue tinted transparent cover makes it easy to detect contamination
- Internal and external color-changing indicators assure proper processing in chemical vapor, steam, or ethylene oxide (EO) gas sterilization
- Medical grade materials
- Latex Free
- Pouch sizes below are calculated as outer dimensions



DYN4460.....	2.25" x 2.75"	.....	Box/200
SPMSP71300.....	7.5" x 13"	.....	Box/200
SPMSP35900.....	3.5" x 9"	.....	Box/200
SPMSP52510.....	5.25" x 10"	.....	Box/200
SPMSP22540.....	2.25" x 4"	.....	Box/200
SPMSP27590.....	2.75" x 9"	.....	Box/200

## BOTTLE COVERS

- Bottle Bags will cover your entire spray or wash bottle, preventing cross contamination
- Transparent, for easy identifying the liquid level in the bottle
- Excellent for keeping your equipment clean and germ free

DYN2164.....	Medium, 6" x 8"	.....	Box/500
DYN2165.....	Large, 6" x 10"	.....	Box/500



## MACHINE BAG

- Protect your machine from cross contamination and splatter
- These covers easily slip over your machine

DYN2162.....	5" x 5"	.....	Box/500
--------------	---------	-------	---------



## DENTAL BIBS

DYN4491 .....	17-3/4" x 12-7/8"	.....	White
DYN4492 .....	17-3/4" x 12-7/8"	.....	Blue
DYN4493 .....	17-3/4" x 12-7/8"	.....	Mint Green
DYN4494 .....	17-3/4" x 12-7/8"	.....	Mauve
DYN4495 .....	17-3/4" x 12-7/8"	.....	Lavender



## DENTAL BARRIER FILM

Eliminates cross contamination while working with handheld equipment. Covers most surface areas.

DYN4488 .....	4" x 6"	.....	Clear
DYN4489 .....	4" x 6"	.....	Blue



## BLACK NITRILE EXAM GLOVES

Extraordinary strength and puncture resistance. Maintaining tactile sensitivity. Contains no allergy causing natural rubber proteins. Commonly used in tattoo, EMS and other first responder markets. Ambidextrous, Non-Sterile, Powder Free. Not made with natural rubber latex. Black.



## DISPOSABLE THONG

Made of soft spunbonded polypropylene material. Each garment is individually wrapped for hygienic purposes. One Size.

DUK900502 ..... Black ..... Cs/250



## DISPOSABLE BIKINI

Made of soft spunbonded polypropylene material. Each garment is individually wrapped for hygienic purposes. One Size.

DUK900506 ..... Black ..... Cs/250



## DISPOSABLE BRA

Backless. Made of soft spunbonded polypropylene material. Each garment is individually wrapped for hygienic purposes. One Size.

DUK900510 ..... Small/Medium ..... Cs/250

DUK900512 ..... Large/X-Large ..... Cs/250



## DISPOSABLE MENS BRIEF

Made of soft spunbonded polypropylene material. Each garment is individually wrapped for hygienic purposes.

DUK900528 ..... Large/X-Large ..... Cs/250



## DISPOSABLE MENS BOXER

- Provides comfort with just the right amount of coverage
- Made of soft spun bonded polypropylene material
- Available in a variety of sizes and styles to accommodate multiple treatments
- Each garment is individually wrapped for hygienic purposes
- Disposable
- Not made with natural rubber latex

DUK900530 ..... Small/Medium ..... Cs/50

DUK900532 ..... Large/X-Large ..... Cs/50



## DISPOSABLE HEADBAND

- Convenient Velcro Closure
- Used to keep hair from face during spa treatments

GRM51823 ..... Pkg/48



## CLEAR ULTRASOUND GEL

- Get "clear results" with our dye free gel
- Wavelength clear ultrasound gel is great for procedures requiring a pigment-free, medium viscosity conductive agent
- Contains no salt, alcohol, fragrance or dye



NTPC202 ..... 250ml ..... Bottle  
NTPC503 ..... 5L ..... Bottle

## COLD PACK GEL CIRCLES

- Reusable soft-sided gel pack, ideal for a variety of cold therapy applications, such as following after injections, blepharoplasty, and breast surgery
- Remains flexible when frozen
- Soft-sided cloth backing for barrier between skin and pack



RAP11331-48 ..... 3" ..... Cs/48  
RAP11351-24 ..... 5" ..... Cs/24

## MAKEUP REMOVER WIPES

UNI151355 ..... Bag/50



## INSULIN SYRINGE

- 3/10cc

BD 320440 ..... 31G x 8mm (5/16").. Box/100  
BB 22-6503 ..... Comfort 31G x 5/16 .. Box/100



## TSK STERIGLIDE™ CANNULA

The biggest advantage of cannulas over needles is the reduced risk of bruising during filler injections. This advantage has resulted in a rapid adoption of the cannula technique worldwide, making cannulas the product of choice for selective procedures, e.g. the tear through region. Other advantages of cannulas include the ability to treat a larger area with only one injection point, better comfort, less trauma, and more acceptability by the patient. After TSK's first cannula introduction with the CSH Cannula, TSK introduced their award winning STERIGLIDE cannula in 2012. The STERIGLIDE outperforms any other cannula available and remains to lead the market as the golden standard.

- Patented Dome shaped tip, Nearest to tip filler delivery
- Proprietary surface treatment
- Clear side-port marking
- Introducer needles included
- ULTRA Thin Wall
- Safety PRC hub

TSK (All Sizes) ..... Box/20



Cannulas / STERIGLIDE™				
Code	Description	Inch	UTW*	Box
PRC-30025ISG	STERIGLIDE Cannula 30G x 25mm	1"	✓	20
PRC-27050ISG	STERIGLIDE Cannula 27G x 50mm	2"	✓	20
PRC-27038ISG	STERIGLIDE Cannula 27G x 38mm	1½"	✓	20
PRC-27025ISG	STERIGLIDE Cannula 27G x 25mm	1"	✓	20
PRC-25050ISG	STERIGLIDE Cannula 25G x 50mm	2"	✓	20
PRC-25038ISG	STERIGLIDE Cannula 25G x 38mm	1½"	✓	20
PRC-22070ISG	STERIGLIDE Cannula 22G x 70mm	2¾"	✓	20
PRC-22050ISG	STERIGLIDE Cannula 22G x 50mm	2"	✓	20

\* UTW = ULTRA Thin Wall

## WOODEN ROUND APPLICATOR STICKS

SSP4204-40-500..... 6" ..... Box/1000



## DISPOSABLE LIP BRUSH

Tightly spaces bristles for easy application.

SSP4017-LIP-25SP ..... Pkg/25



## MICRO APPLICATOR BRUSH

- Tips are lint-free and non-absorbent
- Ultrafine tips
- Dentistry, Beauty, Electronics and craft applications
- 4" handle with narrow bendable neck
- No spill dispenser

SSP4017-MAB-4UF..... Ultrafine, White ..... Box/100



## DISPOSABLE MASCARA BRUSH

- Full, lush head coats lashes evenly
- Nylon bristles

SSP4017-MB-20SP ..... 4" ..... Pkg/20



## DISPOSABLE CLEAR SPATULA

- Smooth Edge, flexible blade
- Used for mixing and applying gel or makeup

SSP4017-SPATULA-20SP ..... Pkg/20



## MANICURE/PEDICURE STICKS

- Square tapered on both ends
- Gently grooms cuticles and cleans under nail tips

SSP4202-60-100SP ..... 7" ..... Bag/144





# INDEX

## A

Abdominal Pads	120
Adhesive Remover	118
Adhesive, Tissue	73
Administration Sets, I.V.	79
AED	38
Air Freshener	45
Airways	113
Alcohol Dispenser	106
Alcohol, Isopropyl	106
Allergy Trays	92
Amniotic Test	88
Aneroids	24
Anesthetic	75-76
Anoscopes	101
Antiseptic Solutions	77,84
Applicators, Cotton Tipped	103
Aprons, Disposable	99
Autoclave	47
Autoclave Cleaner	50
Autoclave Tape	50
Autoclave Pouches & Wraps	49

## B

Bandages	116
Bathroom Tissue	44,45
Batteries	27
Baxedin	77
Benzalkonium Wipes	77,103
Biohazard Bags	89
Biological Indicator	48
Biopsy Kit, Endometrial	88
Biopsy Punches	69
Bleach	46
Blades, Scalpels	68
Blood Collection	81-83
Blood Drawing Chair	90
Blood Lancets	90
Blood Pressure Cuffs	23
Blood Pressure Devices, Digital	20-22
Blood Pressure Devices, Manual	24
Blood Pressure Tubing	23
Bags, Recycling	45
Blunt Fill Needle	95
Bouffant Caps	111,127
Brush, Surgical Instrument	49
Bulbs	27
Butterfly Closures	116

## C

Cap, Surgeons	111,127
Carts	16,17
Cast Padding	125
Casting	124,125
Catheters, I.V.	78,79
Cauteries	13,14
Centrifuge	90
Cervical Dilator	102
Cervical Scraper	102
Chair, Procedure	5
Chair, Transport	18
Chairs	7
Chlorhexadine Pads & Swabs	84
Cleaning Solutions, Instrument	40,41
Cold/Hot Packs	122
Condoms	89

Conform Stretch Gauze	121
Connecting Tube	113
Cotton Balls	107
Cotton Tipped Applicators	103
COVID-19 Antibody	87
COVID-19 Antigen	86,87
CPR Faceshield	114
CPR Pocket Resuscitator	114
Cryosurgery	70
Cuffs, Blood Pressure	23
Cups, Drinking	107
Cups, Medicine	105,107
Curettes	61,62
Cytology Collection Brush	88
Cytospray	89

## D

Defibrillators	38
Dermal Tips	14
Detergent, Laundry	46
Dextrose, 5%	78
Diagnostic Set	24,25
Dental	126-129
Dish Soap	46
Disinfectant Solutions	39,40
Disinfectant Spray/Wipes	39,40
Doppler	31
Drapes, Exam	98
Drapes, Surgical	76
Dressing Tray	74,120
Drinking Cups	107
Drugs of Abuse Test Kit	88

## E

Ear Basin	58
Ear Cleaning	58,61,62
Ear Curettes	57,61,62
ECG	33-35
ECG Electrodes	36
ECG Prep Tape	37
Elastic Bandage	121
Elastic Net Bandage	122
Electrode Skin Prep Pad	37
Electrodes	36
Electrostatic Sprayer	46
Electrosurgical Generator	13,14
Emergency Responder Kit	114
Emesis Basin	105
Endometrial Biopsy	88
Epinephrine Injection	115
EpiPen	115
Exam Tables	4-7
Eye Drops	76
Eye Glasses, Protective	108
Eye Pads	120
Eye Wash Station & Solution	106

## F

Face Mask	108,109
Face Shield	108
Facial Tissue	44
Filter, Breathing	113
Finger Ring Cutter	57
Finger Splints	125
First Aid Kits	120
Flourescein Strips	76
Fridge	15

## G

Garbage Bags	45
Garbage Can	45
Gauze	119
Glasses, Protective	108
Glove Box Holder	64
Gloves	63-67
Glucose Testing	86
Glustitch Tissue Adhesive	73
Gowns, Exam	97
Gowns, Isolation	110

## H

Hand Sanitizer	42
Hand Towels	44
Height Measuring	10,11
Holter Recorders	35
Hot/Cold Packs	122
Hydrogen Peroxide	106
Hyfrecator	13,14

## I

I.V. Access Cannula	84
I.V. Administration Sets	79
I.V. Catheters	82-83
I.V. Connector	84
I.V. Extension Sets	79
I.V. Pole	84
I.V. Solutions	78,80
Indicator Strips, Steam	49
Indicators, Biological	48
Integrators	49
Instrument Stand	17
Instrument Tip Protectors	51
Instrument Cleaning Supplies	62
Instrument Disinfection/Cleaning	40,41
Instrument Lube	62
Instrument Trays	51
Instruments	52-62
Instruments, Sterile	59,60
Insulin Syringes	96
Iodine Cup	105
Irrigation Solutions	80
Isolation Gowns	110
IUD Insertion Kits	74,102
IUD Hook	102
IUD	102

## J

Jars, Sundry	104
--------------	-----

## L

Laceration Trays	74
Lactated Ringers	78
Laundry Detergent	46
Lancets	90
Lights, Exam	11,12
Lubricating Jelly	104

## M

Marcaine	75
Masks, Oxygen/Nebulizer	112
Mayo Stand	17
Medicine Cups	105,107

# INDEX

## N

N95 Fit Test Kit	109
N95 Masks	109
Nasal Cannula	112
Nebulizer Masks	112
Needles & Syringes	91-96
Neuropen	33
Non-Adherent Pads	119
Normal Saline, I.V.	78
Normal Saline, Irrigation	124
Normal Saline, Respiratory	115

## O

Ointment Dressing	123
Opsite Dressing/Spray	123
Ophthalmoscopes	25
Otoscope Speculum	26
Otoscopes	25,26
Oxygen Masks	112
Ovulation Test	88

## P

Packing Strips	120
Paper Towels	44
Patient Monitors	20-22
Peak Flow Meter	112
Pediatric Scales	8
Penlight	60
Pillowcase	99
Pillows	99
Prefilled Syringe	96
Pregnancy Test	87
PRN Luer Lock I.V. Adapters	84
Probe Covers, Thermometer	30
Probe Covers, Ultrasound	32
Proctology Swabs	103
Pulse Oximeters	29
PVP Preps, Swabs & Solutions	77

## Q

Q-Style I.V. Connector	84
Q-Style Vial Adapter	96
Qualitative Fit Test Kit	109

## R

Rapid Test	86,87
Razors	37,73
Refrigerators	15
Resuscitator	113

## S

Sanitary Pads	99
Sanitizer, Hand	42
Scales	8-10
Scalpel Handle	68
Scalpels	68
Scrub Brush, Nylon	77
Sharps Containers	95
Shoe Covers	110,127
Shorts, Exam	99
Sigmoidoscopes	101
Silver Nitrate Applicators	103
Skin Closures	74
Skin Markers	74
Staple Removal	74
Skin Stapler	74

Smoke Shark	14
Soap, Hand	43
Solu IV Swabs	84
Spa & Beauty	131-133
Specimen Containers	89
Speculum, Otoscope	26
Speculum, Vaginal	100-101
Spirometer	33
Sponge Basin	104,105
Sprayer, Disinfection	46
Sterile Instruments	59,60
Sterile Water, Injection	80
Sterile Water, Irrigation	80
Sterilization Pouch Rack	51
Sterilization Pouches & Wraps	49
Sterilizers	47
Stethoscopes	27,28
Stockinette Tubular Bandage	121
Stools	8
Stopcocks	79
Straw, Drinking	107
Strep-A Test Kit	86,87
Stretcher	18
Suction Machines	115
Sundry Jars & Racks	104
Surface Cleaner	39
Surface Disinfectant	39
Surgeons Cap	111,127
Surgical Mask	108,109
Surgilast Bandage	122
Surgitube Bandage	122
Suture Removal Tray	74
Sutures	71,72
Swabs, Proctology	103
Syringe, Prefilled	96
Syringes & Needles	91-96
Syringes, Insulin	96

## T

Table Paper	98
Table, Treatment	4-7
Table, Exam	4-7
Tablets, Ultrasonic Cleaner	48
Tape Measure	10
Tape, Surgical	117,118
Tattoo Supplies	130
Tegaderm Dressing	84,123
Tensor Bandages	121
Thermometer, Fridge	16
Thermometers	30
Tissue, Facial	44
Tissue, Toilet	44,45
Tongue Depressors	103
Topical Anaesthetic	76
Tourniquets	84
Towels, Medical	98
Transport Chair	18
Treatment Carts	16-17
Treatment Table	4-7
Triangular Bandage	120
Tuning Forks	58

## U

Ultrasonic Cleaners	48
Ultrasound Supplies	32

Underpads	99
Urinalysis	85
Urine Analyzer	85
Uterine Dilator	102
Uterine Sound	101

## V

Vaccine Refrigerator	15
Vacutainer Products	95
Vaginal Illumination System	100
Vaginal Speculum	100,101
Veritor Flu Test Kit	86
Vision Screener	26
Vital Signs Monitor	20,22

## W

Wall Systems	19
Warmer, Ultrasound Gel	32
Wart Treatment	70
Water Distilled	50
Wheelchair	18
Wipes, Disinfectant, Surface	39

## X

Xylocaine	75
-----------	----

## Z

Zoll AED	38
----------	----



1.800.232.7286

canadamedical.net | orders@canadamedical.net

**TERMS AND CONDITIONS:**

Net 30 Days on approved accounts.

Taxes extra where applicable.

Free Shipping on most orders (Call us for more details).

Returns require an authorization #

Returns may be subject to restocking fee.

No returns on special order products.



The reproduction of any part of this catalogue is strictly prohibited unless accompanied by written permission from Canada Medical Ltd.  
Errors and Omissions Excepted. Prices are subject to change without notice.