

Imaging Table with Fowler Back and Drop Window



MODEL 80072

- Two-piece top with manual, pneumatic, adjustable backrest and window drop
- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser, cutter and positioning straps included
- Optional side rails & casters available
- 34" wide model available
- Model 80072: 500 lbs load capacity under normal use
- Model 80072-X: 600 lbs load capacity under normal use

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

1



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

80072 SPECIFICATIONS

	length	width	height range
overall	72"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
backrest	38"	27"	
window	10"	10"	

80072-X Extra-Wide 34" width model



Upholstery Colors



OPTIONS

Side Rails (pair)

- 3 height positions
- Plunger height adjustment

098



Power Backrest

- Smooth, power control raises to 80° and lowers to flat

011



Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087

