

# Tabletop Portable Model 300

The Gomco 300 is a heavy duty general use tabletop aspirator for hospitals, clinics and physician's offices. It generates up to -559 mm Hg of vacuum pressure and a 29 lpm flow rate at open flow.

## Key Features

- Maintenance-free and oil-free diaphragm pump ensures maximum reliability
- Similar performance and quality of a mobile aspirator in a tabletop design
- Easy to store, the 300 takes up less than one square foot of space
- Disposable hydrophobic bacteria filter helps prevent fluid and aerosol contamination



**01-22-0300**

## Specifications

Vacuum Range:	-25 to -559 mm Hg
Flow Rate:	29 lpm at open flow
Sound Level:	Less than 60 dBA
Pump Type:	Oil-less diaphragm
Electrical Cord:	6' hospital grade
Dimensions:	13" W x 9" D x 12" H (30 x 23 x 30 cm)
Weight:	19 lbs. (9 kg)
UL/ETL Listed:	No
Warranty:	1 year
Regulatory Approvals:	FDA, HC

Model	Voltage	Glass Bottle	Disposable Canister
01-22-0300	115V/60 Hz	1 - 600 ml	
01-22-3001	115V/60 Hz		1 - 1100 ml
01-23-0300	220V/50 Hz	1 - 600 ml	

*All models come with disposable suction tubing and 3 disposable hydrophobic bacteria filters.*

Catalog #	Description
01-90-2000	Disposable Suction Tubing Kit, one 18" & one 72"
S610100	Disposable Suction Tubing Kit, one 15" & one 72"
S615473	Disposable Suction Tubing Kit, one 15"
S615725	Disposable Suction Tubing Kit, one 72"
01-90-2454	Glass Collection Bottle, with rubber cap, short & long tube assembly, 600 ml
01-90-3100	Disposable Hydrophobic Bacteria Filter, 1/8" NPT x 1/4" hose barb (3/pkg.)
01-90-3695	Disposable Collection Canister, stem inlet & hydrophobic filter, 1100 ml (12/case)

*Disposable collection canisters are shown on page 21 and 22.  
Additional case quantities are available.  
Collection bottles are shown on page 23.  
Suction accessories are shown on page 24 and 25.*