



Baxter

Welch Allyn KleenSpec Single-Use Vaginal Specula and Cordless Illumination System

Our Latest Women's Health Technology to Outfit Your Exam Room.

Welch Allyn KleenSpec Cordless Illumination System and Single-Use Vaginal Specula — in a Variety of Sizes for Your Patient Base.

590 Series KleenSpec Single-Use Vaginal Speculum

Smooth, molded acrylic designed to help prevent pinching or binding of tissue with patient comfort in mind



#59004
Large



#59001
Medium



#59000
Small



#590XS
Extra Small

Look for the distinctive **Welch Allyn** blue, green, white, and yellow coloured handles

HELP REDUCE CROSS-CONTAMINATION

- Using single-use disposable speculum vs. metal speculum helps eliminate the time and expense of reprocessing specula and cross-contamination risk between exams

DESIGNED FOR QUALITY & EFFICIENCY

- Wide, ergonomic handle designed with clear acrylic allows complete visualization of the vaginal canal
- A strong grip and a yoke that's more flexible than low-cost substitutes



KleenSpec Cordless Illumination System

- LED technology offers a white, bright light and true tissue colour during an exam
- No LED bulb replacement necessary
- Features a cordless design so you'll never have to worry about cords getting in your way, breaking, or needing to be cleaned, giving you a full range of motion
- Fits every **Welch Allyn** size speculum
- Uniform spot—no dark or hot spots—for enhanced visualization of the exam area
- LED light source will not generate damaging heat and has a continuous on time of 80 minutes
- A lithium ion battery provides a recharge time of five hours



Ordering Information

KleenSpec Vaginal Specula — Premium 590 Series

590XS	KleenSpec 590 Series Premium Disposable Vaginal Specula, X-Small (24/box; 4 boxes/case)
59000	KleenSpec 590 Series Premium Disposable Vaginal Specula, Small (24/box; 4 boxes/case)
59001	KleenSpec 590 Series Premium Disposable Vaginal Specula, Medium (24/box; 4 boxes/case)
59004	KleenSpec 590 Series Premium Disposable Vaginal Specula, Large (18/box; 4 boxes/case)

KleenSpec Cordless Illumination System

80010	Complete KleenSpec Cordless Illuminator with Charging Station
80000	KleenSpec Cordless Illuminator
74010	Charging Station for Cordless Illuminator
59010	Disposable Sheath for Cordless Illumination System (100/box; 5 boxes/case)

Green Series Exam Light IV — Another LED innovation from **Welch Allyn**

48810	Green Series Exam Light IV with Table/Wall Mount
48950	Mobile Stand for Green Series Exam Light IV
52640	Disposable Sheath for Exam Light III and Green Series Exam Light IV (25/box, 5 boxes/case)

Connect with your local distributor or
Baxter sales representative to learn more.

For the safe and proper use of the devices referenced within,
refer to intended use, warnings, precautions, and complete
instructions for use or the appropriate manual.

Baxter.ca

Baxter, Green Series, KleenSpec and Welch Allyn are
trademarks of Baxter International Inc. or its subsidiaries.
Any other trademarks, product names or brand images
appearing herein are the property of their respective owners.
CA-FLC94-260001 (v1.0) 07/2026