



# CLINITEK Status Connect System Product Specifications

## Transmit a wealth of clinical information with a connected urinalysis solution

[siemens-healthineers.ca/urinalysis](http://siemens-healthineers.ca/urinalysis)



### Semi-quantitative Tests Measured

**Routine urine tests:** albumin, bilirubin, blood, creatinine, glucose, ketone, leukocytes, nitrite, pH, protein, specific gravity, and urobilinogen  
**Kidney disease test:** albumin-to-creatinine ratio (ACR), and protein-to-creatinine ratio (PCR)  
**Pregnancy test:** Human Chorionic Gonadotropin (hCG)

### General Specifications

<b>Depth:</b>	10.7 in. (27.2 cm)
<b>Width:</b>	6.7 in. (17.1 cm)
<b>Height:</b>	7.5 in. (18.5 cm)
<b>Weight:</b>	5.0 lbs (2.3 kg)
<b>Electrical rating:</b>	100V–240V AC, 50–60 Hz (with in-line lead)
<b>Battery operation:</b>	Available for standalone analyzer only
<b>Ambient temperature range:</b>	18–30°C (64–86°F)
<b>Humidity range:</b>	18–80% relative humidity, non-condensing
<b>Compliance:</b>	CLIA-waived, CE-marked, UL, EMC
<b>Calibration:</b>	Automatic, self-calibrating

### Communication Capabilities

**Communication protocol:** HL7, POCT1-A2  
**Data transmission options:** Ethernet or serial connection  
**Interfaces with leading point-of-care data management solutions:** POCcelerator™ Data Management System, RAPIDComm® Data Management System, MAS RALS® Plus, Telcor Quick-Linc®

### Memory

**Patient test results:** 700  
**QC test results:** 200  
**Operator IDs:** 700

### Auto-Checks™ Technology\*

Identify strip type  
Detect humidity exposure†  
Check for common sample interferences‡

### Quality Control Management Features

QC testing and reporting  
QC reminders or mandatory prompt  
QC lock-out

### Barcode Data Entry

Patient ID, operator ID and name, QC control type, lot number, strip lot/expiration, and hCG lot/expiration

### Barcode Symbology Supported

Code 93, Code 39 (with and without check digit), Code 128, Coda Bar, Interleaved 2 of 5 (with and without check digit)

### Operator Security Features

Define up to 700 operators  
Operator lock-out to prevent unauthorized use  
Customize access level for each operator for the following functions:

- Performing patient tests
- Performing QC tests
- Recalling results
- Modifying system settings

\*Available only when using Siemens Healthineers strips with ID or color bands.

†Only available with Siemens Healthineers test strips that have the leukocyte pad.

‡Not available in the U.S.

Note: Current CLINITEK Status analyzers can be upgraded for connectivity. Ask your sales representative for information on upgrade kits.

At Siemens Healthineers, our purpose is to enable healthcare providers to increase value by empowering them on their journey toward expanding precision medicine, transforming care delivery, and improving patient experience, all made possible by digitalizing healthcare.

An estimated 5 million patients globally benefit every day from our innovative technologies and services in the areas of diagnostic and therapeutic imaging, laboratory diagnostics, and molecular medicine, as well as digital health and enterprise services.

We are a leading medical technology company with over 120 years of experience and 18,000 patents globally. Through the dedication of more than 50,000 colleagues in 75 countries, we will continue to innovate and shape the future of healthcare.

Auto-Checks, CLINITEK, CLINITEK Status, POCcelerator, RAPIDComm, and all associated marks are trademarks of Siemens Healthcare Diagnostics Inc., or its affiliates. All other trademarks and brands are the property of their respective owners.

Product availability may vary from country to country and is subject to varying regulatory requirements. Please contact your local representative for availability.

---

**Siemens Healthineers Headquarters**

Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen, Germany  
Phone: +49 9131 84-0  
siemens-healthineers.com

**Legal Manufacturer**

Siemens Healthcare Diagnostics Inc.  
Laboratory Diagnostics  
511 Benedict Avenue  
Tarrytown, NY 10591-5005  
USA  
Phone: +1 914-631-8000

**Local Contact Information**

Siemens Healthcare Limited - Canada  
1577 North Service Road East  
Oakville, ON, L6H 0H6, Canada  
Phone: +1-905-465-8000  
siemens-healthineers.ca