



Here for L.I.F.E.

Leadership. Innovation. Focus. Excellence.

TECHLEM PREMIER SERIES 4500GYN STRETCHER

Part Number: 465-0300-0011

TECHNICAL SPECIFICATIONS

- ✓ Sliding top deck - converts from stretcher to chair
- ✓ Infinite chair positioning - Supine to full chair
- ✓ Pnuematic assist head section..... 0 - 90°
- ✓ Pnuematic assist leg section..... 0 - 105°
- ✓ Trendelenburg/Reverse..... +18/-18°
- ✓ Head, dual side hydraulic height adjustment controls
- ✓ Integrated, adjustable, telescoping heel stirrups
- ✓ Four corner access brake/steer pads with Techlem "Posi-Lock" Braking System
- ✓ Techlem "Accu-Steer" Fifth Wheel Steering System
- ✓ 8" (20cm) All Surface, Tente Four Wheel Locking Caster System
- ✓ Thermoform Base with Built-in Storage Area & Oxygen Tank Holder
- ✓ Built-in Recessed Storage Tray Under Fowler (Head) Section
- ✓ Head End Permanently Telescoping Folding I.V. Pole
- ✓ Head End Push I Pull Handles - Drop Down Design
- ✓ Techlem "Viscolastic" 2 Piece Sectional Mattress
- ✓ 4 IV plus 4 Accessory Attachment Receptacle Wells
- ✓ Total Perimeter Wraparound Bumper with Co-ordinated Accents Available in 14 Colors**

Dimensions :

- ✓ Height Adjustment..... 24" to 36.5" (61 cm to 93 cm)
- ✓ Total Dimensions..... 78" X 32" (198 cm X 81 cm)
- ✓ Side Rails..... 48" X 12" (123 cm X 32 cm)
- ✓ Sectional Mattress Dimensions..... 74" X 26" X 4" (193 cm X 73 cm X 10 cm)

