



Magnifier

UV Light

The ideal light for numerous applications.

Designed for a variety of uses, the Burton UV (Woods) Light is also available in a model with a magnifier and a model with a magnifier and white and UV bulbs. With handheld flexibility, large magnification and a five-year limited warranty, this UV light provides great value for your facility.

Features

- Special fluorescent bulbs produce UVA at approximately 365 nm wavelength
- 4-light magnifier models include optically ground and polished glass 5 diopter magnifier lens with 2.25X total power magnification



burton®

UV Light product dimensions



UV502 & UV503



UV501

Burton's UV lights are depended on worldwide for performance and durability.

Summary of technical data

UV LIGHT

| | |
|------------------------|-------------------------|
| Optics: | 5 Diopter |
| Viewing Area Diameter: | 8" (20 cm) |
| Light Source: | UV or white |
| Rated Life of Bulb: | 6,000-hour |
| Power: | 120v |
| Product Warranty: | 5-year limited warranty |

Ordering Information

UV LIGHT

| Model | Description |
|-------|---|
| UV501 | UV Light, 2 UV Bulbs |
| UV502 | UV Light with Magnifier, 2 UV Bulbs |
| UV503 | UV Light with Magnifier, 2 White and 2 UV Bulbs |

ACCESSORIES

| Model | Description |
|-----------|------------------------------------|
| 1007090 | UV Light Wall Bracket |
| 0001127PK | UV Replacement Bulbs, White, 4/Box |
| 1003073PK | UV Replacement Bulbs, UV, 4/Box |



© 2021 Burton Medical. All rights reserved. Burton Medical reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

burton®

Burton Medical, LLC
2300 W Windsor Ct., Suite C
Addison, IL 60101 USA
Phone (800) 444-9909, Fax (800) 765-1770
Burtonmedical.com