

MICRO-TOUCH® Nitrile Sterile

P e r f o r m a n c e S e r i e s

Nitrile, sterile, powder-free medical exam glove featuring extended cuff, packed as a pair

- **Barrier protection:** Packed as a pair of gloves per pack, these sterile gloves are made from a synthetic, nitrile formulation that provides high elasticity with superior barrier protection and puncture resistance
- **Sterile:** With a minimum 25kGy to provide sterility, they are an ideal solution for a range of sterile medical exam (non-surgical) procedures
- **Extended protection:** The 11.6" extended cuff provides additional protection
- **Secure grip:** Textured fingers provide extra grip and reliability
- **Reduced allergy risk:** Nitrile material prevents latex Type I allergic reactions
- **Specialized protection:** Tested and approved for use with chemotherapy drugs*




*see product packaging for details

Recommended For

- Tested for use with chemotherapy drugs
- Sterile procedures outside of the O.R.
- Extra protection over the wrist and arm
- Sustained, high-risk examination procedures
- Protection from Type I latex allergy in HCW's or patients

FEATURED TECHNOLOGIES



MICRO-TOUCH™ Everyday, Everywhere Protection

MICRO-TOUCH™ premium medical examination gloves are scientifically designed to complement every need of medical professionals. Sophisticated design features include Ansell ergonomic technologies that ensure comfort during extended procedures, wet grip innovations that prevent harmful mishandling, and micro-thin films for maximum sensitivity. Caregivers who are fully prepared and comfortably protected by versatile MICRO-TOUCH styles have the peace of mind to focus on what matters most: the health and protection of themselves and their patients.

MICRO-TOUCH® Nitrile Sterile

Performance Series

PRODUCT INFORMATION

Material	Nitrile	Not Made From Natural Rubber Latex	Yes
Color	Blue	Cuff Length	Extended
Powder Content	Powder-Free	External Glove Surface	Textured Fingers
Freedom from Holes (Inspection level I)	1.5 AQL	Finger Thickness (mm/mil)	0.16/6.3
Palm Thickness (mm/mil)	0.11/4.3	Allergy Prevention	Latex (Type I)
Available Sizes	S (6.5 - 7), M (7.5 - 8), L (8.5 - 9), XL (9.5 - 10)	Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 and US FDA cleared
Sterile	Yes	Antistatic	Not Tested
Glove Length (mm/inches)	290 / 11.6	Product Segmentation	Performance Series

ORDERING INFORMATION

Size	S (6.5 - 7)	M (7.5 - 8)	L (8.5 - 9)	XL (9.5 - 10)
Product Code	6034151	6034152	6034153	6034154

PACKAGING AND STORAGE

Packaging	1 pair per pouch 50 pairs per dispenser 4 dispensers per case 200 pairs per case
Shelf Life	5 Years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.

STANDARDS AND CERTIFICATIONS

ASTM D6319, ASTM D6978, CE 0493, CE 2797, Medical Device Regulation (EU) 2017/745 Class I, EN 16523-1, EN 455 1-4, EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 13485, ISO 9001, ISO 10993-10:2002, UKCA

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 1-844-494-7854
Email:
insidesalesus@ansell.com

Latin America & Caribbean
Region
BR +55 (11) 3356-3100
latam.medical@ansell.com
MX: +52 442 296 2050
latam.medical@ansell.com
CO: +57 1 288 3247
latam.medical@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email:
apac.medical@ansell.com

Australia & New Zealand
Tel: +61 3 9270 7270
Email:
protection@ap.ansell.com