

QUALITY PRODUCTS | VALUE PRICED | TIMELY DELIVERY



2022 PRODUCTS CATALOGUE



www.canadamedical.net



Letter from the President

Dear Friends,

For over 30 years, Canada Medical has partnered with thousands of healthcare professionals and institutions with the goal of increasing the quality of patient care by reducing the cost of their medical supply and equipment requirements without impacting service. Our team's dedication to continuous improvement, quality and service enables Canada Medical to serve our customers year after year.

At our core, we are committed to:

- Distributing branded medical supply merchandise to our loyal customers in our markets
- Servicing our customers' and vendors' needs while exceeding their expectations
- Differentiating ourselves by providing the best customer service to quality customers

Our habits of unrivalled customer service will continue, even as we push to find more ways to

build in quality in a way the customer desires. Relationship building with our customers, our vendors, our teammates and the communities we are a part of has been and will continue to be part of what differentiates us from our competition.

We have also realigned our company structure to a team-based model, working to continuously improve our distribution work flow. Our priority is to continuously improve our service to our customers and vendors while continuing to build on our knowledge and experience.

We strive to be the company that best understands and satisfies the products and services of the medical supply industry.

We look forward to partnering with you and your organization in the future.

Sincerely,
Carla Rodych

Connect with us:



Canada Medical Ltd.

Head Office:
16503 116 Ave NW
Edmonton, AB T5M 3V1
Local: 780.465.2020
Local Fax: 780.468.0991
Toll-Free: 1.800.232.7286
Fax Toll-Free: 1.800.278.8503

orders@canadamedical.net
www.canadamedical.net



*Visit canadamedical.net for more details

PRODUCT CATEGORIES

EXAM ROOM EQUIPMENT & FURNITURE

PAGE 4 - 18

DIAGNOSTICS

PAGE 19 - 38

INFECTION PREVENTION

PAGE 39 - 43

HOUSEKEEPING / JANITORIAL

PAGE 44 - 46

STERILIZATION

PAGE 47 - 51

INSTRUMENTS REUSABLE

PAGE 52 - 57

INSTRUMENTS DISPOSABLE

PAGE 58 - 62

SURGERY AND MINOR PROCEDURES

PAGE 63 - 77

INFUSION & BLOOD COLLECTION

PAGE 78 - 84

LABORATORY

PAGE 85 - 90

NEEDLES & SYRINGES

PAGE 91 - 97

PATIENT CARE & EXAM ROOM SUPPLIES

PAGE 98 - 107

PROTECTIVE APPAREL

PAGE 108 - 111

RESPIRATORY & OXYGEN

PAGE 112 - 115

SKIN & WOUND CARE

PAGE 116 - 125

DENTAL

PAGE 126 - 129

TATTOO

PAGE 130

SPA AND BEAUTY

PAGE 131 - 133



RITTER 224 BARRIER-FREE® EXAMINATION CHAIR




Ritter 224 & 225 Barrier-Free chairs are designed to offer low and high positioning heights so you can provide the best care to all patients, whether treating the elderly, expectant mothers or patients with disabilities. Ritter Barrier-Free Examination Chairs by Midmark adjust in height from 18-37 inches, improving accessibility to the exam chair and reducing the likelihood of staff injury

- Power Base
- 224 Model Air Spring Assist Back
- 225 Model Power Back
- 500lb Patient Weight Capacity
- Foot Control

RIT224-001	Base with Manual Back
RIT224-002	Base with Manual Back, Pelvic Tilt and Drawer Heater
RIT225-002	Base with Power Back
RIT225-003	Base with Power Back, Clean Assist Roller System
RIT002-2009-(Colour Code)	Seamless Top

RITTER® STANDARD COLOURS



RITTER® 204 MANUAL EXAMINATION TABLE



What makes Ritter the most purchased brand of exam tables in the industry? It's our willingness to research, study and listen so we can meet your needs—and the needs of your patients. If you need an exam table with easy access for patients to make them feel secure, three side drawers for optimal supply management, modern styling to coordinate with your office and an additional weight capacity up to 500 pounds, the Ritter 204 Manual Exam Table is the examination table for you.

- One piece seamless drawers
- Stylish and easy to clean
- Removable top system
- Enlarged footstep
- Storage access and volume
- Seamless upholstered top
- Manual back with air spring assist



RIT204-011	Base with Manual Back
RIT002-10138-(Colour Code)	Seamless Top



626 BARRIER-FREE® EXAMINATION CHAIR

The Midmark 626 is the first chair of its kind designed to promote proper patient positioning for a more accurate blood pressure measurement. Connect the IQvitals® device to the Midmark 626 Barrier-Free Exam Chair with the integrated IQscale® to weigh patients quickly and discreetly from the comfort of the exam chair. Then, send data to the EMR. The Midmark IQconnect® framework provides a seamless connection from your devices to your EMR system. In other words, it allows you to operate all of your devices and transfer data to the EMR from one, single pane of glass—your computer screen.

Standard Features:

- Clean Assist™ Roller System The retractable roller base allows you to easily and safely move the chair for cleaning.
- Hand and Foot Control
- Stirrups and Pelvic Tilt
- Drawer Heater
- Active Sensing Technology pauses table movement when disrupting impact is detected
- Quick Exam® and Return to Home Quick Exam elevates the chair to the exam position with one touch of the control, while Return to Home returns the chair to its home position
- Patient Weight Capacity: 650 lb
- 15.5" Low - 37" High



244 BARIATRIC EXAM TABLE

Barrier-Free® Accessibility - Safe and Efficient Care for your Patient

- 850lb Weight Capacity
- 18" Low Height. By lowering the seat to 18" from the floor, the Ritter 244 performs this sometimes difficult task. It becomes a very simple routine, requiring minimal assistance from you and your staff. The patient's dignity is maintained and the likelihood of injury is greatly reduced, providing efficient, effective patient care.
- 34" High Height - The Ritter 244 raises the seat to a height of 34", designed to bring the patient to you while helping to minimize the stress and strain on your back. By having the patient in the right position, the exam can be completed properly and efficiently



630 HUMANFORM™ PROCEDURES CHAIR

When performing in-office procedures, physicians need a table that can achieve proper medical positioning to ensure ready access to the patient and a clear view of the surgical site. The Midmark 630 HUMANFORM™ procedures tables gives physicians the flexibility to position patients where they need them for a wide range of procedures.

The Midmark 630 HUMANFORM™ procedures table is the next generation of procedures tables designed by the people who use them, not the technology behind them. The first of its kind, the Midmark 630 HUMANFORM™ procedures table is the new benchmark in patient and provide centered performance. Designed to move with the human body, this table enhances patient comfort while allowing physicians to remain in a proper working position throughout procedures.

Premium Comfort System utilizes unique contoured upholstery and an adjustable footrest to cradle the patient for a soothing, chair-like feel. Barrier-Free® Design addresses the needs of a changing patient demographic with a best-in-class, 18" low seat height and 650lb weight capacity



HI-LOW POWER TWO SECTION EXAM TABLE

Two Section Table

32" Maximum Height
 17.5" Minimum Height
 27" Width
 78" Length
 39.5" Head Section
 24.5" Seat Section



SOLAYR1008

High-Low Exam Table design provides an 17.5" profile in low position to allow level transfer from standard wheelchair heights, complete with upright seat back.

Pneumatic backrest is adjustable from either side and raises 80 degrees from flat Chair-to-Table design takes patient from consultation to examination level Hands-free low voltage foot control Heavy duty, all welded powder coated steel frame Top sections have welded steel support frames. 2" Firm density foam tops Paper Roll holder (18" to 21") 500 lb weight capacity

HI-LOW POWER THREE SECTION EXAM TABLE

Three Section Table

32"Maximum Height
 17.5" Minimum Height
 27" Width
 78" Length
 39.5"Head Section
 24.5" Seat Section
 13" Foot Section



SOLAYR1009

SOLAYR1009E.....with stirrups

EXAM ROOM BASIC TREATMENT TABLE

- Table Dimensions (24.5" W x 75" L x 28.5" H)
- 350lb weight capacity
- All steel construction with strong 4-leg base
- Nylon floor levelling glide helps prevent scratches on floor

SOL2001



EXAM ROOM BASIC TREATMENT TABLE WITH DRAWERS

- Table Dimensions (24.5" W x 75" L x 30" H)
- 350lb weight capacity
- Mechanical 3-position headrest with adjustable degrees of 25°, 40° & 60°
- 2" thick soft bedded black vinyl top
- Paper roll holder
- Bottom storage shelf
- Two front drawers 20" wide x 18.5" deep x 3.5" high

SOL22401



RITTER SIDE CHAIRS

The Midmark 280 is soft, comfortable seating for your patients in the exam room, waiting room and front office. The wall-saver kickback leg design allows the chair to be close to the wall without damaging it.



COMFORT

Contoured, upholstered seat and back section

STABILITY

Sturdy, powder-coated black steel tube frame with leveling feet

OPTIONS

Available with or without high-impact plastic chair arms

RIT280-001-(Colour Code)*Basic chair without arms

RIT280-002-(Colour Code)*Basic chair with arms

*See Page 4 for colour codes.



RITTER 282 BARIATRIC SIDE CHAIR



The Midmark 282 Side Chair can enhance your clinical workflow while complementing the style of your space with a range of upholstery color options.

COMFORT

Contoured, upholstered seat and back section

STABILITY

Sturdy, powder-coated black steel tube frame with levelling feet, a 30-inch wide seat and a 650-pound weight capacity

RIT282-001-(Colour Code)*Basic chair with arms

*See Page 4 for colour codes.





RITTER® AIR LIFT HEIGHT STOOL, 5 CASTER

These stools offer a black, plastic composite five leg caster base and hood and a single lever height control. The seat height adjusts pneumatically. Black vinyl upholstery is stocked for quick delivery. Other colours available with lead time.

- RIT273-001-(Colour Code)Stool with back
- RIT272-001-(Colour Code)Stool without back

*See Page 4 for colour codes.



STEP ON STOOLS

- AMG116-460Stool
- AMG116-470Stool with handrail
- AMG116-465Replacement rubber tips
- DUK4350Bariatric step stool



DIGITAL PEDIATRIC TRAY SCALE

Scale function that allows user to select and use only pounds or kilograms, reducing the risk of incorrect data entry

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- Display: 0.6" LCD
- Power Source: 9 Volt battery included
- Warranty: 2 Year limited warranty
- Carrying case: optional carrying case (order # 64769)
- Tray Size: 19 1/4" (w) x 13" (d) x 2 5/8" (h)

CON549KLPediatric scale



DIGITAL PEDIATRIC SEAT SCALE

Plastic, easy-to-clean seat limits infant movement while accurately calculating weight measurement

- Capacity: 50 lb / 23 kg
- Resolution: 0 - 20 lb / 0.2 oz; 20 - 50 lb / 0.5 oz; 0 - 9 kg / 5g; 9 - 23 kg / 10 g
- Power Source: 6 AA batteries included, optional 120V adapter (order # ADPT31)
- Connectivity: USB
- Seat Size: 24 1/4" (w) x 14 1/4" (d) x 18 1/2" (h)
- Warranty: 2 Year limited warranty

CON524KLDigital scale with seat

REMOTE DISPLAY DIGITAL FLOOR SCALE

Interfaces with Welch Allyn® Spot Vital Signs® LXi Monitor. The remote display provides versatility with the ability to be mounted on a wall or placed on a table.

- Capacity: 400 lb / 180 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size: 12 1/2" (w) x 12" (d) x 1 7/8" (h)
- Display: 1/2" LCD
- Power Source: 9V battery included, optional 120V power adapter (Order#: ADPT40)
- Connectivity: serial port
- Warranty: 2 Year limited warranty

CON349KLXDigital floor scale



LARGE PLATFORM FLOOR SCALE

- Capacity: 400 lb / 180 kg
- Platform Size: 11" (w) x 12 1/2" (d) x 3" (h)
- Graduation: 1 lb / 1 kg
- Warranty: 1 Year limited warranty

CON160KLPlatform floor scale



MECHANICAL BEAM SCALE

The easy-to-use height rod provides height measurements from 23" to 84".

- Capacity: 500 lb / 200 kg - with included counterweights
- Platform Size: 10 1/2" (w) x 14" (d) x 3 1/4" (h)/267mm (w) x 356mm (d) x 83mm (h)
- Height Rod: 23 5/8" - 84" / 60 cm - 213 cm
- Graduation: 1/4 lb / 100 g
- Warranty: 10 Year limited warranty

CON450KLBeam scale



EYE LEVEL DIGITAL SCALE

Scale has the ability to interface with a computer, EMR software or other electronic device via various Health-o-meter® Professional Connectivity Solutions. The easy-to-use height rod provides height measurements from 30" to 84"/76cm - 214cm

- Capacity: 500 lb/ 220 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size: 13 3/4" (w) x 16 1/2" (d) x 2 3/8" (h)
- Wheels: (2) Wheels for ease of movement
- Power Source: 6 AA batteries included, optional 120V adapter (order # ADPT31)
- Warranty: 2 Year limited warranty - eligible for ScaleSurance extended warranty

CON500KLDigital scale

CON500HB.....Handlebars for 500 series scale

CON500KLHBDigital Scale with Handlebars





HEAVY DUTY EYE LEVEL DIGITAL SCALE W/ SWIVEL DISPLAY

Scale has the ability to interface with a computer, EMR software or other electronic device via various Health-o-meter® Professional Connectivity Solutions. Eye-level rotating display can be read from almost any angle and allows for ease of use and patient privacy. The easy-to-use height rod provides height measurements from 24" to 84"/60cm – 214cm. Scale functions include the ability to calculate the patient's Body Mass Index (BMI).

- Capacity: 600 lb / 272 kg
- Resolution: 0.2 lb / 0.1 kg
- Display: 1" LCD, 350° Swivel
- Connectivity: USB
- Power Source: 110-240V adapter included, 6 C-cell batteries not included
- Warranty: 2 Year limited warranty



CON597KLScale with height rod



DIGITAL PLATFORM SCALE WITH HEIGHT ROD

EMR Connectivity Scale has the ability to interface with a computer, EMR software or other electronic device via various Health-o-meter® Professional Connectivity Solutions. Calculates Body Mass Index (BMI).

- Capacity: 1000 lb / 454 kg
- Resolution: 0.2 lb / 0.1 kg
- Platform Size: 22" (w) x 15 3/4" (d) x 2 1/2" (h)
- Height Rod: Included
- Display: 3/4" LCD, 180° Swivel
- Wheels: (2) Wheels for ease of movement
- Power Source: 110-240V adapter included, 6 AA batteries (not included)
- Connectivity: USB
- Warranty: 2 Year limited warranty



CON1100KLHRScale with height rod



HEIGHT MEASURING TAPE, WALL MOUNT

- The space-saving solution with roll-up mechanism.
- Easy fixing on wall with just one screw.
- Display window in the head piece, it's easy to read the height.
- The durable tape is made of metal

SOLHX400CHeight Measuring Tape



RETRACTABLE TAPE MEASURE

The space-saving solution with roll-up mechanism. Easy fixing on wall with just one screw. Display window in the head piece so it's easy to read the height. The durable tape is made of metal.

AMG116-860Measuring Tape

WALL MOUNTED MECHANICAL MEASURING ROD

The Seca 216 has a variable scale with a range of 1.4" to 90" that is inserted into a 55.3" frame.

- Two different measurement ranges
- Easy read out
- Adjustable head positioner
- Choice of two units of measurement
- Simple assembly

SEC2161814009Measuring rod



TELESCOPIC HEIGHT MEASURING ROD

Mechanical telescopic measuring rod with large measuring range 60-200cm.

- Designed specifically for wall mounting.
- The telescopic measuring rod ensures that the result can be easily read at eye level.
- When not in use, the measuring slice can be folded down.

SOLHX200A.....Measuring Rod



INFANT HEIGHT MEASURING DEVICE

The Seca 210 is a lightweight and space-saving solution for precise length measurement of babies and toddlers while lying down. The Seca 210 is washable and foldable. The fixed head piece and the sliding foot positioner make it simple to use.

SEC2101821004Measuring device



Ritter® 253 LED EXAM LIGHT

Designed to be easy to use, easy to clean, durable and aesthetically pleasing, the Ritter 253 exam light is the simple lighting solution. Because the controls are located on the light head within easy reach, you'll never have to reach around the patient to adjust the lighting, making this the perfect light for lengthy procedures. With table mount and mobile options available, the 253 is a flexible, spacesaving addition perfect for any exam room.

- Life: greater than 50,000 hours
- Flexible arm
- Power cord: 8 ft. (2.5 m) hospital grade plug

RIT253-017Light with wall mount
 RIT253-012Light with table mount
 RIT253-016Light with mobile stand



Ritter® 255 LED PROCEDURE LIGHT

Whether you practice in an office environment or specialized procedure area, your work requires superior, white, shadowfree illumination that is safe and cool to the touch. The Ritter 255 LED Procedure Light by Midmark is designed to deliver both quality and affordability with features that can help you provide the best care to your patients.

- Illumination at 39" (1 m)
- LED Light Module: 16 LEDs arranged in a 4x4 array, 32 Watt rating
- Fully autoclavable handle
- Available in single or dual ceiling mount versions

RIT255-001Single head (8' ceiling mount)
 RIT255-004Dual head (9' ceiling mount)





THE GREEN SERIES EXAM LIGHT IV

Perfect for the OB/GYN office or for specialties that require very focused light in specific areas

- One LED with 3X the typical light output than halogen
- Highly flexible light pipe that is twice as long as conventional pipes with a compact, unobtrusive head that ensures pinpoint accuracy
- Touchless On/Off
- No bulb replacement

WEL48810 GS Exam Light IV with wall mount
 WEL48950 Mobile stand for Green Series Light



GREENSERIES 300 GENERAL EXAM LIGHT

Perfect for the OB/GYN office or for specialties that require very focused light in specific areas

- LED with 3X the typical light output than halogen
- Highly flexible light pipe that is twice as long as conventional pipes with a compact, unobtrusive head that ensures pinpoint accuracy
- Touchless on/off
- No bulb replacement
- 5 Year limited warranty

WEL44410 GS Exam Light 300 with wall mount
 WEL44400 GS Exam Light 300 with mobile stand



600 MINOR PROCEDURE LIGHT

Ideal for minor procedures and exams in any healthcare environment

- Touchless intensity controls provide the output you need, when you need it.
- Features 3 LEDs and a wider head specifically for demanding viewing conditions that require good light.
- 5 year limited warranty

WEL44610 GS600 Light with table/wall mount
 WEL44600 GS600 Light with mobile stand



GS900 PROCEDURE LIGHT

Perfect for minor procedures/procedures in the hospital, ambulatory care center, outpatient facility or physician's office. Features 6 LEDs with intensity controls, outstanding durability, extraordinary maneuverability and unparalleled spot quality. Ceiling mount options gives ability to maximize floor space and direct light wherever it is needed during a procedure. 5 year limited warranty.

WEL44900 with mobile stand
 WEL44900-C with ceiling mount
 WEL44900-W with wall mount

AARON® A942 HYFREACATOR HIGH FREQUENCY DESICCATOR

- 40 watts of power with higher resolution in Bi-Polar mode
- Increase and decrease in 1/10th watt increments from zero to ten watts
- Universal Power with 100-240 VAC
- More User Friendly with features that can save you time and money
- Less to do with fewer steps to accomplish the same procedures
- Less to buy as there are fewer accessories needed for additional functionality
- 3-Button handpiece

AARA942 Aaron® 942 Desiccator



BANTAM|PRO A952 HF DESICCATOR WITH CUT

The Bovie Bantam | PRO, A952 is designed and engineered to be the most reliable and durable generator available today.

- Offers you up to 50 watts in Cut mode, Blended Cut mode, and Pinpoint Coagulation mode. 40 Watts – Fulguration, Bipolar and Micropolar
- The unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat
- Unit cleanup is simplified by its non-absorbing, non-staining, high-impact ABS enclosure
- Standard package includes the unit, power cord, reusable handpiece, non-sterile drapes, 2 sterile drapes, 5 sterile disposable electrodes, 6 disposable sterile dermal tips, 12 disposable non-sterile dermal tips, 2 reusable grounding cord, 5 disposable split grounding pad, wall mount kit and users guide

AARA952 Aaron® 952 Desiccator



DISPOSABLE DERMAL TIPS

Non-Sterile

- AARA804-IMC Sharp Dermal Tip..... Box/100
- AARA806-IMC Blunt Dermal Tip Box/100

Sterile (packaged individually)

- AARA805 Sharp Dermal Tip..... Box/50
- AARA807 Blunt Dermal Tip Box/50



ACCESSORIES

- AARA812 Mobile stand for A941, A942 & A952
- AARA808T Mayo tray for stand
- AARA803 Footswitch for A941, A942 & A952
- AARA827V Bipolar forcep cord for A942 & A952
- AARA825 Bipolar forcep reusable 4.75" straight
- AARESRE Disposable split grounding pads Box/50
- AARA902 Reusable 3-Button handpiece for A941, A942 & A952
- AARA901 Reusable 3-Button handpiece for A940 & A950






Bovie

STERILE HIGH TEMPERATURE CAUTERIES

High temperature cauteries are popular for many types of procedures including pinpoint hemostasis, evacuation of a subungual hematoma, controlling diffuse bleeders, vasectomies and sculpting woven grafts. Individually packaged sterile.

- AARAA01 Fine tip 2200° F (1204° C)
- AARAA03 Loop tip 2200° F (1204° C)
- AARAA21 Elongated tip 1800° F (982° C) vasectomy



Bovie

DISPOSABLE ELECTRODES

- AARES01 Standard blade Box/50
- AARES37 Modified blade Box/25
- AARES20 Standard 3/16" ball Box/25
- AARES02 Standard needle Box/25
- AARES18 Angled blade Box/5
- AARES55 Extended blade 4" Box/25
- AARES03 Extended needle 6" Box/25
- AARES07 Extended 5mm ball Box/5
- AARES60 SuperCut fine needle 3cm Box/5
- AARES62 SuperCut fine needle 3cm Box/5
- AARES24 Short shaft loop 3/8" Box/5



Bovie

SMOKE SHARK™ II EVACUATOR

A modern, light-weight, compact design allowing for simplified use, storage and stackability. Never again order the Smoke Shark™ and still need to order a filter upon receipt, the SE02 is packaged with both the SF35 Extended-Life Filter offers up to 35 hours of performance reducing procedure cost by 38% (+/-) from the previous model, when used at the lowest setting. Covered 3-port filter design accommodates a variety of tubing sizes [1/4" (6.4mm), 3/8" (9.5mm), and 7/8" (22mm)] providing versatility in procedure and accessory setup. The ultra-quiet design offers the same excellent performance with less noise. The SE02 Smoke Shark Smoke Evacuator is -3 and RoHS compliant. Package includes: Unit, Power Cord, User's Guide, Pneumatic Footswitch (FSSEP) and Long Life 35-hour Filter (SF35).

- AARSE02 Smoke Shark™



Bovie

ACCESSORIES

- AAR786T 7/8" x 6' tubing, non-sterile Box/24
- AAR786TS 7/8" x 6' tubing, sterile Box/24
- AARSEPAT Pencil adapter 6 tubing, sterile Case/10
- AARSF35 Replacement 35hr filter for SE02
- AARSF18 Replacement 18hr filter for SE01
- AARSEAS Arm stand to hold tubing

PREMIER PHARMACY/VACCINE UNDERCOUNTER REFRIGERATOR BUILT-IN

- 4.6 Cu. Ft. Capacity
- One solid swing door, right hinged
- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- 2 probes; 1 in air and 1 in glycol bottle
- Adjustable operating temperature range: 36°F - 46°F (2° - 8° C)
- Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable.
- Three shelves (two adjustable/one fixed) with guard rail on back to prevent products from being pushed against the back wall.
- 23 ³/₄" W x 24 ¹/₂" D x 33 ³/₈" H


ABSPH-ABT-HC-UCBI-0404



PREMIER PHARMACY/VACCINE UNDERCOUNTER REFRIGERATOR FREESTANDING

- 2.5 Cu. Ft. Capacity
- One solid swing door, right hinged, not reversible
- Microprocessor temperature controller
- Digital temperature display
- Adjustable operating temperature range: 36°F - 46°F (2°C - 8°C)
- Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- 2 probes; 1 in air and 1 in sample bottle
- Four adjustable shelves (two large, two small) with guard rail on back to prevent product from being pushed against the back wall
- 17 ³/₄" W x 19 ¹/₂" D x 28 ⁷/₈" H

ABSPH-ABT-HC-UCFS-0204



SOLID DOOR PHARMACY PREMIER LABORATORY REFRIGERATOR FREESTANDING

- 10.5 Cu. Ft. Capacity
- Its ample room and six adjustable shelves provide ample space for a variety of contents, while a microprocessor temperature controller makes sure that it's always at a correct, safe temperature conveniently shown on a digital display.
- Its average temperature cycle is 41 degrees Fahrenheit, while its adjustable operating range is between 36 and 46 degrees Fahrenheit.
- High-density, urethane foam insulates the low-temp refrigerator and works in tandem with HFC-free refrigerant (R600a).
- Its door features a magnetic gasket for a positive seal, as well as a keyed door lock to keep contents secure.
- Audible and visual alarms are always on hand to alert nearby users of high or low temperature fluctuations.
- 59" H x 23 ³/₄" W x 26 ¹/₂" D

ABSPH-ABT-HC-10PS





VACCINE PLUS THERMOMETER

Triple display simultaneously shows minimum, maximum and current temperatures. Minimum/maximum monitors high/low readings overnight, on weekends, and for any time period. Alarm is programmable in 1° increments. 10' micro-cable cord.

VWR89094-778



Rubbermaid UTILITY CART, 3 SHELVES

Quiet casters and plastic construction help eliminate noise. Large 4" (10.2 cm) non-marking swivel casters. High style appearance with brushed aluminum uprights.

- 200 LB weight load.
- Black.

UNI139155



STAINLESS STEEL UTILITY CART, 3 SHELVES

- Stainless steel cart
- Strong and elegant
- 4 noiseless smooth casters
- 3 shelves
- 26" Wide x 34.75" High x 17.3" Deep

SOLF20



STAINLESS STEEL EQUIPMENT CART, 1 DRAWER

- Stainless steel cart
- Strong and elegant
- 4 noiseless smooth casters
- 2 shelves, 1 drawer
- 26" Wide x 34.75" High x 17.3" Deep

SOLF21

MULTI PURPOSE CART

- Top surface 25.5" wide x 17" deep
- Stain less steel top rail
- Oxygen tank holder on left side
- 1 x garbage cans on right side
- Built-in I.V. Pole
- Drawer configuration:
- (Includes drawer organiser / slots)
- Drawers 1 and 2, 20" wide x 14.9" deep x 2.5" high
- Flat working tray
- 4-Drawer 9.5" wide x 14.8" deep x 3" high
- 4th Drawers: 21" wide x 14.9" deep x 3" high
- 5th Bottom Large Drawer: 21" wide x 14.9" deep x 5.5" high



SOLF1

EMS CART WITH 5 DRAWERS

- Stainless steel top rail
- Plastic top surface
- Plastic sliding tray under the top surface
- Paper / chart holder (on left side)
- 2 x garbage cans (on right side)
- Built-in I.V. Pole
- CPR Board (hanged at rear)
- Drawer Configuration:
- (Includes drawer organizer / slots)
- Drawers 1 and 2, 20" wide x 14.25" deep x 2.5" high
- 3rd and 4th Drawers: 20" wide x 14.25" deep x 4" high
- 5th Bottom Large Drawer: 20" wide x 14.25" deep x 8.5" high



SOLF45-1

5-LEG MOBILE INSTRUMENT STAND

- Adjustable height (28 - 48") for comfortable use while sitting or standing
- 1" (2.54cm) chrome plated steel tubing
- Stable five-leg 20" base
- Removable stainless steel tray: 12.5" x 19"
- 2" (5.08cm) dual wheel nylon



DUK4365

MAYO STAND WITH TRAY

Mayo Stands with adjustable height, and stainless steel tray.

- DUK4368Double Post
- DUK4352Single Post



STANDARD WHEELCHAIR WITH FIXED OR DESK ARMS

This affordable wheelchair also features dual axle positions that allow for varying seat-to-floor heights, a heavy-duty chrome-plated steel frame, solid snap-on tires, adjustable brakes, sealed bearing assemblies, ergonomic handles, and flame-retardant black naugahyde upholstery.

- Institutional tested durability & proven quality.
- Proprietary contact points are easy to clean & disinfect.
- Dual axle positions allow varying seat heights.
- Maximum weight capacity: 300 lbs (136 Kg).
- Latex free.

DRV700-625



TRANSPORT CHAIR

Anatomically correct armrest height along with heel loops on the easy adjust swing away footrests make it even more secure and comfortable to sit in. Comfort-grip™ wheels provide shock absorption, better brake performance and they're easier to clean.

- 19" Seat
- 300 lb Weight Capacity

DRV700-848



STRETCHERS

- Adjustable Height 76"x22"x16"-28.8" Made of high strength aluminum alloy thick wall pipes.
- Adjustable back
- Foldaway legs are controlled by the handles on both ends of the stretcher.
- The legs can be folded under the stretcher

SOLAG4H





COLLABORATIVE CARE MADE EASIER

Our Plus devices let you move from an optical exam to digital in a few clicks, expanding your capabilities. Through our free iExaminer Pro App, save and share eye and ear images and videos for tracking and trending patient exams and consultations with specialists.

Explore our new systems with Plus models or update your device heads which are backwards compatible.



#777-PM3WAS



#71-PM3LXE



#71-PM3LDE

For more information on how Hillrom is bringing the future of care into view, please contact your Hillrom representative.

hillrom.ca

6950 Creditview Rd, Unit 4, Mississauga, ON, L5N 0A6

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

© 2021 Welch Allyn, Inc. ALL RIGHTS RESERVED. APR225603 Rev 1 21-MAY-20201 ENG-CA

1.800.232.7286 | canadamedical.net | orders@canadamedical.net

Welch Allyn Connex® Vital Signs Monitor 6000 Series

Experience accurate, automatic spot-check vitals and continuous monitoring sent to your EMR from the bedside to help support improved patient safety and clinical decisions.

The Connex Vital Signs Monitor represents the next generation in vital signs monitoring with an easy-to-use color touchscreen and advanced vital signs technology to monitor vitals, respiration and motion. Plus, EMR integration gives clinicians immediate access to the information they need—where and when they need it.



Connect to the future of vital signs today

Designed specifically for busy floors like yours to increase efficiency

- > Captures both spot-check vitals and continuous monitoring
- > EMR integration for accurate electronic vitals documentation
- > Vivid color touchscreen simplifies training and operation
- > Comprehensive point-of-care documentation including customizable data fields
- > Advanced Patient ID including two-way name confirmation
- > Patient data protection for HIPAA compliance
 - Protect patient identifiable information with customizable options like abbreviations
 - Allow automated deletion of patient data after sending to EMR
 - Apply rules to check record completeness before saving
- > Easy-to-understand error messages for simple troubleshooting

Industry-leading parameters for accurate patient information

- > SureBP® blood pressure technology
 - Adult, pediatric and neonatal blood pressure modes
 - Customizable blood pressure monitoring programs help nurses follow hospital protocol
 - 15-second readings for greater patient comfort
- > SureTemp® Plus and/or Braun ThermoScan® PRO 6000 ear thermometry (optional—choose one or both)
- > Masimo® or Nellcor® SpO₂ algorithms (optional)—both sensors and signal processing
- > Masimo noninvasive continuous monitoring of total hemoglobin (optional)
- > Manual measurements entry (weight, height, respiration rate, pain) and documentation of modifiers like body position and O₂ therapy
- > Customizable documentation fields to fully document at the point of care

Configurations Available

CONNEX VSM	6700	6800
Parameters		
SureBP NIBP	●	●
SureTemp Plus Thermometry	○	○
Braun PRO 6000 Ear Thermometry	○	○
Masimo or Nellcor SpO ₂	○	○
EarlySense Contact-Free Motion	⬆	⬆
Covidien Microstream EtCO ₂	⬆	⬆
Masimo RRA	⬆	⬆
Masimo SpHb	⬆	⬆
Communications		
4 USB Ports	●	●
Ethernet	●	●
802.11 a/b/g Wireless Radio	⬆	●
Accessories		
Barcode Scanner	○	○
Weight Scale	○	○
EMR Connectivity	○	○
Printer	○	○

● Included ○ Optional ⬆ Upgradeable

Features:

- SureBP technology can provide BP in only 15 seconds BP averaging to help ensure accurate hypertension detection and diagnosis
- SureTemp Plus Oral / Axillary Thermometry
- Option for Nonin® SpO2 technology
- 7" color touchscreen for a simplified workflow Options for documenting patient observations (i.e. Respiration Rate, BMI)
- EMR connectivity helps minimize documentation errors
- 8-hour Lithium ion battery to maximize time before recharging
- Each device ships with one adult blood pressure cuff



- WEL44XT Spot 4400 with NIBP and SureTemp® Thermometer
- WEL44WT Spot 4400 with NIBP, SureTemp Thermometer and Nonin® SpO2
- WEL4400-DST Spot 4400 Desk Stand
- WEL4400-MBS Spot 4400 Mobile Stand
- WEL77794-M4400 Green Series™ 777 Integrated Wall System for Spot 4400
- WEL77794-2M4400 Green Series 777 Integrated Wall System with PanOptic™ Ophthalmoscope for Spot 4400
- WELBATT22 Connex® Spot, 2-Cell LI-Ion Battery
- WEL106275 Micro USB Cable
- WELREUSE-11 Cuff Reusable Adult
- WELREUSE-12 Cuff Reusable Large Adult
- WEL4500-35 Fast BP Hose with FlexiPort®, 10 ft

The Welch Allyn ProBP 2400 is not a typical electronic blood pressure device. It is an affordable, professional device that you can trust to accurately measure your patients' blood pressures despite many common underlying circumstances. The ProBP 2400 is specifically designed for practices like yours. It gives you the quality, durability and features you don't get from home blood pressure devices that are mainly designed for general household use, which some practices resort to because of cost. It provides you consistent technique that's not assured with manual blood pressure gauges. And it does so at a budget-friendly price that other professional blood pressure monitors with the same features don't offer.

Features

- More affordable than similarly featured professional blood pressure monitors
- 3-Reading Average Mode can give a broader view of patient's true blood pressure
- Irregular heartbeat detection provides data to help you make testing decisions
- Manual Mode allows it to serve as a digital gauge when auscultation is indicated
- Microprocessor-controlled deflation rate provides consistent technique
- Built-in folding handle allows you to carry it safely and comfortably
- Optional mobile stand makes transporting easy for maximum location flexibility
- Optional wall mount bracket ensures device is exactly where and when you need it



- WEL2400 ProBP 2400 BP Unit with Adult & Large Adult Cuffs
- WEL2400-MS Mobile Stand for ProBP 2400
- WEL2400-WMB Wall Bracket for ProBP 2400
- WELREUSE-09-2400 Blood Pressure Cuff, Child (14 - 22cm)
- WELREUSE-11-2400 Blood Pressure Cuff, Adult (22 - 32cm)
- WELREUSE-12-2400 Blood Pressure Cuff, Large Adult (32 - 42cm)



WelchAllyn  Hillrom. **Connex[®] SPOT MONITOR**

The Connex[®] Spot Monitor features an easy-to-use, vivid touchscreen display and provides accurate vital signs measurement including blood pressure averaging, spot checking, interval monitoring and custom scoring across patient populations. The device connects wirelessly to your EMR and sends vitals to the patient chart from the point of care to help improve efficiency. Its upgradeable design and custom configurations include optional wall unit integration, industry-leading thermometry and pulse oximetry choices, connectivity and other mounting/mobility options. In addition, the Accessory Power Management Stand with WhisperDrive™ technology features additional power backup, providing up to 17 hours of on-time, alleviating the need to charge the monitor as frequently. Ships standard with one Adult and one Large Adult cuff.

Features

- SureBP[®] technology can provide BP in only 15 seconds
- Connex Scoring App calculator for custom protocols (e.g., NEWS, MEWS)
- BP averaging to help ensure accurate hypertension diagnosis
- Supports resident list workflows for long-term care settings
- SureTemp[®] Plus or Braun ThermoScan[®] PRO 6000
- Thermometry Partnerships with more than 90 leading EMRs including Epic[®] and Cerner[®]
- Choice of leading SpO2 technologies: Masimo[®], Nellcor[®] or Nonin[®]
- Add to integrated wall system for complete patient exam workflow
- Scale connectivity to Detecto, Health-o-meter[®] and SECA scales
- Mobile stand option includes work surface with integrated battery for extended usage

WEL71XX-B	Digital BP
WEL71WX-B	Digital BP, Oximeter
WEL71WT-B	Digital BP, Oximeter & SureTemp Thermometer
WEL7000-MS3	Mobile Stand
WEL7000-GCX	Wall Mount
WEL7000-DST	Desk Stand



Affordable and easy to use, Welch Allyn FlexiPort Reusable Blood Pressure Cuffs offer all the advantages of FlexiPort technology—including simplified connectivity—in a dependable, durable design.

Features

- FlexiPort technology enables every cuff to be used as a one- or two-tube cuff
- Tested to withstand demanding clinical environments
- Latex-free material minimizes risk of allergic reactions
- Rotatable port, improves patient comfort and cuff durability
- Folded edge reduces risk of cuts and scrapes, maximizing patient comfort
- Meets all latest clinical guidelines for proper fit from AAMI and AHA
- Primary sizes are color-coded to make it easy to find the proper cuff size
- Three-year warranty

FlexiPort™ Cuffs (No Tubing)

- WELREUSE-06 Small Infant
- WELREUSE-07 Infant
- WELREUSE-08 Small Child
- WELREUSE-09 Child
- WELREUSE-10 Small Adult
- WELREUSE-11 Adult
- WELREUSE-12 Large Adult
- WELREUSE-13 Thigh
- WELREUSE-MLT 8-Cuff Kit

Complete Inflation System (includes Cuff, Tubing, Bulb and Valve)

- WELREUSE-09-2BV Child
- WELREUSE-11-2BV Adult
- WELREUSE-12-2BV Large Adult
- WELREUSE-13-2BV Thigh



WelchAllyn® INFLATION BULB & VALVE COMBINATION

- WEL5088-01 Bulb & Valve Combination
- WEL5087-01 Air Valve Only
- WEL5086-01 Bulb Only



WelchAllyn® COIL TUBING with TWIST FIT CONNECTOR

- WEL5089-13 8 "



WelchAllyn® 767 SERIES ANEROIDS



The versatile 767 Aneroid is the smart choice for accurate and reliable readings throughout your facility. Improved ergonomics and a new modern design make for an easy-to-use blood pressure aneroid, developed with staff and facility workflow efficiency in mind. Recessed, laser-engraved dial face with high-contrast pointer provides certified accuracy to reduce errors and improve patient care.

- WEL7670-01 Aneroid with adult cuff, wall mount
- WEL7670-03 Aneroid with adult cuff, mobile stand

WelchAllyn® DS58 HAND ANEROID



With gear-free DuraShock™ technology, you can count on the DS58 aneroid to remain in calibration longer than traditional aneroids, providing more accurate readings over time and a lower lifetime cost, helping to improve practice efficiencies and your bottom line.

- Gear-free DuraShock™ technology Certified accuracy to +/- 3 mmHg
- Laser-engraved dial face for accuracy
- Jeweled movement for long life
- Lifetime calibration warranty
- Latex-free for safety
- Thumbscrew bleed valve and ergonomic inflation bulb
- Unsurpassed reliability in a traditional, handheld design
- Designed to enhance comfort
- Withstands a 60" drop, exceeds AAMI's shock resistant requirement

- WELDS58-11 With adult cuff and zippered case
- WELDS58-MC With 4 sized cuffs and zippered case

WelchAllyn® DIAGNOSTIC SET, COAXIAL



Welch Allyn 3.5V Diagnostic Set with PanOptic Basic LED Ophthalmoscope, MacroView Basic LED Otoscope, one Lithium Ion Plus Rechargeable Power Handle, and Hard Case

- WEL71-PM2LXUDiagnostic set

WelchAllyn® 3.5V NICKEL-CADMIUM RECHARGEABLE HANDLE

Includes #72300 battery. A choice of convenient rechargeable power sources that fit all the WelchAllyn 3.5 V instrument heads. Smooth, non abrasive knurled finish ensures a no-slip grip. Chrome-plated brass construction.

WEL71000-A



WelchAllyn® 3.5V LITHIUM ION RECHARGEABLE HANDLE

Twice the on-time, half the weight and size. Fits all Welch Allyn 3.5 V Instrument Heads. 120 minutes of on-time (compared to 60 minutes for standard handles). Low battery indicator, no "memory" effect. Includes AC charging module.

WEL71900



WelchAllyn® COAXIAL OPHTHALMOSCOPE HEAD

- Halogen HPX lamp illumination for true tissue color.
- 18 unique aperture/filter combinations for greater versatility
- Polarizing filter virtually eliminates corneal blue filter
- Compatible with all existing Welch Allyn 3.5 V power sources
- Look for the traditional "flat-top" design you've come to trust

WEL11720-L Ophthalmoscope head only



WelchAllyn® PANOPTIC™ BASIC LED OPHTHALMOSCOPE

- Welch Allyn PanOptic Basic LED Ophthalmoscope; power handle not included
- 5X larger view of the fundus vs. standard ophthalmoscopes in an undilated eye
- Digitally capture, store and share fundus images by combining with iExaminer
- Greater working distance improves comfort for both practitioner and patient
- Dynamic focusing wheel allows continuous, smooth action and more precise control
- Hallogen HPX® lamp provides bright, white light
- Compatible with all existing Welch Allyn 3.5 V power sources

WEL118-2 Head only





WelchAllyn® MACROVIEW™ OTOSCOPE HEAD

- A nearly complete view of the tympanic membrane
- Better clarity and definition of landmarks
- The ability to adjust focus for variable ear canal length or farsighted eyes
- Lip grip for secure fastening and ejection of ear specula
- Fiber optics produce cool light with no reflections, no obstruction
- Twice the field of view and 30% greater magnification than traditional otoscopes

WEL238-2 MacroView Basic LED Otoscope
 WEL238-3 MacroView Plus LED Otoscope



WelchAllyn® DISPOSABLE OTOSCOPE SPECULUM

WEL52432-U Pediatric (2.75mm) Bag / 850
 WEL52434-U Adult (4.25mm) Bag / 850



WelchAllyn® SPOT™ VISION SCREENER

The Spot Vision Screener is a handheld, portable device designed to help users quickly and easily detect vision issues on patients from 6 months of age through adult. Spot screens both eyes at once from a nonthreatening 3-foot distance. The touch-screen display allows for one-touch activation, simple management of patient data entry, and easy configuration for both vision screening and autorefracton applications.

Features:

- Supports AAP screening guidelines for early detection of amblyopic risk factors
- Automated screening provides thorough and easy-to-understand results
- Easy to use, easy to implement with minimal user training required
- Adequately powered for mass screening events
- Captures readings 97% of the time helping to screen otherwise difficult patients
- Wireless printing and WiFi-enabled for easy export of data
- Manual exports available via USB port
- Uses lights and sounds to help engage children

WELVS100S-B



WelchAllyn® RETEVAL-DR™ ELECTRORETINOGRAPH

The Welch Allyn RETeval-DR Electroretinograph is designed for healthcare providers who want to improve management of their patients with diabetes. In just three minutes, you can objectively screen for vision-threatening diabetic retinopathy (VTDR) wherever patients receive care. Welch Allyn RETeval-DR Electroretinograph enables efficient, comfortable, non-mydratic screening for VTDR that improves patient compliance to help preserve vision in your patients with diabetes.

Features:

- No dilation required, enhancing patient comfort
- Short test time—only three minutes per patient
- No interpretation and minimal training necessary
- Immediate results provided on the device and in a PDF report
- Portable, handheld device goes wherever patients with diabetes receive care

WELRETEVAL-001

WelchAllyn® REPLACEMENT BULBS

- WEL03000-U6For Ophthalmoscope (11710)
- WEL03100-U6For Otoscope
- WEL04200-U6For Exam Light (48600 & 48700)
- WEL04900-U6For Ophthalmoscope (11720)
- WEL06300-U6For Exam Light (44100 & 44110)
- WEL06500-U6For MacroView™ (23810 & 23820)
- WEL07800-U6For Vaginal Illuminator (78000)
- WEL08800-U6For Vaginal Illuminator (78810)



WelchAllyn® RECHARGEABLE BATTERIES

- WEL72200Fits 71000 & 71670 Handles (Black)
- WEL72300Fits 71000-A Handle (Orange)



3M LITTMANN™ CARDIOLOGY IV STETHOSCOPE

- Outstanding acoustics with better audibility of high-frequency sounds (adult side, firm pressure).
- Tunable diaphragm; Hear high or low frequency sounds by slightly adjusting pressure on the chestpiece.
- Dual-lumen tubing: Two sound baths in one tube eliminates the rubbing noise of traditional double tubes.

- | | |
|-----------------------|-----------------------------|
| 3M 6152.....Black | 3M 6153 Burgundy |
| 3M 6154.....Navy Blue | 3M 6155..... Hunter Green |
| 3M 6156.....Plum | 3M 6157..... Carribean Blue |
| 3M 6158.....Raspberry | 3M 6161 Rose Pink |



PROSCOPE™ DUAL HEAD STETHOSCOPE

The Proscope™ Dual Head combines lightweight construction with truly affordable pricing. Combination diaphragm/bell chestpiece for wide frequency response Adjustable aluminum binaurals with PVC eartips. Overall length 32.25" (81cm).

- ADC670BD Burgundy
- ADC670BK Black
- ADC670N Navy
- ADC670R Red
- ADC670RB Royal Blue
- ADC670TL Teal
- ADC670V Purple



LITTMANN™ CLASSIC III STETHOSCOPE

The 3M™ Littmann® Classic III™ Stethoscope offers high acoustic sensitivity for exceptional performance when doing general physical assessments. It features dual tunable diaphragms and an updated design that is easier to clean and maintain.

The 3M™ Littmann® Classic III™ Stethoscope is the latest version of the stethoscope that helps millions of medical professionals achieve their best. The Classic III stethoscope offers high acoustic sensitivity for exceptional performance, plus a versatile two-sided chestpiece with tunable diaphragms. The Classic III stethoscope is used by students and medical professionals alike to identify, listen to, and study heart, lung, and other body sounds for physical assessment, patient monitoring, and diagnosis.

Tunable diaphragm technology—a 3M invention—lets clinicians hear different frequency sounds by simply adjusting the pressure on the chestpiece. Hold the chestpiece with light pressure to hear low frequency sounds; press a little more to hear higher frequency sounds. This time-saving feature allows the healthcare provider to focus on the patient, while the non-chill diaphragm provides patient comfort.

Headset tension is easily adjusted for head size and comfort by squeezing together or pulling apart the eartubes. Snap-tight, soft-sealing eartips conform to individual ears for an excellent acoustic seal and comfortable fit. Eartips snap firmly onto the ends of the eartubes and, for safety, require a significant effort to remove.

Resilient next-generation tubing retains its shape and flexibility even after folding tightly into a pocket. It provides longer stethoscope life due to improved resistance to skin oils and alcohol, and is less likely to pick up stains. Next-generation tubing is not made with natural rubber latex, a benefit to sensitive medical professionals and their patients. And, to play its part in the protection of human health and the environment, the next-generation tubing does not contain phthalate plasticizers.

The Classic III stethoscope is made in the USA for consistent quality. It comes complete with large and small soft-sealing eartips, non-chill bell sleeve, and instructions.

Features:

- Two-sided chestpiece with tunable diaphragms on both the adult and pediatric sides.
- Single-piece tunable diaphragm is easy to attach, and easier to clean because its surface is smooth without crevices.
- Pediatric side converts to a traditional open bell by replacing the single-piece diaphragm with a non-chill rim.
- Next-generation tubing provides longer life due to improved resistance to skin oils and alcohol; less likely to pick up stains.
- Small tunable diaphragm is useful for pediatric, small, or thin patients; around bandages; and for carotid assessment.
- Stainless steel chestpiece is precision-machined into an aesthetically pleasing, less angular shape. The stem features open side indicator.
- Snap-tight eartips have a soft, smooth surface providing a comfortable acoustic seal and comfortable fit.

3M 5620..... Black
 3M 5626..... Raspberry
 3M 5832 Lavender
 3M 5839 Lemon Lime
 3M 5639..... Rose Pink
 3M 5621..... Grey
 3M 5629..... Orange
 3M 5829 Lime Green

3M 5622..... Navy Blue
 3M 5627..... Burgundy
 3M 5835..... Turquoise
 3M 5840..... Emerald
 3M 5624..... Hunter Green
 3M 5623..... Caribbean Blue
 3M 5630..... Ceil Blue
 3M 5831..... Plum



DIGITAL FINGER TIP PULSE OXIMETER

Moisture proof, one-button operation, requires two AAA batteries, one-year limited warranty.

- Checks heart rate & respiratory function. Suitable for finger thickness: .37" - .86" (9.5mm - 21.9mm)
- Not made with natural rubber

DYN7088Adult
 DYN7090Pediatric



DIGITAL FINGER TIP PULSE OXIMETER

- Dual colour LED displays SpO2 (blood-oxygen saturation) and pulse rate
- Suitable for adult and pediatric patients
- Adjustable brightness
- Automatically powers off; low power consumption
- Includes 2 AAA batteries, instruction manual, carrying case.



ZENGG6500

PULSE OXIMETER, HAND HELD

The Edan H100B pulse oximeter with powerful memory (300 hours) and configurable alarms provides an easy to use tool to measure oxygen levels (SpO2) and heart rate information. The Edan H100B has a wide variety of sensors options (neonate, pediatric, adult) to fit any patient sizes. Memory from the Edan H100B can be downloaded and printed using the optional software and USB cable.

Edan H100B Hand Held Pulse Oximeter Includes:

- Operator's Manual & Quick Reference Card
- Adult Finger Clip Sensor
- 4 "AA" alkaline batteries
- Operator's Manual
- Deluxe Carrying Case

SOLH100BPulse Oximeter
 SOL02.01.210121.....Pediatric Sensor
 SOL01.57.040199.....Infant Disposable Sensor



PULSE OXIMETER, HAND HELD

Affordable, dependable and easy to use, the SPECTRO2™ 10 Pulse Oximeter is a low cost, spot-check device displaying SpO2, pulse rate, signal strength and Pulse amplitude Index (PI). This device delivers reliable performance, meeting the needs of health care providers in primary care environments. Its unique optional accessories allow clinicians to customize a solution to meet the needs of virtually any clinical environment.

- Large, bright, easy-to-read LED for at-a-glance assessment
- Unique on-board cradle to store and protect finger sensor
- Compatible with Docking Station, Attachable Printer, and BCI® spot-check sensor

BCIWW1000ENPulse Oximeter
 BCI1303Infant disposable sensor 10/BX
 BCI1301Pediatric disposable sensor 10/Bx





WelchAllyn  Hillrom

THERMOSCAN® PRO 6000 TYMPANIC EAR THERMOMETER

- Innovative PerfectTemp™ technology adjusts for variability in probe placement
- ExacTemp™ technology detects stability of the probe during measurement
- Designed to take quick and accurate readings, to improve patient satisfaction
- 60 second pulse timer assists with manual measurement of pulse rate, respiration
- Electronic security features help prevent theft, loss
- Small or large cradles for probe cover storage that pop-up for easy attachment
- Optional charging station stores 200 probe covers, with rechargeable battery
- 6- or 9-foot security tether options keep device attached to cradle

WEL06000-200 Braun PRO 6000 with Small Cradle (holds 1-pc box)

WEL06000-005 Braun PRO 6000 Probe Covers Box/200



WelchAllyn  Hillrom

CARETEMP® TOUCH-FREE THERMOMETER

- Infrared sensor reads temperature 4-6 cm from the center of forehead
- 60 second pulse timer, memory recall button, C/F scale conversion button
- No patient skin contact required, helps minimize the risk of cross-contamination
- Measurement feedback system detects proper operating distance for accuracy
- Large, backlit LCD display that's intuitive and easy to read
- Takes approximately 3,000 readings on two AAA batteries
- 2 year device warranty

WEL105801



WelchAllyn  Hillrom

SURE TEMP® PLUS 690 ORAL THERMOMETER

- Fast: quickly captures accurate readings on all patient types
- Large LCD displays temperatures in Fahrenheit or Celsius
- Convenient storage for 25 probe covers
- Last temperature recall button
- Optional security wall mount
- Takes approximately 6,000 readings on 3 AA batteries
- Interchangeable, removable wells to help reduce the risk of cross-contamination

WEL01690-200 Thermometer

WEL05031-105 Probe Covers Box/250



THERMOMETER NON-CONTACT

- AMSINO® Infrared, 1 year warranty, 2x AAA batteries included.
- Designed to be a screening device to measure body and forehead temperature without contact.

AMG016-662

EDAN® SONOTRAX DOPPLER

- A high performance Fetal Doppler with heart rate
- Compact design
- Built-in audio recorder
- Large LCD display & backlit
- Accurate FHR detection with clear sound
- Automatic / Average / Manual fetal heart rate counting
- Highly sensitive interchangeable probes
- 9V alkaline battery



SOLMS91046613mHz OB Probe with LCD Display

IMEX ELITE 100 DOPPLER

Imex Elite 100 Fetal Dopplers have superior performance and broad beam to find the fastest and easiest OB Heartbeat and signal location as well as arterial and venous blood flow. The comfortable shape of the Imex Elite probe is easy to hold and features broad beam technology, making it easy to find and lock onto the signal of interest. Obstetrical applications include detecting pregnancy, reassuring the mother, and determining multiple pregnancies

- Imex Elite 100 Fetal Dopplers feature a large built-in speaker for superior sound
- New interchangeable probes incorporate Static Suppression Doppler (SSD) technology for a clear, static-free signal
- Imex Elite 100 Fetal Dopplers operate on a standard 9V alkaline or optional rechargeable batteries
- Imex Fetal Doppler probes can be gas sterilized or sheathed
- 5 Year Limited Imex Parts Warranty



IMXEN30 with 3mHz OB Probe

IMXEN50 with 5mHz Vascular Probe

IMXEN80 with 8mHz Vascular Probe

IMEX ELITE 200 DOPPLER WITH DISPLAY

Imex Elite 200 Fetal Dopplers are Nicolet's top of the line doppler offering a very dependable and easy to read digital display as well as clean audio of the fetal heart. The Imex Elite 200 Fetal Doppler comes in an easy to use handheld size and features a comfortable easy to hold probe shape. The Elite 200 Doppler automatically calculates fetal heart or pulse rate based on probe application in order to assure accurate results.

Other features of the Imex Elite 200 Fetal Doppler are as follows:

- Digital heart rate display to easily determine fetal heart rate
- Obstetrical applications include detecting pregnancy, reassuring the mother, and determining multiple pregnancies
- New interchangeable probes incorporate Static Suppression Doppler (SSD) technology for a clear, static-free signal
- Imex Elite 200 Fetal Dopplers operate on a Standard 9V alkaline or optional rechargeable batteries
- The comfortable shape of the Elite probe is easy to hold and features broad beam technology, making it easy to find and lock onto the signal of interest
- Elite 200 Fetal Dopplers feature a large built-in speaker for superior sound
- Imex Fetal Doppler probes can be gas sterilized or sheathed
- 5 Year Limited Imex Parts Warranty



IMXED30 with 3mHz OB Probe

IMXED50 with 5mHz Vascular Probe

IMXED80 with 8mHz Vascular Probe



AQUA-SONIC® 100 ULTRASOUND GEL

Blue Gel

- PAR01-01 20mL Foil Package - Sterile Box/48
- PAR01-08 250mL
- PAR01-50 5L

Clear Gel

- PAR03-08 250mL
- PAR03-50 5L



AQUAFLEX® GEL PAD

- PAR04-02 2cm x 9cm Box/6



ECLIPSE® ULTRASOUND PROBE COVERS

- PAR38-01 Eclipse Pre-lubed Probe Cover Box/100
- PAR38-03 Eclipse 3D Pre-lubed Probe Cover Box/100



THERMASONIC ULTRASOUND GEL WARMER

- PAR82-01 Single Bottle Warmer
- PAR82-03 Three Bottle Warmer



TRANSEPTIC® CLEANSING SOLUTION

Multi-purpose spray for professional use in cleansing ultrasound transducer/probe surfaces between procedures.

- PAR09-25 250 ml Spray Bottle

DERMLITE®



The versatile DermLite has a multitude of applications. Viewing pigmented lesions, varicose veins, hair transplants, skin irritations, or skin lesions is just the beginning. Together with the DermLite's compact size and lightweight design (only 3.5 oz/99g), the possibilities are virtually infinite. To make sure you will not run out of power, the DermLite's endurance measures over 1200 lesion examinations per battery, due to its diode lighting technology and long-life lithium power.

DERDL100DL-100 Unit

NEUROPEN® SENSATION SCREENING DEVICE



Includes: Neuropen, 1 x 10g Monofilament, 1 x Neurotip

The Neuropen® is a reusable device that features a calibrated spring mechanism to deliver a 40g sharpness test using Neurotips® single use neurological examination pins, to assess reduced sensation to sharpness/ pain in small nerve fibers and a 10g pressure test using Neuropen Monofilaments® to assess protective touch/pressure sensation in large nerve fibres.

MDPNT0100Neuropen Device
 MDPNT010410g Monofilaments Box/5
 MDPNT540540g Neurotips Box/100

SE-1200 Express 12-channel ECG with Interpretation



- Large high resolution color touch screen
- Alphanumeric keyboard
- One-touch operation
- Enhanced weak signal detection
- Comprehensive filters & anti-noise technology
- Real-time ECG waveforms freezing and review
- Advanced Glasgow algorithm for automatic measurements and interpretation
- Built-in high resolution thermal printer with external USB printer (optional) support
- Extended printing
- Large internal memory which can be enlarged by USB flash disk
- Built-in data management with reports review
- LAN/Wi-Fi (optional) connections
- Bi-directional communication with SE-1515 Data Management System
- Multi-format reports export (optional)
- Barcode scanner support (optional)

SOLSE1200-EXBSE-1200 ECG Basic
 SOLSE1200-EXPRESSS..... SE-1200 ECG with touch screen



SPIROMETER DISPOSABLES

WEL703419Disposable Flow Transducer Pkg/100
 WEL100680.....Spirometry Noseclips Pkg/25





Welch Allyn CP 150
Transfer information
via USB or Ethernet

Welch Allyn
CardioPerfect
Workstation
Software



Record test data and
patient information to
individual patient records
in your EHR for analysis
and permanent storage



CP 150™ ELECTROCARDIOGRAPH

The CP 150 combines the best traditional box ECG features (e.g. full-size printouts and clinical decision support with MEANS interpretation) with connectivity enhancements designed especially for today's primary care practices.

- ECG print previews can help save time and reduce paper costs
- Seven-inch color, touchscreen display for easy entry of patient data
- Optional wireless connectivity helps increase device mobility and utilization
- Optional spirometry lets you perform two tests with one device
- Optional DICOM compatibility enhances interoperability
- Optional adult and pediatric MEANS ECG interpretation software
- Three-, six-, or full 12-lead ECG preview options
- Battery operation lets you easily transport your device to the patient
- Simplified EMR connectivity without a complicated interface
- Instant-on powers the device up quickly to take an ECG at a moment's notice

WELCP150-1ENB ECG Non-interpretive

WELCP150A-1ENB..... ECG Interpretive

WEL105353 CP100/150/200 Recording paper..... Box/200

WEL105341 CP150 - Office Cart

WELCH ALLYN® DIAGNOSTIC CARDIOLOGY SUITE



The Welch Allyn Diagnostic Cardiology Suite now includes ECG and SPIROMETRY capabilities, helping you shorten the path to diagnosis for improved patient outcomes.

Features:

- Captures best 10 seconds of ECG data from memory, which helps reduce need for repeat ECGs
- ECG complies with ACC, AHA and HRS recommendations for adult/pediatric patients
- Easily access prior spirometry exams to monitor treatment plans over time
- Utilize one device for spirometry and ECG exams for a common workflow and setup
- ECG and spirometry applications launch directly from the EMR
- Instantly verify spirometry exam quality and test performance



WELCC-RXX-WAXX..... Connex Cardio PC based 12-Lead Multi-Channel Resting ECG - software, Connex Cardio WAM wireless acquisition module, AHA 10 wire, 12-lead lead set with banana plugs, CardioSens Ultra II resting tab electrodes (pouch of 100), Clear Choice adapter clips (set of 10), Connex Cardio UTK

WELCC-RXX-WADXConnex Cardio PC based 12-Lead Multi-Channel Resting ECG - software, Connex Cardio WAM wireless acquisition module, AHA 10 wire, 12-lead lead set with banana plugs, CardioSens Ultra II resting tab electrodes (pouch of 100), Clear Choice adapter clips (set of 10), Connex Cardio UTK, DICOM

WELCC-RXX-AAXX.....Connex Cardio PC based 12-Lead Multi-Channel Resting ECG - software, Connex Cardio AM12 wired acquisition module, AHA 10 wire, 12-lead lead set with banana plugs, CardioSens Ultra II resting tab electrodes (pouch of 100), Clear Choice adapter clips (set of 10)

WELCC-RXX-AADXConnex Cardio PC based 12-Lead Multi-Channel Resting ECG - software, Connex Cardio AM12 wired acquisition module, AHA 10 wire, 12-lead lead set with banana plugs, CardioSens Ultra II resting tab electrodes (pouch of 100), Clear Choice adapter clips (set of 10), DICOM

WELDCSS-CXXDiagnostic Cardiology Suite Spirometry with Calibration Syringe

WELDCSS-NXX.....Diagnostic Cardiology Suite Spirometry without Calibration syringe

WEL105341Office Cart

IQECG™ DIGITAL ECG

The Midmark IQecg makes your job easier with a range of reviewing and editing tools to help increase productivity and efficiency. At the touch of a button, you can capture an abnormal ECG as it appears on screen, enhance the diagnosis with tools like zoom and electronic calipers, quickly edit statements and perform serial comparison.

Features include:

- Lightweight, portable design, weighing less than 11 oz.
- View and monitor live, 12-lead ECG in selectable formats on screen
- Interpret and analyze data for neonate, pediatric and adult patients via the Midmark 12-lead Resting Electrocardiogram Analysis Program
- Review and edit ECG results on-screen with free text or customizable quick edit statements
- Perform serial comparison with overlay of up to 4 ECG tests
- Save time and reduce data entry errors by populating patient information directly from your EHR/EMR software
- Facilitate paperless office strategy-eliminate costs of thermal paper and the hassles of scanning and filing patient charts

RIT4-000-0062IQECG™ Digital ECG with Interpretive



IQSPIRO™ DIGITAL SPIROMETRY SYSTEM

The Midmark IQspiro Digital Spirometry System combines a digital spirometer with innovative software solutions that make it fast and easy to capture, interpret and retrieve patient information. Our digital spirometer also features correct patient positioning that can optimize testing and an eject lever that saves time and alleviates handling of used mouthpieces.

- View real-time Flow/Volume Loop and Volume/Time graphs with fully customizable parameters.
- Pre and post-bronchodilator testing with automatic interpretation.
- New eject lever alleviates handling of used mouthpiece
- Accurate and precise; meets or exceeds ATS 2005 Standards
- Lightweight, portable design weighs only 9 ounces
- Correct, ergonomic patient testing orientation
- Easy grasp, new T-shaped handle design for patient support and stability during testing
- Improved mouthpiece stability and protection with new easy-latch door design
- Easy to print full page reports on standard paper

RIT4-000-0028IQ Spiro System with calibration syringe

RIT4-100-1200IQ Diagnostic Workstation Software

RIT2-100-1205Disposable Mouthpieces Box/25

RIT2-100-1206Disposable Mouthpieces Case/100



**Call your Canada Medical account manager
to order the appropriate Midmark
software solution**



3M RED DOT™ RESTING ECG TAB ELECTRODE

3M2330E Pkg/100



WelchAllyn® RESTING TAB ELECTRODE

WEL108071 Box/500



3M RED DOT™ MONITORING ELECTRODE

3M 2231.....With Softcloth® Tape and Solid Gel Bag/50



3M RED DOT™ MONITORING ELECTRODE

3M 2670-5.....Flexible Non-Woven and Sticky Gel..... Case/1000



3M RED DOT™ DIAPHORETIC ELECTRODE

3M 2271-50With Softcloth® Tape and Solid Gel Bag/50

ELECTRODE SKIN PREP PAD

Formulated to remove excess skin oil enhancing ECG electrode to skin contact.

DYN1508Skin Prep Pad.....Box/100



3M RED DOT™ ECG TRACE PREP TAPE

Designed to remove the non-conductive skin layer.

3M 2236..... 18mm x 5m Roll



SAFETY PREP RAZOR

PER75-4008Single Blade.....Box/100



PERSONNA PREP RAZOR (PREMIUM)

- Specially designed comb (rake)
- Wrapped blade reduces cuts and nicks
- Polymer Coated Stainless Steel Blade
- Heavy weight handle design
- Fixed head
- Blades precision ground for comfort.

PER75-4004Single Edge.....Box/100



TWIN BLADE RAZORS

- Long-lasting, Superior Quality Blade
- Clear plastic guard
- Disposable
- Unisex

DYN4250Twin Blade.....Box/50





AECPLUS®

SEMI-AUTOMATIC DEFIBRILLATOR



The AED Pro is designed for both basic and advanced life support. While it offers guidance to the BLS rescuer, it also provides more advanced capabilities to the professional rescuer. The AED Pro not only relays advanced CPR feedback, it's also a vital signs patient monitor that can be configured to the way you work.

- A three-lead cable allows basic monitoring using ECG electrodes.
- Flexible operating modes let the professional rescuer apply manual override to access more advanced capabilities, including heart rate displays, and to decide if and when to deliver treatment.
- A high-resolution auto-gaining LCD display keeps the ECG trace within a specific area, always displaying the gain, so the rescuer can easily see the patient's rhythm.

ZOL22600010102011060



AECPLUS® METAL WALL CABINET WITH ALARM

ZOL8000-0855



AECPLUS® ELECTRODES CPR-D PADZ®

ZOL8900-0810-01..... Pediatric

ZOL8900-0800-01..... Adult



AECPLUS® LITHIUM BATTERIES

ZOL8000-0807 Pkg/10



ZOLL AED 3 DEFIBRILLATOR



Enhanced Real CPR Help, featuring a full-color display with vivid rescue images, CPR cycle timer, and a large color bar gauge that shows CPR compression depth. RapidShock Analysis™ enables the industry's shortest rhythm analysis for more continuous life-saving CPR. CPR Uni-padz™ electrodes can be used on both adult and pediatric patients, so there is no need for two separate sets of pads. Cloud connectivity that can automatically report AED readiness. 5-year shelf life for pads and batteries increases readiness and decreases frequency of maintenance. Available in semi-automatic and fully automatic models.

ZOL8501-001102-06

CAVICIDE1® SURFACE DISINFECTANT

CaviCide1, a multi-purpose disinfectant/decontaminant cleaner, can be used on hard, non-porous surfaces. CaviCide1 is effective in killing bacteria and fungi in 1 minute, including TB, Norovirus, HIV-1, HBV, HCV and many more. When used as directed, fragrance-free CaviCide1 will effectively clean and disinfect surfaces and can help reduce the risk of cross-contamination.

- MET11-5024 CaviCide1, 24oz spray
- MET11-5001 CaviCide1, 1 Gallon
- MET11-5100 CaviWipes1, Small, 6" x 7"
- MET11-5150 CaviWipes1 X-Large, 9" x 12"



Metrex
INFECTION PREVENTION

INFECTION PREVENTION

CLOROX® PROFESSIONAL DISINFECTING BLEACH WIPES

- Kills 58 Microorganisms in 3 minutes or less.
- 30 second kill times on 32 bacteria
- 30 seconds to 1 minute on 19 viruses
- 1 minute on norovirus
- 3 minutes on C. difficile sporesm TB and fungi

- CLX01308 6.75" x 9" Wipes
- CLX01309 12" x 12" Wipes
- CLX01310 12" x 12" Wipe Refill Pack



CLOROX
HEALTHCARE™

CLOROX® PROFESSIONAL HYDROGEN PEROXIDE WIPES

Speed for Hard Surface Disinfection. Fast kill times on more than 40 microorganisms. If your wipe doesn't keep large surface areas wet long enough, it can lead to costly re-wiping. Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes keep larger surfaces wet longer, enabling you to use fewer wipes for better value, less waste and less worry.

- 30 second to 1 minute contact times on most bacteria and viruses.

- CLX01456 6.75" x 9" Wipes
- CLX01458 12" x 11" Wipes
- CLX01459 12" x 11" Wipes Bucket Refill



SANI-CLOTH® PLUS SURFACE DISINFECTANT WIPES

- Germicidal Disposable Cloth
- Bactericidal / Virucidal / Cleans and Disinfects
- Ready-to-use disposable cloths are pre-moistened with bactericidal disinfectant
- For disinfection and cleaning of instruments and all hard surfaces such as counter-tops and equipment
- Effectively disinfects against: Influenza A2/HK, Hepatitis C, Respiratory Syncytial Virus (HRS), Methicillin resistant Staphylococcus aureus (MRSA), Vancomycin resistant Enterococcus faecalis (VREF), Escherichia coli (E.Coli), Pseudomonas aeruginosa, and HIV-1 (AIDS) virus, Herpes Simplex Type II and Salmonella
- Ideal for hospitals, health care facilities, first aid rooms, ambulances and other areas where avoidance of cross-contamination is essential
- Packaged in a convenient dispenser for easy access
- 6" x 6.75" wipe, 2 minute contact time is effective against enveloped / non-enveloped viruses, bacteria, TB, fungi, multi-drug resistant organisms and bloodborne pathogens2
- Designed to be compatible with hard nonporous surfaces and equipment.

- PDIQ55172 Super Sani-Cloth Wipes 160/12ea/cs





GLUTE-OUT NEUTRALIZER FOR GLUTARALDEHYDE & OPA

Disposing of used disinfectants can be extremely dangerous as glutaraldehyde and OPA products remain highly toxic and biocidal even after 14 days of re-use. Improper disposal of disinfectants by pouring down a drain has proven to cause high levels of exposure for the clinician and is considered illegal in many regions. Glute-Out Neutralizer offers a solution to dangerous disinfectant disposal problems. Glute-Out deactivates glutaraldehydes and OPA by using a glycine powder base that neutralizes solutions in just five minutes.

CV-610-1045 2oz bottle (Neutralizes 1 Quart)
 CV-610-1099 5oz bottle (Neutralizes 1 Gallon)



CIDEX OPA SOLUTION, DEVICE DISINFECTANT SOLUTION

- 14-Day High Level Medical Device Disinfectant Solution
- Sterilization in 10 hours at 25°C
- Activation required

JOH20394 3.8 Litre
 JOH20393 Test Strips, Bottle / 15 Box/2



METRICIDE 28® DISINFECTANT SOLUTION

MetriCide 28 must maintain the minimum effective concentration (MEC) level of 1.8% glutaraldehyde or greater during the 28-day use life to effectively sterilize and high-level disinfect. It is recommended to test the glutaraldehyde solution prior to each use to help prevent healthcare-associated infections (HAIs).

- Recommended for immersible and heat-sensitive equipment and instruments
- Contains surfactants

MET11-2801 3.8 Litre
 MET10-304 Test Strips (1.8%) Bottle/60



EASYTASK CLEANING SYSTEM 200WIPES/BOX

10"x12" wipe, add any cleaning or disinfecting solutions to make a durable wipe. Resealable bag keeps wipes moist. Holds 1/2 gallon of liquid, low lint.

CFSN-ETA110CZGW.....Dry Wipe Box/200 wipes



BIOESQUE DISINFECTING SOLUTION

No Harmful Chemicals, No PPE needed. Surfaces must remain wet for 1 minute to eliminate viruses*,

2 minutes to eliminate bacteria, 3 minutes to eliminate fungi & M.bovis BCG (TB), and to disinfect Norovirus 4min. DIN 02486858

CML910216.....3.78L
 CML910208.....946ml

EMPOWER® DUAL ENZYMATIC DETERGENT

EmPower® is used to pre-clean critical or semi-critical medical devices prior to terminal sterilization or high level disinfection. It contains two protease enzymes that offer the broadest cleaning action on a variety of protein soils. Inclusion of surfactants provides additional cleaning on carbohydrates, lipid and protein soils. Works fast to remove blood, tissue, mucous and other protein-rich bodily fluids from instruments.

- Dual-enzymatic detergent cleaner
- Gentle on instruments
- Less corrosive than distilled water
- Low-foaming formula is ideal for use in automated machines
- Fragrance free option for respiratory-sensitive users

MET10-4100 3.8 Litre



EMPOWER® FOAM SPRAY

EmPower® Foam is used to pre-clean critical or semi-critical medical devices prior to terminal sterilization or high level disinfection. It is a ready-to-use, foaming dual-enzymatic spray that starts the pre-cleaning process by breaking down blood, tissue, mucus and other protein-rich bodily fluids from instruments. The foaming spray delivery system is ideal for use in the operating room, sterile processing, and other departments where instrument cleaning is not immediately available.

- Ideal for instruments with the presence of moderate to heavy organic debris
- Covers instruments, reducing splashes and spills
- Ready-to-use in a convenient sprayer
- Covers instruments to rapidly begin the cleaning process

MET10-4224 24 ounce



ENZYCLEAN PROTEASE INSTRUMENT DETERGENT

Enzyclean Protease Enzyme Detergent is specially formulated for use as a presoak or manual detergent.

The protease and amylase enzymes break down organic and inorganic soils, residues, films and deposits and ensures complete cleaning.

WEIEBL1HC 3.78 Litre



CUTS-IT GEL INSTRUMENT PRE-SOAK

Clings to and dissolves blood, fat, tissue and other body fluids and keeps bioburden from drying on instrumentation.

Neutral pH, Non-corrosive, Phosphate-free, Biodegradable. Perfect for a transport gel.

WEIM60005HC 710 ml





Purell ES4 MANUAL HAND SANITIZER DISPENSER

ES4 Hand Sanitizer Dispensers

- GOJ5020-01..... Wall Dispenser, White
- GOJ5024-01..... Wall Dispenser, Black/Chrome

Purell Advanced Hand Sanitizer

- GOJ05060-02..... Gel for ES4 Manual Dispenser, refill



Purell ES6 TOUCH FREE HAND SANITIZER SYSTEM

ES6 Touch Free Hand Sanitizer Dispensers

- GOJ6420-01..... Wall Dispenser, White
- GOJ6424-01..... Wall Dispenser, Black

Purell Advanced Hand Sanitizer

- GOJ6460-02..... Gel for ES6 Touch Free Dispenser, 1200mL Refill



Purell GEL HAND SANITIZER

Kills 99.9% of most common germs that may cause illness. Specially formulated with moisturizers to leave hands feeling soft and refreshed.

Gel

- GOJ3770-12..... 354mL Pump Bottle
- GOJ9674-12..... 8oz Bottle with Aloe



HAND SANITIZING GEL 540ML



Ecolab's waterless hand sanitizer formulated to promote skin health and reduce the spread of infections. Contains moisturizers to condition and protect the skin and contains 70% ethyl alcohol.

- ECO6000325 Quikcare Gel Hand Sanitizer 540ml

Purell ES4 MANUAL SOAP DISPENSER

ES4 Hand Soap Dispensers

- GOJ5030-01.....Wall Dispenser, White/White
- GOJ5034-01.....Wall Dispenser, Black

Healthcare Foam Soap

- GOJ5072-02.....1200mL Refill



Purell ES6 TOUCH FREE SOAP DISPENSER

ES6 Touch Free Soap Dispensers

- GOJ6430-01.....Wall Dispenser, White
- GOJ6434-01.....Wall Dispenser, Black

Health Care Foam Soap

- GOJ6472-02.....1200ml Refill



ECOLAB BACTI-STAT® LIQUID ANTIMICROBIAL HANDWASH

Ecolab's Bacti-Stat is formulated to provide outstanding efficiency and exceptional mildness. Ecolab's hand soaps are designed specifically for the handwashing needs of healthcare workers and help reduce the potential for transmission of healthcare acquired infections, as well as promote skin health.

- ECO6017285..... 4 Litre



ECOLAB BACTISTAT HAND SOAP 540ML

- 0.3% Triclosan formula meets FDA requirements for Healthcare Personnel Handwash.
- Lathers quickly, rinses freely which minimizes residual lather than can cause irritation.
- Contains Aloe Vera which helps protect against irritation and dry, cracked skin.

- ECO6048722Bactistat Hand Soap 540ml



GREEN FOAM HAND SOAP 222ML

A mild, biodegradable foam soap formula certified by Green Seal™ and EcoLogo® to meet their joint standard for hand cleaners and hand soaps. The certifications mean that the product has been thoroughly evaluated by these unbiased, third-party. Unscented.

- GOJ5715-06.....Green Foam Hand Soap 222ml





TORK® SINGLE-FOLD TOWELS

CFSSK1850A Pack/250 Case/16 Packs



TORK® PREMIUM HAND TOWELS

CFS101298 Pack/100 Case/21 Packs
 CFS552120 Dispenser



TORK PREMIUM HAND TOWEL ROLLS

CF5774820 Dispenser
 CFS8634550 Roll Towel Case/12Rolls



WHITE SWAN® 2 PLY HOUSEHOLD PAPER TOWELS

CFS01890 90 Sheets/Roll Case/21 Rolls



WHITE SWAN® 2 PLY QUICK WIPES

CFS08500 80 Wipes/Box Case/135 Boxes



WHITE SWAN® 2 PLY FACIAL TISSUE

CFS08301 100 Tissues/Box Case/30 Boxes



2 PLY BATHROOM TISSUE

CFSTM1601A 500 Sheets/Roll Case/48 Rolls



TORK® ADVANCED MINI-JUMBO BATH TISSUE

CFS12024402..... 2 Ply, 751' Roll Case/12
CFS555020A..... Dispenser



GARBAGE CAN, WHITE

- Steel square outer can and lid with removable plastic inner liner.
- Enamel Finish Size.
- 32 quart (8 gallon) capacity.
- 11.5"D x 12.25"W x 18.25"H (29.21cm x 31.12cm x 46.36cm)

BRE35266



GARBAGE/RECYCLING BAGS

CFSTT-20W..... Garbage Bag, Regular Strength, White (20" x 22") Cs/500
 CFSTT-24..... Garbage Bag, Regular Strength, Black (22" x 24") Cs/500
 CFSTT-24W..... Garbage Bag, Regular Strength, White (22" x 24") Cs/500
 CFSTT-32..... Garbage Bag, Regular Strength, Black (26" x 36") Cs/250
 CFSTT-52..... Garbage Bag, Regular Strength, Black (30" x 38") Cs/250
 CFSTT-32TB..... Recycling Bag, Regular Strength, Blue (26" x 36") Cs/250



METRIMIST® ODOUR ELIMINATOR

MET10-1158..... 240mL, 8 oz.



LYSOL® DISINFECTANT SPRAY

Lysol® Disinfectant Spray kills 99.9% of viruses and bacteria on hard surfaces and eliminates odours by killing the bacteria that cause them. Also effective on mold and mildew.

UNI34052..... 350g, Crisp Linen Scent





TURBOPRO™ HANDHELD ELECTROSTATIC SPRAYER

Get your entire space Clorox® clean in minutes*. The Clorox® TurboPro™ Electrostatic Sprayer is engineered to clean high-touch areas such as schools, offices, athletic facilities and healthcare settings more quickly and efficiently with the convenience of touchless coverage. Clorox® electrostatic sprayers are paired with EPA-registered chemistries that are designed especially for sprayer devices. When used with Clorox® Turbo Pro™ Disinfectant Cleaner (sold separately), this cordless electrostatic sprayer disinfects hard-to-reach places in airports, daycare centers, hotels, laundry rooms and long-term care facilities and more.

CLX29561TurboPro™ Handheld Electrostatic Sprayer



CLOROX® TURBO PRO™ DISINFECTANT CLEANER FOR SPRAYER DEVICES

Safe for sprayer devices: Tested for use through hand held, backpack and cart sprayer devices without compromising safety, compatibility and performance
 Kills SARS-CoV-2: EPA-approved to kill the virus that causes COVID-19* in 2 minutes and 99.9% of bacteria in 5 seconds on hard, nonporous surfaces
 Kills 51 Pathogens: Kills germs, including cold and flu viruses, MRSA and norovirus in 2 minutes

CLX60091Clorox® Turbo Pro™ Disinfectant Cleaner for Sprayer Devices



CLOROX® GERMICIDAL BLEACH 3.5L 3/CS

- Concentrated 7.4% bleach
- Health Canada registered to kill a total of 35 microorganisms, including C. difficile spores in 5 minutes
- DIN# 02408961 Concentrated 7.4% bleach

CLX01383 Clorox® Germicidal Bleach 3.5L 3/CS



DAWN ULTRA DISH SOAP 532ML

UNI74836 Dawn Ultra Dish soap 532ml



TIDE LIQUID HE LAUNDRY DETERGENT 2.95L

UNI0888564Loads per Jug

Ritter® M9D STERILIZER

The Ritter M9D AutoClave Sterilizer, an economical option to the Ritter M9, offers easy, effective sterilization. One button starts the cycle and you must manually open the door within 30 minutes of completion to activate the dry cycle. This extra step is the only difference between the Ritter M9D and its fully automatic counterpart (shown below).

RITM9D-042..... M9D Sterilizer



Ritter® M9 & M11 AUTOMATIC STERILIZER

Today, more than ever, you need products that are reliable and simple to use. The Ritter M9 & M11 couldn't be easier and is time tested for reliability. You load it... set it... and forget it. Select your cycle and the time is set for you. Press "Start" and the Ritter M9 fills with the correct amount of water. Once the sterilization process is complete, the door opens automatically and quietly to dissipate the steam to give you fast drying time for your instruments. It's real convenience!

The Ritter M9 & M11 sterilizer also has easy front access for filling and draining and controls that are easy to understand. Programmable cycle features allow custom creation of different cycle parameters for special sterilization needs.

Programmable buttons allow you to change the following:

- Time (3-90 minutes), Temp (23-275 °F), Dry Time (0-60 minutes), Vent (slow or fast and two program buttons allow you to store the custom cycles.
- Pre-programmed sterilization cycle controls: Unwrapped (132°C or 270°F for 3 min), Pouched (132° C or 270° F for 5 minutes), and Packs (121°C or 250°F for 30 mins).
- The display Indicates cycle selected, cycle temperature and exposure time for selected cycle. During the cycle, it shows messages describing status of cycle. When the cycle enters sterilization mode and dry mode, remaining cycle time is displayed as well as temperature and pressure.

RITM9-042 M9 Sterilizer 9" x 15" chamber
RITM11-042..... M11 Sterilizer 11" x 18" chamber
RIT9A599001 Thermal Printer for M9 & M11
RIT9A682001 Data Logger for M9 & M11
RIT060-0016-00 Thermal Printer Paper for M9 & M11(042 Model)



Tuttnauer® EZ11PLUS STERILIZER

Your Sterilization & Infection Control Partners

Features:

- 11" diameter X 19.8" deep chamber and flexible loading options.
- Printer Included
- USB Port enables users to easily transfer cycle data via a flash drive to a PC
- EZGlide door features a sleek modern look and effortless opening and closing.
- Closed door drying, meaning the door remains closed through the dry cycle to maintain sterility
- Comprehensive 2 year parts & labor warranty and 10 year chamber warranty

TutEZ11Plus-P





QUICKCLEAN™ ULTRASONIC CLEANERS

Table top unit, basket is included. QuickClean™ washes instruments for sterilization, allowing caregivers and staff to spend their time with patients, rather than washing instruments.



POWERFUL CLEANING. CONSISTENT RESULTS. SAFER FOR EVERYONE.

RITQC1-01	1.2 Gallon
RITQC3-01	3.3 Gallon
RITQC6-01	6.6 Gallon

MINIBIO AUTO READER

Tuttnauer MiniBio will allow 20 Minute & 1 hour process times simultaneously, holds 3 indicators with different incubation times, at the same temperature, USB connection with free software download, thermal printer

TUTWTL198-0057



ULTRA-RAPID BIOLOGICAL INDICATORS 50/BOX

Tuttnauer Ultra-Rapid BI's are to be used on the Bionova or MiniBio units. 20 minute read out, Can be used in gravity or pre/post vacuum autoclaves.

TUTWTL198-0072



EZTEST® STERILIZATION MONITORING SYSTEM

Complete in office sterilization monitoring system - incubator, EZTest (25), English log book. Everything you need for effective monitoring of steriliation cycles in your office. Visual results in 24 hours.

SSP4145-SK/25B



DELTA STEAM BIOLOGICAL INDICATOR

- EZTest is a self-contained biological indicator (BI) for monitoring sterilization
- Steam, ethylene oxide, propylene oxide, hydrogen peroxide
- Bacterial spores on a carrier material packaged within a small, thermoplastic culture tube.
- Inside is a sealed-glass ampoule of specially formulated soybean casein digest culture medium containing a color indicator which turns a dramatic yellow when growth occurs.
- Easy-to use in-office; no laboratory transfer needed
- Dramatic visual results within 24 hours (48 hours for ethylene oxide)

SSP4145-BV5-100 Box/100



3M COMPLY™ STERIGAGE™ CHEMICAL INTEGRATOR

Chemical integrators are used for pack control which validates that the sterilant has penetrated to the point of placement in the pack and confirms that sufficient exposure conditions have been met. 3M™ Comply™ Steam Chemical Integrators consist of a paper wick and steam and temperature sensitive chemical pellet contained in a paper/film and foil laminate envelope of laminate materials.

- Moving front ink technology provides easy to read Accept or Reject result
- Can be used in all 250-275°F / 121-135°C steam sterilization cycles
- Response correlates with a Biological Indicator at 3 time/temperature relationships under ideal steam sterilization conditions



3M 1243B 100/Bag

COMPLY™ STEAM CHEMICAL INDICATOR STRIPS

3M™ Comply™ Indicator Strip enables pack control by validating that the sterilant has penetrated to the point of placement in the pack and confirming that sufficient exposure conditions have been met. Use in steam sterilization cycles.



3M 1250 Class 4 Box/480

STERILIZATION WRAP SHEETS (CSR)

MAXGEM211212"x12" CS/1000
 MAXGEM211515"x15" CS/1000
 MAXGEM211818"x18" CS/1000
 MAXGEM212020"x20" CS/500
 MAXGEM212424"x24" CS/500
 MAXGEM213030"x30" CS/500



STERILIZATION POUCHES, SELF SEALING

SPMSP22540 2.25" x 4" **Box/200**
 GSM442100 3.5" x 6.25" **Box/200**
 SPMSP35900 3.5" x 9" **Box/200**
 SPMSP52510 5.25" x 10" **Box/200**
 SPMSP71300 7.5" x 13" **Box/200**
 GSM222555 12" x 15" **Box/200**
 SPMSP27590 2.75"x9" **Box/200**



INSTRUMENT AUTOCLAVABLE BRUSHES

SRG2872550-50 Nylon
 SRG2972550-50 Stainless Steel



SINGLE END STAINLESS STEEL



SINGLE END NYLON

3M COMPLY™ AUTOCLAVE TAPE

3M 1322-18MM..... 3/4" x 60 Yards
3M 1322-24MM..... 1" x 60 Yards



DISTILLED WATER

TRA5-2100 4 Litre



SPEED CLEAN® AUTOCLAVE CLEANER

Speed Clean is an autoclave cleaner solution that can be used in all brands of autoclaves, to clean and remove chlorides.

RIT002-0396-05..... 16oz



CHAMBER BRITE® AUTOCLAVE CLEANER

Chamber Brite powder is simply the easiest way to keep autoclaves clean and free of water deposit build-up. Designed specifically for steam sterilizers, Chamber Brite effectively cleans and descales the reservoir, chamber and tubing in minutes.

TUTCB0010 Box/10



INSTRUMENT TRAYS, STAINLESS STEEL

- AMG020-142 10.5" x 6.5" x 2.5"
- AMG020-140 Cover For 020-142
- AMG020-162 12.75" x 7" x 2.5"
- AMG020-160 Cover For 020-162




FLAT TRAY, STAINLESS STEEL

- AMG020-527 10" x 6.5" x 3/4"
- AMG020-525 15" x 10.5" x 3/4"



STERILIZER POUCH RACK

- RIT002-2108-00..... Rack




INSTRUMENT PROTECTORS

Instrument Protectors will protect instrument tips and edges from damage during handling and sterilization.

Variations Include:

- perforated protectors (holes) which reduce vacuum seal & moisture build-up
- steotome protectors
- double protectors
- retractor protectors
- validated for pre-vacuum sterilization
- available in a wide array of sizes for almost every application
- steam and gas permeable
- non sterile
- latex free

Assorted Sizes and Colours Available. Call for details.



STERILIZATION MARKERS

- CML 13601 Non-Sterile, withstands temperatures up to 500°F, Black



THUMB DRESSING FORCEP

ALMM06-0120..... 5"



TISSUE FORCEP

ALMM08-0120 5", 1 x 2 teeth
 ALMM08-0140 5.5", 1 x 2 teeth
 ALMM08-0160..... 6", 1 x 2 teeth



ADSON TISSUE FORCEP

ALMM08-0420 4.75", 1 x 2 teeth



ADSON DRESSING FORCEP

ALMM06-0520..... 4.75", Serrated



SPLINTER FORCEP

ALMM10-0120 4"



WILDE NASAL & EAR FORCEP

ALMM06-0710..... 6"
 ALMM06-0720..... 7"



HARTMANN ALLIGATOR FORCEP

ALMM06-0800..... 3.5"
 ALMM06-0830..... 6.5"



KELLY HEMOSTAT FORCEP

ALMM18-0320 5.5" Straight
 ALMM18-0340 5.5" Curved



HALSTED MOSQUITO FORCEP

ALMM18-0200 5" Straight
 ALMM18-0220 5" Curved



BOZEMAN UTERINE DRESSING FORCEPS

ALMM14-0300 10" Straight
 ALMM14-0320 10" Curved



FOERSTER SPONGE FORCEP

ALMM20-0180 7" Straight Serrated
 ALMM20-0200 7" Curved Serrated
 ALMM20-0220 7" Straight Smooth
 ALMM20-0240 7" Curved Smooth

ALMM20-0260 9.5" Straight Serrated
 ALMM20-0280 9.5" Curved Serrated
 ALMM20-0300 9.5" Straight Smooth
 ALMM20-0320 9.5" Curved Smooth



SIMS UTERINE SOUND

ALMM14-0500 12.5"



SCHROEDER TENACULUM FORCEP

ALMM14-0340 9.5" Straight
 ALMM14-0360 9.5" Curved



OPERATING SCISSORS

ALMM04-0220 5.5" Straight, Blunt/Blunt
 ALMM04-0240 5.5" Straight, Sharp/Blunt
 ALMM04-0260..... 5.5" Straight, Sharp/Sharp



ALMM04-0500 5.5" Curved, Blunt/Blunt
 ALMM04-0520 5.5" Curved, Sharp/Blunt
 ALMM04-0540 5.5" Curved, Sharp/Sharp



IRIS SCISSOR

ALMM04-1200 3.5" Straight
 ALMM04-1220 3.5" Curved
 ALMM04-1240 4.5" Straight
 ALMM04-1260..... 4.5" Curved



MAYO SCISSORS

ALMM04-0800 5.5" Straight
 ALMM04-0840 6.75" Straight
 ALMM04-0870 9.0" Straight

ALMM04-0820 5.5" Curved
 ALMM04-0830 6.0" Curved
 ALMM04-0872 9.0" Curved



DOUBLE ENDED PROBE

ALMM34-0130 6.0"



GROOVED DIRECTOR


ALMM34-0240 5.5"



METZENBAUM SCISSORS

ALMM04-0910 5.5" Straight
 ALMM04-0920 7.0" Straight
 ALMM04-0960..... 9.0" Straight

ALMM04-0912 5.5" Curved
 ALMM04-0940 7.0" Curved
 ALMM04-0980 9.0" Curved



SPENCER STITCH SCISSORS

ALMM04-1100 3.5"
 ALMM04-1120 4.5"



UNIVERSAL BANDAGE SCISSORS

Available in Black, Blue, Orange or Red

ALMM02-0500 6.0"



LISTER BANDAGE SCISSORS

ALMM02-0140 5.5"
 ALMM02-0160..... 7.0"



CRILE WOOD NEEDLE HOLDER

ALMM12-0420 6"



DERF NEEDLE HOLDER

ALMM12-0120 4.75"



HALSEY NEEDLE HOLDER, SERRATED

ALMM12-0160..... 5"



MAYO-HEGAR NEEDLE HOLDER

ALMM12-0320 5"
ALMM12-0340 5.5"



OLSEN-HEGAR NEEDLE HOLDER

ALMM12-0280 5.5"
ALMM12-0300 6.5"
ALMM12-0310 7.0"



BAUMGARTNER NEEDLE HOLDER

ALMM12-0260..... 5.5"



NAIL NIPPER (WITH BARREL SPRING)

ALMM29-0720 5.25"



NAIL NIPPER (WITH DOUBLE SPRING)

ALMM29-0744 5.5"



TAYLOR PERCUSSION HAMMER

ALMM24-0140 8"



QUEEN SQUARE HAMMER

ALMM24-0170 9"
ALMM24-0180 15.5"



FINGER RING SAW

ALMM30-0920 with Blade Protector
ALMM30-0940 Blade Only



CURETTES - STAINLESS STEEL - AUTOCLAVABLE

Straight Buck Curettes


ALMM22-0300 Size 00
ALMM22-0320 Size 0
ALMM22-0340 Size 1
ALMM22-0360 Size 2
ALMM22-0380 Size 3
ALMM22-0400 Size 4
ALMM22-0420 Size 5

Angled Buck Curettes

ALMM22-0500 Size 00
ALMM22-0510 Size 0
ALMM22-0520 Size 1
ALMM22-0530 Size 2
ALMM22-0540 Size 3
ALMM22-0550 Size 4
ALMM22-0560 Size 5

Billeau Ear Curettes

ALMM22-0600 Size 1
ALMM22-0620 Size 1
ALMM22-0640 Size 1 Each



FOX DERMAL CURETTE - AUTOCLAVABLE

ALMM30-0300 3mm
ALMM30-0320 4mm
ALMM30-0340 5mm
ALMM30-0360 6mm



EAR WASH SYSTEM

- Eliminates messy outflow by automatically dislodging, removing and disposing of cerumen
- Transparent tip seals tightly and directs outflow back to sink and away from you and your patient.
- You control and monitor suction, water temperature and pressure
- Safe - reduces typical side effects like vertigo, nausea and tympanic membrane damage
- Turn cerumen removal into a clean, quick, safe procedure with the first system to simultaneously irrigate and suction discharge using Hydrovac Action.

WEL29350 Ear Wash System
 WEL29360 Disposable Eartips Pkg/25
 WEL29330 Replacement Tubing (10')



OTOCLEAR® EAR WASH SYSTEM

Otoclear spray wash kit, manual, intermittent lavage

- Delivers greater volume compared to syringes or other manual devices
- Efficient design allows for two-handed procedure
- Built-in temperature gauge on bottle for water temperature

BNX7290 Ear Wash System w/ ear basin
 BNX7200 Otoclear Irrigation Tips Box/40 - Disposable



EAR SYRINGE, METAL

MILV919-385 4oz



EAR BASIN, PLASTIC

GFI1280 Autoclavable



TUNING FORKS

MIL19-102 C128
 MIL19-104 C256
 MIL19-106 C512



ELEPHANT EAR WASH SYSTEM

Designed by a physician, the spraybottle provides enough water pressure to quickly and thoroughly clean ears, and holds enough solution to clean both ears without refilling. The unique hose design allows for easy positioning of the solution stream.

STA555-EW Ear Washer with 3 Tips
 STA555-TW Tips Only Pkg/20 - Disposable
 STA555-SHIELD Shield Only



STERILE SINGLE-USE INSTRUMENTS


Our single-use sterile instrument trays are designed to give you a good floor quality stainless steel instrument in a peel top plastic package that is sterile and ready for use.



IRIS SCISSOR, STERILE

Stainless Steel instrument in sterile peel top plastic package


- SUPSS.390.11..... 4.5" Straight Scissor
- SUPSS.393.11..... 4.5" Curved Scissor



SPENCER STITCH SCISSOR, STERILE

Stainless Steel instrument in sterile peel top plastic package


- MAXDYND04037 4.5" Straight Scissor



O.R. SCISSORS, STERILE

Stainless Steel instrument in sterile peel top plastic package


- SUPSS.030.14..... 5.5" Sharp/Sharp Scissor
- SKL96-2521..... 5.5" Sharp/Blunt Scissor



METZENBAUM SCISSORS, STERILE

Stainless Steel instrument in sterile peel top plastic package


- SUPSS.120.18..... 7" Straight Scissor
- SUPSS.121.18..... 7" Curved Scissor



LISTER BANDAGE SCISSORS, STERILE

Stainless Steel instrument in sterile peel top plastic package


- SUPPL.014.18 5.5" Scissor



ADSON FORCEP, STERILE

Stainless Steel instrument in sterile peel top plastic package


- SUPDT.318.12 4.75" Serrated
- SUPDT.320.12 4.75" 1x2 Teeth



HALSTEAD MOSQUITO FORCEP, STERILE

Stainless Steel instrument in sterile peel top plastic package


- SUPAF.060.13 5" Straight
- SUPAF.061.13 5" Curved



KELLY FORCEP, STERILE

Stainless Steel instrument in sterile peel top plastic package


- MAXDYNJ04048 5.5" Straight
- SUPDT.320.12 5.5" Curved



BAUMGARTNER NEEDLE HOLDER, STERILE

Stainless Steel instrument in sterile peel top plastic package


- MAXDYND04043 5.5"



SCHROEDER TENACULUM FORCEP, STERILE

Stainless Steel instrument in sterile peel top plastic package

- SKL96-2624 9.5" Straight



PENLIGHT

- AAR66RN Disposable Pkg/6



LIGHTED EAR CURETTES WITH MAGNIFICATION

Standard Box (50 per box & 1 Light Source & 1 Magnifying Lens)

- BNX2210Lighted AngleLoop® Box
- BNX2220Lighted MicroLoop® Box
- BNX2230Lighted WaveCurette™ Box
- BNX2240Lighted FlexLoop® Box
- BNX2250Lighted CeraSpoon® Box
- BNX2260Lighted InfantScoop® Box
- BNX2270.....Lighted VersaLoop® Box
- BNX2245Variety Pack (Includes 7 tip styles) Box

Clinic Packs (200 per box & 4 Light Sources & 4 Magnifying Lens)

- BNX2209 Variety Pack (Includes 7 tip styles) Box
- BNX2200 Light Source Only Box



LIGHTED PLACEMENT TOOL

Gives Audiologists the ability to safely and accurately place foam blocks in the ear.

- BNX2850Lighted Placement Tool Box/15



LIGHTED SUCTION FOR CERUMEN REMOVAL


This combination provides the ability to quickly remove cerumen or foreign bodies with the utmost ease and safety. Suction unit sold separately or works with any existing suction unit.

- BNX2625Lighted Suction Box/25



LIGHTED FORCEPS FOR FOREIGN BODY REMOVAL

- BNX2750 Forceps with Light Source & Magnifying Lens Box/10
- BNX2775 Forceps Only Box/10



LIGHTED ARTICULATING EAR CURETTE

The easy to maneuver tip articulates up to 90 degrees to assist in the removal of cerumen from the ear canal and foreign bodies from the ear and nose. The Lighted Articulating Ear Curette is the only lighted curette that articulates to assist in the removal of cerumen from the ear canal and works because of the unique hinged design.

- BNX2511..... Lighted Articulating Curette Box/50



MILTEX DISPOSABLE DERMAL CURETTES

- MIL33-522mm Box/50
- MIL33-533mm Box/50
- MIL33-544mm Box/50
- MIL33-555mm Box/50
- MIL33-577mm Box/50



-
- BNX9555..... Safe Ear Currettes™ - White FlexLoop® 4mm **Box/50**
 - BNX4888..... Safe Ear Currettes™ - Blue InfantScoop® 2mm..... **Box/50**
 - BNX1222..... Safe Ear Currettes™ - Green MicroLoop® 3mm **Box/50**
 - BNX6333..... Safe Ear Currettes™ - Yellow CeraSpoon® 4mm **Box/50**
 - BNX3444..... Safe Ear Currettes™ - Red AngleLoop® 4mm..... **Box/50**
 - BNX4111..... Safe Ear Currettes™ - Purple VersaLoop® 3mm **Box/50**
 - BNX2999..... Safe Ear Currettes™ - Orange ControLoop® 4mm..... **Box/50**
 - BNX5777..... Safe Ear Currettes™ - Variety Pack (10/each style & 15/white) **Box/75**

INSTRUMENT CLEANING SUPPLIES

- MIL3-740 Surgical Instrument Stain Remover, 3oz
- MIL3-720 Surgical Instrument Cleaner, 8oz
- SRG2872550-50 Surgical Instrument Brush, Nylon
- SRG2972550-50 Surgical Instrument Brush, Stainless Steel



METRILUBE® INSTRUMENT LUBRICANT

MetriLube is an oil-based instrument lubricant that coats metal instruments to prevent spotting, staining, and rusting of surgical instruments and helps keep hinges working freely.

- MET-10-3425 24 oz



INSTRUMENT LUBRICANT MILK

Instrument Lubricants inhibit rust and prolong the life of surgical instruments. Weiman® lubricants contain dual-purpose corrosion inhibitors that protect all types of metal instruments, including carbon steel. Do not interfere with steam sterilizing.

- WEIT2 1 Gallon





The PrimaTouch® Medical Examination and Procedure Glove line offers a high quality, comprehensive collection of gloves for every procedure requiring medical-grade hand protection.

As a focused single-use glove supplier, priMED is committed to staying on top of the latest glove trends and providing clinicians with the most cutting edge, single-use hand protection available.

FIT NITRILE 200/BX

- PriMed PrimaTouch® Fit A latex-free alternative,
- Fully micro-textured with enhanced finger texturing,
- Tested for use with chemotherapy drugs,
- 3mm thickness, Powder-Free
- Purple color



- PRIPM6-2005.....Extra Small
- PRIPM6-2006.....Small
- PRIPM6-2007.....Medium
- PRIPM6-2008.....Large
- PRIPM6-2009.....Extra Large (180/BX)

MOVE NITRILE GLOVE 150/BX

PriMed Move nitrile exam glove dons easily and conforms to every movement of the wearer's hand. This 100% latex-free glove fits, feels, and protects like a latex glove without the risk of Type 1 latex allergies.



- PRIPM6-1505.....X-Small
- PRIPM6-1506.....Small
- PRIPM6-1507.....Medium
- PRIPM6-1508.....Large
- PRIPM6-1509.....X-Large (130/BX)



OMNITRUST #212 SERIES NITRILE

Powder Free Examination Glove micro-texture fingertips and matte finish for added grip control without loss of tactile sensitivity. The #212 series is also tested for safe use with chemotherapy drugs.



OmniTrust Nitrile Small 100/BX

- OMI212-21..... Small
- OMI212-22..... Medium
- OMI212-23..... Large
- OMI212-24..... Extra Large

NOAH STRETCH NITRILE - BOX OF 100

- Finger Textured
- Nitrile
- Powder free
- 3.0mil Thickness

SPM88896Small
 SPM88897Medium
 SPM88898Large



COPOLYMER GLOVE

AMG014-200One Size Fits All Box/100



GLOVE BOX HOLDER WIRE

BOWGL011-0613Single Glove Wire
 BOWGL022-0613Double Glove
 BOWGL033-0613Triple Glove Wire
 BoWGL044-0613.....Quad Glove Wire



GLOVE BOX HOLDER PLASTIC

BOWGP-013Single Glove Plastic
 BOWGP-014Double Glove Plastic
 BOWGP-015Triple Glove Plastic
 BOWGP-016Quad Glove Plastic



TRIPLE GLOVE BOX HOLDER - WIRE

Holds three boxes of standard examination gloves.
 Dimensions: 17.1 x 4.25 x 7.5 inch (43.7 x 10.8 x 19 cm) White Wire
 ** No hardware with product**

AMG011-910 Triple Glove Box Holder - Wire



VIBRANT LATEX - BOX OF 100

- Powder free
- Beaded cuff eliminates roll-up
- Micro Textured
- Non-Sterile
- 4.5 mil Thickness

SPM98225	Extra Small
SPM98226	Small
SPM98227	Medium
SPM98228	Large
SPM98229	Extra Large

AURELIA
Sign of Comfort™



DELIGHT VINYL - BOX OF 100

- Powder free
- Smooth finish
- Latex-free
- hypoallergenic
- Protein-free
- Non-sterile
- clear white

SPM38225	Extra Small
SPM38226	Small
SPM38227	Medium
SPM38228	Large
SPM38229	Extra Large



BLUSH® NITRILE PINK EXAMINATION GLOVES - BOX OF 200

- Nitrile (100% Latex-Free)
- Non Sterile
- Glove Design & Feature: Ambidextrous and Finger Textured
- Cuff: Beaded
- Color: Blush Pink
- Thickness: 2.5 mil
- Packaging: 200 pieces per box / 10 boxes per case

SPM78885	Extra Small
SPM78886	Small
SPM78887	Medium
SPM78888	Large



BOLD NITRILE BLACK - BOX OF 100

- Latex-Free
- Non-Sterile
- Fully Textured
- Black Color
- Ambidextrous

SPM73995	X-Small
SPM73996	Small
SPM73997	Medium
SPM73998	Large
SPM73999	X-large



NITRIDERM® NITRILE STERILE EXAM GLOVES

Nitrile synthetic contains no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort. Pair packaged, inner wrapper in ozone resistant multi-vac packaging.

IHC106100.....	Small
IHC106200.....	Medium
IHC106300.....	Large
IHC106350.....	X-Large
50 Pairs/Box	



PRESTIGE® DHD LATEX SURGICAL GLOVES

Damp don gel interior provides effortless donning. Curved finger design increases comfort. Powder-Free to reduce risk of powder induced allergies, adhesion, and granulomas.

IHC139550.....	Size 5.5
IHC139600.....	Size 6.0
IHC139650.....	Size 6.5
IHC139700.....	Size 7.0
IHC139750.....	Size 7.5
IHC139800.....	Size 8.0
IHC139850.....	Size 8.5
IHC139900.....	Size 9.0
50 Pairs/Box	



NITRIDERM® NITRILE STERILE EXAM GLOVES SINGLE

Nitrile synthetic contains no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort. Pair packaged, inner wrapper in ozone resistant multi-vac packaging.


IHC109100.....	Small
IHC109200.....	Medium
IHC109300.....	Large
IHC109350.....	E-Large
100/BX	



NITRIDERM® NITRILE STERILE EXAM GLOVE

NitriDerm Nitrile Surgical Gloves, These nitrile synthetic surgical gloves contain no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort, fit, and protection, while offering exceptional tactile sensitivity.

IHC135260.....	Size 6.0
IHC135265.....	Size 6.5
IHC135270.....	Size 7.0
IHC135275.....	Size 7.5
IHC135280.....	Size 8.0
IHC135285.....	Size 8.5
IHC135290.....	Size 9.0
50 Pairs/Box	



PRESTIGE® CHLOROPRENE SURGICAL GLOVES

Made from non-latex, polychloroprene synthetic polymer, these new surgical gloves offer all of the benefits of natural rubber latex, with none of the allergen problems associated with latex. The unique formulation allows them to be worn for extended periods without hand fatigue. Slightly roughened finish provides the ideal combination of gripping ability and tactility

IHC134550.....	Size 5.5
IHC134600.....	Size 6.0
IHC134650.....	Size 6.5
IHC134700.....	Size 7.0
IHC134750.....	Size 7.5
IHC134800.....	Size 8.0
IHC134850.....	Size 8.5
IHC134900.....	Size 9.0
Box/25 Pairs	



DERMASSIST® LATEX SURGICAL GLOVES

IHC133550.....	Size 5.5
IHC133600.....	Size 6.0
IHC133650.....	Size 6.5
IHC133700.....	Size 7.0
IHC133750.....	Size 7.5
IHC133800.....	Size 8.0
IHC133850.....	Size 8.5
IHC133900.....	Size 9.0
Box/50 Pairs	



LATEX SURGICAL GLOVES POWDER FREE

- One pair (one right, one left) per paper wallet
- Maximum strength with no loss of tactile sensitivity
- Textured surface in palm and finger tips
- Powder Free
- Sterile

DYN2360.....	Size 6.0
DYN2365.....	Size 6.5
DYN2370.....	Size 7.0
DYN2375.....	Size 7.5
DYN2380.....	Size 8.0
DYN2385.....	Size 8.5
DYN2390.....	Size 9.0
Box/50 Pairs	



MAGNA DISPOSABLE SCALPELS

These scalpels are individually packaged and are attached to a plastic handle with a plastic blade cover for protection.

- ALMM92-10 #10
- ALMM92-11 #11
- ALMM92-12 #12
- ALMM92-15 #15
- ALMM92-20 #20
- ALMM92-21 #21
- ALMM92-22 #22
- ALMM92-23 #23

Box/10



MEDI-CUT® DISPOSABLE SAFETY SCALPELS

Blade retracts safely with one simple motion, and can be permanently locked and covered for safe disposal.

- DYN4160..... #10
- DYN4161..... #11
- DYN4165..... #15

Box/10



MAGNA STAINLESS STEEL SURGICAL BLADES

- ALMM90-10#10
- ALMM90-11#11
- ALMM90-12#12
- ALMM90-15#15
- ALMM90-20#20
- ALMM90-21#21
- ALMM90-22#22
- ALMM90-23#23
- ALMM90-24#24
- ALMM90-25#25

Box/100



SCALPEL HANDLE STAINLESS STEEL

- ALMM36-0100..... #3 (Fits 10, 11, 12, 15 blades)
- ALMM36-0120..... #4 (Fits 20 - 25 blades)




#3



#4

MILTEX® DISPOSABLE BIOPSY PUNCHES

MIL33-31AA.....	1.0mm
MIL33-31A.....	1.5mm
MIL33-31.....	2.0mm
MIL33-31B.....	2.5mm
MIL33-32.....	3.0mm
MIL33-33.....	3.5mm
MIL33-34.....	4.0mm
MIL33-35.....	5.0mm
MIL33-36.....	6.0mm
MIL33-37.....	7.0mm



DYNAREX DISPOSABLE BIOPSY PUNCHES

DYN4086.....	1.5mm
DYN4087.....	2.0mm
DYN4088.....	2.5mm
DYN4089.....	3.0mm
DYN4090.....	3.5mm
DYN4091.....	4.0mm
DYN4093.....	5.0mm
DYN4094.....	6.0mm
DYN4096.....	8.0mm



BIOPBLADE

Used for skin/tissue saucerization and shave biopsy, sterile, disposable. Flexible design allows the blade to be positioned at various angles.

MIL33-100



BiopBlade™

SUTUREOUT KIT

SutureOut is a sterile, disposable product that is intended for removing multiple stitches from a single patient. Kit includes SutureOut tool, plastic forcep and 1ea 4x4 gauze.

CHS85-4337



BRYMILL LIQUID NITROGEN GUN

Units come with four spray tips with apertures of various sizes (.04 - .016 inches) plus a 20-gauge straight spray and 20-gauge bent spray for more accurate, controlled treatments

- CYGCRY-AC..... 500mL
- CYGCRY-AC-3..... 350mL
- CYGCRY-BABY.....150ml "BABY"



LIQUID NITROGEN TANKS AND ACCESSORIES

- CYGTANK XC2020 Litre Dewar Tank
- CYGD12.....Dipper Cup for Liquid Nitrogen
- CYGCRYO GLOVESLiquid Nitrogen Safety Gloves

HISTOFREEZER® PORTABLE CRYOSURGICAL SYSTEM

- On the spot treatment for common warts, skin tags, benign lesions and plantar warts
- 2 year shelf life
- Freezing Temp: -55°C
- Each kit includes 2 x 80ml cannister. 24ea 2mm foam tips for small lesions, 36ea 5mm foam tips for large lesions
- 40 treatments per kit
- Instruction manual



PAL2879100580

CANTHARONE® WART TREATMENT SOLUTIONS

- For topical use of removal of warts.

- DOR9001-975M..... Cantharone Wart Treatment 7.5ml
- DOR9002-975M..... Catharone Plus Wart Treatment 7.5ml





LOOK™ SUTURES

Nylon

LOO917B 5-0 Nylon with C17 Needle Box/12
 LOO921B 4-0 Nylon with C17 Needle Box/12
 LOO916B 6-0 Nylon with C17 Needle Box/12



SURGERY & MINOR PROCEDURES

Reverse Cutting Needles

	Look	Ethicon
11mm	3/8 Circle C1	
12mm	C17	C-2
14mm	•	•
16mm	C22	FS-3
18mm	C6	FS-2
24mm	C7	FS-1
26mm	C8	FS
1/2 Circle	C31	X1

Precision Reverse Cutting Needles

	Look	Ethicon
7mm	•	G-6
8mm	•	P-6
11mm	C1	P-1
13mm	C3	P-3
16mm	PC34	PS-3
18mm	PC31	PS-2
23mm	PC33	PS-1

Canada Medical carries the complete line of Look® and Ethicon Sutures.
 Contact us for item codes and pricing.

ETHICON

a *Johnson & Johnson* company

ETHICON™ SUTURES

Chromic Gut

- ETH634G..... 5-0 Chromic Gut with FS-2 Needle..... Box/12
- ETH635H..... 4-0 Chromic Gut with FS-2 Needle..... Box/36
- ETH636H..... 3-0 Chromic Gut with FS-2 Needle..... Box/36

Plain Gut

- ETHH820G..... 5-0 Plain Gut with FS-2 Needle Box/12
- ETHH821H..... 4-0 Plain Gut with FS-2 Needle Box/36
- ETHH822H..... 3-0 Plain Gut with FS-2 Needle Box/36

Ethilon

- ETH668G..... 5-0 Nylon with C-2 Needle..... Box/12
- ETH661G..... 5-0 Nylon with FS-2 Needle..... Box/12
- ETH662G..... 4-0 Nylon with FS-2 Needle..... Box/12

Silk

- ETH682G..... 5-0 Silk with FS-2 Needle..... Box/12
- ETH683G..... 4-0 Silk with FS-2 Needle..... Box/12
- ETH684G..... 3-0 Silk with FS-1 Needle..... Box/12

Prolene

- ETH8661G..... 5-0 Prolene with FS-2 Needle Box/12
- ETH8683G 4-0 Prolene with FS-2 Needle Box/12
- ETH8665G..... 3-0 Prolene with FS-2 Needle Box/12



Reverse Cutting Needles

	Look	Ethicon	Syneture
	3/8 Circle	C1	C-1
		C17	C-2
		•	•
		C22	FS-3
		C6	FS-2
		C7	FS-1
		C8	FS
	1/2 Circle	C31	X1

Precision Reverse Cutting Needles

	Look	Ethicon	Syneture
	3/8 Circle	•	G-6
		•	P-6
		C1	P-1
		C3	P-3
		PC34	PS-3
		PC31	PS-2
		PC33	PS-1

Canada Medical carries the complete line of Ethicon™ Sutures.

Contact us for item codes and pricing.

3M STERI-STRIP™ SKIN CLOSURES

3M R1540 1/8" x 3" 5/Env Box/50 Env
 3M R1541 1/4" x 3" 3/Env Box/50 Env
 3M R1542 1/4" x 1-1/2" 6/Env Box/50 Env
 3M R1547 1/2" x 4" 6/Env Box/50 Env



PRIME3 PRIMAFLEX SKIN CLOSURES

PRIPM3-71831... 1/8" x 3" 5/Env Box/50 Env
 PRIPM3-71431.. 1/4" x 3" 3/Env Box/50 Env



GLUSTITCH TISSUE ADHESIVE

BLK200..... Violet 0.2mL Box/12
 BLK201..... Clear 0.2mL Box/12



TWIN BLADE RAZOR

- Long-lasting, Superior Quality Blade
- Clear plastic guard
- Disposable
- Unisex

DYN4250 Box/50



PREP RAZOR

- Double Edge Premium Razor

PER75-4004 Personna prep razor Box/100



SAFETY PREP RAZOR

DYN4251 Gallant Safety Prep Razor Box/50



SURGICAL SKIN MARKERS WITH RULER, STERILE

SRGS-3000-R0 Multi Tip Box/50



3M MULTI-SHOT SKIN STAPLERS

3M DS-15 15 Shot, Disposable Box/12



SUTURE REMOVAL TRAY, STERILE

Suture Removal Trays help you run your procedure more efficiently and help reduce the risk of cross-contamination.

DYN4521-IMCSuture Removal Tray



DYN4521-IMC

SKIN STAPLE REMOVER, DISPOSABLE, STERILE

BUS716-IMCSkin Staple Remover



BASIC DRESSING TRAY, STERILE

This tray contains:

- 1 Three Compartment Tray
- 1 Plastic Transfer Forceps
- 2 Metal Forceps
- 4 10 x 10 cm Gauze Sponges
- 4 5 x 5 cm Gauze Sponges
- 1 Waste Bag with Tie
- 1 Safety Pin
- 2 Towel/Underpad

TRA85-2000TY



IUD INSERTION KIT, STERILE

Kit Includes:

- 1ea Scissors Metal, 7" Curved
- 1ea Tenaculum Metal, 9.5" Straight
- 1ea Uterine Sound Metal
- 10ea 4x4" Gauze
- 30x30" Sterile Wrap

TRASMA30IUOPB .



LACERATION TRAYS, DISPOSABLE, STERILE

- One Tray
- One 5" floor-grade, Baumgartner (serrated jaw) needle holder
- One 4-1/2" sharp/sharp scissors
- One 5" metal tissue forceps
- Two absorbent towels
- Five 4" x 4" gauze sponges
- One polylined sterile field
- One polylined fenestrated drape

BUS751




XYLOCAINE®

- APNVR89896CA1Xylocaine 1% Plain 10ml 50/BX
- APNVR44571CA1Xylocaine 1% Plain 20ml
- APNVR89892CA1Xylocaine 1% Plain 2ml Poly Ampoule 50/BX
- APNVR44569CA1Xylocaine 1% Plain 50ml
- APNVR89893CA1Xylocaine 1% Plain 5ml Poly Ampoule 50/BX Box/50 vials
- APNVR44553CA1Xylocaine 1% with Epinephrine 20ml
- APNVR89897CA1Xylocaine 2% Plain 10ml 50/BX
- APNVR44567CA1Xylocaine 2% Plain 20ml
- APNVR89894CA1Xylocaine 2% Plain 2ml Poly Ampoule 50/BX
- APNVR44565CA1Xylocaine 2% Plain 50ml
- APNVR89895CA1Xylocaine 2% Plain 5ml 50/BX
- APNVR44551CA1Xylocaine 2% with Epinephrine 20ml
- APNVR89868CA1Xylocaine Endotracheal Spray 10mg/50ml
- APNVR44900CA1Xylocaine Jelly 2% 30ml Tube
- APNVR44927CA1Xylocaine Ointment 5% 35gm Tube
- APNVR44677CA1Xylocard 100mg 5ml 10amp/box



MARCAINE®

- WIN02073020... 20mL, 0.25% Plain, Single Use
- WIN02083020... 20mL, 0.5% with epi, Single Use Box/10
- WIN02080020... 20mL, 0.25% with epi, Multi-Dose Box/10
- WIN1587050 50mL, 0.25% Plain, Multi-Dose
- WIN1610050 50mL, 0.5% Plain, Multi-Dose



LIDOCAINE

- TEL0122AG01Lidocaine 1% Plain 20ml
- TEL0122AJ01Lidocaine 1% Plain 50ml
- TEL0126AF01Lidocaine 2% Plain 10ml Tray/20
- TEL0127AG01Lidocaine 2% Plain 20ml
- TEL0127AJ01Lidocaine 2% Plain 50ml
- TEL0126AD01Lidocaine 2% Plain 5ml 20/bx Box/20
- TEL0128AG01Lidocaine 2% with Epinephrine 20ml



PAIN EASE® INSTANT TOPICAL ANESTHETIC

Gebauer's Pain Ease® instant topical anesthetic skin refrigerant (vapocoolant) is FDA cleared to temporarily control the pain associated with needle procedures and minor surgical procedures such as immunizations, allergy injections, foreign body removal, suturing and suture removal, IV starts and blood draws. Pain Ease can be used prior to, or in some cases in place of lidocaine injections.

SEK0008-023.5oz Mist Spray
 SEK0008-033.5oz Medium Spray



EMLA® CREAM

- Topical Anesthetic

APNVR44812CA130g tube



MINIMS® TETRACAINE HYDROCHLORIDE DROPS

BL 10666CA0.5% Pkg/20
 BL 10667CA1.0% Pkg/20



MINIMS® FLOURESCEIN SODIUM DROPS

BL 10692CA0.5mL Drop, 2% Pkg/20



DISPOSABLE TOWEL DRAPES, STERILE

- For use in minor surgical or therapeutic procedures where a sterile field or drape is required
- Poly-interlined, three-ply construction - tissue/poly/tissue - provides absorbency while preventing possible moisture or bacteria migration
- Eliminates all labor costs associated with cloth towels or drapes
- Available in both fenestrated and plain

DYN4409-IMCFenestrated, 18"x26" Box/50
 DYN4410-IMCPlain, 18"x26" Box/50



PROVIDINE IODINE SOLUTION 10%

Providine Iodine solution is an antiseptic cleanser for wound cleansing. Effective in destroying certain bacteria and removing impurities to provide antiseptic cleansing. Air dries quickly requiring no wiping. For external use only.

LER109-08500mL
 LER109-09115mL



BAXEDIN® PRE-OP ANTISEPTIC SOLUTION

- 0.5% Chlorhexidine Gluconate (CHG) with 70% Isopropyl Alcohol
- For skin asepsis prior to surgery
- Disinfectant

OMGL000003.....500mL Bottle, Tinted



STANHEXIDINE® LOW SUDSING AQUEOUS SOLUTION

- 2% Chlorhexidine Gluconate (CHG) with 4% Isopropyl Alcohol
- Anti-bacterial for the skin
- Treatment for superficial skin infections

OMGL000009450mL Bottle



BENZALKONIUM CHLORIDE WIPES

For use wherever antiseptic cleansing/refreshing is needed for sting free sensitive skin areas. Saturated with benzalkonium chloride 0.13%. Skin antiseptic.

LER126-01Wipes Box/100



PVP PREPS & SWABS

PDIB51200Povidone Iodine 10% Prep Pad Box/100
 PDIS42050.....Povidone Iodine 10% Swabstix Box/50



NYLON SCRUB BRUSH

AMG118-140.....



E-Z SCRUB SPONGE

Proprietary polyethylene soft bristles enhance mechanical action while reducing skin irritation. Soft sponge provides gentle cleaning on irritated areas or when firm mechanical action is not required.

BD 3774794% Chlorhexidine Box/30



Baxter NORMAL SALINE (0.9% SODIUM CHLORIDE) I.V. SOLUTIONS

- BAXJB1322P250mL
- BAXJB1323500mL
- BAXJB13241000mL



Baxter LACTATED RINGERS, I.V. SOLUTIONS

- BAXJB2322P250mL
- BAXJB2323500mL
- BAXJB23241000mL



Baxter DEXTROSE 5%, I.V. SOLUTIONS

- BAXJB0062P250mL
- BAXJB0063500mL
- BAXJB00641000mL



Baxter DEXTROSE 3.3% & 0.3% SALINE, I.V. SOLUTIONS

- BAXJB1033500mL
- BAXJB10341000mL



B|BRAUN SHARING EXPERTISE **STERILE WATER, I.V. SOLUTIONS**

- TRABBL8502250mL
- TRABBL8501-01500mL



Baxter MINI BAG® I.V. SOLUTIONS

Dextrose 5% Water Solution

- BAXJB0081P 50mL per bag, 4 bags/pack
- BAXJB0082P 100mL Case/96

Sodium Chloride 0.9% Solution


- BAXJB1301PK 50mL per Bag, 4 bags per pack Pack
- BAXJB1302PK 100mL per Bag, 4 bags per pack Pack



Baxter 4 WAY STOPCOCKS

Stopcock 4-way large bore, with rotating male luer lock adapter.

BAX2C6204 Case/50



Baxter CLEARLINK® NEEDLE-FREE I.V. SETS

IV Needleless administration set with Clearlink® luer activated valve and male luer lock adapter. 76" Length, 10 Drops/mL

BAX2C8401 Case/48



IV Needleless Continu-Flo solution set with 2 Clearlink luer activated valve and male luer lock adapter. 104" Length, 10 Drops/mL

BAXJC8519 Case/48



Baxter SOLUTION SET (10 DROPS / ML)

BAXJC8109..... with Luer Lock Adapter (101")



Baxter I.V. EXTENSION SETS

IV Microbore catheter, 8" extension set with Interlink® injection site and luer lock adapter.

BAX2N3374 Case/200



STERILE WATER FOR IRRIGATION

BAXJF7113 500mL Pour Bottle
 BAXJF7114 1000mL Pour Bottle



STERILE WATER FOR INJECTION

Bacteriostatic
 WIN03977030 30mL vial

No Preservative
 WIN04887010 10mL vial
 WIN04887020 20mL vial Cs/25



NORMAL SALINE 0.9% (NaCl) FOR IRRIGATION

BAXJF7123 500mL Pour Bottle
 BAXJF7124 1000mL Pour Bottle



NORMAL SALINE 0.9% (SODIUM CHLORIDE) FOR INJECTION

No Preservative
 WIN04888010 10mL vial
 WIN04888020 20mL vial Case/25

Bacteriostatic
 WIN01966010 10mL vial Case/25
 WIN01966020 20mL vial
 WIN01966030 30mL vial



Baxter MINI BAG® SOLUTION SETS

Dextrose 5% Water Solution
 BAXJB0081P 50mL per bag, 4 bags/pack
 BAXJB0082P 100mL Case/96


Sodium Chloride 0.9% Solution
 BAXJB1301PK 50mL per Bag, 4 bags per pack Pack
 BAXJB1302PK 100mL per Bag, 4 bags per pack Pack



BD INSYTE™ AUTOGUARD™ BC PRO WITH BLOOD CONTROL TECHNOLOGY

New blood-control technology proven to reduce the risk of blood exposure by 95%

- No clean up of blood spills
- No need to worry about stopping blood
- No add-on costs for clean-up supplies
- Push-button shielding technology proven to reduce needlestick injuries by 95%
- BD Vialon™ Biomaterial—softens up to 70% in the vessel, resulting in longer dwell times; reduces the chance of mechanical phlebitis by up to 50%
- No risk of exposure to blood-borne pathogens
- Due to dramatically less blood leakage at insertion,
- BD Insyte™ Autoguard™ BC can mean less potential cross contamination; no spills means there's no blood to get on gloves, clothing, shoes, bedding, floors, or equipment.



BD 381012	24g x 0.75" Yellow
BD 381023	22g x 1" Blue
BD 381033	20g x 1" Pink
BD 381034	20g x 1.16" Pink
BD 381037	20g x 1.88" Pink
BD 381044	18g x 1.16" Green
BD 381047	18g x 1.88" Green
BD 381054	16g x 1.16" Gray
BD 381057	16g x 1.77" Gray

BD **INSYTE™ I.V. CATHETER**

New Standard in I.V. therapy. When faced with the everyday challenges of infusion therapy, Insyte™ I.V. catheters deliver.

- Increases the ease of penetration
- Reduces the risks of infiltration
- Exhibits greater kink resistance
- Reduces the risks of phlebitis up to 50%
- Reduces I.V. restarts due to reduced complications
- Reduces the total I.V. therapy cost to use due to longer indwelling times

BD 38121224g x 0.75" Yellow
BD 38122322g x 1.00" Blue
BD 38123320g x 1.00" Pink
BD 38123420g x 1.16" Pink
BD 38123720g x 1.88" Pink
BD 38124418g x 1.16" Green
BD 38124718g x 1.88" Green
BD 38125416g x 1.16" Gray
BD 38125716g x 1.77" Gray
BD 38126714g x 1.75" Orange

Box/50



BD **INSYTE™ AUTOGUARD™ WINGED I.V. CATHETER**

Protect from blood leakage

The BD Insyte™ Autoguard™ BD shielded IV catheter with blood control technology is clinically demonstrated to reduce the risk of blood exposure by 95%. BD blood control technology minimizes blood exposure, cleanup time and use of supplies associated with blood spills. BD Instaflash™ needle technology on the 20-G, 22-G and 24-G sizes incorporates a notched needle to improve first-stick success and reduce painful hit-and-miss insertions.

- No venous compression required during insertion
- Less mess and minimizes healthcare worker risk resulting from blood leakage
- Less blood leakage means less disruption during the patient's care

With Wings

BD 381944 18g x 1.16"
BD 381933 20g x 1"
BD 381934 20g x 1.16"
BD 381923 22g x 1"
BD 381912 24g x .75"

Box/50



BD **VACUTAINER® SAFETY-LOK™ BLOOD COLLECTION SETS**

The BD Vacutainer® Safety-Lok™ Blood Collection Set is simple, easy to use, and is safety-engineered. Single-handed safety activation with irreversible locking mechanism, confirmed with an audible click.

BD 367296 21g x 0.75", 12" tubing
BD 367297 23g x 0.75", 12" tubing
BD 367298 25g x 0.75", 12" tubing


Box/50



TERUMO® SURFLO® WINGED INFUSION SET

The soft, interlocking wings of these infusion sets readily conform to the body's contour and are angled for easy positioning of the needle. The precisely honed, ultra-thin needle design assures ease of penetration and continued access with a single insertion. Coil-resistant tubing reduces recoil during use and is both latex and DEHP free.

TERSV.21BLK.....	21g x 0.75", 12" tubing	Box/50
TERSV.22BLK.....	22g x 0.75", 12" tubing	
TERSV.23BLK.....	23g x 0.75", 12" tubing	
TERSV.25BLK.....	25g x 0.75", 12" tubing	




TERUMO® SURSHIELD® SAFETY WINGED INFUSION SET

Terumo's Surshield Safety Winged Needle Sets with One-Handed Activation provides a solution with a cost-effective, easy to use design. Press shield over the needle until it locks into place.

- Unique needle locking mechanism helps reduce the risk of needlestick injuries.
- One handed activation for safety and ease of use
- No change in needle insertion technique required

TERSV.S19BL	19g x 3/4", 12" tubing	Box/50
TERSV.S19BLS	19g x 3/4", 3.5" tubing	
TERSV.S21BL	21g x 3/4", 12" tubing	
TERSV.S21BLS	21g x 3/4", 3.5" tubing	
TERSV.S23BL	23g x 3/4", 12" tubing	
TERSV.S23BLS	23g x 3/4", 3.5" tubing	
TERSV.S25BL	25g x 3/4", 12" tubing	
TERSV.S25BLS	25g x 3/4", 3.5" tubing	



TERUMO® SURSHIELD® SAFETY I.V. CATHETER

Terumo's innovative catheter is the only one with Surflash, the exclusive catheter placement indicator, plus a passive safety feature.

- Clear visual confirmation of catheter placement, for increased confidence and first-stick success
- Reduced risk of needlestick injury - without adding extra steps to your procedures
- Ultra-sharp, double-bevel needle and medical-grade lubricant for smoother venipuncture

TERSRSFA1832A	18g x 1.25"	Box/50
TERSRSFA2032A	20g x 1.25"	
TERSRSFA2225A	22g x 1"	
TERSRSFA2419A	24g x 0.75"	



TOURNIQUETS, DISPOSABLE

Latex Free, 1x18"

PRIPM5-1018 Box/25



3M TEGADERM® TRANSPARENT FILM DRESSING

3M 1622W.....4.4cm x 4.4cm Box/100

3M 1624W.....6cm x 7cm Box/100



3M SoluPrep™ ANTISEPTIC SWABS

Chlorhexidine gluconate and 70% v/v alcohol

3M 101062% CHG & Alcohol Pad Box/200

3M 102020.5% CHG & Alcohol Swabstick Box/50

3M 102032% CHG & Alcohol Swabstick Box/50



BD I.V. ACCESS CANNULA

The BD Plastic Cannula is a truly universal cannula and can be used with all types of split-septum IV injection sites and vials designed for needleless access. It provides simplified access to the Interlink injection site with smoother insertion, reduced penetration force and greater fluid dispersion.

BD 303345 Blunt Plastic Cannula Box/100



BD PRN LUER LOCK IV ADAPTERS

Does not contain natural rubber latex. 0.1mL/min. Sterile.

BD 385111Short 0.75" Injection site, sterile, latex free Box/50

BD 385112Mid-Size 0.0625" Injection site, sterile, latex free Box/50



BD Q-SYTE™ IV NEEDLELESS CONNECTOR, LUER ACCESS

The BD Q-Syte™ Luer Access Split Septum combines the convenience of a Luer access with the benefits of a split septum design to provide closed access for IV therapy. The benefits of this unique split septum design include reliable performance, high flow rates, and a visible fluid path.

BD 385100Split Septum Box/50



DELUXE I.V. POLE, 4 LEG

This I.V. has a four point base, twin wheeled swivel casters, and a circular "Ram's Horn" design. As well, the Economy 2 Hook I.V. Stand features a positive friction grip, chrome plated steel pole, stable base design and friction grip lock adjustment.

AMG775-751I.V. Pole



CLARITY AUTOMATIC URINE ANALYZER

- Large LCD Display
- Auto calibration and self-check
- Flagging of abnormal patient results
- Internal Memory of last 2000 Tests
- RS 232 port or USB for connectivity



CLACD-MCPROMO5 Clarity Urocheck Urine Analyzerch
 Includes: Analyzer, 4 Bottles of 10SG, and 1 Bottle of Microalbumin
 CLACD-UCTL30 Clarity Urocheck Urine Control 10ml 3/BXch
 For use with Clarity Urine Strips and Analyzers only, REFRIGERATE UPON ARRIVAL.
 Control solution must be used with each new box of strips.
 CLADTG-10SG..... Clarity Urocheck urine strips 10SG 100/Bt Bottle
 CLADTG-UAPPR..... Clarity Analyzer paper rolls 4/PK Pack

LABORATORY

BIOSTRIP® URINE TEST STRIPS

URINE TEST STRIP	BLOOD	LEUKOCYTES	UROBILINOGEN	BILIRUBIN	NITRATE	SPECIFIC GRAVITY	GLUCOSE	PROTEIN	KETONES	pH
4192-2P										
4192-U7										
4192-U10										

SSP4192-2P.....Urinalysis for protein and glucose Bottle/100
 SSP4192-U7Urinalysis reagent strip, 7 parameters Bottle/100
 SSP4192-U10.....Urinalysis reagent strip, 10 parameters Bottle/100



SIEMENS® URINE TEST STRIPS

URINE TEST STRIP	BLOOD	LEUKOCYTES	UROBILINOGEN	BILIRUBIN	NITRATE	SPECIFIC GRAVITY	GLUCOSE	PROTEIN	KETONES	pH
Labstix										
Multistix 5										
Multistix 8 SG										
Multistix 10 SG										
Uristix										

SEM2300AMultistix 10 SG Bottle/100
 SEM2304AMultistix 8 SG Bottle/100
 SEM2309Multistix 5 Bottle/100
 SEM2855UUristix Bottle/100



CHEMSTRIP™ URINE TEST STRIPS

URINE TEST STRIP	SG	pH	Leukocytes	Nitrite	Protein	Glucose	Ketone	Urobilinogen	Bilirubin	Blood
Chemstrip 10										
Chemstrip 9										
Chemstrip 7										
Chemstrip 5										
Chemstrip GP										

BMC1120347911 Chemstrip 10 Bottle/100
 BMC11008552119 Chemstrip 7 Bottle/100
 BMC11893467119 Chemstrip 5 Bottle/100
 BMC1020074311 Chemstrip GP Bottle/100



VERITOR® STREP & FLU TEST KIT

For Rapid Detection of Flu A+B & Strep

- Technology-verified Performance-Incorporates state of the art technologies that deliver results you can trust.
- Multiple Sample Types-Validated to be used with nasopharyngeal washes/aspirates, and nasopharyngeal swab in transport
- Convenient Kit Configuration-Full customized kit to test and delivery results in minutes
- Objective Results-Easy to read digital display eliminates subjective interpretation



Convenient Workflow Delivers Results in Minutes



Easy sample processing
Prefilled unitized tube facilitates workflow



Ready in minutes
Test device is ready to insert into reader 10 minutes after sample is added



Insert and read
Simple one-touch button readies the reader for test device insertion



Objective results
Once the test device is inserted in the reader, an objective, digital test result is displayed in 10 seconds

BD 256055	Veritor System Reader
BD 256040	Veritor Strep Test Kits Box/30
BD 256045	Veritor Flu A+B Test Kits Box/30
BD 256089	Veritor SARS COVID-2 TEST Box/30

MEDISURE MULTI USER BLOOD GLUCOSE MONITOR

- Results in six seconds or less
- Stores 960 test results
- Tiny 0.6 uL blood sample size
- Large LCD display. Large buttons with simple functions

MDIDG001MM	Glucose Meter
MDIDG001S	Test Strips 100/BX
MDIDG001C	Control Solution



STREP-A TEST KIT

SSP4010-IST-501Box/25



PREGNANCY TEST STRIPS

SSP4001-422-50.....Box/50



PREGNANCY TEST, CASSETTE

SSP4001-441-25.....Box/25



COVID-19 ANTIGEN RAPID TEST 25 TESTS/BOX

Our new Rapid Response® COVID-19 Antigen Rapid Test Device allows for an efficient and cost-friendly method of detection. In just 15 minutes, the test can detect the SARS-CoV-2 viral nucleoprotein antigen through visual interpretation. Non Returnable

GLMCOV-19C25Box/25



COVID-19 ANTIBODY RAPID TEST 25 TEST/BOX

The Rapid Response® COVID-19 IgG/IgM Rapid Test Device with just a prick of blood the test device can rapidly & qualitatively test for the presence of anti-SARS-CoV-2 IgM and anti-SARS-CoV-2 IgG antibodies. Results in 15 minutes. FINAL SALE - No Returns

GLMCOV-13C25Box/25



BIOSTRIP® LH OVULATION PREDICTION TEST

SSP4220-13 Box/6 Tests



DRUGS OF ABUSE TEST KIT

SSP4192-D700-25D.A.T. 7 Drug Adulteration Test Box/25
 SSP4215-DUA-184-125Quickcheck 8 Panel Urine Drug Test 25/KT
 SSP4215-DUA-1127A-125 EACH.....Quickcheck 12 cup Urine Drug Screen Test

* 10 Panels: AMP/BAR/BZD/BUP/COC/MTD/MET/OPI-II/PCP/THC



AMNICATOR® AMNIOTIC FLUID INDICATOR SWABS

- A rapid and simple method of detecting ruptured amniotic membranes during pregnancy.
- Highly sensitive indicator.
- Non threatening to patients
- Each to read-color change visable in seconds

SSP4080-MWAM1-10 Box/10



ROVERS® CERVEX-BRUSH® COMBI WITH REMOVER TUBE

SSP4008-38010100RT Bag/25 Brushes



CYTOLOGY BRUSH


SSP4201-CB-8B8" Non-Sterile Bag/100



ENDOMETRIAL BIOPSY SAMPLER

Bioteque's 3.0 mm Preferred Curette is used to obtain samples of mucosal lining for histological study in cancer screening, endometrial dating, or bacterial culturing.

MDG022721MEDGYN 3mm Currette Sampler
 BTQM0014BIOTEQUE 3mm Currette Sampler



MEDGYN ENDO SAMPLER STERILE KIT

MDG0227203mm with 10cc Syringe



URINE SPECIMEN CONTAINERS, STERILE

- Individually Wrapped

DYN4253-IMC..... 4oz Bag/100



URINE COLLECTOR

- Fits securely on standard toilets

DYN4232



CONDOMS

ANB5711 Non-Lubricated Box/144

ANB1511 Lubricated Box/144



SPECIMEN BIO-HAZARD TRANSPORT BAGS

VWR11215-684 6" x 9" Cs/500



CYTOSPRAY FIXATIVE SPRAY

Ensures cells adhere to the cytology slide during staining and that cell morphology is maintained.

SSP4240-CYT-017..... 100mL



ZIP-IQ TT TEST TUBE CENTRIFUGE

- The LW Scientific Zip-IQ TT centrifuge is an impressive 6-place test tube centrifuge for spinning blood, urine, and other fluids in 1.0ml to 15ml test tubes and microtubes, motor spins at 5,000rpm (2,500g) producing serum in as little as 5 minute

LWS06AD-15T3



BLOOD DRAWING CHAIR (PHLEBOTOMY CHAIR)

- Contoured padded seat and padded, arched backrest
- Padded flip arm and armrests
- Height adjustable armrest
- Comfortable, 2" padding throughout (5cm)
- Easy-clean, durable vinyl
- 1-1/4" square, heavy duty, all-welded, tubular steel frame
- Adjustable height and depth flip arm
- Flip arm positions in place when in use, then flips out of the way for easy exit
- Neutral Gray, powder-coated frame
- Adjustable foot levelers
- 400lbs. load capacity under normal use

SOLS573




LORIS SAFETY LANCETS, STERILE

- Tri-bevel tip for virtually painless sampling
- Needle automatically retracts
- For capillary blood samples
- No lancing device required
- Pressure activated lancet ready to use - no loading needed
- Latex free
- 1.4 mm puncture depth

LER201-01 Cs/100


Directions for Use



1 Twist off protective cap.



2 Gently push the safety lancet down against test site to activate.



3 Apply pressure near the puncture site to obtain the required blood volume.

PRICK LANCETTER ALLERGY LANCETS

1mm point length to optimize skin test response. Sterile and individually wrapped. Disposable. One lancet used per puncture to eliminate antigen carryover.

OMGLH8375ZA Box/100



BD ECLIPSE LUER LOCK SYRINGE & NEEDLE


Needle Only

BD 305757	30g x 1/2"
BD 305758	27g x 1/2"
BD 305767	25g x 1-1/2"
BD 305761	25g x 1"
BD 305759	25g x 5/8"
BD 305762	23g x 1"
BD 305763	22g x 1-1/2"
BD 305768	22g x 1"
BD 305765	21g x 1-1/2"
BD 305764	21g x 1"
BD 305766.....	18g x 1-1/2"

Syringe & Needle Combo

BD 305778	1ml 30g x 1/2"	Box/50
BD 305789	1ml 27g x 1/2"	Box/50
BD 305780	1ml 25g x 5/8"	Box/50
BD 305787	3ml 25g x 1"	Box/50
BD 305781	3ml 25g x 5/8"	Box/50
BD 305782	3ml 23g x 1"	Box/50
BD 305783	3ml 22g x 1-1/2"	Box/50
BD 305779	3ml 21g x 1"	Box/50
BD 305785	5ml 22g x 1-1/2"	Box/50
BD 305786	10ml 22g x 1-1/2"	Box/50

Safety Engineered Needles



NEEDLES AND SYRINGES

BD SYRINGES ONLY

BD 309659	1ml Slip tip	Box/200
BD 309628	1ml Luer lock	Box/100
BD 309657	3ml Luer lock	Box/200
BD 309646.....	5ml Luer lock	Box/125
BD 302995	10ml Luer lock	Box/200
BD 302830	20ml Luer lock	Box/48
BD 302832	30ml Luer lock	Box/58
BD 309653	60ml Luer lock	Box/40
BD 305819	Sterile Caps	Box/200



BD PRECISION GLIDE™ NEEDLES ONLY

BD 305128	30g x 1"	Box/100
BD 305106	30g x 1/2"	Box/100
BD 305136	27g x 1-1/4"	Box/100
BD 305109	27g x 1/2"	Box/100
BD 305111	26g x 1/2"	Box/100
BD 305122	25g x 5/8"	Box/100
BD 305125	25g x 1"	Box/100
BD 305127	25g x 1-1/2"	Box/100
BD 305145	23g x 1"	Box/100
BD 305155	22g x 1"	Box/100
BD 305156	22g x 1-1/2"	Box/100
BD 305167	21g x 1-1/2"	Box/100
BD 305195	18g x 1"	Box/100
BD 305196	18g x 1-1/2"	Box/100
BD 305198	16g x 1-1/2"	Box/100



BD ALLERGY TRAYS

BD 305535	1/2cc 27g x 1/2"
BD 305540	1cc 27g x 1/2"
BD 305541	1cc 27g x 3/8"



BD SYRINGE & NEEDLE COMBINATIONS

BD 309623	1ml 27g x 1/2"	Box/100
BD 309625	1ml 26g x 3/8"	Box/100
BD 309626	1ml 26g x 3/8"	Box/100
BD 309570	3ml 25g x 5/8"	Box/100
BD 309571	3ml 23g x 1"	Box/100
BD 309572	3ml 22g x 1"	Box/100
BD 309574	3ml 22g x 1-1/2"	Box/100
BD 309575	3ml 21g x 1"	Box/100
BD 309577	3ml 21g x 1-1/2"	Box/100
BD 309578	3ml 20g x 1"	Box/100
BD 309579	3ml 20g x 1-1/2"	Box/100
BD 309581	3ml 25g x 1"	Box/100
BD 309582	3ml 25g x 1-1/2"	Box/100



SOL-CARE™ SAFETY NEEDLE

- Easy, one-handed technique of activating the protective arm.
- Audible click indicates safety mechanism activation.
- Patented locking mechanism ensures needle is secure and completely locked in protective arm.
- Two options for activation of the protective arm.
- Compatible with all standard luer lock and slip tip syringes.
- 100/Box



SN3012.....SOL-CARE™ Safety Needle 30Gx.5"
SN2712.....SOL-CARE™ Safety Needle 27Gx.5"
SN2515.....SOL-CARE™ Safety Needle 25Gx1.5"
SN2510.....SOL-CARE™ Safety Needle 25Gx1"
SN2558.....SOL-CARE™ Safety Needle 25Gx5/8"
SN2315.....SOL-CARE™ Safety Needle 23Gx1.5"
SN2310.....SOL-CARE™ Safety Needle 23Gx1"

SN2215.....SOL-CARE™ Safety Needle 22Gx1.5"
SN2210.....SOL-CARE™ Safety Needle 22Gx1"
SN2115.....SOL-CARE™ Safety Needle 21Gx1.5"
SN2110.....SOL-CARE™ Safety Needle 21Gx1"
SN2015.....SOL-CARE™ Safety Needle 20Gx1.5"
SN2010.....SOL-CARE™ Safety Needle 20Gx1"
SN1815.....SOL-CARE™ Safety Needle 18Gx1.5"

SOL-CARE™ SAFETY NEEDLE WITH SYRINGE

12558SNSOL-CARE 1ml Luer Lock Syringe w/Safety Needle 25Gx5/8
12705SNSOL-CARE 1ml Luer Lock Syringe w/Safety Needle 27Gx.5
13005SNSOL-CARE 1ml Luer Lock Syringe w/Safety Needle 30Gx.5
32110SNSOL-CARE 3ml Luer Lock Syringe w/Safety Needle 21Gx1
32115SNSOL-CARE 3ml Luer Lock Syringe w/Safety Needle 21Gx1.5
32210SNSOL-CARE 3ml Luer Lock Syringe w/Safety Needle 22Gx1
32215SNSOL-CARE 3ml Luer Lock Syringe w/Safety Needle 22Gx1.5
32310SNSOL-CARE 3ml Luer Lock Syringe w/Safety Needle 23Gx1
32510SNSOL-CARE 3ml Luer Lock Syringe w/Safety Needle 25Gx1
32558SNSOL-CARE 3ml Luer Lock Syringe w/Safety Needle 25Gx5/8

SOL-M LUER LOCK SYRINGE WITHOUT NEEDLE

SOM180001PP.....SOL-CARE 1cc Luer Lock Syringe 100/BX
SOM180003.....SOL-CARE 3CC Luer Lock Syringe 100/BX
SOM180005.....SOL-CARE 5CC Luer Lock Syringe 100/bx
SOM180010.....SOL-CARE 10CC Luer Lock Syringe 100/BX
SOM180020.....SOL-CARE 20cc Luer Lock Syringe 100/BX
SOM180030.....SOL-CARE 30CC Luer Lock Syringe 50/BX
SOM180060.....SOL-CARE 60CC Luer Lock Syringe 30/BX

SOL-M LUER SLIP SYRINGE WITHOUT NEEDLE

SOM180011.....SOL-CARE 1cc slip tip syringe 100/bx

SOL-M NEEDLE ONLY, 100/BOX

SOM1181018g x 1"	SOM1231023g x 1"
SOM1181518g x 1-1/2"	SOM1231523g x 1-1/2"
SOM1191019g x 1"	SOM1250625g x 5/8"
SOM1191519g x 1-1/2"	SOM1251025g x 1"
SOM1201020g x 1"	SOM1251525g x 1-1/2"
SOM1201520g x 1-1/2"	SOM1261226g x 1/2"
SOM1211021g x 1"	SOM1263826g x 3/8"
SOM1211521g x 1-1/2"	SOM1265826g x 5/8"
SOM1221022g x 1"	SOM1270527g x 1/2"
SOM1221522g x 1-1/2"	SOM1300530g x 1/2"

SOL-M BLUNT FILL NEEDLE

SOM110022.....18g x 1-1/2

TERUMO[®] SYRINGE ONLY

TERSS-01T	1cc TB Syringe only, Slip tip	Box/100
TERSS-03L.....	3cc TB Syringe only, Luer lock tip	Box/100
TERSS-05L.....	5cc TB Syringe only, Luer lock tip	Box/100
TERSS-10L.....	10cc TB Syringe only, Luer lock tip	Box/100
TERSS-20L2	20cc TB Syringe only, Luer lock tip	Box/50
TERSS-30L.....	30cc TB Syringe only, Luer lock tip	Box/25
TERSS-60L.....	60cc TB Syringe only, Luer lock tip	Box/25



TERUMO[®] NEEDLE ONLY

TERNN-1825R.....	18g x 1"
TERNN-1838R.....	18g x 1-1/2"
TERNN-2025R	20g x 1"
TERNN-2038R.....	20g x 1-1/2"
TERNN-2125R.....	21g x 1"
TERNN-2138R.....	21g x 1-1/2"
TERNN-2225R.....	22g x 1"
TERNN-2238R.....	22g x 1-1/2"
TERNN-2325R.....	23g x 1"
TERNN-2516R.....	25g x 5/8"
TERNN-2525R.....	25g x 1"
TERNN-2713R	27g x 1/2"



TERUMO[®] SYRINGE WITH NEEDLE

TERSS-01T2516.....	1cc TB Syringe w/25g x 5/8" needle	Box/100
TERSS-01T2713.....	1cc TB Syringe w/27g x 1/2" needle	Box/100
TERSS-03L2025	3cc Syringe w/20g x 1" needle	Box/100
TERSS-03L2038	3cc Syringe w/20g x 1-1/2" needle	Box/100
TERSS-03L2125	3cc Syringe w/21g x 1" needle	Box/100
TERSS-03L2138	3cc Syringe w/21g x 1-1/2" needle	Box/100
TERSS-03L2238	3cc Syringe w/22g x 1-1/2" needle	Box/100
TERSS-03L2325	3cc Syringe w/23g x 1" needle	Box/100
TERSS-03L2525	3cc Syringe w/25g x 1" needle	Box/100



BD BLUNT FILL NEEDLE

- BD 305180 Blunt fill needle only, 18g x 1-1/2" Box/100
 BD 305211 Blunt fill filter needle only, 18g x 1-1/2" Box/100
 BD 305060 3cc syringe w/Blunt full needle, 18g x 1-1/2" Box/100

For Preparation and Filling


BD™ Blunt Fill Needle – designed for vial access

BD™ Blunt Filter Needle – designed for safe ampoule medication preparation

BD™ Blunt Fill Needle	
BD Catalog #	Description
305180	18 G x 1 1/2"
305181	18 G x 1"
305183	20 G x 1"

BD™ Blunt Fill Needle, with BD Luer-Lok™ Tip	
BD Catalog #	Description
305060	18 G x 1 1/2" 3 mL Syringe
305062	18 G x 1 1/2" 5 mL Syringe
305064	18 G x 1 1/2" 10 mL Syringe

BD™ Blunt Filter Needle	
BD Catalog #	Description
305211	18 G x 1 1/2"




BD VACUTAINER® PRODUCTS

BD provides a full offering of BD Vacutainer® Blood Collection Tubes in both glass and plastic. BD Vacutainer® Plus Plastic Tubes offer a safe method of blood collection. They also reduce the potential for tube breakage and specimen spoilage, thereby reducing the potential for exposure to bloodborne pathogens.

The BD Vacutainer® Push Button Collection Set is the only blood collection system with push button technology. The BD Vacutainer® Push Button Collection Set offers split second protection for that single moment which could potentially change your life. The push button safety mechanism instantly helps protect you against needle-stick injury.



- BD 367290 Vacutainer Sterile multi sample luer adapter, latex free Box/100
 BD 367899 Vacutainer Plus collection plastic tube, 13mm x 100mm, 6ml, Pink closure cap, sterile Box/100
 BD 367925 Vacutainer Plus collection plastic tube, 13mm x 100mm, 6ml, Grey closure cap, sterile Box/100
 BD 367815 Vacutainer Plus collection plastic tube, 13mm x 100mm, 6ml, Red closure cap, sterile Box/100
 BD 367861 Vacutainer Plus collection plastic tube, 13mm x 75mm, 6ml, Lavender closure cap, sterile Box/100
 BD 367986 Vacutainer Plus collection plastic tube, 13mm x 100mm, 5ml, Gold closure cap, sterile Box/100
 BD 363083 Vacutainer Plus Citrate plastic tube, 13mm x 7mm, 2.7ml, Light Blue closure cap, sterile Box/100
 BD 363080 Vacutainer Plus Citrate plastic tube, 13mm x 7mm, 1.8ml, Light Blue closure cap, sterile Box/100
 BD 367352 Vacutainer Push button collection set with pre-attached holder, 21g x 0.75", 12" tubing, sterile Cs/100
 BD 368656 Vacutainer Push button collection set with pre-attached holder, 23g x 0.75", 12" tubing, sterile Cs/100
 BD 368607 Vacutainer Eclipse safety plastic blood collection needle only, 21g x 1.25", sterile Box/48
 BD 368608 Vacutainer Eclipse safety plastic blood collection needle only, 22g x 1.25", sterile Box/48
 BD 368650 Vacutainer Eclipse safety blood collection needle with pre-attached tube holder, 21g x 1.25", sterile Cs/100
 BD 368651 Vacutainer Eclipse safety blood collection needle with pre-attached tube holder, 22g x 1.25", sterile Cs/100
 BD 368652 Vacutainer Safety-Lok blood collection needle with pre-attached tube holder, 21g x 1.25", sterile Box/25
 BD 368653 Vacutainer Safety-Lok blood collection needle with pre-attached tube holder, 23g x 1.25", sterile Box/25

BD POSIFLUSH™ PREFILLED SYRINGES

BD PosiFlush™ is a range of ready to use sterile pre filled flush syringes containing sodium chloride (NaCl) 0.9%, specifically designed to enhance best clinical practice.

BD 306573 Sterile Pathway (SP) 3ml Box/30
 BD 306574 Sterile Pathway (SP) 5ml Box/30
 BD 306575 Sterile Pathway (SP) 10ml Box/30



BD Q-SYTE™ VIAL ACCESS ADAPTER

The BD Q-Syte™ Vial Access Adaptor Split Septum is a true split septum device with a design that can make a significant impact reducing CRBSIs in your organization. A split septum needleless access system has 64%-70% lower CR-BSI rates than mechanical valves. Purposefully simple in design and function, split septum devices eliminate the complexities of mechanical valves, and with them, the places that may harbor bacteria.



BD 385108 BD Q-Syte™ Vial Access Adaptor Box/25

BD INSULIN SYRINGES

- Ease of penetrating the skin
- Ease of dispensing the total dosage
- Ease of measuring to the correct dose line
- Ease of reading the scale on the syringe
- Consistent quality from syringe to syringe
- Comfort of injection into the skin
- Smooth action of the plunger when drawing up insulin



BD 329420Micro-fine IV 1ml syringe with 28g x 1/2" regular bevel sterile needle Box/100
 BD 329424Micro-fine 1ml syringe with 28g x 1/2" regular bevel sterile needle Box/100
 BD 329461Micro-fine 0.5ml Lo-Dose syringe with 28g x 1/2" sterile needle Box/100
 BD 320440Micro-fine 0.3ml syringe with 31g x 8mm thin wall sterile needle Box/100

SURE COMFORT INSULIN SYRINGE

BB 22-65033/10CC Syringe with 31g x 5/16 Box/100
 BB 22-60051/2cc syringe with 30g x 5/16 Box/100



BD SHARPS CONTAINERS

- BD 300460 1.4 Litre Sharps Container
- BD 305970 Non Locking Tray Style Bracket
- BD 300450 3 Litre Sharps Container
- BD 300974 5.1 Litre Sharps Container
- BD 305447 Locking Wall Bracket for 5.1L
- BD 300440 7.6 Litre Sharps Container
- BD 305485 Locking Wire Bracket
- BD 300452 10.3 Litre Sharps Container



COVIDIEN PHLEBOTOMY SHARPS CONTAINER

Portable phlebotomy containers that fit into present blood drawing trays for the efficient disposal of blood, hypodermic and IV needles. These containers facilitate temporary and final closures in a variety of sizes and colors.

- COV8906..... 1 Quart, Yellow



COVIDIEN MULTI-PURPOSE SHARPS CONTAINER

Covidien offers a complete line of brackets and holders to ensure your safety and compliance with all OSHA regulations.

- COV8970 7.6L with clear rotor opening, red
- COV8975 Wall bracket for 7.6L container
- COV8963..... Locking wall bracket for 7.6L container



COVIDIEN SHARPS CONTAINER, HORIZONTAL DOOR

Covidien offers a complete line of brackets and holders to ensure your safety and compliance with all OSHA regulations.

- COV8967..... 7.6L Container w/clear horizontal drop lid, Yellow



COVIDIEN SHARPS CONTAINER, SHARPSTAR™ LID

Covidien offers a complete line of brackets and holders to ensure your safety and compliance with all OSHA regulations.

- COV8506Y 5L, Yellow



TABLE PAPER

- BAN170912IMC Smooth, 18" x 225' Cs/12
- BAN170914IMC Smooth, 21" x 225' Cs/12
- BAN171002IMC Crepe, 18" x 125' Cs/12
- BAN171004IMC Crepe, 21" x 125' Cs/12



EXAM DRAPES, DISPOSABLE

- BAN1710824IMC..... 40" x 48" Cs/100



EXAM GOWNS, DISPOSABLE

- BAN1710847IMC..... 30" x 42" Cs/50
- BAN910521..... 32" x 42" (Extra Large, Blue) Cs/50



GOWN (HALF WITH TIES), DISPOSABLE

- BAN1710857IMC..... 30" x 21" Cs/100



MEDICAL TOWELS, 3-PLY TISSUE, 13" x 18"

- BAN918101..... White Cs/500



PILLOWCASE

BAN179365IMC...21" x 30", Disposable, White Cs/100
GFI385721" x 27", Reusable, White, Vinyl



PILLOWS (FLUID PROOF)

These pillows feature a soft, vinyl ticking that is inherently flame-retardant, antistatic and fluid proof. Product has antimicrobial properties built in to protect the product. The product does not protect users or others against bacteria, viruses or germs.

MAXMDT219684 18" x 24"



EXAM SHORTS, DISPOSABLE

Disposable exam shorts have multi-layer design for comfort and fluid resistancy.

GRM10000Small/Medium, Blue Cs/50
GRM10001Large/X-Large, Blue Cs/50



APRONS, DISPOSABLE

AMG018-400Plastic, White Box/100



UNDERPAD, DISPOSABLE, FLUFF FILLED

DYN134017" x 24" Bag/100



SANITARY PADS

CFSMT-4 Individually Wrapped Box/250



WelchAllyn® KLEENSPEC® CORDLESS VAGINAL ILLUMINATION SYSTEM COMPLETE WITH CHARGER

- Superior Visualization - advanced LED light source provides better visualization of the exam area by supplying brighter, whiter light, right where it is needed.
- Cordless and Rechargeable - no cord to get in the way or break during procedures, continuous on time of 80 minutes, recharge time of six hours, lithium ion battery technology.
- Automatic on/off - no more hassle with turning the unit on and off.

WEL80000.....Illuminator Only
 WEL80010.....Complete Illumination System



WelchAllyn® KLEENSPEC® CORDED LIGHTING SYSTEMS

- Unique integrated light source transmits superior illumination for a more accurate diagnosis.
- Eliminates positioning of an external light source
- Integrated on/off switch close to the task area for operational convenience.
- New ergonomic design allows for easy cleaning
- New lamp design with improved lamp life

WEL78900.....Speculum Lighting System



WelchAllyn® DISPOSABLE VAGINAL SPECULUM

Designed to eliminate the time and expense of reprocessing metal specula and to reduce the risk of cross-contamination. Single-use speculum is manufactured to deliver durability and reliability. Use with the KleanSpec Illumination System or alone with an external light source.

WEL-590XS.....KleanSpec 590 Series, Extra-Small..... Box/24
 WEL-59000.....KleanSpec 590 Series, Small..... Box/24
 WEL-59001.....KleanSpec 590 Series, Medium Box/24
 WEL-59004.....KleanSpec 590 Series, Large Box/18



DISPOSABLE VAGINAL SPECULUM

GAL8309.....Small Cs/100
 GAL8310.....Medium Cs/100
 GAL8311.....Large Cs/100



DISPOSABLE VAGINAL SPECULUM W/LIGHT SOURCE CAPABILITY

DYN4914.....Small Box/25
 DYN4915.....Medium Box/25
 DYN4916.....Large Box/25



PEDERSON STAINLESS STEEL VAGINAL SPECULUM

- ALMM14-0180Small
- ALMM14-0200Medium
- ALMM14-0220Large



GRAVES STAINLESS STEEL VAGINAL SPECULUM

- ALMM14-0100Virgin
- ALMM14-0120Small
- ALMM14-0140Medium
- ALMM14-0160.....Large



ANOSCOPES, DISPOSABLE

- WEL53110.....Welch Allyn Box/25
- GAL82420.....Sani-Scope Box/25



SIGMOIDOSCOPES, DISPOSABLE

- WEL53130.....KleenSpec Disposable Sigmoidoscope Box/25
individually wrapped



Miltex UTERINE SOUND, DISPOSABLE

- MIL30-6000.....Uterine Sounds, Sterile Pkg/25



Miltex UTERINE SOUND, STAINLESS STEEL

- MIL30-650SS.....13" Stainless Steel



IUD EXTRACTOR/HOOK

MIL30-940.....10.25"



IUD THREAD RETRIEVER DISPOSABLE

SMBRETRIEVER.....Sterile, Single Use Box/5



CERVICAL SCRAPER DUE END

SSP51-200.....7" Ayres Cervical Scraper, Non-Sterile Box/200



HEGAR UTERINE DILATOR

ALMM14-0610.....1 - 2mm, Double Ended

ALMM14-0620.....3 - 4mm, Double Ended

ALMM14-0630.....5 - 6mm, Double Ended



IUD INSERTION KIT, STERILE

Kit Includes:

- 1ea Scissors Metal, 7" Curved
- 1ea Tenaculum Metal, 9.5" Straight
- 1ea Uterine Sound Metal
- 10ea 4x4" Gauze
- 30x30" Sterile Wrap

TRASMA30IUOPB



CERVICAL DILATOR, SINGLE USE

Cervical Dilators are especially useful in locating and dilating the external OS, the cervical canal and the internal OS.

- Length: 215mm
- Taper: 23mm (1.75mm Tip to 3.8mm OD)
- Disposable
- Used on post-menopausal patients with cervical stenosis

BTQOS265BA-1



LORIS™ BZK TOWELETTES

The LORIS™ 0.13% BZK Wipe, is saturated with 0.13% Benzalkonium Chloride and is used as an antiseptic topical skin cleanser.

- Will not sting
- Ideal for perineal and maternity care

LER126-01Individually Packaged Box/100



SILVER NITRATE APPLICATORS

AMG118-3956", 75% Silver and 25% Nitrate Pack/100



RAYON TOP PROCTOLOGY SWABS

DYN43218" Proctology Swab, Rayon Tipped Box/100

HPC81616" Proctology Swab, Rayon Tipped Box/50



COTTON TIPPED APPLICATORS

AMD8033" Non-Sterile Pkg/100

AMD562006" Non-Sterile Pkg/100

AMD4225-IMC6" Sterile (2/Pkg) Box/100 Pkg




TONGUE DEPRESSORS

AMD50000Junior Non-Sterile Box/500

AMD4205-IMCSenior Non-Sterile Box/500

AMG118-421Senior Sterile Box/100



AQUAGEL® LUBRICATING JELLY

Convenient flip cap lid. Non-runny, viscous gel. Bacteriostatic. Clear, water-soluble, non-irritating lubricant. Compatible with natural and synthetic substances.

PAR57-055oz



LUBRICATING JELLY

Water soluble and latex free lube is intended for medical purposes to lubricate body orifices and facilitate entry of diagnostic or therapeutic devices. Comes in 5 gram individual packets.

DYN12515 grams, Foil Packs Box/72



SUNDRY JARS

- DUK4010Glass Sundry Jars, Labelled, Set of 5 Set
- DUK4011Plastic Sundry Jars, Labelled, Set of 5 Set
- DUK4020Wall Mount 5 Jar Rack
- DUK4023Wall Mount 3 Jar Rack



RECTANGULAR WASH BASIN, 7L, NOT AUTOCLAVABLE

DYN4227Plastic, Mauve



RECTANGULAR WASH BASIN, 7L, AUTOCLAVABLE

AMG025-184Plastic, Blue



SPONGE BASIN, STAINLESS STEEL

AMG020-5657.5" x 3"



MEDICINE CUP, GRADUATED, STAINLESS STEEL

AMG90-62002oz



IODINE CUP, STAINLESS STEEL

ALM90-60002oz (60cc)
 ALM90-60053oz (90cc)
 ALM90-60104oz (120cc)
 ALM90-60206oz



EMESIS BASIN, SINGLE USE

Made in plastic, this emesis basin has a 16 oz. (0.47 L) capacity. Non-autoclavable.

AMG760-680plastic, Beige



EMESIS BASIN, STAINLESS STEEL

AMG020-5068", 13oz
 AMG020-50510", 26oz



EYE WASH STATION AND SOLUTION

BIO MED WASH® is an all-natural, sterile first aid wash for cuts and lacerations, is the first product that replaces saline solution for washing abrasions and lacerations (saline is still preferred for deep intramuscular cuts). The product is invaluable in the initial washing of skin that has come in contact with irritating chemicals such as gasoline, pepper spray, tear gas or plant toxins; and the washing of eyes irritated by foreign particles, heavy smoke, dust or perspiration. BIO MED WASH contains 100% Bio-Logic Aqua® tissue-culture grade water, proven to have the ideal pH balance, hypotonicity and mild mineral nutrient content to maximize moisture absorption and promote healing.

RBW47085Eyewash Station, 2 x 7oz. Cannisters, Mirror
 RBW400077oz. Cannisters



EYE WASH SOLUTION

Effective against redness, smarting and mild congestion of the eyes and eyelids due to light wind, dust, smoke or other pollutants.

SFA04077Isotonic Solution Sterile, 32oz / 1L



HYDROGEN PEROXIDE

TRASHP500500mL, 3%



ISOPROPYL ALCOHOL

LER111-01500mL, 70%
 TRASCA599500mL, 99%



ALCOHOL DISPENSER

DUK40249oz. Dispenser



COTTON BALLS

DYN3170Medium Bag/2000



PLASTIC MEDICINE CUPS, GRADUATED

DYN42581oz Pkg/100



PAPER MEDICINE / PORTION CUPS

DYN42440.5oz Pkg/250



PAPER DRINKING CUPS

SPMC000044oz. Poly Coated Sleeve/100
 SPMC000055oz. Poly Coated Sleeve/100



DRINKING CUPS, PLASTIC

DYN4236.....5oz., White, Plastic, Disposable Pkg/50
 DYN42375oz., Blue, Plastic, Disposable Pkg/50
 DYN42385oz., Mint Green, Plastic, Disposable Pkg/50
 DYN42395oz., Mauve, Plastic, Disposable Pkg/50
 DYN42405oz., Lavender, Plastic, Disposable Pkg/50



FLEXIBLE DRINKING STRAW

DYN4260.....Individually Wrapped Box/400



PROTECTIVE EYE GLASSES



FACE SHIELD, FULL, DISPOSABLE

Pre-assembled, one-piece, face shield that includes an optical-grade lens, foam headband, and adjustable elastic strap to ensure a comfortable fit. Latex-free.



EARLOOP FACE MASK

Earloop masks in Level 1, 2, and 3



EARLOOP FACE MASK, WITH FACE SHIELD

Available with optically clear plastic wraparound face shield which provides protection for the face and eyes against pathogens contained in body fluids. Mask has fluid resistant shield. Latex-free.



PRIMED CHILD SIZE FACE MASK

Child Size Face Mask with earloops is soft, non-irritating and easy to breathe through. The dye free mask is not made with natural rubber latex and is comfortable to wear for pediatric size patients.



3M QUALITATIVE FIT TEST KIT

3M Respirator Fit Test Kits can be used with any particulate respirator to test face mask to face seal fit. The 3M Fit Test Kit can help satisfy regulatory requirements for face mask fit testing. Available in either a bitter or sweet test solution kit. Simply spray the solution into the hood environment. If the wearer does not detect the test agent's taste, an acceptable respirator fit has been demonstrated. The 3M Respirator Fit Test Kit can fit test approximately 150 employees with the included supplies.

Kit Includes:

- 55mL Sensitivity Solution
- 55mL Fit Test Solution
- Two nebulizers (FT-13)
- Laminated instruction book (English/Spanish/French/Portuguese)
- One collar assembly (FT-15)
- Two sets replacement nebulizer inserts
- One hood (FT-14)

3M FT-10..... Test Kit with Sweet Sensitivity Solution

3M FT-30..... Test Kit with Bitter Sensitivity Solution



SENSITIVITY SOLUTION, SWEET, 55ML



3M FT-11.....Box/6

TEST SOLUTION, SWEET, 55ML



3M FT-12.....Box/6

SENSITIVITY SOLUTION, BITTER, 55ML



3M FT-31.....Box/6

TEST SOLUTION, BITTER, 55ML



3M FT-32.....Box/6

3M N95 FACE MASK, FLAT FOLD



3M 9210PLUS Box/20

3M N95 FACE MASK, MOLDED



Regular Size

3M 8210..... Box/20

Small Size

3M 8110S Box/20

3M N95 FACE MASK, FLAT FOLD, FLUID RESISTANT



3M 1870PLUS Box/20

3M N95 FACE MASK, MOLDED FLUID RESISTANT



Regular Size

3M 1860..... Box/20

Small Size

3M 1960S..... Box/20

ISOLATION GOWN 50/CS

- Universal
- Offers economical, comfortable and reliable protection
- Full length for maximum coverage
- Elastic cuffs
- Sewn seams provide superior strength
- Latex free
- Extra length on waist ties allows gown to be easily secured



DYN2145White, Impervious
 DYN2141Yellow, Fluid Resistant

ISOLATION GOWN

- Level 3 single, lightweight layer of protective polyethylene film with an open backdesign, Thumbbloops, Fluid Resistance.
- Not chemotherapy approved.

PRIPM4-88335protective gown 20 per bag



ISOLATION GOWN - CHEMOTHERAPY APPROVED

- High quality PE-coated SMS material provides maximum protection against fluids and penetration of harsh chemotherapy drugs. ASTM F739-12, AAMI 2

PRIPM4-99450chemotherapy approved 10 per bag



ISOLATION GOWN - REUSABLE

- Designed with a densely woven fabric made from a non-linting continuous filament polyester to provide excellent fluid resistance and protection.
- 99% polyester / 1% carbon, 3" stretch knit cuffs, Gown Length 46".
- Machine washable and dryer safe, Level 3

EDEISO40MFisolation gown 12 per box



SHOE COVERS

- Non-Skid, Non-Conductive.
- Will cover shoe size from S to L shoes (SPM6007) and from M to XL (SPM6006) (gender non-specific).
- Manufactured of spunbond polypropylene
- Sewn seams, not glued - offers greater strength. Blue in color.
- Elastic opening provides a snug comfortable fit.

SPM6007 Universal Bag/100
 SPM6006 XLarge Bag/100



SHOE COVERS PLASTIC/WATERPROOF

Resistant to water, snow and slush, AMG plastic shoe covers are convenient, cost effective and safe. Made from embossed plastic to prevent slipping, these shoe covers are waterproof and feature an elastic band opening for universal fitting. Blue.

AMG018-285 Universal Size Cs/250 Pairs



BOUFFANT CAP

- Lightweight and breathable spunbond

PRIPM4-6644.....Blue Small 100 per bag



SURGEONS CAPS, DISPOSABLE

PRIPM4-800 Surgeon Cap, Premium SMS...Box/100

DYN2121 Surgeon CapBox/100



COVER YOUR COUGH COMPLIANCE KIT

Cover Your Cough (CYC) Compliance Kits hold an assortment of personal protection equipment in a convenient and attractive dispenser. Kit 101 holds two boxes of earloop face masks, one - two boxes of facial tissue, and one hand sanitizer bottle or hand sanitizing wipe canister (consumables not included). The sign holder is vertical and a variety of signs can be downloaded and printed from www.bowmandispensers.com. Proper and consistent use of the products in this dispenser reduces the chance of spreading infectious bacteria and viruses, such as those found in influenza and the common cold. CYC Compliance Kits are ideal for use in public waiting areas, offices, and lobbies. 18" W x 59" H

Kit 101 includes the following products:

- BOWMAN® Respiratory Hygiene Station (RS001 0512)
- BOWMAN® Floor Stand (KS010-0412)
- BOWMAN® 8.5 x 11" Sign Holder - Vertical (MP-070)

BOWRS001-0212..... Hygiene Station Only

BOWKS010-0412 Floor Stand Only, Steel

BOWMP-070..... 8.5" x 11" Sign Holder, Vertical, Clear Plastic



INFECTION PROTECTION SYSTEM

Holds flat pack or launderable gowns, three boxes of gloves, one box of earloop face masks and one box of face shields. Keyholes for wall mounting or may be used with a wall or door hanger. Quartz ABS plastic.

BOWLD-070



PPE CART II PROTECTION SYSTEM MOBILE

Easy roll 3" casters with toe-locking brakes, steel doors with magnetic closures, safety yellow powder-coated aluminum and steel

BOWCT030-0000Mobile PPE Cart



PEAK FLOW METER

International, independent testing has proven that both the Mini-Wright and the Mini-Wright ATS are the most accurate peak flow meters available. All meters are individually calibrated by hand to ensure consistent accuracy and reproducibility.

- MCA3103085 Adult
- MCA3104050 Pediatric



DISPOSABLE MOUTHPIECES

Our disposable one-way-valved and single patient use mouthpieces are designed to provide protection against possible cross infection in common respiratory procedures.

- MCA3122200 Adult Bag/100
- MCA3123200 Pediatric Bag/100



OXYGEN MASKS & NASAL CANNULA

Cannula

- HUD1822 Adult Cannula with Tubing
- HUD1101 Child Cannula with Tubing

Masks-Medium Concentration

- HUD1041 Adult Mask with Tubing
- HUD1042 Child Mask with Tubing

Oxygen Tubing

- HUD1115 7 foot
- HUD1118 14 foot
- HUD1119 25 foot



NEBULIZER MASKS

- HUD1707 Child Mask & Cup
- HUD1705 Adult Mask & Cup
- HUD1732 Mouthpiece & Cup



CONNECTING TUBE, STERILE

CHS70-6072S..... 3/16" x 72"



DAR® BARRIERBAC "S" VIRAL BREATHING FILTER

COV350/5879 Box/25



MANUAL RESUSCITATOR

Controlled ventilation. Disposable.

ZEN4050 Adult, 66 lb +
 ZEN4010 Pediatric, 22 - 66 lb
 ZEN4025 Infant, 0-22 lb



GUEDEL AIRWAY KIT

Single use, colour coded airways. Eight sizes (40mm, 50mm, 60mm, 70mm, 80mm, 90mm, 100mm, 110mm). Comes in clear poly case with airway tray. Soft, disposable.

ZEN2200-1 Guedel Airway Kit Set/8



GUEDEL AIRWAY NON STERILE

Individual easy open packaging.

ZEN2104 40mm - Pink
 ZEN2105 50mm - Blue
 ZEN2106 60mm - Black
 ZEN2107 70mm - White
 ZEN2108 80mm - Green
 ZEN2109 90mm - Yellow
 ZEN2110 100mm - Red
 ZEN2111 110mm - Orange



CPR FACESHIELD WITH KEYCHAIN CASE

CPR Protective Barrier with Belt/Keyring 11 x 11cm Pouch (Red), One Pair of Latex Free Gloves.

ZEN3052R



CPR POCKET RESUSCITATOR

CPR Pocket Ventilator, Head Strap, Oxygen Inlet, Filter, Gloves and Complete with hard plastic clamshell carry case.

ZEN3100-10



EMERGENCY RESPONDER KIT

- Includes Oxygen
- 1 Guedal Airway Kit
- 1 Universal Face Mask, BVM Reservoir, Oxygen Tubing
- 1 Cylinder "D" with 6LPM Regulator
- 1 Adult Mask 3 Valves with Tubing Non-Breathing
- 1 Pediatric Mask 3 Valves with Tubing Non-Breathing

ZEN911..... Complete Kit



PORTEX® SOFT SEAL LMA MASK, DISPOSABLE

- BCI100/222/200..... Size 2.0 Infant 10kg-20kg
- BCI100/222/250..... Size 2.5 Infant 20kg-30kg
- BCI100/222/300..... Size 3.0 Infant 30kg-60kg
- BCI100/222/400..... Size 4.0 Infant 50kg-70kg
- BCI100/222/500..... Size 5.0 Infant 70kg-100kg



ASPIRATOR (MANUAL VACUUM/SUCTION PUMP)

State of the art, high vacuum, light, portable emergency suction device. Built in 3M HEPA Filter for bacterial and viral protection. Adult and Child settings, Built-in overflow protection seal, 100 degree swivel. Disposable collection jar and catheter.



ZEN7000 Suction Machine
 ZEN7000-AD Disposable Jar and Catheters

OPTIVAC® AC/DC PORTABLE SUCTION UNIT

Long life rechargeable battery. Compact design. Heavy Duty. 3 year warranty, 1 year on battery.



AHPG180 Suction Machine

EPINEPHRINE INJECTION, SINGLE DOSE AMPOULE

WIN07241L0..... 1mL Ampule

*DIN00721891



EPIPEN® EPINEPHRINE AUTO INJECTOR

WIN00121 EpiPen®
 WIN00111 EpiPen® Jr.



NORMAL SALINE, 0.9% SODIUM CHLORIDE

For respiratory therapy.

CHS116..... 15mL Box/24



PLATIC STRIP ADHESIVE BANDAGES

QMDASP19/76 3/4" x 3" Box/100



FABRIC STRIP ADHESIVE BANDAGES

AMDAF0334 3/4" x 3" Box/100



SPOT ADHESIVE BANDAGES

DUK4930-IMC 7/8" Box/100



FABRIC KNUCKLE BANDAGES

QMDASFK3876 1.5" x 3" Box/100



BUTTERFLY WOUND CLOSURES

DYN3615 Sterile, 3/8" x 1 - 3/16" Box/100



KIDS BANDAGES

DIS912816 Spiderman box/20

DIS912824 Frozen box/20

DUK15606 Emoji box/100



LEUKOPLAST STRONG FABRIC BANDAGES

BSN7646007 2.2cm x 7.2cm, Fabric, Strip Box/100

BSN7646005 3.8cm x 2.2cm, Fabric, Island Box/100

BSN7646009 3.8cm x 3.8cm, Fabric, Island Box/100

BSN7646013 5.0cm x 7.2cm, Fabric, Island Box/100

BSN7646012 5.0cm x 4.5cm, Fabric, SM Finger Tip Box/50

BSN7646011 5cm x 6.3cm, Fabric, Finger Tip Box/100

BSN7646010 7.5cm x 3.8cm, Fabric, Anchor Box/100



COVERPLAST® CLASSIC FABRIC DRESSING

- BSN7115104..... 4cm x 5m Box
- BSN7115103..... 6cm x 5m Box
- BSN7115102..... 8cm x 5m Box



PAPER SURGICAL ADHESIVE TAPE

- QMDSUP1X10 1" x 10 Yards Box/12



TRANSPARENT SURGICAL ADHESIVE TAPE

- QMDSUR1X10 1" x 10 Yards Box/12



CLOTH SURGICAL ADHESIVE TAPE

- QMDST1X10 1" x 10 Yards Box/12
- QMDST2X10 2" x 10 Yards Box/12



HYPAFIX® TAPE

A cross-elastic non-woven tape that provides light compression. Excellent for areas requiring flexibility and freedom from constraint while providing superior holding power. Air and exudate permeable. Hypoallergenic.

- BSN7144300..... 2.5cm x 10m
- BSN7144301..... 5.0cm x 10m
- BSN7144302..... 10 cm x 10m
- BSN7144303..... 15 cm x 10m
- BSN7144304..... 20 cm x 10m



3M MICROPORE™ PAPER TAPE

A latex free, hypoallergenic paper tape that is gentle to the skin yet adheres well and leaves minimal adhesive residue upon removal. An economical, general purpose, breathable surgical tape.

- 3M 1530-0..... 1/2" x 10 Yards Box/24
- 3M 1530-1..... 1" x 10 Yards Box/12
- 3M 1530-2..... 2" x 10 Yards Box/6



3M TRANSPORE™ TRANSPARENT TAPE

A latex free, hypoallergenic, transparent and perforated plastic tape that offers strong adhesion. Transpore tape is easy to tear into very thin strips, and is easy to use even with gloves.

- 3M 1527-0..... 1/2" x 10 Yards Box/24
- 3M 1527-1..... 1" x 10 Yards Box/12
- 3M 1527-2..... 2" x 10 Yards Box/6



3M BLENDERM™ TRANSPARENT TAPE

A latex free, hypoallergenic, occlusive and transparent plastic tape. This flexible medical tape is designed to provide a barrier from external fluids and contamination.

- 3M 1525-1..... 1" x 5 Yards Box/12
- 3M 1525-2..... 2" x 5 Yards Box/6



TENOPLAST® ELASTIC ADHESIVE TAPE

- BSN7206902..... 2.5cm x 4.5m
- BSN7205218..... 5.0cm x 4.5m Box/6 RL
- BSN7205219..... 7.5cm x 4.5m Box/4 RL



ADHESIVE REMOVER

- Effectively removes tape and residue from skin
- Reduces patient discomfort
- Contains C11-12 Isoparaffin, Isopropyl Alcohol,
- Petroleum Naphtha, Isopropyl Palmitate, Isopropyl
- Myristate, Isopropyl Stearate

- PDIB16400Remover Pad Box/100



WOVEN GAUZE SPONGES, NON-STERILE

8-Ply Non-Sterile

- AMDD3001 2" x 2" Pkg/200
- DYN3232 3" x 3" Pkg/200
- AMDD3005 4" x 4" Pkg/200

12-Ply Non-Sterile

- DYN3223 2" x 2" Pkg/200
- DYN3233 3" x 3" Pkg/200
- DYN3243 4" x 4" Pkg/200



WOVEN GAUZE SPONGES, STERILE 2'S

8-Ply Sterile 2's

- DYN3322 2" x 2" Box/50 Pkg of 2
- DYN3342 4" x 4" Box/25 Pkg of 2



NON-WOVEN GAUZE SPONGES 4-PLY, NON-STERILE

- AMDA2101-CH 2" x 2" Pkg/200
- AMD2100-IMC 4" x 4" Pkg/200



NON-WOVEN GAUZE SPONGES 4-PLY, STERILE 2'S

- AMDA2240 2" x 2" Box/50 Pkg of 2
- DYN3363 3" x 3" Box/50 Pkg of 2
- AMD4440-IMC 4" x 4" Box/50 Pkg of 2



NON-ADHERENT PADS

- AMDA1923 2" x 3" Box/100
- COVZ6017 2" x 3" w/ Adhesive Box/100
- AMDA1934 3" x 4" Box/100
- DYN3438 3" x 8" Box/50



ABDOMINAL PADS

DYN3501 Sterile, 5" x 9" Tray/20
 DYN3503 Sterile, 8" x 10" Tray/24



EYE PADS, STERILE

DYN3360..... 1-5/8" x 2-5/8" Box/50



PACKING STRIPS, STERILE, PLAIN & WITH IODOFORM

PLAIN, 1 per Bottle

DYN3401-IMC..... 1/4" x 5 Yards
 DYN3402-IMC..... 1/2" x 5 Yards
 DYN3403-IMC..... 1" x 5 Yards

WITH IODOFORM, 1 per Bottle

DYN3411-IMC..... 1/4" x 5 Yards
 DYN3412-IMC..... 1/2" x 5 Yards
 DYN3413-IMC..... 1" x 5 Yards



DRESSING TRAY, STERILE

- 1 Three Compartment Tray
- 1 Plastic Transfer Forceps
- 2 Metal Forceps
- 4 - 10 x 10cm Gauze Sponges
- 4 - 5 x 5cm Gauze Sponges
- 1 Waste Bag with Tie
- 1 Safety Pin
- 2 Towel/Underpad



TRA85-2000TYEach
 TRA85-2000..... Cs/30

FIRST AID KITS

SFA50100B..... Alberta No. 1 Kit



CONFORM STRETCH GAUZE, NON-STERILE

- QMDQ134050.....2" x 5 Yards Box/12
- QMDQ134075.....3" x 5 Yards Box/12
- QMDQ134100.....4" x 5 Yards Box/12
- QMDQ134150.....6" x 5 Yards Box/6



CONFORM STRETCH GAUZE, STERILE

- QMDS134050.....2" x 5 Yards Box/12
- QMDS134075.....3" x 5 Yards Box/12
- QMDS134100.....4" x 5 Yards Box/12



ELASTIC TENSOR BANDAGES

- QMDQ133SP0502" x 5 Yards
- QMDQ133V75.....3" x 5 Yards
- QMDQ133SP1004" x 5 Yards
- QMDQ133V1506" x 5 Yards



STOCKINETTE TUBULAR 100% COTTON BANDAGE

- QMDR10-0502" x 25 Yards
- QMDR10-0753" x 25 Yards
- QMDR10-1004" x 25 Yards
- QMDR10-1506" x 25 Yards



TRIANGULAR BANDAGE

- QMDCT4056



SURGILAST® ELASTIC NET BANDAGE, LATEX FREE, 25 Yards

- ZENGL-701 Size 1 Roll
- ZENGL-702 Size 2 Roll
- ZENGL-703 Size 3 Roll
- ZENGL-704 Size 4 Roll
- ZENGL-705 Size 5 Roll
- ZENGL-706..... Size 5.5 Roll
- ZENGL-707 Size 6 Roll
- ZENGL-708 Size 7 Roll
- ZENGL-709 Size 8 Roll



SURGITUBE® TUBULAR GAUZE BANDAGE FOR APPLICATOR, 50 Yards

- ZENGL-219 Size 1, Small Fingers & Toes Roll
- ZENGL-220 Size 2, Large Fingers & Toes Roll
- ZENGL-221 Size 3, Hands, Wrist & Feet Roll
- ZENGL-222 Size 4, Arms & Lower Legs Roll
- ZENGL-223 Size 5, Legs, Thighs & Head Roll



SURGITUBE® APPLICATORS

- ZENGL-230 Size 0
- ZENGL-231 Size 1
- ZENGL-232 Size 2
- ZENGL-234 Size 3
- ZENGL-236..... Size 4



HOT & COLD PACKS, REUSABLE, GEL

- RAP12246-24..... 4" x 6"
- RAP12259-24 5" x 9"



HOT & COLD PACKS, INSTANT, DISPOSABLE

- DYN4512Cold 5"x9 CS/24
- DYN4516.....Hot 5"x9



CIRCULAR COOL GEL PACKS

Reusable.

- RAP11331-48 3" Cs/48
- RAP11351-24 5" Cs/24



POLYSPORIN® OINTMENT

IDC75620..... 15g Tube
 IDC75647..... 30g Tube



BACITRACIN ZINC OINTMENT

Helps prevent infection in minor cuts, scrapes, and burns

DYN1171 0.9g Foil Pack Box/144



CUTICELL® DRESSING

Cuticell® is soothing and low-adherent minimizing pain at dressing change. The primary wound contact layer is made of smooth acetate fabric impregnated with non-medicated ointment on Eucerit® basis. The result: an ointment dressing that ensures that the wound is not disturbed and helps maintain a moist environment that supports wound healing.

BSN7253901..... 3" x 3" Box/10



JELONET® DRESSING

SMN7403..... 5cm x 5cm Box/50



TEGADERM™ TRANSPARENT FILM DRESSING

3M 1624W..... 6cm x 7cm Box/100
 3M 1626W..... 10cm x 12cm Box/50



OPSITE® TRANSPARENT FILM DRESSING

Flexifix Film

SMN66000040 5cm x 10m
 SMN66000041 10cm x 10m

Flexigrid Dressing

SMN4628 6cm x 7cm Box/100
 SMN4630 10cm x 12cm Box/50



OPSITE® SPRAY

SMN66004978 100mL Spray Bottle



SALJET® SINGLE DOSE NORMAL SALINE FOR IRRIGATION

Single use 30mL containers. Saves time, wastage and assures sterility.

WLB10-18815100009-08 30mL Box/40



IRRIGATION SOLUTIONS

Normal Saline

BAXJF7633P 500mL
 BAXJF7634 1000mL

Sterile Water

BAXJF7623P 500mL
 BAXJF7624 1000mL



STERILE WATER FOR IRRIGATION

MAXDYND40570 110mL cup Cs/48



NORMAL SALINE 0.9% (SODIUM CHLORIDE) FOR IRRIGATION

MAXDYND40540 100mL



DYNACAST® PRELUDE SPLINT SYSTEM (ROLL FORMAT)

Fiberglass.

BSN7141800 2.5cm x 4.6m Roll
 BSN7142000 5.0cm x 4.6m Roll
 BSN7142100 7.5cm x 4.6m Roll
 BSN7141500 10cm x 4.6m Roll
 BSN7141600 12.5cm x 4.6m Roll
 BSN7141700 15cm x 4.6m Roll



DELTA-LITE® PLUS CAST TAPE

Fiberglass.

BSN7345800.....	2.5cm x 1.8m, White	Box/10
BSN7345801.....	5.0cm x 1.8m, White	Box/10
BSN7345802.....	7.5cm x 3.6m, White	Box/10
BSN7345803.....	10cm x 3.6m, White	Box/10
BSN7345804.....	12.5cm x 3.6m, White	Box/10

Other colours available



SOFFBAN CAST PADDING - SYNTHETIC

BSN7146401.....	5cm x 2.7m	Bag/12
BSN7146701.....	7.5cm x 2.7m	Bag/12
BSN7148605.....	10cm x 2.7m	Bag/12
BSN7146501.....	15cm x 2.7m	Bag/12



GYPSONA® PLASTER OF PARIS CAST BANDAGE

BSNLPL23	5cm x 2.7m	Box/24
BSNLPL33	7.5cm x 2.7m	Box/12
BSNLPL43	10cm x 2.7m	Box/12
BSNLPL63	15cm x 2.7m	Box/12
BSNLPL65	15cm x 4.5m	Box/12
BSNLPL83	20cm x 2.7m	Box/12
BSNLPL85	20cm x 4.5m	Box/12



GYPSONA® PLASTER OF PARIS CAST SLABS

BSNL5330	7.5cm x 75cm	Box/50
BSNL5430	10cm x 10cm	Box/50
BSNL5530	12.5cm x 75cm	Box/50



ZIP STICK CAST REMOVAL AID

BSN7204629.....	2.5cm X 48cm
-----------------	--------------



ALUMAFOAM FINGER SPLINTS

CCO6212	1/2" x 18"	Box/12
CCO6234	3/4" x 18"	Box/12
CCO6210.....	1" x 18"	Box/12



DENTAL BIBS

- DYN449117-3/4" x 12-7/8", White
- DYN449217-3/4" x 12-7/8", Blue
- DYN449317-3/4" x 12-7/8", Mint Green
- DYN449417-3/4" x 12-7/8", Mauve
- DYN449517-3/4" x 12-7/8", Lavender



DENTAL BARRIER FILM

Eliminates cross contamination while working with handheld equipment. Covers most surface areas.

- DYN4488 4" x 6", Clear Box/1200
- DYN4489 4" x 6", Blue Box/1200



IMPRESSION TRAYS

- Disposable perforated plastic trays
- Consistent impressions
- Small, medium and large upper and lower trays, quadrant and anterior
- 12 pieces per bag, 10 bags per case
- Not made with natural rubber latex

- DYN2220 No. 1, Large Upper..... Bag/12
- DYN2221 No. 2, Large Lower..... Bag/12
- DYN2222 No. 3, Medium Upper..... Bag/12
- DYN2223 No. 4, Medium Lower..... Bag/12
- DYN2224 No. 5, Small Upper Bag/12
- DYN2225 No. 6, Small Lower Bag/12
- DYN2234 No. 7, Quadrant Tray, UL/LR..... Bag/12
- DYN2235 No. 8, Quadrant Tray, UR/LL..... Bag/12
- DYN2238 No. 9, Anterior Tray, Upper Bag/12 (Bag Size 25x17)
- DYN2239 No. 10, Anterior Tray, Lower..... Bag/12 (Bag Size 25x17)



EARLOOP FACE MASK

Procedure Masks offer a flexible nose piece that helps form a closure for better protection against exposure. Pleated earloop procedure Masks with glass free filter. Not made with natural rubber latex.

- Earloop Masks in level 1, 2, and 3 Box/50



BOUFFANT CAP

- Lightweight and breathable spunbond

PRIPM4-6644.....Blue Small.....100 per bag



SHOE COVERS

- Non-Skid, Non-Conductive.
- Will cover shoe size from S to L shoes (SPM6007) and from M to XL (SPM6006) (gender non-specific).
- Manufactured of spunbond polypropylene
- Sewn seams, not glued - offers greater strength. Blue in color.
- Elastic opening provides a snug comfortable fit.



SPM6007..... UniversalBag/100

SPM6006..... XLarge..... Bag/100

SHOE COVERS PLASTIC/WATERPROOF

Resistant to water, snow and slush, AMG plastic shoe covers are convenient, cost effective and safe. Made from embossed plastic to prevent slipping, these shoe covers are waterproof and feature an elastic band opening for universal fitting. Blue.



AMG018-285 Universal SizeCs/250 Pairs

SURGEONS CAPS, DISPOSABLE

PRIPM4-800 Premium SMS..... Box/100

DYN2121 Box/100



ISOLATION GOWN TEAL UNIVERSAL 5/30bag/cs

- PriMed Overhead Fullback
- Thumbloop with Elastic Cuff
- 34gsm Teal SMS
- AAMI Level 3

PRIPM4-22430.....5/BG



NITRIDERM® NITRILE STERILE EXAM GLOVES

Nitrile synthetic contains no natural rubber latex proteins. Soft, highly elastic formulation provides the maximum in comfort. Pair packaged, inner wrapper in ozone resistant multi-vac packaging.

- IHC106100Small
- IHC106200Medium
- IHC106300Large
- IHC106350X-Large

50 Pairs/Box



NOAH STRETCH NITRILE - BOX OF 100

- Finger Textured
- Nitrile
- Powder free
- 3.0mil Thickness

- SPM88896Small
- SPM88897Medium
- SPM88898Large



BLUSH® NITRILE PINK EXAMINATION GLOVES - BOX OF 200

- Nitrile (100% Latex-Free)
- Non Sterile
- Glove Design & Feature: Ambidextrous and Finger Textured
- Cuff: Beaded
- Color: Blush Pink
- Thickness: 2.5 mil
- Packaging: 200 pieces per box / 10 boxes per case

- SPM78885 Extra Small
- SPM78886 Small
- SPM78887 Medium
- SPM78888 Large



PRIMED PRIMATOUGH® FIT NITRILE GLOVES (NON-POWDERED) - BOX OF 200

- PRIPM6-2005Extra Small
- PRIPM6-2006Small
- PRIPM6-2007Medium
- PRIPM6-2008Large
- PRIPM6-2009Extra Large



See more glove options on pages 63-67,
or check out the website for all available gloves.

RUBBER BANDS

- Made of high quality rubber
- Fits tattoo machines that use size #12 & #13
- Helps keep your tattoo needle in place while tattooing
- A must-have for all standard tattoo machines
- One of the most commonly replaced parts on a tattoo machine



DYN2166..... 1l..... b Bag Bag

GROMMET NIPPLES

- Gives the needle a snug fit and keeps it in place
- Easily slip your needle over the top of the nipple, secures needle tightly on the A-bar
- Rubber nipples designed specifically for tattoo machine use
- Packaged in bags of 100



DYN2179 BlackBag/100

INK CUPS

- Used to hold pigment or ink
- Enables you to place your ink cup directly on the work surface without the need for an ink cup stand
- Transparent to allow a better view of the ink color
- Disposable
- Non-Sterile
- Available in Small, Medium, Large & Extra Large.



DYN2175 Medium, 14mm..... Bag/1000

SOAKER PADS

- Absorbent center is wrapped and sealed on all four sides by a polyethylene film with perforations on one side
- Pads allow tattoo to breathe and not dry out
- Non-sterile
- Non-stick for easy removal
- Absorbs excess ink after tattoo application



DYN2168..... 4" x 7", Black..... Box/500

DYN2169..... 6.25" x 7", Black..... Box/500

STERILIZATION POUCHES, SELF SEALING

- Self sealing, no tape or heat sealing needed
- Blue tinted transparent cover makes it easy to detect contamination
- Internal and external color-changing indicators assure proper processing in chemical vapor, steam, or ethylene oxide (EO) gas sterilization.
- Medical grade materials
- Latex Free
- Pouch sizes below are calculated as outer dimensions.



DYN4460.....	2.25" x 2.75"	Box/200
SPMSP71300.....	7.5" x 13"	Box/200
SPMSP35900.....	3.5" x 9"	Box/200
SPMSP52510.....	5.25" x 10"	Box/200
SPMSP22540.....	2.25" x 4"	Box/200
SPMSP27590.....	2.75" x 9"	Box/200

BOTTLE COVERS

- Bottle Bags will cover your entire spray or wash bottle, preventing cross contamination
- Transparent, for easy identifying the liquid level in the bottle
- Excellent for keeping your equipment clean and germ free



DYN2164.....	Medium, 6" x 8"	Box/500
DYN2165.....	Large, 6" x 10"	Box/500

MACHINE BAG

- Protect your machine from cross contamination and splatter
- These covers easily slip over your machine



DYN2162.....	5" x 5"	Box/500
--------------	---------------	---------

DENTAL BIBS

DYN4491	17-3/4" x 12-7/8", White
DYN4492	17-3/4" x 12-7/8", Blue
DYN4493	17-3/4" x 12-7/8", Mint Green
DYN4494	17-3/4" x 12-7/8", Mauve
DYN4495	17-3/4" x 12-7/8", Lavender



DENTAL BARRIER FILM

Eliminates cross contamination while working with handheld equipment. Covers most surface areas.

DYN4488	4" x 6", Clear
DYN4489	4" x 6", Blue



BLACK NITRILE EXAM GLOVES

Extraordinary strength and puncture resistance. Maintaining tactile sensitivity. Contains no allergy causing natural rubber proteins. Commonly used in tattoo, EMS and other first responder markets. Ambidextrous, Non-Sterile, Powder Free. Not made with natural rubber latex. Black.



DISPOSABLE THONG

Made of soft spunbonded polypropylene material. Each garment is individually wrapped for hygienic purposes. One Size.

DUK900502 Black..... Cs/250



DISPOSABLE BIKINI

Made of soft spunbonded polypropylene material. Each garment is individually wrapped for hygienic purposes. One Size.

DUK900506 Black..... Cs/250



DISPOSABLE BRA

Backless. Made of soft spunbonded polypropylene material. Each garment is individually wrapped for hygienic purposes. One Size.

DUK900510 Small/Medium..... Cs/250

DUK900512 Large/XLarge..... Cs/250



DISPOSABLE MENS BRIEF

Made of soft spunbonded polypropylene material. Each garment is individually wrapped for hygienic purposes.

DUK900528 Large/XLarge..... Cs/250



DISPOSABLE MENS BOXER

- Provides comfort with just the right amount of coverage
- Made of soft spun bonded polypropylene material
- Available in a variety of sizes and styles to accommodate multiple treatments
- Each garment is individually wrapped for hygienic purposes
- Disposable
- Not made with natural rubber latex

DUK900530 Small/Medium..... Cs/50


DUK900532 Large/XLarge..... Cs/50



DISPOSABLE HEADBAND

- Convenient Velcro Closure
- Used to keep hair from face during spa treatments.

GRM51823 Pkg/48



CLEAR ULTRASOUND GEL

- Get "clear results" with our dye free gel.
- Wavelength clear ultrasound gel is great for procedures requiring a pigment-free, medium viscosity conductive agent.
- Contains no salt, alcohol, fragrance or dye.



NTPC202 250mL Bottle
 NTPC503 5L Bottle

COLD PACK GEL CIRCLES

- Reusable soft-sided gel pack, ideal for a variety of cold therapy applications, such as following after injections, blepharoplasty, and breast surgery.
- Remains flexible when frozen
- Soft-sided cloth backing for barrier between skin and pack



RAP11331-48 3" Cs/48
 RAP11351-24 5" Cs/24

MAKEUP REMOVER WIPES

UNI151355 Bag/50



INSULIN SYRINGE

- 3/10cc

BD 320440 31G x 8mm (5/16").. Box/100
 BB 22-6503 Comfort 31g x 5/16 ... Box/100



TSK STERiGLIDE™ CANNULA

The biggest advantage of cannulas over needles is the reduced risk of bruising during filler injections. This advantage has resulted in a rapid adoption of the cannula technique worldwide, making cannulas the product of choice for selective procedures, e.g. the tear through region. Other advantages of cannulas include the ability to treat a larger area with only one injection point, better comfort, less trauma, and more acceptability by the patient. After TSK's first cannula introduction with the CSH Cannula, TSK introduced their award winning STERiGLIDE cannula in 2012. The STERiGLIDE outperforms any other cannula available and remains to lead the market as the golden standard.

- Patented Dome shaped tip, Nearest to tip filler delivery
- Proprietary surface treatment
- Clear side-port marking
- Introducer needles included
- ULTRA Thin Wall
- Safety PRC hub

TSK (All Sizes) Box/20



Cannulas / STERiGLIDE™					
Code	Description	Inch	UTW*	Box	
PRC-30025ISG	STERiGLIDE Cannula 30G x 25mm	1"	✓	20	
PRC-27050ISG	STERiGLIDE Cannula 27G x 50mm	2"	✓	20	
PRC-27038ISG	STERiGLIDE Cannula 27G x 38mm	1½"	✓	20	
PRC-27025ISG	STERiGLIDE Cannula 27G x 25mm	1"	✓	20	
PRC-25050ISG	STERiGLIDE Cannula 25G x 50mm	2"	✓	20	
PRC-25038ISG	STERiGLIDE Cannula 25G x 38mm	1½"	✓	20	
PRC-22070ISG	STERiGLIDE Cannula 22G x 70mm	2¾"	✓	20	
PRC-22050ISG	STERiGLIDE Cannula 22G x 50mm	2"	✓	20	

* UTW = ULTRA Thin Wall

WOODEN ROUND APPLICATOR STICKS

SSP4204-40-500..... 6" Box/1000



DISPOSABLE LIP BRUSH

Tightly spaces bristles for easy application

SSP4017-LIP-25SP Pkg/25



MICRO APPLICATOR BRUSH

- Tips are lint-free and non-absorbent
- Ultrafine tips
- Dentistry, Beauty, Electronics and craft applications
- 4 inch handle with narrow bendable neck
- No spill dispenser

SSP4017-MAB-4UF..... Ultrafine, White Box/100



DISPOSABLE MASCARA BRUSH

- Full, lush head coats lashes evenly.
- Nylon bristles

SSP4017-MB-20SP 4" Pkg/20



DISPOSABLE CLEAR SPATULA

- Smooth Edge, flexible blade.
- Used for mixing and applying gel or makeup.

SSP4017-SPATULA-20SP Pkg/20



MANICURE/PEDICURE STICKS

- Square tapered on both ends.
- Gently grooms cuticles and cleans under nail tips.

SSP4202-60-100SP 7" Bag/144



INDEX

A		
Abdominal Pads	120	
Adhesive Remover	118	
Adhesive, Tissue	73	
Administration Sets, I.V.	79	
AED	38	
Air Freshener	45	
Airways	113	
Alcohol Dispenser	106	
Alcohol, Isopropyl	106	
Allergy Trays	92	
Amniotic Test	88	
Aneroids	24	
Anesthetic	75-76	
Anoscopes	101	
Antiseptic Solutions	77,84	
Applicators, Cotton Tipped	103	
Aprons, Disposable	99	
Autoclave	47	
Autoclave Cleaner	50	
Autoclave Tape	50	
Autoclave Pouches & Wraps	49	
B		
Bandages	116	
Bathroom Tissue	44,45	
Batteries	27	
Baxedin	77	
Benzalkonium Wipes	77,103	
Biohazard Bags	89	
Biological Indicator	48	
Biopsy Kit, Endometrial	88	
Biopsy Punches	69	
Bleach	46	
Blades, Scalpels	68	
Blood Collection	81-83	
Blood Drawing Chair	90	
Blood Lancets	90	
Blood Pressure Cuffs	23	
Blood Pressure Devices, Digital	20-22	
Blood Pressure Devices, Manual	24	
Blood Pressure Tubing	23	
Bags, Recycling	45	
Blunt Fill Needle	95	
Bouffant Caps	111,127	
Brush, Surgical Instrument	49	
Bulbs	27	
Butterfly Closures	116	
C		
Cap, Surgeons	111,127	
Carts	16,17	
Cast Padding	125	
Casting	124,125	
Catheters, I.V.	78,79	
Cauteries	13,14	
Centrifuge	90	
Cervical Dilator	102	
Cervical Scraper	102	
Chair, Procedure	5	
Chair, Transport	18	
Chairs	7	
Chlorhexadine Pads & Swabs	84	
Cleaning Solutions, Instrument	40,41	
Cold/Hot Packs	122	
Condoms	89	
Conform Stretch Gauze	121	
Connecting Tube	113	
Cotton Balls	107	
Cotton Tipped Applicators	103	
COVID-19 Antibody	87	
COVID-19 Antigen	86,87	
CPR Faceshield	114	
CPR Pocket Resuscitator	114	
Cryosurgery	70	
Cuffs, Blood Pressure	23	
Cups, Drinking	107	
Cups, Medicine	105,107	
Curettes	61,62	
Cytology Collection Brush	88	
Cytospray	89	
D		
Defibrillators	38	
Dermal Tips	14	
Detergent, Laundry	46	
Dextrose, 5%	78	
Diagnostic Set	24,25	
Dental	126-129	
Dish Soap	46	
Disinfectant Solutions	39,40	
Disinfectant Spray/Wipes	39,40	
Doppler	31	
Drapes, Exam	98	
Drapes, Surgical	76	
Dressing Tray	74,120	
Drinking Cups	107	
Drugs of Abuse Test Kit	88	
E		
Ear Basin	58	
Ear Cleaning	58,61,62	
Ear Curettes	57,61,62	
ECG	33-35	
ECG Electrodes	36	
ECG Prep Tape	37	
Elastic Bandage	121	
Elastic Net Bandage	122	
Electrode Skin Prep Pad	37	
Electrodes	36	
Electrostatic Sprayer	46	
Electrosurgical Generator	13,14	
Emergency Responder Kit	114	
Emesis Basin	105	
Endometrial Biopsy	88	
Epinephrine Injection	115	
EpiPen	115	
Exam Tables	4-7	
Eye Drops	76	
Eye Glasses, Protective	108	
Eye Pads	120	
Eye Wash Station & Solution	106	
F		
Face Mask	108,109	
Face Shield	108	
Facial Tissue	44	
Filter, Breathing	113	
Finger Ring Cutter	57	
Finger Splints	125	
First Aid Kits	120	
Flourescein Strips	76	
Fridge	15	
G		
Garbage Bags	45	
Garbage Can	45	
Gauze	119	
Glasses, Protective	108	
Glove Box Holder	64	
Gloves	63-67	
Glucose Testing	86	
Glustitch Tissue Adhesive	73	
Gowns, Exam	97	
Gowns, Isolation	110	
H		
Hand Sanitizer	42	
Hand Towels	44	
Height Measuring	10,11	
Holter Recorders	35	
Hot/Cold Packs	122	
Hydrogen Peroxide	106	
Hyfrecator	13,14	
I		
I.V. Access Cannula	84	
I.V. Administration Sets	79	
I.V. Catheters	82-83	
I.V. Connector	84	
I.V. Extension Sets	79	
I.V. Pole	84	
I.V. Solutions	78,80	
Indicator Strips, Steam	49	
Indicators, Biological	48	
Integrators	49	
Instrument Stand	17	
Instrument Tip Protectors	51	
Instrument Cleaning Supplies	62	
Instrument Disinfection/Cleaning	40,41	
Instrument Lube	62	
Instrument Trays	51	
Instruments	52-62	
Instruments, Sterile	59,60	
Insulin Syringes	96	
Iodine Cup	105	
Irrigation Solutions	80	
Isolation Gowns	110	
IUD Insertion Kits	74,102	
IUD Hook	102	
IUD	102	
J		
Jars, Sundry	104	
L		
Laceration Trays	74	
Lactated Ringers	78	
Laundry Detergent	46	
Lancets	90	
Lights, Exam	11,12	
Lubricating Jelly	104	
M		
Marcaine	75	
Masks, Oxygen/Nebulizer	112	
Mayo Stand	17	
Medicine Cups	105,107	

INDEX

N

N95 Fit Test Kit	109
N95 Masks	109
Nasal Cannula	112
Nebulizer Masks	112
Needles & Syringes	91-96
Neuropen	33
Non-Adherent Pads	119
Normal Saline, I.V.	78
Normal Saline, Irrigation	124
Normal Saline, Respiratory	115

O

Ointment Dressing	123
Opsite Dressing/Spray	123
Ophthalmoscopes	25
Otoscope Speculum	26
Otoscopes	25,26
Oxygen Masks	112
Ovulation Test	88

P

Packing Strips	120
Paper Towels	44
Patient Monitors	20-22
Peak Flow Meter	112
Pediatric Scales	8
Penlight	60
Pillowcase	99
Pillows	99
Prefilled Syringe	96
Pregnancy Test	87
PRN Luer Lock I.V. Adapters	84
Probe Covers, Thermometer	30
Probe Covers, Ultrasound	32
Proctology Swabs	103
Pulse Oximeters	29
PVP Preps, Swabs & Solutions	77

Q

Q-Style I.V. Connector	84
Q-Style Vial Adapter	96
Qualitative Fit Test Kit	109

R

Rapid Test	86,87
Razors	37,73
Refrigerators	15
Resuscitator	113

S

Sanitary Pads	99
Sanitizer, Hand	42
Scales	8-10
Scalpel Handle	68
Scalpels	68
Scrub Brush, Nylon	77
Sharps Containers	95
Shoe Covers	110,127
Shorts, Exam	99
Sigmoidoscopes	101
Silver Nitrate Applicators	103
Skin Closures	74
Skin Markers	74
Staple Removal	74
Skin Stapler	74

Smoke Shark	14
Soap, Hand	43
Solu IV Swabs	84
Spa & Beauty	131-133
Specimen Containers	89
Speculum, Otoscope	26
Speculum, Vaginal	100-101
Spirometer	33
Sponge Basin	104,105
Sprayer, Disinfection	46
Sterile Instruments	59,60
Sterile Water, Injection	80
Sterile Water, Irrigation	80
Sterilization Pouch Rack	51
Sterilization Pouches & Wraps	49
Sterilizers	47
Stethoscopes	27,28
Stockinette Tubular Bandage	121
Stools	8
Stopcocks	79
Straw, Drinking	107
Strep-A Test Kit	86,87
Stretcher	18
Suction Machines	115
Sundry Jars & Racks	104
Surface Cleaner	39
Surface Disinfectant	39
Surgeons Cap	111,127
Surgical Mask	108,109
Surgilast Bandage	122
Surgitube Bandage	122
Suture Removal Tray	74
Sutures	71,72
Swabs, Proctology	103
Syringe, Prefilled	96
Syringes & Needles	91-96
Syringes, Insulin	96

T

Table Paper	98
Table, Treatment	4-7
Table, Exam	4-7
Tablets, Ultrasonic Cleaner	48
Tape Measure	10
Tape, Surgical	117,118
Tattoo Supplies	130
Tegaderm Dressing	84,123
Tensor Bandages	121
Thermometer, Fridge	16
Thermometers	30
Tissue, Facial	44
Tissue, Toilet	44,45
Tongue Depressors	103
Topical Anaesthetic	76
Tourniquets	84
Towels, Medical	98
Transport Chair	18
Treatment Carts	16-17
Treatment Table	4-7
Triangular Bandage	120
Tuning Forks	58

U

Ultrasonic Cleaners	48
Ultrasound Supplies	32

Underpads	99
Urinalysis	85
Urine Analyzer	85
Uterine Dilator	102
Uterine Sound	101

V

Vaccine Refrigerator	15
Vacutainer Products	95
Vaginal Illumination System	100
Vaginal Speculum	100,101
Veritor Flu Test Kit	86
Vision Screener	26
Vital Signs Monitor	20,22

W

Wall Systems	19
Warmer, Ultrasound Gel	32
Wart Treatment	70
Water Distilled	50
Wheelchair	18
Wipes, Disinfectant, Surface	39

X

Xylocaine	75
-----------	----

Z

Zoll AED	38
----------	----



1.800.232.7286

canadamedical.net | orders@canadamedical.net

TERMS AND CONDITIONS:

Net 30 Days on approved accounts.

Taxes extra where applicable.

Free Shipping on most orders (Call us for more details).

Returns require an authorization #

Returns may be subject to restocking fee.

No returns on special order products.



The reproduction of any part of this catalogue is strictly prohibited unless accompanied by written permission from Canada Medical Ltd.
Errors and Omissions Excepted. Prices are subject to change without notice.